DSR-20MD/20MDP

SERVICE MANUAL





Photo: DSR-20MD

US Model Canadian Model DSR-20MD

AEP Model Australian Model New Zealand Model DSR-20MDP

E MECHANISM

SPECIFICATIONS

System

Recording format DVCAM format

Video signal

DSR-20MD:

EIA STANDARD, NTSC color.

system

DSR-20MDP: CCIR STANDARD, PAL colour

system

Usable cassettes Recording time

Standard-DVCAM cassettes and

Mini-DVCAM cassettes

184 minutes (when using the PDV-184ME cassette)

40 minutes (when using the

PDVM-40ME cassette)

Clock

Quartz locked

DSR-20MD:

12-hour cycle display DSR-20MDP: 24-hour cycle display

Power back-up

Built-in self-charging capacitor

Back-up duration: up to about 100

hours

(After 8-hour charges)

Inputs and outputs

Video input

BNC connector

Input signal: 1 Vp-p

(75 ohms unbalanced) BNC connector

Video output

Output signal: 1 Vp-p

S video input

(75 ohms unbalanced) Mini DIN 4-pin

Luminance signal: 1 Vp-p

(75 ohms unbalanced) Chrominance signal:

0.286 Vp-p (DSR-20MD) 0.3 Vp-p (DSR-20MDP) (75 ohms unbalanced)

S video output

Audio input

Mini DIN 4-pin Luminance signal: 1 Vp-p

(75 ohms unbalanced)

Chrominance signal;

0.286 Vp-p (DSR-20MD) 0.3 Vp-p (DSR-20MDP)

(75 ohms unbalanced)

Phono jack (L, R)

Input level: 2 Vrms (full bit) Input impedance: more than

47 kohms

Phono jack (L, R) Audio output

Output level: 2 Vrms (full bit)

Output impedance: less than

10 kohms

BNC connector Monitor output

Output signal: 1 Vp-p (75 ohms unbalanced)

Control S input Stereo minijack (1)

For the optional DSRM-10 Remote

Control Unit

Foot switch input Stereo minijack (1)

For the optional Foot Switch

LANC input/output

Stereo mini-mini jack (1)

For the optional RM-95 Remote

Commander

- Continued on next page -

DVCAM |

DIGITAL VIDEO CASSETTE RECORDER



SONY

RS-232C input/output

D-sub 9-pin connector (1) Output: 3 kilohms at load Typ ±9V

Input: 5 kilohms at load High level 5 to 15V

Low level -5 to -15V Headphones output

Stereo minijack (1)

DV input/output 4-pin jack (1)

General

Power requirements

120 V AC, 60 Hz DSR-20MD:

12 V DC, 2.0 A (4.0 A at the peak) DSR-20MDP: 220 – 240 V AC, 50 Hz

12 V DC, 2.0 A (4.0 A at the peak)

Power consumption

0.45 A at 77°F, 120 V AC, 60 Hz DSR-20MD:

(during playback)
DSR-20MDP: 0.35 A at 25°C, 220 – 240 V AC,

50 Hz (during playback)

Operating temperature

5°C to 40°C (41°F to 104°F)

Storage and transport temperature

-20°C to +60°C (-4°F to +140°F)

Operating humidity


20% to 80%

Storage and transport humidity

20% to 80%

Dimensions

Approx. 212 × 98 × 395 mm (8 1/4 × 3 1/4 × 15 1/4 inches) (w/h/d, including projecting parts and controls)



Mass

Approx. 5 kg (11 lb.)

Supplied accessories

Remote commander (1) Size AA (R6) batteries (2) AC power cord (1) Cleaning cassette (1) Instructions for Use (1)

Design and specifications are subject to change without notice.

SAFETY-RELATED COMPONENT WARNING!!

COMPONENTS IDENTIFIED BY MARK A OR DOTTED LINE WITH MARK A ON THE SCHEMATIC DIAGRAMS AND IN THE PARTS LIST ARE CRITICAL TO SAFE OPERATION. REPLACE THESE COMPONENTS WITH SONY PARTS WHOSE PART NUMBERS APPEAR AS SHOWN IN THIS MANUAL OR IN SUPPLEMENTS PUB-LISHED BY SONY.

ATTENTION AU COMPOSANT AYANT RAPPORT À LA SÉCURITÉ!


LES COMPOSANTS IDENTIFIÉS PAR UNE MARQUE A SUR LES DIAGRAMMES SCHÉMATIQUES ET LA LISTE DES PIÈCES SONT CRITIQUES POUR LA SÉCURITÉ DE FONCTIONNEMENT. NE REMPLACER CES COM-POSANTS QUE PAR DES PIÈCES SONY DONT LES NUMÉROS SONT DONNÉS DANS CE MANUEL OU DANS LES SUPPLÉMENTS PUBLIÉS PAR SONY.

SAFETY CHECK-OUT

(US Model only)

After correcting the original service problem, perform the following safety checks before releasing the set to the customer:

- Check the area of your repair for unsoldered or poorly-soldered connections. Check the entire board surface for solder splashes and bridges.
- 2. Check the interboard wiring to ensure that no wires are "pinched" or contact high-wattage resistors.
- Look for unauthorized replacement parts, particularly transistors, that were installed during a previous repair. Point them out to the customer and recommend their replacement.
- Look for parts which, though functioning, show obvious signs
 of deterioration. Point them out to the customer and recommend their replacement.
- 5. Check the line cord for cracks and abrasion. Recommend the replacement of any such line cord to the customer.
- 6. Check the B+ voltage to see it is at the values specified.
- Check the antenna terminals, metal trim, "metallized" knobs, screws, and all other exposed metal parts for AC leakage. Check leakage as described below.



LEAKAGE TEST


The AC leakage from any exposed metal part to earth ground and from all exposed metal parts to any exposed metal part having a return to chassis, must not exceed 0.5 mA (500 microamperes). Leakage current can be measured by any one of three methods.

- A commercial leakage tester, such as the Simpson 229 or RCA WT-540A. Follow the manufacturers' instructions to use these instruments.
- A battery-operated AC milliammeter. The Data Precision 245 digital multimeter is suitable for this job.
- 3. Measuring the voltage drop across a resistor by means of a VOM or battery-operated AC voltmeter. The "limit" indication is 0.75V, so analog meters must have an accurate lowvoltage scale. The Simpson 250 and Sanwa SH-63Trd are examples of a passive VOM that is suitable. Nearly all battery operated digital multimeters that have a 2V AC range are suitable. (See Fig. A)

SERVICE NOTE

Note for Repair

Make sure that the flat cable and flexible board are not cracked or bent at the terminal. Do not insert the cable insufficiently nor crookedly.



Cut and remove the part of gilt which comes off at the point. (Take care that there are some pieces of gilt left inside)




TABLE OF CONTENTS

Sect	ion <u>Title</u>	<u>Page</u>	Section	<u>Title</u> <u>F</u>	Page
SER\	/ICE NOTE	4	•	VA-102 (AUDIO) Schematic Diagram 4	-83
			•	VA-102 (HI MICOM) Schematic Diagram 4	1-86
	•		•	RS-78 Printed Wiring Board 4	-89
1.	GENERAL		•	RS-78, VA-102 (RS MICOM)	
				Schematic Diagram 4	1-89
	ures		•	MD-63, MD-64, MD-65, FP-406	
	s on Video Cassettes			Printed Wiring Boards 4	1-94
	s on Recording /Playback		•	, ,,	
	tion and Function of Parts			Schematic Diagram	
•	pack		•	HP-100 Printed Wiring Board4	
	ections for Playback		•	HP-100 Schematic Diagram4	
	ngs for Playback		•	CM-56 Printed Wiring Board4	
	pack Procedure		•	CM-56 Schematic Diagram4	
	back Functions		•	FR-136 Printed Wiring Board4	
	ording		•	FR-136 Schematic Diagram 4	
	ections for Recording		•	HG-1 Printed Wiring Board4	
	ngs for Recording		•	HG-1 Schematic Diagram4	
Reco	rding Procedure	1-7	•	U-1 Printed Wiring Board 4	
	nging Menu Settings		•	U-1 Schematic Diagram 4	
	ging the SET UP MENU Settings		•	U-2 Printed Wiring Board4	
	Contents		•	U-2 Schematic Diagram 4	1-12:
	n Messages				
	bleshooting		- A	DUCTMENTO	
Note	s on Use	1-9	5. A	DJUSTMENTS	
			5-1. M	IECHANICAL SECTION ADJUSTMENTS 5	5-1
2.	DISASSEMBLY		5-1-1.	Information 5	5-1
				ow to Search Reference Pages for Removal 5	
2-1.	Removal of Upper Case	2-1		hase Adjustment Mark " PH- " 5	
2-2.	Removal of Front Panel Assembly			Preparation for Mechanical Check,	
2-3.	Removal of VA-102 Board			Adjustment and Maintenance 5	5-2
2-4.	Removal of Power Block		2-1.	FL Block Assembly 5	
2-5.	Removal of MD Block Assembly		2-2.	Cassette Positions 5	
2-6.	Removal of CM-56 Board		2-3.	Loading/Unloading 5	
2-7.	Removal of JC-19 Board		2-4.	Manual Up/Down the FL Block 5	
2-8.	Removal of RP-228 Board		2-5.	Service Jigs List 5	
2-9.	Circuit Boards Location		5-1-3.		
L-0.	Official Education	_ 0	3-1.	Phase Adjustment	
			- 1:	(Loading/Unloading Driving Section) 5	5-6
3.	BLOCK DIAGRAMS		3-2.	Phase Adjustment	
Э.	BLOOK DIAGNAMO			(S/L Cassette Selection Section) 5	6-8
3-1.	Overall Block Diagram 1	2.1	3-3.	Phase Adjustment	
3-1.				(Mechanism Chassis Upper Surface Parts) 5	5-10
3-3.	Overall Block Diagram 2 Overall Block Diagram 3		5-1-4.		
J-J.	Overall Block Diagram 3	3-9	4-1.	Cleaning of Rotary Drum Assembly 5	
			4-2.	Cleaning of Tape Path System 5	-12
	PRINTER WIRING POARRO AND		4-3.	Periodic Checks	
4.	PRINTED WIRING BOARDS AND		5-1-5.		
	SCHEMATIC DIAGRAMS		5-1.	Tape Stopper, Motor FPC Assembly and	
			0 1.	Elastic Connector	-14
4-1.	Frame Schematic Diagram	4-1	5-2.	Drum Assembly and Drum Base 5	
4-2.	Printed Wiring Boards and Schematic Diagrams		5-3.	LM Cover, LM Worm Wheel, LM Holder and	
	RP-228 Printed Wiring Board		0 0.	LM Motor Assembly 5	-15
	RP-228 Schematic Diagram		5-4.	TG3/4 Catcher Block Assembly,	. 10
	JC-19 Printed Wiring Board	4-14	J- 4 .	Pinch Driving Gear and TC Arm Assembly 5	-15
	JC-19 (AD/DA CONVERTER)		5-5.	Pinch Arm Assembly, Pinch Limiter and	-13
	Schematic Diagram		5-5.	Tension Coil Spring (Pinch)	-16
	JC-19 (S1 AFC) Schematic Diagram	4-25	5-6.	HC Arm, HC Roller Assembly,	- 10
	JC-19 (U1) Schematic Diagram	4-29	5-0.	Pinch Retainer, Pinch Cam Gear and	
	JC-19 (D1) Schematic Diagram	4-34		TG5/6 Catcher Block Assembly 5	16
	JC-19 (C1 SPCON) Schematic Diagram	4-37	5-7.	RL Arm and RL Link	
	JC-19 (MODE) Schematic Diagram	4-41	5-7. 5-8.	Gooseneck Guard and Gooseneck Arm Assembly 5	
	JC-19 (DV IN/OUT) Schematic Diagram	4-45			-17
	JC-19 (AUDIO CORE) Schematic Diagram	4-50	5-9.	Tension Coil Spring (TG2), Spring Adjustor,	
	JC-19 (DIGITAL AUDIO) Schematic Diagram	4-53		TG2 Spring Hook, TG2 Selection Arm and Damper	-10
	JC-19 (AUDIO D/A, A/D CONVERTER)		E. 10	Sheet	-10
	Schematic Diagram		5-10.	Tension Coil Spring (TG7), Spring Adjustor and	10
	RE-32 Printed Wiring Boards	4-61	E 44	TG7 Spring Hook	
	RE-32 Schematic Diagram		5-11.	TG2 Retainer, TG2 Arm Assembly (TG2 Plate Spring	
	VA-102 Printed Wiring Board			and ET Magnet), S Tension Regulator Band Assemb	
	VA-102 (IF) Schematic Diagram		E 10	and TG2 Load Arm Assembly 5-	
	VA-102 (VIDEO IN) Schematic Diagram		5-12.	TG7 Retainer, TG7 Arm Assembly (TG7 Plate Spring	y Ju
	VA-102 (VIDEO OUT) Schematic Diagram	4-75	100	and ET Magnet), T Tension Regulator Band Assembly	
	• VA-102 (UVIC, DV IN/OUT) Schematic Diagram	4-79	E 40	and TG7 Load Arm Assembly 5-	
	A VA 100 (MONITOR OUT) Cohematic Diagram	4.04	5-13.	S Reel Table Block Assembly 5-	-21

Section	<u>Title</u>	<u>Page</u>	<u>Section</u>	<u>Title</u>	<u>Page</u>
5-14.	T Reel Holder and T Reel Table Block Assembly		1.	Power Supply Voltage Check, Power Block (U-2 Board)	5-62
5-15.	S Reel Plate Assembly		0	Video/Audio Block Power Supply Voltage Check,	. 5-02
5-16.	T Reel Plate Assembly		2.	Power Block (U-2 Board)	5-62
5-17.	TG1/8 Base Assembly, S Brake and T Ratchet		22 64	stem Control System Adjustment	
5-18	TG3/6 Roller Assembly and TG3/6 Lock Spring	3-20		Initializing the C, D, E Page Data	
5-19.	FL Joint Gear, TG5/6 Retainer and	E 06	1. 2.	Input of C Page Initial Data	
	Capstan Motor	5-26		Input of D Page Initial Data	
5-20.	FL Selection Arm, FL Relay Gear and FL Joint Arm		3.		
	Assembly		4.	Input of E Page Initial Data	5-63
5-21.	Rotary Switch, TC Gear and Relay Gear		5. 0	Page C Address List	
5-22.	GL Arm Retainer and GL Arm	5-28	6.		
5-23.	M Slider and M Slider Arm	5-28	7.	Page D Address List	
5-24.	TG7 Selection Arm, TG7 Cam Gear and T Cam Gear	5-28	8.	Page E Address List	
5-25.	Main Cam, TG2 SL Arm Assembly and		3-4. Se	ervo System Adjustments	. 5-66
	Tension Coil Spring (TG2 SL)	5-29	1.	Switching Position Adjustment (CM-56 Board)	. 5-66
5-26.	TG3/4 Arm Block Assembly (TG3/4 Arm Assembly	,	2.	Capstan FG Duty Adjustment (CM-56 Board)	
	TG3/4 Limiter Spring and TG3/4 Gear), TG3/4 Bas		3-5. Vid	deo System Adjustments	. 5-67
	Block Assembly (TG3/4 Base Assembly)			RP-228 Board Adjustments	
5-27.	TG5/6 Arm Block Assembly (TG5/6 Arm Assembly		1.	Recording Current Adjustment (RP-228 Board)	. 5-67
0 2.7.	TG5/6 Limiter Spring and TG5/6 Gear), TG5/6	•	2.	PLL fo Adjustment (RP-228 Board)	
	Base Block Assembly (TG5/6 Base Assembly)	5-32	3.	CLK Delay Adjustment (RP-228 Board)	
5-28.	Reel Motor		4.	AGC Center Level Adjustment (RP-228 Board)	
5-20. 5-29.	RS Arm Assembly		5.	AEQ Adjustment (RP-228 Board)	
	RS Gear Assembly, Mic Press Spring and	0-04	6.	PLL Capture Range Adjustment	
5-30.		5-24	Ų.	(RP-228 Board)	5-70
E 04	Mic Lever	J-J4	7.	IC774 41.85 MHz VCO Check (RP-228 Board)	
5-31.	Rack Joint Gear, Rack Holder, Mic Holder,	E 25		JC-19 Board Adjustments	
c 00	Rack (LC) and Rack (SC)	5-35 E 36		A/D Converter Reference Voltage Adjustment 1	
5-32.	Plate Link Assembly		1.	(JC-19 Board)	5.71
5-33.	Roller Shaft Assembly and Roller Belt		•	A/D Converter Reference Voltage Adjustment 2	. 5-71
5-34.	Lid Opener		2.		5-71
5-35.	C Door			(JC-19 Board)	5-71
5-36.	Damper Arm and Tension Spring (DB)		3.	Y Signal Clamp Reference Voltage Adjustment	E 71
5-37.	Gear (A), Gear (B), and C Worm	5-40		(JC-19 Board)	5-71
5-38.	Tension Coil Spring (HS), Tension (DB), Shift Plate		4.	CR Signal Clamp Reference Voltage Adjustment	E 74
	Spring and C Sloat Block Assembly		_	(JC-19 Board)	. 5-/1
5-1-6.			5.	CB Signal Clamp Reference Voltage Adjustment	F 74
6-1.	Adjustment Position	5-42		(JC-19 Board)	5-/1
6-2.	Adjustment Order	5-43	6.	Playback Y Signal Level Adjustment	F 70
6-3.	Adjustment and Checking Method	5-44		(JC-19 Board)	5-72
6-3-1.			7.	Playback CR Signal Level Adjustment	
6-3-2.	TG2/TG7 Height Adjustment	5-45		(JC-19 Board)	5-72
6-3-3.		5-45	8.	Playback CB Signal Level Adjustment	
6-3-4.		5-46		(JC-19 Board)	
6-3-5.	FWD/RVS Back Tension Adjustment	5-48	9.	IC422 27MHz XTAL fo Adjustment (JC-19 Board)	
	Adjustment Preparations and		10.	AFC Preliminary Adjustment (JC-19 Board)	
	RF Waveform Check	5-49	11.	AFC Picture Frame Adjustment (JC-19 Board)	5-74
6-3-7.	Tracking Adjustment		12.	AFC Adjustment (JC-19 Board)	5-74
6-3-8.			3-5-3.	VA-102 Board Adjustments	
	CUE and REV Check	5-52	1.	AGC Adjustment (VA-102 Board)	. 5-75
	. Rising Check		2.	Analog E-E VIDEO Signal Output Level	
	. Tape Path Check			Adjustment (VA-102 Board)	. 5-75
0-0-11	. Tape I dill Official III.	• • •	3.	Analog E-E Y Signal Output Level Check	
5-2. SI	ERVICE MODE	5-55	<u> </u>	(VA-102 Board)	5-76
	Adjusting Remote Commander		4.	Analog E-E Chroma Signal Output Level Check	
	Used Adjustment Remote Commander	5-55		(VA-102 Board)	5-76
1.		0-00	5.	Decoder VXO Freerunning Frequency Adjustmen	
2.	Precautions Upon Using The Adjusting Remote	5-55	J.	(VA-102 Board)	
	Commander		6	Video Input Y/C Separation Adjustment	0
5-2-2.	Data Processing		6.	(VA-102 Board)	5,77
5-2-3.			-	Decoder HUE Adjustment (VA-102 Board)	5.79
1.	Emergence Memory Address		7.	Decoder HUE Adjustment (VA-102 Board)	5-70 5-70
1-1.	EMG Code (Emergency Code)	5-57	8.	REC Y Level Adjustment (VA-102 Board)	5-79 E 70
			9.	REC CR Level Adjustment (VA-102 Board)	5-79
5-3. VI	DEO SECTION ADJUSTMENTS	5-58	10.	REC CB Level Adjustment (VA-102 Board)	5-80
3-1. Pi	reparations Before Adjustment		11.	Encoder Freerunning Frequency Adjustment	. 5.00
3-1-1.		5-58		(VA-102 Board)	
	Connection of Equipment		12.	Playback Y Level Check (VA-102 Board)	5-80
	Adjusting Connectors (RP 228 Board CN775)		13.	Playback Chroma Level Adjustment	
	Checking the Input Signals	5-60		(VA-102 Board)	
1.	S Video Input		14.	Playback Burst Level Adjustment (VA-102 Board)	
2.	Video Input			(PAL Model Only)	5-81
	Alignment Tapes		3-5-4.	General Adjustment	5-82
	Input/Output Level and Impedance		.1.	Y Output Level Adjustment (JC-19 Board)	. 5-82
3-2. Po	ower Supply System Adjustment	5-62			

<u>Section</u>	<u>Title</u>	<u>Page</u>
2.	Encoder R-Y Input Level Adjustment (JC-19 Board)	5-82
3.	Encoder B-Y Input Level Adjustment	
	(JC-19 Board)	5-83
4.	Decoder HUE Input Adjustment (JC-19 Board)	
5.	Battery Down Adjustment and Confirmation	
3-5-5.	BIST Check	5-85
1.	Playback System Check	
	(JC-19, RP-228 Boards)	5-85
2.	Record System Check	5-86
3-6. Au	dio System Adjustments	
1.	Playback Level/Indicator Check	
2.	E-E Level Check	5-87
3.	Recording/Playback Check (Audio Lock Mode)	
3-7. Arı	rangement Diagram for Adjustment Parts	5-88
6. RI	EPAIR PARTS LIST	
6-1. Ex	ploded Views	6-1
6-1-1.	Overall Assembly	6-1
	Chassis Assembly	
	MD Block Assembly	
6-1-4.	FL Cassette Compartment Assembly	6-4
6-1-5.	Mechanism Chassis Assembly (1)	
	(Top Side View (1))	6-5
6-1-6.	Mechanism Chassis Assembly (2)	
	(Top Side View (2))	6-6
6-1-7.	Mechanism Chassis Assembly (3)	
	(Top Side View (3))	6-7
6-1-8.	Mechanism Chassis Assembly (4)	
	(Bottom Side View (1))	6-8
6-1-9.	Mechanism Chassis Assembly (5)	
	(Bottom Side View (2))	6-9
6-2. Ele	ectrical Parts List	
Hardwar	e l ist	6-42



Features

The DSR-20MD/20MDP is a 1/+ inch digital video The DSR-20MD/20MDP is a \(\frac{1}{\sinch}\) inch digital video cassette recorder that uses the DVCAM digital recording format. This system achieves stable, superbicture quality by digitally processing video signals that are separated into color difference signals and luminance signals (component video). The unit is equipped with a full-fledged analog interface to support hybrid systems that combine conventional analog equipment with digital equipment,

The DSR-20MD/20MDP's main features are described

DVCAM Format

DVCAM is based on the consumer DV format, which uses the 4:1:1 component digital format (DSR-20MD) or the 4:2:0 format (DSR-20MDP), and provides a '4-inch digital recording format for professional use.

High picture quality, high stability

Video signals are separated into color difference Video signals are separated into color distrence signals and huminance signals, which are encoded and compressed to one-fifth size before being recorded to ensure stable and superb pliciture quality. Because the recording is digital, multi-generation dubbing can be performed with virtually no deterioration of quality.

Wide track pitch

The recording track pitch is 15 µm, fully 50 percent wider than the DV format's 10 µm track pitch. Thanks to this feature, the DVCAM format sufficiently meets the reliability and precision requirements of professional editing.

High-quality PCM digital audio

PCM recording makes for a wide dynamic range and a high signal-to-noise ratio, thereby enhancing sound quality.
There are two recording modes: 2-channel mode

nere are two recording modes: 2-channel mode (48 kHz sampling and 16-bit linear code), which offers sound quality equivalent to the DAT (Digital Audio Tape) format, or 4-channel mode (32 kHz sampling and 12-bit nonlinear code).

Playback compatibility with DV format

A DV cassette recorded on a DV-format VCR can be played back on this unit. (Casacites recorded in LP mode cannot be played back.)

Choice of two cassette sizes

- The unit can use both standard-size and mini-size DVCAM cassettes.

 •According to cassette size, it automatically changes the position of the real drive plate.

 •The maximum recording/playback times are 184 minutes for standard size cassettes and 40 minutes for

Other Features

Compact size

The unit achieves compact size suitable for using on a demonstration. The unit is also equipped with basic functions that are needed for videocassette recorders and players used in professional digital video editing

- DC IN connector

Menu system for functionality and operation settings

The unit provides a menu system to make its various functions easier to use and set up its operation conditions.

Superimposition function

Time code, menus, error messages, and other text data can be superimposed and output in analog composite video signals.

Remote control

The unit can be operated by remote control from an editing controller that supports the RS-232C interface or from a SIRCS¹⁰-aystem remote controller or foot switch such as the optional DSRM-10 or SVRM-100A.

1) SIRCS
SIRCS-system has the same function as CONTROL S-system

star 1 Ovandaw 36

Features

Notes on Video Cassettes Usable cassettes

Uso Standard-DVCAM cassettes or Mini-DVCAM cassettes with this VCR. PDV-184ME can record programs for 184 minutes and PDVM-40ME can record for 40 minutes. You can get the highest quality pictures with this digital video cassette recorder using DVCAM cassettes, You may not be able to get as good quality with other cassettes. We recommend using DVCAM cassettes so that you can record your one-time events in highest quality.





Cassette memory

Castette memory is an optional feature that is mounted on some Standard DVCAM cassettes and Mini DVCAM cassettes. When you record a program, the recording date and time, and the program's position on the tape are stored in the cassette memory so that you can quickly locate the program later on. (H16K indicates that you can use the cassettes 16 kbits of date can be stored on. On this VCR, you can use the cassettes up to 16 kbits of data can be stored on. On this VCR,

To save a recording

To prevent accidental erasure of a recording, slide in the safety switch on the cassette so that the red portion becomes visible. To record on a tape, slide out the switch so that the red portion is hidden.




Note

[DVCAM], DY, "DY and CI'l are trademarks.

Checking the tape for slack

Using a paper cilp or a similar object, turn the reel gently in the direction shown by the arrow. If the reel does not move, there is no slack, insert the cassette into the cassette compartment, and after about 10 seconds take it



Notes on Recording / Playback

Copyright precautions

On recording

You cannot record any software having copyright protection signals on this VCR. If you start recording protected video and audio signals, a warning message appears on the monitor screen and the VCR stops recording.

recording.
On playback
When you play back software having copyright protection signals on this
VCR, you may not be able to copy it onto other equipment.

Limitations caused by the difference in format

This VCR can record, play back and edit the tapes recorded in DVCAM format. It can also play back the tapes recorded in DV format (SP mode). However, due to the difference in format, you may not be able to record or edit some tapes affected by recording conditions of the tape (e.g., A tape originally recorded in DV format is dubbed in DVCAM format). For details, refer to "Compatibility of DVCAM and DV format" on page 42,

No compensation for contents of the recording

Contents of the recording cannot be compensated for if recording or playback is not made due to a malfunction of the VCR, video tape, etc.

Note

You cannot play back a DVCAM tape recorded in other color systems on this VCR.

O Cassette compartment

Accepts standard-size or mini-size DVCAM digital videocassettes. When using a mini-size cassette, insert it into the middle of the compartment. For details of usable cassettes, see page 4.

♦ © ON/STANDBY switch and ON/STANDBY lamp Press this switch to turn on the power, and the ON/ STANDBY lamp lights in green. Press it again to turn to standby mode, and the lamp goes off.

When the REMOTE/LOCAL switch is set to REMOTE, you cannot turn the unit to standby mode.

during recording or playback.

The audio signal you want to monitor can be selected


with the AUDIO MONITOR selector inside of the door ([3]).

Ø → PHONE LEVEL control knob Controls the volume of the headphones connected to the \(\) PHONES jack.

• COUNTER RESET button

Press this button to reset the tape counter in the display window to "0:00:00 (0:00:000)." This button does not work when displaying the time code or the remaining time.

Press this button to eject a cassette. 1 Display section



O Audio tevel meter

Indicates the recording level during recording or EB mode" and the playback level during playback. When the audio level exceeds 0 dB, the red indicator lights.

If you play back the tape whose audio was only recorded on channel 2, the audio level meter for CH2/4 may not function.

- Time counter display Indicates the following:

 Time data: count value of the time counter, time code and remaining time
 Alarm messages (see page 35)
 Messages for solf-diagnosis function (see page 39)

- For DSR-20MDP: Time code is set to the non drop
- frame mode only.
 Time code is indicated as follows:
 Drop frame: "00:00.00:00" ("00:00,00:00" on the monitor) (DSR-20MD only) Non drop frame: "00:00:00:00"

Video input signals indicators

Indicates the currently selected video input signals.
INPUT VIDEO, INPUT S VIDEO or INPUT DV

O NS (Non Standard) AUDIO indicator Unit (Non Standard) AUDIO indicator
Lights when the VCR plays back a tape whose audio recording was made in the unlock mode, or when unlock mode signals are input through the DV Lick.

For details of unlock mode, see page 42.

 Audio mode indicators Indicates the audio mode during playback or recording, or while in BE mode,

During playback it indicates the audio mode in which

the tape was recorded.

During recording or while in BE mode, it indicates
the currently aclected audio recording mode. You can
select audio recording mode by setting "AUDIO
MODE" menu (see page 33).

F\$32k: Lights when playing the tapes recorded in
4-channel mode, or recording a tape in 4-channel
mode.

P\$48k: Lights when playing the tapes recorded in
2-channel mode, or recording a tape in 2-channel
mode. the tane was recorded.

Note

When recording in 4-channel mode on this VCR, audio signals are recorded only in channels 1/2.

@ TIMER indicator Lights when setting the TIMER switch to REPEAT or REC.

O DVCAM indicator Lights except playing back the DV-formatted tapes.

1) EE mode

EE mode
"HE" stands for "Electric to Electric". When in this mode, the video and audio signals that are input to the VCR's
recording circuitry do not pass through any magnetic conversion circuits but instead are output via electric circuits only.
This mode is used to check input signals and adjust input levels.

Chapter 1 Overview 7^{as}


Location and Function of Parts

O Cassette indicator

Lights when inserting a digital video cassette available for this VCR. It flashes when ejecting a cassette,

Mini indicator
Lights when inserting mini-size digital video cassette,

2 Tape transport control section



• Tape Transport indicators

2 DUP (duplicate) button

Use this button to make a work tape having the same time code as the source tape.

For details of duplicate, see page 30.

REC (record) button

When you press the PLAY button while holding down this button, the indicator lights and recording begins. To set the VCR to recording pause mode, press the II PAUSE button while holding down this button.

II PAUSE button When you press this button, the indicator lights and the VCR is set to pause mode.

♦ ►► FWD (fast forward) button When you press this button, the indicator lights and the tape is fast forwarded. During fast forward, the picture does not appear on the monitor (you can see the picture of the EE mode during fast forward).
To search forward, hold this button down during fast

PLAY button

When you press this button, the indicator lights and playback begins. If you press this button while holding down the ◀ REW button during stop, the tape is rewound to its beginning and starts playing automatically (during rewind, the REW indicator lights and the PLAY indicator flashes).

● **REW** (rewind) button


Q ◆ REW (rewind) button when you press this button, the indicator lights and the tage starts rewinding. During rewind, the picture does not appear on the monitor (you can see the picture of the EE mode during rewind). To search backward, hold this button down during

If you press the PLAY button while holding down It you pless the PLA1 button winte nothing at this button during stop, the tape is rewound to its beginning and starts playing automatically (during rewind, the REW indicator lights and the PLAY

G # STOP button

button to stop the current tape transport

3 Inside of the door



O INPUT SELECT button

Select video input signals. Each press of this button cycles through three video signal selection options: video, S-video, and DV input. When you select one of these options, the corresponding indicator in the display lights up.

O COUNTER SELECT selector

Select the type of time data in the time counter display.
TC: Time code
COUNTER: Count value of the time counter
TAPE REMAIN; Remaining time

 AUDIO MONITOR selector Use to select the audio track you want to hear when playing back a tape recorded in 4-channel mode (Fs32k).

PS32k), CH-1/2: Channels 1/2 only MIX: Channels 1/2 and channels 3/4 (mix) CH-3/4: Channels 3/4 only

O TIMER selector

▼TIMER selector
Uses to select timer recording or outo repeat using an external AC timer (not supplied).
REPEAT: When the power is supplied to this VCR, a tupe rewinds to its beginning automatically and playback starts. The VCR repeats the playback from the beginning to the first index (if there is no index on the tape, to the unrecorded portion; if no unrecorded portion, to the tape end). Autor repeat also functions if you set this selector to REPEAT during playback. during playback. OFF: Timer is released.

REC: When the power is supplied to this VCR, recording starts.

⊘ REMOTE/LOCAL switch

Selects whether the unit is operated from its front panel or from external (remote) equipment.

REMOTE: The unit is operated from an editing controller connected to the RS-232C 🗷 conne

controller connected to the RS-232C 22 connector
on the rear panel. No operation on the front panel
works except sliding the switch or selectors.
LOCAL: The unit is operated from its front panel,
from an external equipment connected to the
LANC e jack on the rear panel, or from a SIRCSsystem remote controller connected to the
REMOTE 22 CONTROL S jack on the rear panel.

@ MENII button

→ MENU button

Press this button to display the menu on the monitor
screen. Press it again to return from the menu display
to the usual display.

Note

If you set the REMOTE/LOCAL switch to REMOTE while the menu display is on the monitor, it returns to the usual display.

On how to use the menu, see Chapter 3 "Menu Sestings."

SET button

Press this button to save selected menu items to the unit's memory.

② Arrow (分 号) buttons
Use these buttons to move around the menu items.

Chapter 1 Overview 9 68

8^{aa} Chapter 1 Overview

OUTPUT @ connector O INPUT -O connectors

• INPUT -o connectors
Input video and audio signals. To connect a VCR equipped with the S-video output jack, use the S VIDEO jack on this VCR.

OUTPUT or connectors

Output video and audio signals. To connect a VCR equipped with the S-video input jack, use the S VIDEO jack on this VCR.

MONITOR & connectors
 Output video and audio signals for monitoring.

O SYNC switch

Selects the reference signal. The video signal is locked to V-sync or H-sync, but not locked to sub-career. The sync phase is not adjusted. The video signal is not

to v-sync or H-sync, but not locked to sub-career. The sync phase is not adjusted. The video signal is not locked to DV input. INT: Selects the playback signal on this VCR as the reference signal. EXT: Selects the input video signal from the external equipment connected to this VCR as the reference signal.

Notes

- NAMES

 The picture and the sound may be distorted if:

 You set the SYNC switch during playback.

 The analog signal is input from the INPUT ocnnectors during playback with the SYNC switch act to EXT.
- If the SYNC switch is set to EXT during playback, the INPUT SELECT button does not work.

10⁶⁸ Chapter I Overview

@ REMOTE 122 CONTROL Sinck

Connect a SIRCS-system remote controller. When controlling this VCR from a remote controller such as the DSRM-10 or SVRM-100A (not supplied), connect the unit to the editing controller via this jack.

Note

SIRCS-system has the same function as CONTROL S-system.

© REMOTE № FOOT SWITCH ≥ jack
Connect the optional Foot Switch to control this VCR.

The Foot Switch must be conformed with Standard UL2601-1/EN60601-1,

For details on the Foot Switch, consult with authorized Sony

O LANC & lack

When you connect the LANC 9 jacks on this VCR and the other VCR, you can control this VCR (player) from the other VCR. The LANC connection transmits signals such as control signals, time code and time counter data and status data. You can control this VCR by connecting the optional RM-95 Remote Commander to this jack.

- The other VCR (recorder) receives the time code data from the LANC o Jack only when this VCR (player) is set to show the time code indications.

 If the REMOTE/LOCAL switch is set to REMOTE,

⊕ RS-232C
 □ connector (D-sub 9-pin)
 Connect an editing controller or a personal computer with the RS-232C interface for remote-control of this VCR.

DV I Jack
 The DV I jack is i.I.RNK compatible. Use when the equipment connected to the VCR has a DV I jack. If you connect the VCR and the other equipment using DV I jacks, you can minimize deterioration of picture quality during dubbing, editing or capating setsill pictures into a personal computer by digital processing. For details, refer to the instruction manual of the equipment you use.

is a trademark of Sony Corporation and indicates

that this product is in agreement with IEEE1394-1995 specifications and their revisions.

@ ~ AC IN connector

Connect to an AC power outlet using the supplied power cord.

D == DC IN connector Connect the optional BP-90A Ni-Cd Battery Pack with the battery adaptor and DC cable.

- If the voltage of the Battery Pack falls less than 11 V, a beep sound is output (when BEEP in the menu is set to OFF, it is not output) and "delo" appears in the display window. Replace the battery by a charged one or remove it to use the AC power outlet.

 If the voltage of the Battery Pack falls less than 10.5 V, a beep sound is output (when BEEP in the menu is set to OFF, it is not output) and the VCR is set to the standby mode. As you cannot turn on the VCR at this moment, replace the battery by a charged one or remove it to use the AC power outlet.

② Equipotential ground terminal

Used to connect to the equipotential plug to bring the various parts of a system to the same potential. Refer to "Important asfeguards/ordices for use in the medical environment" on page ii.

Chapter 1 Overview 11 0

Location and Function of Parts

Supplied Remote Commander 0 **⊕ EJECT** butto O On/at/ O COUNTER RESET 35 DISPLAY bullon O DATA CODE button SEARCH SELECT 0:0 Buttone for playing a various eposds 88 00 REW PLAY PE O >> FF button © ■ STOP button

- O COUNTER RESET button
- 6 INPUT SELECT button
- ♦ SEARCH SELECT buttons
 Press these buttons to search for sonnes using the index riess these buttons to search for scenes using the in function.

 For details, see "Searching using the index function" on page 19.

12^{GB} Chapter 1 Overview

- 6 Buttons for playing at various speeds
- B buttons
- × 1/10 button × 1/5 button × 1/5 button

FRAME - 11/11 buttons
For details, see "Ploying at various speeds" on page 18.

O II PAUSE button

O > PLAY button Q - REW button

O On/standby switch

© DISPLAY buttonPress this button to see the indications, such as tape counter, on the monitor screen.

O DATA CODE button Press this button to see tape information on the monitor screen.
For details, see "Displaying tape information" on page 22.

@ Buttons for menu operation MENU button

SET button

REC buttons
 When you press these buttons at the same time, the indicator lights and recording begins.

►► FF button

@ STOP button

When using the supplied remote commander, set REMOTE CONTROL in the menu to VTR4 (see page 34). Otherwise, you cannot operate this VCR with the supplied remote commander.

Battery installation

1 Push and slide the lid to open.



2 Install the two size AA (R6) batteries (supplied) ith the correct polarity



3 Replace the lid.



- Notes on batteries

 Make sure that the battery orientation is correct when inserting batteries.

 Do not mix an old battery with a new one, or
- different types of batteries.

 If you will not use the Remote Commander for a long time, remove the batteries to avoid damage from battery leakage. If batteries have leaked, remove them, wipes the battery compartment dry and replace the batteries with new ones.

Chapter 1 Overview 13⁶⁹


Playback

This section describes the necessary connections, settings and operations to perform playback on this unit. The same settings and operations apply whether you are using the unit as part of an editing system, for dubbing, or as a stand-alone videocassette player,

Connections for Playback

To digital video equipment with DV jack

The video and audio signals are sent with hardly any degradation, enabling high-quality editing. The signal flow is automatically detected so you need not make separate connections for input and output.



- Set DV RE OUT in the menu to OFF (see page 24).


 *Audio signals are not output during playing at various speeds.

 *With the DV connection, the sound is recorded in the same audio recording mode as that of the source tage. To exect all a different audio recording mode from the source tage, to sected in a different audio recording mode from the source tage, use the INPUT-© connectors
- instead.

 With the DV connection, tape information (recording date, camcorder data, sic.) recorded on the source tape is transmitted from this VCR (player). As a result, when you play back a recorded tape and press the DATA CODE button, the same tape information recorded on the source tape is displayed on the monitor screen. However, contents of the cassette memory are not transmitted. In addition, the time code is newly recorded on the tape on the other VCR, except when copying a tape in Duplicate mode.
- mode.
 As for the LANC connection, see "Notes for LANC connection" on the

Chapter 2 Playback and Reco.

To video equipment without DV jack



- When you connect output lacks of the recorder to input jacks of this VCR, select the input correctly to prevent a humming noise.

 Distorted signals (e.g., when played back at a speed other than normal) will not be recorded properly.

 The indication displayed on the monitor screen are output only via the MONITOR G- connector.

Notes for LANC connection

- With the LANC connection, refer to the instruction for use supplied with the recorder VCR.

- the recorder VCR.

 *The LANC connection transmits signals such as control signals, time code, time counter data and status data.

 *If the other VCR has a LANC @ jack of 5-pin DIN type, connect with the VK-8 lo Control L connecting cable (not supplied).

 *The Jacks labeled CONTROL L have the same function as LANC @ jacks. The jacks labeled REMOTE on other equipmont may also have the same.

 *The other VCR (recorder) receives the time code data from the LANC @ jack only when this VCR (player) is set to show the time code indications.

 *With the LANC connection, this VCR only works as a slave unit.

Chapter 2 Playback and Recording 15^{as}

Playback

Settings for Playback

Preparation on the player (this VCR)

- 1 Power on the video monitor, then set the monitor's input according to the input signals from the recorder.
- 2 Set up the recorder.

 For details, see "Preparation on the recorder" below.
- 3 Power on this unit by pressing the ON/STANDBY switch.

enter entre en une deputer en entre en La companya de la co

The ON/STANDBY lamp lights in green.

- 4 If the other equipment that controls this VCR has the time code function, set the COUNTER SELECT selector to TC (see page 9).
- When you play back a tape recorded in 4-channel mode (Fe 32k), set the AUDIO MONTTOR selector to MIX (see page 9). Then select the precise belance between the tracks with the AUDIO MIX BALANCE in the menu (see page 3).

- With the DV connection, the playback VCR's AUDIO MONITOR (sound selection) and AUDIO MIX BALANCE (audio balance adjustment) do not function on the source audio output through the DV I.
- jock.

 You cannot change the input signal selection during playback or playback pause mode.

Preparation on the recorder

- Insert a tape for recording.
 Select the formats of video and audio input signal to be recorded.
 Set the LANC mode to M.

- Editing is not possible with a tape that is copyright protected.

 You cannot use the video equipment that has no LANC mode switch as a


Playback Procedure

- When controlling this unit from an editing controller or a personal computer, set the REMOTE/LOCAL switch to REMOTE. When not, set the switch to LOCAL (See page 9).

 Do not insert the cassette forcibly. The VCR may be damaged.
- After checking the tape for slack, hold the cassette so that the tape window is facing upward, then insert it into this unit as illustrated below.

 For details of checking the tape for slack, see page 5.

The cassette is automatically drawn into the unit.



2 Press - PLAY.

This starts the playback operation.

Searching using the index function

Three kinds of search are available on this VCR:

PHOTO SEARCH.

INDEX BEARON

Three kinds of search are available on this VCR:

Searching for the beginnings of recordings: Index search

Searching for a point on the tape where the recorded date changes:

Date search

Searching for scenes recorded in the photo mode with a digital

camcorder: Photo search

Searching with the cassette memory
If the tape has a cassette memory, the recordings are listed in the
chronologically in the order they were made. You can search using this enronological nat.

If the tape does not have a cassette memory, you cannot search for scenes in the chronological order.

1 Press SEARCH SELECT to select the search type: INDEX, DATE or

The VCR starts searching and when it locates the recording, begins playback. During Photo search, the VCR pauses.

The chronological list appears on the monitor screen.

2 Press + or > to select a recording.

About the cassette memory

If you use a tape with the mark, the cassette memory stores up to 135 index signals. (The number changes depending on the data size combination of index, date, and photo data stored on a tape.) This VCR is capable of storing and retrieving up to 16 kbits of cassette memory.

To locate recordings whose signals are disabled to be stored in the cassette memory, or to locate recordings in order of their position on the tape, set CASSETIE MEMORY SHARCH in the menu to OFF (see page 34). You can use the same procedure to search for a recording on a tape without cassette memory.

Notes

**Each program is indexed at its beginning. If you record another program over the beginning of the first program, you will not be able to locate the original program.

_	Inde	calge	al la	
Α		В	С	If D is recorded over the beginning of B
	ŧ		ganno	be searched
A	D	В	C	

You cannot add indexes after recording.
To add indexes only for Auto Repeat, start recording from the point you want to start indexing.
You cannot erase indexes after recording.
To delete indexes for Auto Repeat, set INDEX WRITE in the mens to OFF (see page 34). Then record over the index signal you want to crase.
Searching may not be done correctly if the signals were not recorded on a Sony-brand digital video equipment.

Playback Functions

Playing at various speeds

You can enjoy playback functions using supplied remote com-

Send for the end of the send of the end of t

Playback options	Operation
Play at 1/10 of normal speed	Press x 1/10 during playback
Play at 1/5 of normal speed	Presa x 1/6 during playback
Play at normal speed	Press x 1 during playback
Play at twice the normal speed	Press x 2 during playback
Play frame by frame	Press FRAME -11/11> during pause.

To hear the sound during playing at various speeds If you want to hear the sound during playing at various speed WITH SOUND in the menu to ON (see page 33). eeds, set JOG

18^{GB} Chapter 2 Playback and Re-

Searching without cassette memory
When you use a tape without a cassette memory, the VCR searches in the
order of the actual positions of the recordings, regardless of the setting of
CASSETTE MEMORY SEARCH in the meau.

When you use a tape with a cassotte memory, set CASSBTTE MEMORY SEARCH in the menu to OFF (see page 34).

1 Press SEARCH SELECT to select the search type.

INDEX 0

2 Press or >> repeatedly to locate the recording you want.

The VCR starts searching backwards or forwards until the index number comes to zeto, then plays back the recording. During Photo search, the VCR pauses.

How signals are recorded

The VCR marks the tape when TREC button is pressed.

There are three different signals for each search method, The type of signal recorded and where it is recorded (on the tape or in the cassette memory) depends on the video equipment used for recording, Please note that if the signals for certain search type are not recorded, you cannot do that type of

When you record with a Sony digital camcorder

Signals for	In cassette memory	On tape
index search*	No	No
Date search	Yes	Yes
Photo search	Yes	Yes

When you record on this VCR

Signals for	in cassette memory	On taps	
Index search*	Yes	Yes	
Date search	No	Yes	
Photo search	No	No	


The signals for Index search are recorded when you start recording in stop mode.

Chapter 2 Playback and Recording 21^{db}

Displaying tape information

If you record on a tape using a Sony digital camcorder DSR-200/200P/ 200A/200A/PD100/PD100/PD100A/PD

Press DATA CODE during playback.
Each time you press DATA CODE, the display changes as follows.



- •When the information was not recorded, "---" appears instead,
 •The carneorder data displayed on the monitor screen by this VCR are
 partially different from those shown by the digital carneorder.

Auto Repeat

This VCR can repeat the playback of all or a part of the tape.

1 Set the TIMER selector on the front panel to REPEAT.


The TIMER indicator on the front panel lights.

- 2 Press REW to rewind the tape to its beginning.
- 3 Press > PLAY.

Playback starts automatically. The VCR repeats the playback from the beginning to the first index (if there is no index on the tape, to the surrecorded portion; if no unrecorded portion, to the tape end).

Auto Repeat using an external AC timer
If you connect an external AC timer (not supplied) to this VCR, you can
repeat playback automatically at the preset time.

1 Connect an external AC timer (not supplied) to this VCR.



2 Set the TIMER selector on the front panel to REPBAT.

The TIMER indicator in the display window lights,

3 Set the timer-on time on the external AC timer.

At the preset time, the power turns on, and Auto Repeat playback starts automatically within one minute. The VCR repeats the playback from the beginning to the first index (if there is no index on the tape, to the unrecorded portion; if no unrecorded portion, to the tape end).

22^{GP} Chapter 2 Playback and Recording

Playback

- NAGES

 The VCR cannot search for an index or unrecorded portion within 20 seconds from the beginning of the tape.

 While a tape is running, do not turn off the power using an AC timer. The VCR and a tape may be damaged. When turning off the power of the VCR, make sure to press the STOP button on this VCR first to stop the tape transport, then turn off the power.

To stop Auto Repeat Press the STOP button.

To release Auto Repeat mode Set the TIMER selector to OFF.

Chapter 2 Playback and Recording 23⁰⁸


Recording

Ints section describes the necessary connections, actings and operations to perform recording on this unit. The same settings and operations apply whether you are using the unit as part of an editing system, for dubbing, or as a stand-alone recorder.

Connections for Recording

To digital video equipment with DV jack

The video and audio signals are sent with hardly any degradation, enthigh-quality editing. The signal flow is automatically detected so you not make separate connections for input and output. ted so you need



- Audio signals are not output during playing at various speeds.
 With the DV connection, the sound is recorded in the same audio recording made as that of the source tage. To record in a different audio recording mode from the source tage, use the INPUT o connectors.
- Instead.

 With the DV connection, tape information (recording date, camcorder data, etc.) recorded on the source tape is transmitted from the other VCR (player). As a result, when you piley back a recorded tape and press the DATA CODE button, the same tape information recorded on the source tape is displayed on the monitor screen. However, contents of the cassette memory are not transmitted. In addition, the time code is newly recorded on the tape on this VCR, except when copying a tape in Duplicate mode.

 As for the LANC connection, see "Notes for LANC connection" on the next page.

Chapter 2 Recording and Playback 25^{cs}

Settings for Recording

Preparation on the recorder (this VCR)

- **Before recording, set the clock on the VCR so that the recording time can be written into the index signal. You can set the clock by setting the CLOCK SET menu (see page 34).

 When controlling this unit from an editing controller or a personal computer connected to the RS-232C C2 connector, set the REMOTE/LOCAL switch to REMOTE. When not, set the switch to LOCAL (See
- page 9).

 Editing is not possible with a tape that is copyright protected.
- Power on the video monitor, then set the monitor's input according to the input signals from this unit.
- 2 Set up the player to play back a tapo.

 For details, see "Preparation on the player" on the next page.
- 3 Power on this unit by pressing the ONSTANDBY switch.

The ON/STANDBY lamp lights in green.

4 Use the COUNTER SELECT selector to select the type of time data to

Type of time data	Set the selector to
Count value of the time counter	COUNTER
Time code	TC

5 Select the video and audio input signals to be recorded.


Press INPUT SELECT to select the desired signal. Each press of this button cycles through three video signal selection options: video, S-video, and DV input. Each selection is shown by a lit indicator in the display window.

Note

Once you have started recording, you cannot change the input signal selection (except during recording pause mode).

Chapter 2 Recording and Playback 27^{as}

To video equipment without DV jack



- When recording the analog input signals, this VCR can digitally output the signals from the DV & jack for backup. Set DV EB OUT in the menu to ON (see page 34).

 When you connect output jacks of this VCR to input jacks of the player, select the input correctly to prevent a humming noise.

 Poistorted signals (e.g., when played back at a speed other than normal) will not be recorded properly.

 The indications displayed on the monitor screen are output only via the MONITOR © connector.

Notes for LANC connection

- With the LANC connection, refer to the instruction for use supplied with
- the player VCR.

 The LANC connection transmits signals such as control signals, time

- *The LANC connection transmits signals such as control signals, time code and time counter data and status data.

 If the other VCR has a LANC ♥ jack of 5-pin DIN type, connect with the VK.910 Control L connecting cable (not supplied).

 The jack shelded CONTROL L has the same function as LANC ♥ jacks. The jacks labeled REMOTE on other equipment may also have the same.

 *This VCR (recorder) receives the time code data from the LANC ♥ jack only when the other VCR (player) is set to show the time code indications.

· With the LANC connection, this VCR only works as a slave unit,

Recording

26⁶⁸ Chapter 2 Recording and Pfs

6 When using the line connections (INPUT

connectors), select the audio mode.

Select the desired mode by setting the AUDIO MODE menu.

Audio mode	Set the menu to
2-channel mode	F848k
4-channel mode	Fs32k

On how to use the menu, see Chapter 3 "Menu Settings,"

- In the DVCAM format, there are two audio recording modes, with either two channels at 48 kHz or four channels at 32 kHz. It is not possible to select other modes (for example with four channels at
- possible to select other modes (for example with four channels at 48 kHz).

 When recording in 4-channel mode on this VCR, audio signals are recorded only in channels I/2.

 Once you have started recording, you cannot change the audio mode selection.
- 7 Use the AUDIO INPUT LEVEL control knobs to adjust audio

Watching the audio level meter (see page 7), adjust the level so that the meter does not indicate higher values than 0 dB when the audio signal is at its maximum. When the level exceeds 0 dB, sound distortion occurs.

With the DV connection, the recorder VCR's AUDIO MODE (sound selection) and AUDIO INPUT LEVEL (audio balance adjustment) do not function.

Preparation on the player

- Insert a source tape.
 If the player VCR has an EDIT switch, set it to ON,
 Turn off the on-screen display.
 Set the LANC mode to M.

With the DV connection, the playback VCR's AUDIO MONITOR (some selection) and AUDIO MIX BALANCE (audio balance adjustment) do a function on the source audio output through the DV L jack.

Recording Procedure

Note
When controlling this unit from an editing controller or a personal
computer connected to the RS-232C 22 connector, set the REMOTE
LOCAL switch to REMOTH. When not, set the switch to LOCAL (See page 9).

1 After checking that the cassette's anfety switch is set to write enabled position and the tape for slack, hold the cassette so that the tape window is facing upward, then insert it into this unit.

For details of the cassene's safety switch, see page 4. For details of checking the tape for slack, see page 5.

The cassette is automatically drawn into the unit and the tape is wound round the head dram. The tape is stationary while the head dram

2 Press the playback button on the player.

This starts the player's playback operation

3 Press and hold ● REC, and press ► PLAY.

This starts the recorder's recording operation.

To stop recording Press the STOP button.

To record using the optional Foot Switch

1 Press the pedal of the Foot Switch when the VCR is in stop mode.

The VCR starts recording.

2 Press the pedal again.

The recording stops and the VCR goes into recording pause mode.

To stop recording
Press the STOP button on the VCR.

- Vou should set the REMOTE/LOCAL switch to REMOTE to prevent concurrent use of the foot switch and the front panel controls.

 The beginning of the recording (for about two seconds) cannot be made. If you immediately start recording, press the pedal twice to go into recording pause mode, then start recording.

 The Foot Switch must be conformed with Standard UL/2601-1/EN60601-1.

Chapter 2 Recording and Playback 29 as

28^{dir} Chapter 2 Recording and Playback

• The recording pause mode will be automatically released after five minutes to protect the tape, and the VCR goes into stop mode.
• The foot switch operation works even if the VCR is in any operation mode. To provent accidental erasure of a recording, you should slide in the safety switch on the cassette so that the red portion becomes visible before you insert the cassette into this VCR.

Duplicate

If you copy a source taps, using the DUP (duplicate) button on this VCR, you can copy the time code recorded on the source tape as they are. You can easily make a work tape having the same time codes as the source

tispe,
The duplicate function on this VCR works only when using a source tape recorded in DVCAM format and making DV connections.

- 1 Connect this VCR and the other (playback) VCR, using an i.LINK cable (DV connecting cable) (not supplied) and select DV with the INPUT SELECT selector on this VCR.
- 2 Locate the points where you want to start playback and recording.
- 3 Press STOP on this VCR to atop the tape transport operation.
- 4 Press and hold DUP on this VCR, and press > PLAY.

The DUP indicator flashes and this VCR enters into duplicate-standby

- of the other (playback) VCR has already started playback, the DUP indicator lights and duplicate starts immediately.

 If the other (playback) VCR is in the playback pause mode, duplicate starts immediately and this VCR continues to record a still picture and a certain time code.
- 5 Press the play button on the other VCR to start playback.

The DUP indicator lights and duplicate starts.

To adjust the point where duplicate starts
In step 4 above, press and hold the DUP button instead of the ▶ FLAY
button, and press the \$1 PAUSE button. This VCR remains recording
standby mode until you press the \$1 PAUSE button again.

After the other VCR starts playback, press the \$1 PAUSE button at the
point where you want to start duplicate.

To stop duplicate Press the M STOP button.

30th Chapter 2 - Recording and Pisylo


oter 2 Recording and Playback 31 de

Recording

AC timer recording

By connecting this VCR to an external AC timer (not supplied), you can start recording at a preset time.

1 Connect this VCR to an external AC timer (not supplied).



- 2 Insert a tape for recording.
- 3 Press INPUT SELECT to select the recording source.
- 4 Set the timer-on time on the conn

At the preset time, the power of this VCR and the recording source turn on automatically and recording starts about soveral to 10 seconds later. Set the timer allowing a margin for the recording to start.

5 Set the TIMER selector at the front to REC.

You need not press ● REC.

If the tape ends before the recording source stops operation. The tape stops without rewinding. If you set AUTO REWIND in the menu to ON, the tape rewinds to its beginning automatically (see page 34).

To stop recording during the timer recording Press the **E** STOP button.

To release AC timer recording Set the TIMER selector to OFF.

While a tape is running, do not turn off the power using an AC timer. The YCR and a tape may be demaged. When turning off the power of the YCR, make sure to press the STOP button on this VCR first to stop the tape, then turn off the power.

Changing Menu Settings

This VCR has various functions available, and you can set and check them on the monitor acreen. Before operation, set the clock by setting the CLOCK SET menu.

You can change the menu settings on the SET UP MENU screen.

If necessary, change the settings manually during editing, etc.

Changing the SET UP MENU Settings

Follow the instructions below to change the settings.

Notes

*During duplicate, do not change the speed of the player's tape or set it to pause mode. Otherwise, the time code of the recorded tape becomes out of sequence and you cannot use it for editing.

*During duplicate, time counter does not appear. Check it in the other VCR.

VCR.
When you start duplicating, the first part of the source tape may be dropped on the copied tape. Play back the source tape from the preceding point. You cannot completely copy the tape if the source tape is recorded

from its beginning point.

- You may not be able to copy the first part or an unrecorded portion of the source tape. Locate the recorded portion on the source tape, then start

source tape. Locate the recorded portion on the source tape, then start copying.

The recording does not stop the moment you press the **W** STOP button to stop editing. The source picture may be recorded a little longer than you expected.

If you duplicate a tape by using two DSR-20MD/20MDPs, set DV EB OUT in the menu of the player to OFF (see page 34).

The index signals are not recorded when the duplicate starts.

If you set the REMOTE/LOCAL switch to REMOTE during duplicate, the tage starts.

The SET UP MENU appears on the monitor screen. To cancel the menu settings, press MENU again.



2 Press $^{4}/^{4}$ to select the option you want to change, and press SET.

Each menu option appears on the monitor screen (see the table below).

3 Press ऐ/₺ to change the setting, and press SBT.

The menu disappears from the monitor screen,

Menu Contents

Initial settings are indicated in bold letters.

Menu options	Set this option to	Description of settings
AUDIO MIX BALANCE		If you set the AUDIO MONITOR selector to MIX, you can select the precise balance between channels 1/2 and channels 3/4 by five steps
AUDIO MODE	Fa48k Fa32k	 To set the sudio mode to 2-channel mode (19bit mode). This mode uses the whole sudio teste to record one stereo track. You can get higher sound quality. To set the sudio mode to 4-channel mode (12bit mode). This mode separates the sudio area to to two parts. You can record two kinds of audio, stereo 1 and stereo 2. When recording on this VCR, sudio signals are recorded only in glannels 1f2.
JOG WITH SOUND	OFF	To listen to the sound when playing a tape in various speads. To turn off the sound when playing a tape in various speads.

32⁶⁸ Chapter 2 Recording and Playb

330

Set this option to Description of settings Menu options REMOTE CONTROL Description of settings Set the command mode (VTR1 to 6, INST) on this YCR. Change this setting when using infrared remote commander or external (remote) equipment to remotely control the unit. When using the supplied remote commander, setel YTR4 (finite setting). When using the remote controller such as the optional DSRM-10 or SVRM-100A, setect INST, When selecting OFF, you cannot remotely control the unit. -To sat the baud rate with an editing controller that supports RS-232C interactor to 9800bps. -To set the baud rate to 19200bps. BS232C BAUD BATE 19200bps To display the tape counter in the center of the monitor screen. To display the tape counter in the lower right of the monitor screen. To display the tape counter in the lower right of the monitor screen. CENTER LOWER R To display the alarm message on the monitor screen. Not to display the alarm message. To output a base sound when an illegical operation is made. To deadivate it. CAUTION DISPLAY ON OFF To depot cluste alignals when recording begins. Not to record index alignals. To espect recordings with the cassette memory. If the tape does not have a cassette memory, but offer the signals recorded on the tape fisc INDEY WOITE CASSETTE MEMORY SEARCH ALL DATA INDEX DATA DATE DATA PHOTO DATA CASSETTE MEMORY EDASE NOID When using the cassette whose mamory can store over 16 kbits of data, you can only select ALL DATA. You cannot erase index data on the tape To set the time ends to the same one as skeady recorded on the TIME CODE (DSR-20MD only) AUTO · To set the time code to the same one as already tape. • To set the time code to Non Drop Frame. • To set the time code to Drop Frame. NOF To set the line code to Drop Frame. To set the line code to Drop Frame. It you use AUTO and start recording at the beginning of the lape, the line code is set to Non Drop Frame. If the VCR automatically if there is no operation for an hour To turn off the VCR automaticelly if there is no operation for an hour during stop mode (Auto Off). To deactivate Auto Off. AUTO OFF ON OFF to usewheat AUIO-UII. To rewind the tape to lie beginning automatically if the tape reaches a med (Auto Rewind). To deadlivate Auto Rewind. To prevent the picture from blurring when playing a tape recorded in phole mode. AUTO REWIND РНОТО РВ photo mode. photo mode clear plotter when playing a still plotters. Date clear plotter when playing a still plotters. When using FRAME, the plotter recorded in photo mode may blur. Set the cleck on this VCR so that the recording time can be written in the ridde signal. Using 9/9 and SET buttons, set the date and time. The still hours moler keeps curvalative coults of the head dram. The still hours moler keeps curvalative coults of the head dram. FRAME CLOCK SET HOURS METER The digital hours meter keeps cumulative counts of the head drum rotation lime and the number of unthreading operations. These counts can be displayed on the monitor screen and are unresettable. The cumulative total hours of dum rotation with tape threaded is displayed in 10-hour increments. The cumulative number of tape unthreading operation is displayed in 10-operation increments. DRUM ROTATI THREADING DV EE OUT

To output the selected line input signals from the DV & jack.
 To output only playback video and audio signals from the DV & jack

34 Chapter 3 Monu Settings

Alarm Messages

Various messages appear on the monitor screen ("Err" appears in the display window). Check them with the

Message	Meaning / Remedy
PLEASE CONFIRM THE SAFETY SWITCH OF THE CASSETTE	Check that the protect tab is slid in so that the red portion visible -> Slide back the safety switch (see page 4).
NO CASSETTE MEMORY	You try to erase cassette memory when there is no cassette memory.
VCR IS RECORDING	You press a certain operation button during recording or editing.
PLEASE INSERT A NEW CASSETTE	Though no cassette is inserted in the cassette compartment, you press ➤ PLAY, etc. → Insert a cassette.
THE TAPE IS REWOUND	You press ◄◄ REW at the beginning of the tape.
PLEASE REWIND OR INSERT A NEW CASSETTE	You try to start playback or recording at the tape end. → Rewind the tape or insert a new cassette.
PLEASE SET THE CLOCK	When turning on the power, the clock has not been set, Set the clock in the menu (see page 34).
THIS PROGRAM IS COPYRIGHT PROTECTED	You try to dub the tape on which copyright protect signals are recorded.
CASSETTE MEMORY IS TOO LARGE TO ERASE	You try to erase "INDEX DATA," "DATE DATA," or "PHOTO DATA" on a tape having more than 18 Kbits memory capacity. → Erase "ALL DATA" on the tape (see page 34).
WRITING ON CASSETTE MEMORY, PLEASE WAIT	You do certain operation white the VCR is writing on casestle memory. — Operate after writing on cassatte memory is complete.
VCR IS IN DUP MODE	You press a certain operation button during duplicate.

Chapter 4 Meintenance and Troybleshootho 356

Troubleshooting

If the VCR does not function or functions incorrectly, check the following.

ON

Symptom	Cause / Remedy
The power cannot be turned on.	The power plug is disconnected Connect the plug.
The unit will not operate even if the power has been turned on.	 The REMOTE/LOCAL switch is set to REMOTE. → Set it to LOCAL (See page θ).
	Moisture condensation occurs. Turn of the power and disconnect the power plup. After about one minute, connect the plug and turn on the power. Wait for about one hour with the power turned on.
	 The cassette is not inserted straight. → insert it straight.
The unit cannot be controlled using buttons on the unit,	The REMOTE/LOCAL switch is set to REMOTE. → Set it to LOCAL (See page 9).
The cassette cannot be ejected.	The REMOTE/LOCAL switch is set to REMOTE. → Set it to LOCAL (See page 9).
The cassette cannot be inserted, or it is ejected promptly.	There is moisture condensation on the head drum, → Walt for about an hour. The cassette is not inserted straight. → insert it straight.
No picture,	The video heads are dirty, Clear the video heads using the cleaning cassette.
Noise appears on the screen.	A damaged cassette is inserted. → Insert other cassette. The video heads are dirty. → Clear the video heads using cleaning cassette.
No picture via the DV jack.	Reconnect an i.LINK cable (DV connecting cable) (not supplied).
The audio is noisy.	A damaged cassette is inserted, Insert other cassette.
The playback automatically starts when the power is turned on.	The TIMER selector is set to REPEAT Set it to OFF (See page 9).
The recording automatically starts when the power is turned on.	The TIMER selector is set to REC. → Set it to OFF (See page 9).
The remote commander does not function,	The batteries are dead. → Replace the batteries, Something is blocking the interest rays. → Remove the obstacle. The command mode is wrong. → Set up REMOTE CONTROL in the SET UP MENU (See page 34).
The menu does not appear.	Connect the video monitor to the MONITOR G- connector.

Notes on Use

Notes on the video cassette recorder

Do not install the unit in a place subject to direct smallght or heat sources If you do, its cabinet, mechanical parts, etc., may be

Do not install the unit in an extremely bot place If the unit is left in a car parked with its windows closed (especially in summer), its cabinet may be damaged or it may not work correctly.

If the unit is brought directly from a cold to a warm

location

Moisture may condense inside the unit and cause
damage to the video head and tape. If you use the unit
in a place subject to direct cold currents from an air
conditioner, moisture may also condense inside the

Do not place a heavy objects on the unit The cabinet may be damaged, or the VCR may not work correctly.

Do not handle the recorder roughly Avoid rough handling or mechanical shock.

To avoid damaging the cabluet finish
Plastic is often used for the surface finishing of the
recorder. Do not spray a volatile solvent such as an
insecticide toward the cablnet or place rubber or vinyl products on the cabinet for a long time, 11 you as, ... finish of the cabinet may be damaged or the coating ets on the cabinet for a long time. If you do, the

Do not clean the cabinet with thinner or benzing The cabinet may be damaged or its coating may come off. When you use a chemical-impregnated cloth, use it according to its directions.

Clean the cabinet with soft dry cloth
When the cabinet is very dirty, clean it with a soft dry
cloth lightly moistened with a mild detergent solution
and finish it with dry cloth,

Do not put magnetic objects close to the unit Magnetic fields may damage the recording.

Checking the video heads every 1000 hours A VCR is a high-precision piece of equipment that records and plays back the picture on a magnetic tape. In particular, the video head and other mechanical in particular, the video head and other mechanics parts become dirty or worn. To maintain a clean picture, we recommend maintenance every 1000 hours, though the using condition may differ depending on temperature, humidity, dust, etc.

Cleaning of the video heads

If the video heads are contaminated, the pictores cannot be recorded properly or the playback pictures become noisy. If the following phenomena occur, use the cleaning caseate PDVM-12CL (supplied) or PDV-12CL (not supplied) to clean the heads.

ed noise appears on the playback picture.

 A part of the playback picture does not move.
 The playback picture does not appear on the screen. cloms caused by contaminated video heads







To use the cleaning cassette Refer to your cleaning cassette's operating instructions.

After prolonged use, the video heads may become Arter protonged use, the video heads may become worn out. If Optimum picture quality is not restored even after you have cleaned the video heads with the cleaning cassette, the video heads may have worn out. In that case, you have to replace the video heads with new ones. Picase consult your Sony dealer.

Notes on the video cassettes

Cleaning the terminal
If the terminal of the Standard-DVCAM or MiniDVCAM cassette gets dirty, or dust sticks to the
terminal, the VCR may not work correctly.
Clean the terminal with the swab once every ten times you eject a cassette.



When affixing a label on the cassette
Be sure to affix a label on only the correct location so
as not to cause malfunction of the VCR.

After using a cassette
After use, please be sure to rewind the tape completely
(to prevent picture and sound distortion). Return it to
its case and store in upright position.

About moisture condensation

If the unit or tape is brought directly from a cold to a warm location, moisture may condense inside or outside the unit or tape. If you use the tape or video heads in this condition, the tape may adhere to the head drum, and the video heads or the tape may be damaged, or malfunction may occur.

Moisture condensation is likely to occur under the

- following conditions:
 The unit is brought from the cold outdoors to a warm indoor location.

 The unit is brought from the air-conditioned indoors

• The unit is brought from the air-conditioned indoors to the hot outdoors.
• The unit is used in a place subject to cold currents from an air conditioner.
When bringing the unit from a cold place to a warm place or vice versa, put it in a plastic bag and seal the bag tightly. After bringing it into the new place, leave the bag on for about an hour, and remove the bag when the air temperature inside it has reached the temperature surrounding it.

If moisture condensation occurred

At mossure consensation occurred
You cannot operate the unit except to press
EIECT.
If you insert a cassente, it is ejected automatically. If
this occurs, turn on the power, wait about an hour for
the moisture to evaporate.

Digital hours meter

The digital hours moter keeps cumulative counts of the head drum rotation time and the number of unthreading operations. These counts can be displayed on the monitor screen. Use them as guidelines for scheduling maintenance.

In general, consult your Sony dealer about necessary periodic maintenance checks.

The digital hours meter has the following two display modes and you can check them in the HOURS METER menu (see page 34).

- DRUM ROTATION mode
 The cumulative total hours of drum rotation with tape threaded is displayed in 10-hour increments.

 THREADING mode
 The cumulative number of tape unitreading operation is displayed in 10-operation increments.

Self-diagnosis function

The unit is equipped with the self-diagnosis function that works to prevent the VCR from malfunctioning. A two-digit service number appears in the display window. In this case, check the following table.

Message	Symptom	Remedy
22	The video heads are dirty.	Clear the heads. (See page 38)
32	To prevent the unit from malfunctioning, the self-diagnosis function has worked.	Disconnect the power cord. After reinstalling the power source, operate the unit, Remove the cassette or turn on/off the unit,
21	Moisture condensation has occurred.	Remove the cassette and leave the unit for al least one hour.


If you are mable to resolve the problem, contact your Sony dealer or local authorized Sony service facility and inform them of the number.

SECTION 2 DISASSEMBLY


Note: Follow the disassembly procedure in the numerical order given.

2-1. REMOVAL OF UPPER CASE


2-3. REMOVAL OF VA-102 BOARD




2-2. REMOVAL OF FRONT PANEL ASSEMBLY



2-4. REMOVAL OF POWER BLOCK




2-5. REMOVAL OF MD BLOCK ASSEMBLY 2-7. REMOVAL OF JC-19 BOARD




2-6. REMOVAL OF CM-56 BOARD

(CN003) (a) Two claws (b) Screw (CN005, 007, 008) (c) Two flat cables (c) Two flat c

2-8. REMOVAL OF RP-228 BOARD




2-9. CIRCUIT BOARDS LOCATION







SECTION 3 BLOCK DIAGRAMS


3-1. OVERALL BLOCK DIAGRAM 1







3-2. OVERALL BLOCK DIAGRAM 2





3-3. OVERALL BLOCK DIAGRAM 3







SECTION 4
PRINTED WIRING BOARDS AND SCHEMATIC DIAGRAMS

4-1

SP1268 / Druck 7



4-2



FRAME

DSR-20MD/20MDP

4-2. PRINTED WIRING BOARDS AND SCHEMATIC DIAGRAMS

RP-228 BOARD (SIDE A)

A-4 B-6 A-5

A-1 B-5 C-5

C-1 D-3 B-1 B-3 C-5

C-2 B-3 A-1 A-1 B-2 B-2 A-3 B-3

CN101 A-2

CN102 CN771

CN775

D772 D774

IC770 IC771 IC772

IC775 IC777


Q105

Q109 Q774 Q775 Q776 Q777

THIS NOTE IS COMMON FOR PRINTED WIRING **BOARDS AND SCHEMATIC DIAGRAMS.** (In addition to this, the necessary note is printed in each block)


For printed wiring boards:

- : Pattern from the side which enables seeing. (The other layers' pattern are not indicated)
- Circled numbers refer to waveforms.
- Through hole is omitted.
- · There are few cases that the part printed on diagram isn't mounted in this model.
- Chip parts.



For schematic Diagram:

- All capacitors are in μF unless otherwise noted. pF: μμF 50V or less are not indicated except for electrolytics and
- Chip resistors are $^{1}/_{10}W$ unless otherwise noted. $k\Omega$: $1000\Omega,~M\Omega$: $1000k\Omega.$
- Caution when replacing chip parts.
- New parts must be attached after removal of chip. Be careful not to heat the minus side of tantalum capacitor,
- because it is damaged by the heat.
- Some chip part will be indicated as follows.



- Constants of resistors, capasitors, ICs and etc with XX indicate that they are not used. In such cases, the unused circuits may be indicated.
- Parts with ★ differ according to the model/destination. Refer to the mount table for each function.
- All variable and adjustable resistors have characteristic curve B, unless otherwise noted.
- Signal name

XEDIT → EDIT PB/XREC → PB/REC : nonflammable resistor.

- : fusible resistor.
- : panel designation. : B+ Line.*
- --- : B- Line.* : IN/OUT direction of B line (+, -).*
- : adjustment for repair.*
- Circled numbers refer to waveforms.*

Measuring conditions voltege and waveform:

- · Voltages and waveforms are measured between the measurement points and graound when color bar signal input. They are reference values and reference waveforms.* (VOM of DC 10 M Ω input impedance is used)
- Voltage values change depending upon input impedance of VOM used.
- Indicated by the color red.

Note:	Note:
The components identified by mark △ or dotted line with mark △ are critical for safety. Replace only with part number specified.	Les composants identifiés par une marque ∆ sont critiques pour la sécurité. Ne les remplacer que par une piéce portant le numéro spécifié.

number, please include the board

RP-228 (REC/PB AMP) PRINTED WIRING BOARD

- Ref. No.: RP-228 board; 3,000 series -

• For Printed Wiring


• RP-228 board is si

of layers 2 to 5 hav • There are few case

is printed on this d • Chip transistor



RP-228 BOARD (SIDE A)




• For Printed Wiring Board.

- RP-228 board is six-layer print board. However, the patterns
- of layers 2 to 5 have not been included in the diagram.

 There are few cases that the part isn't mounted in this model is printed on this diagram.
- Chip transistor






RP-228 BOARD


CN103 D-5


IC773 IC774

Q701 Q702 Q772 Q773 Q778


RP-228 BOARD (SIDE B)








4-9



The same accounts to the same account of the s



RP-228 BOARD 3.4 Vp-p 0.04 μsec IC772 ② – ⑧ PB 3.3 msec IC777 3 REC 2 41.8 MHz 6.7 msec IC777 ②, ③, ⑥, ⑦ REC IC772 12 PB 8 1 1.6 Vp-p 3.3 msec IC775 @ PB IC774 (9) REC/PB 4 • 0.5 Vp-p 3.6 Vp-p 8.38 MHz 3.3 msec IC774 ⑤ REC/PB IC775 @ PB IC774 @ PB 6 1 3.4 Vp-p 0.6 Vp-p 0.25 μsec 6.7 msec IC777 (9 PB IC774 ① - ④, ᠖ REC 6 B 3.4 Vp-p 3.0 Vp-p 0.25 μsec 6.7 msec IC777 2 REC/PB IC774 59 - 63 PB 0 6.0 Vp-p 3.3 msec IC777 @ REC

4-13


(TAPE DETECT)

JC-19 (AD/DA CONVERTER, S1 AFC, U1, D1, C1 SPCON, MODE, DV IN/OUT, AUDIO CORE, DIGITAL AUDIO, AUDIO D/A, A/D CONVERTER) PRINTED WIRING BOARD JC-19 BOARD (SIDE A) - Ref. No.: JC-19 board; 2,000 series -A-6 A-3 A-9 A-10 D-5 D-10 CN103 CN103 CN104 CN411 CN412 CN501 CN831 C-2 C-2 C-1 B-4 C-11 E-2 E-1 D-4 D-4 B-6 B-10 B-8 F-9 E-10 E-7 F-6 E-4 D-8 E-9 D-9 D-11 IC013 IC014 IC018 IC019 IC200 IC209 IC210 IC211 IC212 IC213 IC214 IC410 IC410 IC410 IC701 IC702 IC801 IC805 IC807 IC807 IC803 IC803 IC803 JC-19 BOARD (SIDE A) Q039 Q040 Q041 Q042 Q043 Q044 Q045 Q050 Q051 Q052 Q053 Q200 Q201 Q801 B-3 B-3 B-3 B-2 B-2 A-3 B-4 A-2 B-1 B-1 C-3 D-7 (31)/ IF, VIDEO IN/OUT, UVIC, (RS-232C IF) DV IN/OUT, MONITOR OUT, AUDIO, HI MICOM, RS MICOM (POWER 2) RP-228 (REC/PB AMP) RE-32 MD-65 (V/A IN/OUT) (TAPE DETECT) MD-63 (TAPE DETECT) OVER DISCHARGE PROTECT HP-100 AD/DA CONVERTER, S1 AFC, U1, HEADPHONE D1, C1 SPCON, MODE, DV IN/OUT, & REC VR AUDIO CORE, DIGITAL AUDIO, AUDIO D/A, A/D CONVERTER FR-136 (DISPLAY CONTROL) (SERVO, SYSTEM CONTROL)


AD/DA CONVERTER, S1 AFC, U1, D1, C1 SPCON, MODE, DV IN/OUT, AUDIO CORE, DIGITAL AUDIO, AUDIO D/A, A/D CONVERTER
SP1268 / Druck 14


- For Printed Wiring Board.
 JC-19 board is six-layer print board. However, the patterns of layers 2 to 5 have not been included in the diagram.
 There are few cases that the part isn't mounted in this model is printed on this diagram.
 Chip transistor







JC-19 BOARD (SIDE B)






JC-19 BOARD (SIDE B)







10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 | 3 2-12 IC011 R139 1k C168 CXD23000-T4 Y CLAMP REF REG IC013 C150 0.01u IC018 R113 ≱ R138 10k 0048 UN5213-TX IC018 TC74HC4066AFT(EL) --€∑ 3-10 ≪ 3-14 ≪ 3-14 ≪ 3-15 TO JC-19 BOARD(2/10) (SEE PAGE 4-25) IC010 C123 IC016 CXD2300Q-T4 IC014 R149 1800 C149 0.01u Q043 2SA1576A-T106-R 2 9 R012 0 C141 C155 10u 10u 63V 6.3V C151 ± T C160 10u T T 0 01u IC009 IC012 C133 1 CB CLAMP REF REG IC015 IC015 22) TO JC-19 BOARD(5/10) (SEE PAGE 4-37,39) C148 C163 0.01u 0.01u SIGNAL PATH VIDEO SIGNAL CHROMA Y Y/CHROMA REC \Rightarrow


TO VA-102 BOARD(1/8 CN101 (SEE PAGE 4-69)








JC-19 BOARD (2/10) 3.8 Vp-p 13.5 MHz IC200 99, 99 REC/PB IC205 ⑦ REC Ø 0.15 μsec IC200 (11)- (18) REC/PB IC205 ® REC/PB 1 B 3.0 Vp-p IC205 18, 29 REC/PB IC200 @- @ REC/PB • Н IC205 28 REC/PB IC200 [®] PB IC205 49 REC/PB IC200 (6) PB 1 IC200 (6) PB























4-37








JC-19 BOARD (6/10)


20 MHz
IC501 @ REC/PB


32.768 KHz


IC501 @ REC/PB



4-45







3 MHz
IC804 ® REC/PB

3 MHz
IC804 ® REC/PB

3.6 Vp-p

12.3 MHz
IC805 ® REC/PB

18.16 MHz
IC805 ® REC/PB


1.1 msec
IC805 ® REC


3.1 Vp-p


1.1 msec
IC805 ® PB


JC-19 BOARD (8/10)

IC801 @ REC/PB










AUDIO D/A, A/D CONVERTER
SP1268 / Druck 37 JC-19 (10/10)



- For Printed Wiring Board.
 There are few cases that the part isn't mounted in this model is printed on this diagram.

RE-32 (V/A IN/OUT) PRINTED WIRING BOARD

- Ref. No.: RE-32 board; 7,000 series -


RE-32 BOARD RE-32 BOARD Α V/A IN/OUT -REF.NO.:7,000 SERIES-XX MARK:NO MOUNT В D C CH-2 EXIT Ε RS-78 (RS-232C IF) (IF, VIDEO IN/OUT, UVIC, DV IN/OUT, MONITOR OUT, AUDIO, HI MICOM, RS MICOM) U-2 (POWER 2) RP-228 (REC/PB AMP) G RE-32 MD-65 (V/A IN/OUT) (TAPE DETECT) MD-63 (TAPE DETECT) 69-385-12 Н D101 MA3075WA-(TX) OVER DISCHARGE PROTECT U-1 (POWER 1) / AD/DA CONVERTER, S1 AFC, U1, \ D1, C1 SPCON, MODE, DV IN/OUT, (HEADPHONE & REC VR) AUDIO CORE, DIGITAL AUDIO, AUDIO D/A, A/D CONVERTER

(SERVO, SYSTEM CONTROL)

(TAPE DETECT)

(DISPLAY CONTROL)

SCHEMATIC DIAGRAM



VA-102 RS-78 / IF, VIDEO IN/OUT, UVIC, (RS-232C IF) DV IN/OUT, MONITOR OUT, U-2 AUDIO, HI MICOM, RS MICOM (POWER 2) RP-228 VA-102 (IF, VIDEO IN/OUT, UVIC, DV IN/OUT, MONITOR OUT, • For Printed Wiri (REC/PB AMP) AUDIO, HI MICOM, RS MICOM) PRINTED WIRING BOARD • VA-102 board is RE-32 of layers 2 to 5 ha - Ref. No.: VA-102 board; 1,000 series -(V/A IN/OUT) MD-65 There are few ca (TAPE DETECT) is printed on this MD-63 Chip transistor (TAPE DETECT) VA-102 BOARD (SIDE A) CN051 CN101 A-1 E-5 E-6 B-5 A-2 D-8 E-2 B-7 E-3 OVER DISCHARGE CN101 CN102 CN401 CN602 CN701 CN702 LANC PROTECT VA-102 BOARD(SIDE A) U-1 (POWER 1) AD/DA CONVERTER, S1 AFC, U1, / HEADPHONE D1, C1 SPCON, MODE, DV IN/OUT, & REC VR E-4 C-5 C-5 A-4 A-4 B-4 D-1 C-1 D401 AUDIO CORE, DIGITAL AUDIO, D403 D404 D851 D852 D853 AUDIO D/A, A/D CONVERTER FR-136 (DISPLAY CONTROL) MD-64 (SERVO, SYSTEM CONTROL) (TAPE DETECT) D866 D867 000000 IC051 B-2 B-3 E-2 E-2 E-2 E-2 E-3 D-3 B-6 C-6 C-7 C-7 B-7 E-1 B-2 B-4 E-8 B-1 IC052 IC101 IC102 IC103 IC104 IC105 IC106 IC202 IC401 IC402 IC652 IC701 IC702 IC703 IC706 IC715 IC715 IC715 В IC852 IC853 IC856 IC862 IC863 IC865 D-6 C-7 C-7 E-6 C-6 E-6 C-5 E-6 C-5 E-6 C-5 E-6 C-5 E-6 C-5 E-6 E-7 D-1 A-1 Q214 Q215 Q217 Q218 Q221 Q222 Q223 Q224 Q25 Q231 Q32 Q404 Q407 Q410 Q411 Q412 Q415 Q417 Q420 Q852 Q855 Q855 30 D 05 2 3 5

T,

T)

• For Printed Wiring Board.


- VA-102 board is six-layer print board. However, the patterns of layers 2 to 5 have not been included in the diagram.
- There are few cases that the part isn't mounted in this model is printed on this diagram.

• Chip transistor




HARGE)


U1, OUT, O,





VA-102 BOARD (SIDE B)


	30, 11 12 (0		
CN601 CN852	A-6 E-1	Q103 Q107 Q108 Q201 Q202 Q203 Q204 Q205 Q206 Q207 Q208 Q209 Q210 Q211 Q212 Q213 Q216 Q226 Q227 Q228 Q233 Q403 Q406 Q416 Q419 Q604 Q604 Q604 Q604 Q605 Q607 Q605 Q607 Q605 Q607 Q607 Q608 Q709 Q709 Q709 Q709 Q710 Q711 Q712 Q713 Q853	
D051 D052 D201 D202 D855 D856 D858 D859 D861 D862 D863 D864	B-6 B-6 D-2 D-3 C-8 C-8 E-1 E-1 A-6 E-1 C-8		
ICO53 IC201 IC203 IC204 IC205 IC206 IC207 IC208 IC210 IC213 IC404 IC405 IC406 IC407 IC602 IC651 IC652 IC704 IC705 IC708 IC708 IC709 IC712 IC713 IC712 IC713 IC714 IC717 IC718 IC851 IC856 IC865 IC866 IC867 IC866	B-7 B-2 B-2 C-2 D-3 D-4 E-5 B-4 D-7 C-5 C-1 D-1 D-1 D-1 C-1 C-1 C-1 B-1 B-5 C-8 B-2 C-1 B-1 C-1 C-1 C-1 C-1 C-1 C-1 C-1 C-1 C-1 C		
Q101 Q102	E-3 E-4	Q854	C-















TEMPHOS BUSH S














VIDEO OUT VA-102 (3/8)




































RS-78 (RS-232C IF) PRINTED WIRING BOARD


- Ref. No.: RS-78 board; 6,000 series -


- For Printed Wiring Board.
 There are few cases that the part isn't mounted in this model is printed on this diagram.


RS-78 BOARD(SIDE A)





RS-78 BOARD(SIDE B)













MD-63, MD-64, MD-65 (TAPE DETECT), FP-406 (TAPE SENSOR) PRINTED WIRING BOARDS AND SCHEMATIC DIAGRAM


- Ref. No.: MD-63 board; 6,000/MD-64 board; 7,000/MD-65 board; 5,000/FP-406 boad; 5,000 series -





- For Printed Wiring Board.
- There are few cases that the part isn't mounted in this modis printed on this diagram.




- For Printed Wiring Board.
- There are few cases that the part isn't mounted in this model is printed on this diagram.











SIGNAL PATH

SIGNAL PATH			
	REC	REC/PB	PB
Drum speed servo			
Drum phase servo			
Drum servo (speed and phase)			
Capstan speed servo			
Capstan phase servo			
Capstan servo (speed and phase)			
Ref. signal			


HP-100 (HEADPHONE & REC VR) PRINTED WIRING BOARD


- Ref. No.: HP-100 board; 5,000 series -





- For Printed Wiring Board.
- HP-100 board is six-layer print board. However, the patterns of layers 2 to 5 have not been included in the diagram.
- There are few cases that the part isn't mounted in this model is printed on this diagram.






SP1268 / Druck 50

DSR-20MD/20MDP

CM-56 (SERVO, SYSTEM CONTROL) PRINTED WIRING BOARD

- Ref. No.: CM-56 board; 4,000 series -



• For Printed Wiring Board.

• There are few cases that the part isn't mounted in this model is printed on this diagram.

Chip transistor



RS-78 (RS-232C IF) (VA-102 IF, VIDEO IN/OUT, UVICED IN/OUT, MONITOR AUDIO, HI MICOM, RS U-2 (POWER 2) RP-228 (REC/PB AMP) RE-32 MD-65 (TAPE DETECT) MD-63 (TAPE DETECT) OVER DI PROTEC U-1 (POWER 1) HP-100 AD/DA CONVERTER, S1 AI D1, C1 SPCON, MODE, DV AUDIO CORE, DIGITAL AU (HEADPHONE & REC VR AUDIO D/A, A/D CONVERT (DISPLAY CONTROL) (SERVO, SYSTEM CONTROL) (TAPE DETECT)

D004 D011 D501 D502 E-2 C-6 D-2 B-5 IC003 IC005 IC008 IC009 IC014 IC017 IC021


F-8 A-7 A-9 A-2 A-6

C-7 B-9 C-4 B-7 F-6 F-7 A-6

Q002 Q012 Q500 Q502 Q503 F-8 E-2 E-5 D-2 E-9

SERVO, SYSTEM CONTROL SP1268 / Druck 51

t mounted in this model is



RS-78
(RS-232C IF)
(IF, VIDEO IN/OUT, UVIC, DV IN/OUT, MONITOR OUT, AUDIO, HI MICOM, RS MICOM)

RE-32
(V/A IN/OUT)

RE-32
(V/A IN/OUT)

HG-1
(OVER DISCHARGE PROTECT)

U-1 (POWER 1)

JC-19

AD/DA CONVERTER, S1 AFC, U1, D1, C1 SPCON, MODE, DV IN/OUT, AUDIO CORE, DIGITAL AUDIO, AUDIO D/A, A/D CONVERTER

CM-56
(SERVO, SYSTEM CONTROL)


DETECT)


4-105


4-107

SP1268 / Druck 53

4-108











FR-136 (DISPLAY CONTROL) PRINTED WIRING BOARD

- Ref. No.: FR-136 board; 5,000 series -










FR-136 (DISPLAY CONTROL) SCHEMATIC DIAGRAM




HG-1 (OVER DISCHARGE PROTECT) PRINTED WIRING BOARD AND SCHEMATIC DIAGRAM


- Ref. No.: HG-1 board; 8,000 series -


For Printed Wiring Board.

- : Pattern from the side which enables seeing.
- : Pattern on the rear side.
- There are few cases that the part isn't mounted in this model is printed on this diagram.

HG-1 BOARD








U-1 (POWER 1) PRINTED WIRING BOARD


- Ref. No.: U-1 board; 10,000 series -


- For Printed Wiring Board.
- :: Pattern from the side which enables seeing.
- : Pattern on the rear side.
- There are few cases that the part isn't mounted in this model is printed on this diagram.

U-1 BOARD		
CN1	A-4	
CN3	A-1	
D1	B-4	
D4	B-2	
D5	B-3	
D6	B-2	
D8	B-2	
D9	B-3	
D10	B-1	
D11	A-2	
D12	A-2	
Q1	B-3	
Q3	A-2	
Z1	A-3	
Z2	A-2	
Z3	A-2	



DSR-20MD: 1-468-441-DSR-20MDP: 1-468-442-





The components identified by mark \triangle or dotted line with mark \triangle are critical for safety. Replace only with part number specified.


Les composants identifiés par une marque \triangle sont critiques pour la sécurité. Ne les remplacer que par une piéce portant le numéro spécifié.

DSR-20MD/20MDP


U-2 (POWER 2) PRINTED WIRING BOARD

- Ref. No.: U-2 board; 20,000 series -

U-2 BOARD				
CN1	B-1			
CN2	B-1			
CN10	A-5			
CN11	A-4			
CN12	A-1			
CN13	B-1			
D1	B-2			
D2	B-3			
D3	B-3			
D4	A-2			
D5	B-4			
D6	B-5			
D7	A-2			
D8	A-3			
Q1	B-1			
Q2	A-2			
Q3	A-4			
Q4	A-5			
Z1	B-2			
Z2	A-1			
Z3	B-3			
Z4	A-3			
Z5	B-4			
Z6	A-2			
Z7	B-4			
Z8	B-5			
Z10	A-5			
Z11	A-3			




- For Printed Wiring Board.
 : Pattern from the side which enables seeing.
- : Pattern on the rear side.
- There are few cases that the part isn't mounted in this model is printed on this diagram.



DSR-20MDP: 1-468-442-

U-2 (POWER 2) SCHEMATIC DIAGRAM




SECTION 5 **ADJUSTMENTS**

5-1. MECHANICAL SECTION ADJUSTMENTS


5-1-1. INFORMATION

1-1. HOW TO SEARCH REFERENCE PAGES FOR REMOVAL




1-2. PHASE ADJUSTMENT MARK "PH-"

Numerous phase adjustments must be performed for removing and attaching parts (replacing parts) of the E mechanism. When removing and attaching parts, be sure to check the phase adjustment of corresponding parts. Parts that need phase adjustment are indicated with PH- mark. When replacing parts indicated with PH- mark, check their positions and phases so that the parts are attached smoothly in later.




In case of the above figure, refer to (B) and (E) of "5-1-3. PHASE ADJUSTMENTS"




5-1-2. PREPARATION FOR MECHANICAL CHECK, ADJUSTMENT AND MAINTENANCE

2-1. FL BLOCK ASSEMBLY


· Removing/Attaching




2-2. CASSETTE POSITIONS



2-3. LOADING/UNLOADING



2-4. MANUAL UP/DOWN THE FL BLOCK



2-5. SERVICE JIGS LIST

Ref. No.	Name	Part No.	Fixtur No.	Usage, Others Application, etc
J-1	Cleaning fluid	Y-2031-001-0		For cleaning drum assembly and tape guide
J-2	Wiping cloth	7-741-900-53	· · ·	For cleaning drum assembly
J-3	Super fine applicator (Made by NIPPON APPLICATOR (P752D))		·	For cleaning tape guide
J-4	Mirror (Small oval type)	J-6080-840-A	GD-2038	Tape path
J-5	Tracking tape (XH2-1AST) Standard cassette Tracking tape (XH2-1ASE) Standard cassette Tracking tape (XH2-1A1) Mini cassette	8-967-999-01 8-967-999-06 8-967-999-03	0.0 2000	Tape path (for tape top checking) Tape path (for tape end checking)
J-6	Mini DV torque cassette	J-6082-360-A		Tape path (for checking) For adjusting FWD/RVS back tension
J-7	Cassette standard plate (D/E mechanism)	J-6082-330-A		
J-8	Reel standard plate (D/E mechanism)	J-6082-331-A		For adjusting tape guide and reel table For adjusting reel table
J-9	TG2/7 preset plate	J-6082-358-A		For adjusting tape guide
J-10	Screwdriver for tape path	J-6082-026-A		For adjusting tape guide
J-11	Adjusting remote commander (RM-95 remodeled partly) Note 1	J-6082-053-B		Tape path
J-12	Torque driver	J-9049-330-A		Machanian de la
J-13	Tension regulator adjustment board	J-6082-359-A		Mechanism check and replacement
J-14	CPC 8-jig	J-6082-388-A		Electric tension regulator adjustment Tape path


Other equipment used


- Oscilloscope
- DC power supply
- · Digital voltage meter

Note 1: If the micro processor IC in the adjusting remote commander is not the new micro processor (UPD7503G-C56-12), the pages cannot be switched.

In this case, replace with the new micro processor (8-759-148-35).

- 1. Make a checking and adjustment at the tape top using the XH1-1AST tape.
- 2. Then, make a checking with the XH2-1ASE (for tape end) and XH2-1A1 (Mini cassette for tape top and end).
- 3. Again make a checking with the XH2-1AST.






5-1-3. PHASE ADJUSTMENTS


• This section classifies the phase adjustment into three blocks for clarity. The attaching order of each part is not described here. For details of the attaching order, refer to "5-1-5. MECHANISM SECTION CHECKS AND REPLACEMENTS".

3-1. PHASE ADJUSTMENT (Loading/Unloading Driving Section)


Note 1: Adjust it at the **UNLOADING** position unless otherwise specified. Note 2: (a) to (a) shown below are the orders for the phase adjustment.




PHASE ADJUSTMENT (A)




PHASE ADJUSTMENT ®




PHASE ADJUSTMENT ©




PHASE ADJUSTMENT ®




PHASE ADJUSTMENT (E)



PHASE ADJUSTMENT (F)




PHASE ADJUSTMENT @




3-2. PHASE ADJUSTMENT (S/L Cassette Selection Section)

Note 1: Adjust if at the S/L cassette position unless otherwise specified.


Note 2: (H) to (M) shown below are the orders for the phase adjustment.




PHASE ADJUSTMENT (H)




PHASE ADJUSTMENT ①




PHASE ADJUSTMENT (6)



PHASE ADJUSTMENT (L)




PHASE ADJUSTMENT M




3-3. PHASE ADJUSTMENT (Mechanism Chassis Upper Surface Parts)

Note: Adjust if at the **UNLOADING** position unless otherwise specified.


PHASE ADJUSTMENT N




PHASE ADJUSTMENT (P)




PHASE ADJUSTMENT @




PHASE ADJUSTMENT ®




PHASE ADJUSTMENT ®




PHASE ADJUSTMENT ①




PHASE ADJUSTMENT (1)




PHASE ADJUSTMENT (V)



PHASE ADJUSTMENT W



PHASE ADJUSTMENT ®



5-1-4. PERIODIC CHECK AND MAINTENANCE

 Carry out the following maintenance and periodic checks not only to fully display the functions and performance of the set, but also for the equipment and tape. After repairing, service the set as follows, regardless of the length of use.

4-1. CLEANING OF ROTARY DRUM ASSEMBLY

1) Press a wiping cloth (Ref No. J-2) moistened with cleaning fluid (Ref No. J-1) against the rotary drum assembly gently, and clean it while rotating the upper rotary drum assembly slowly with your finger in the counterclockwise direction.

Note: Do not rotate the motor on power or rotate the upper rotary drum assembly in the clockwise direction with your finger. The head tip will also be damaged if the wiping cloth is moved perpendicularly against it. Therefore, be sure to follow the above instructions when

Therefore, be sure to follow the above instructions when cleaning the rotary drum assembly.

4-2. CLEANING OF TAPE PATH SYSTEM (See Fig. 1.)

1) In the EJECT mode, clean the tape path systems (TG-1, 2, 3, 4, 5, 6, 7, 8, capstan) and the lower drum using a superfine applicator (Ref No. J-3) moistened with cleaning fluid.

Note 1: Make sure that no oil or grease of the link mechanisms sticks to the superfine applicator. (Ref No. J-3)

Note 2: Do not use a applicator moistened with alcohol to the other guide cleaning. But clean the pinch roller using alcohol




Fig. 1.

4-3. PERIODIC CHECKS

Location of Maintenance and Check		Hours of Use (H)										Remarks
		500	1000	1500	2000	2500	3000	3500	4000	4500	5000	Hemaiks
	Cleaning of tape path surface	0	0	0	0	0	0	0	0	0	0	Take care not to adhere the oil.
	Cleaning and degaussing of rotary drum assembly	0	0	0	0	0	0	0	0	0	0	
Driving System	Capstan shaft (Bearing)	_	☆	_	☆	_	☆	_	☆	_	☆	Make sure that no oil gets on the tape path surface.
	Loading motor	-	☆	_	☆	-	☆	_	☆	-	☆	A-7026-007-A
Performance Confirmation	Abnormal noise	☆	☆	☆	☆	☆	☆	☆	☆	☆	☆	
	Back tension measurement	_	☆	_	☆	-	☆	-	☆	-	. ಭ	
	Brake system	_	☆	_	☆	-	☆	-	☆	_	☆	
	FWD RVS Torque measurement	_	☆	-	☆	-	☆	_	☆	-	☆	

O: Cleaning ☆: Confirmation

Note: When overhauling, refer to the checks above and replace parts.

Note: Grease


Be sure to use the specified the grease. (The SG-055G is used all in the E mechanism)
 Check the quantity of grease when installing the parts which is needed to apply the grease.

• FLOIL (SG-055G): Part No. 7-651-000-09

5-1-5. MECHANISM SECTION CHECKS AND REPLACEMENTS


5-1. TAPE STOPPER, MOTOR FPC ASSEMBLY AND ELASTIC CONNECTOR

· Removing/Attaching



5-2. DRUM ASSEMBLY AND DRUM BASE


· Removing/Attaching



(PH-G) : Page 5-7 (PH-S) : Page 5-11

5-3. LM COVER, LM WORM WHEEL, LM HOLDER AND LM MOTOR ASSEMBLY

· Removing/Attaching



5-4. TG3/4 CATCHER BLOCK ASSEMBLY, PINCH DRIVING GEAR AND TC ARM ASSEMBLY


· Removing/Attaching Removing: After removing the LM motor assembly (Refer to 5-3.), remove each part. Attaching: After attaching each part and the LM motor block assembly, adjust them according to the flowchart (START-3) on page 5-43. (Only when the TG3/4 catcher block assembly is removed) Screw (M1.4 × 2.5) Never turn this screw as it has been adjusted. TG3/4 catcher block assembly PH-G Pinch driving gear COMMY. PH-S TC arm assembly Claw Fixing torque: 0.0686 Nem (0.7 kgecm)

5-5. PINCH ARM ASSEMBLY, PINCH LIMITER AND TENSION COIL SPRING (PINCH)


1. Removing

①. Set the **UNLOADING** position. (Refer to page 5-3)

2. Pinch arm block assembly.



3. Pinch arm assembly and pinch limiter.



2. Attaching

1. Attach the parts in the order of $\bigcirc \rightarrow \bigcirc \rightarrow \bigcirc$.

2. Adjust them according to the flowchart (START-3) on page 5-43.


5-6. HC ARM, HC ROLLER ASSEMBLY, PINCH RETAINER, PINCH CAM GEAR AND TG5/6 CATCHER BLOCK ASSEMBLY

1. Removing


1. Set the **UNLOADING** position. (Refer to page 5-3)

2. Pinch arm block assembly. (Refer to 5-5)


3. Flexible board.



4. HC arm, HC roller assembly and pinch retainer.



(5). Pinch cam gear and TG5/6 catcher block assembly.




2. Attaching

1. Attach the parts in the order of $\textcircled{1} \rightarrow \textcircled{5} \rightarrow \textcircled{4} \rightarrow \textcircled{3} \rightarrow \textcircled{2}$.

2. Adjust them according to the flowchart (START-3) on page


5-7. RL ARM AND RL LINK

• Removing/Attaching (L cassette) position. (Refer to page 5-2))



5-8. GOOSENECK GUARD AND GOOSENECK ARM ASSEMBLY

• Removing/Attaching (L cassette) position. (Refer to page 5-2))




5-9. TENSION COIL SPRING (TG2), SPRING ADJUSTOR, TG2 SPRING HOOK, TG2 SELECTION ARM AND DAMPER SHEET


1. Removing

①. Set the **UNLOADING** position. (Refer to page 5-3)


2. Tension coil spring (TG2) and spring adjustor.



3. LM motor block assembly.



4. TG2 spring hook and TG2 selection arm.



2. Attaching

1. Attach the parts in the order of $\bigcirc \rightarrow \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc$.


2. Adjust them according to the flowchart (START-2) on page 5-43

5-10. TENSION COIL SPRING (TG7), SPRING ADJUSTOR AND TG7 SPRING HOOK


1. Removing

①. Set the **UNLOADING** position. (Refer to page 5-3)

2. Tension coil spring (TG7) and spring adjustor.



3. TG7 spring hook.



2. Attaching


1. Attach the parts in the order of \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc .

2. Adjust them according to the flowchart (START-2) on page 5.43


5-11. TG2 RETAINER, TG2 ARM ASSEMBLY (TG2 PLATE SPRING AND ET MAGNET), S TENSION REGULATOR BAND ASSEMBLY AND TG2 LOAD ARM ASSEMBLY

1. Removing


- ①. Tension coil spring (TG2), spring adjustor, LM motor block assembly and TG2 spring hook. (Refer to 5-9)
- 2. TG2 retainer.




- 3. Set the **LOADING** position. (Refer to page 5-3)
- 4. TG2 arm assembly and S tension regulator band assembly.



S tension regulator band assembly, TG2 plate spring and ET magnet.



6. TG2 load arm assembly.




- 1. Set the **UNLOADING** position. (Refer to page 5-3)
- 2. Attach the parts in the order of $(6) \rightarrow (3) \rightarrow (5) \rightarrow (4) \rightarrow (2) \rightarrow (1)$
- 3. Operation check: **LOADING**/**UNLOADING**. (Refer to page 5-3)
- **4.** Adjust them according to the flowchart (START-2) on page 5-43.


5-12. TG7 RETAINER, TG7 ARM ASSEMBLY (TG7 PLATE SPRING AND ET MAGNET), T TENSION REGULATOR BAND ASSEMBLY AND TG7 LOAD ARM ASSEMBLY

1. Removing


①. Tension coil spring (TG7), spring adjustor.




②. TG7 spring hook.




- 3. Set the **LOADING** position. (Refer to page 5-3)
- 4. TG7 arm assembly and T tension regulator band assembly.



(5). TG7 plate spring, ET magnet and T tension regulator band assembly.



6. TG7 load arm assembly.



2. Attaching

1. Set the UNLOADING position. (Refer to page 5-3)

2. Attach the parts in the order of $\textcircled{6} \rightarrow \textcircled{3} \rightarrow \textcircled{5} \rightarrow \textcircled{4} \rightarrow \textcircled{2} \rightarrow \textcircled{1}$.

3. Operation check: **LOADING**/**UNLOADING**. (Refer to page 5-3.)


4. Adjust them according to the flowchart (START-2) on page 5-43.

5-13. S REEL TABLE BLOCK ASSEMBLY


1. Removing

1. Set the (L cassette) position. (Refer to page 5-2)


2. RL arm




3. Gooseneck guard.



4. S band adjustor.



⑤. S reel table block assembly.




- 2. Adjust them according to the flowchart (START-1) on page 5-43.

5-14. T REEL HOLDER AND T REEL TABLE BLOCK ASSEMBLY


1. Removing

1. Set the L cassette position. (Refer to page 5-2)

2. T band adjustor.



3. T reel holder and T reel table block assembly.



2. Attaching

1. Attach the parts in the order of $\bigcirc \rightarrow \bigcirc \rightarrow \bigcirc$.


2. Adjust them according to the flowchart (START-1) on page 5-43.

5-15. S REEL PLATE ASSEMBLY

1. Removing


1). Set the L cassette position. (Refer to page 5-2)

②. S band adjustor.




3. Set the **LOADING** position. (Refer to page 5-3)


4. S brake.




(5). RL arm and RL link.




6. Gooseneck guard.



7. Plate spring



8. S reel plate assembly.




- **1.** Attach the parts in the order of $(1) \rightarrow (8) \rightarrow (7) \rightarrow (4) \rightarrow (3) \rightarrow (2) \rightarrow (6) \rightarrow (5)$.
- Adjust them according to the flowchart (START-1) on page 5-43.

5-16. T REEL PLATE ASSEMBLY

1. Removing


1. Set the Leassette position. (Refer to page 5-2)

2. T band adjustor.




3. Set the **LOADING** position. (Refer to page 5-3)


4. T ratchet.



⑤. Plate spring.



6. T reel plate assembly.




- **1.** Attach the parts in the order of $\textcircled{1} \rightarrow \textcircled{6} \rightarrow \textcircled{5} \rightarrow \textcircled{3} \rightarrow \textcircled{4} \rightarrow \textcircled{2}$.
- 2. Adjust them according to the flowchart (START-1) on page 5-43.

5-17. TG1/8 BASE ASSEMBLY, S BRAKE AND TRATCHET


1. Removing

1. Set the **LOADING** / **L cassette** positions. (Refer to pages 5-2 to 5-3)

②. TG1/8 base assembly.




3. S brake and T ratchet.



- **1**. Attach the parts in the order of \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc .
- 2. Adjust them according to the flowchart (START-2) on page 5-43.

5-18, TG3/6 ROLLER ASSEMBLY AND TG3/6 LOCK SPRING


· Removing/Attaching



Note: After attaching each part, adjust them according to the flowchart (START-3) on page 5-43.

5-19. FL JOINT GEAR, TG5/6 RETAINER AND CAPSTAN MOTOR


• Removing/Attaching



PH-A : Page 5-6 PH-D , PH-G : Page 5-7

5-20. FL SELECTION ARM, FL RELAY GEAR AND FL JOINT ARM ASSEMBLY

· Removing/Attaching




5-21. ROTARY SWITCH, TC GEAR AND RELAY GEAR

1. Removing


1. Set the **UNLOADING** position. (Refer to page 5-3)

2. Rotary switch, TC gear and relay gear.



2. Attaching

 Remove the LM motor block assembly and LM worm wheel. (To synchronize phase of the pinch driving system (front side) and the loading driving system (back side))




Attach the TC gear, relay gear and rotary switch.

3. Attach the LM worm wheel and LM motor block assembly.

5-22. GL ARM RETAINER AND GL ARM

• Removing/Attaching (UNLOADING) position. (Refer to page 5-3))




5-23. M SLIDER AND M SLIDER ARM

1. Removing

1. Set the **UNLOADING** position. (Refer to page 5-3)

②. GL arm retainer and GL arm. (Refer to 5-22)

3. M slider and M slider arm.



2. Attaching

• Attach the parts in the order of $\bigcirc \rightarrow \bigcirc \rightarrow \bigcirc$.

PH-B : Page 5-6 PH-R : Page 5-10 PH-D , PH-E , PH-F : Page 5-7 PH-V : Page 5-11

5-24. TG7 SELECTION ARM, TG7 CAM GEAR AND T CAM GEAR


1. Removing

1. Set the UNLOADING position. (Refer to page 5-3)

②. GL arm retainer and GL arm. (Refer to 5-22)

3. M slider and M slider arm. (Refer to 5-23)

4. TG7 selection arm, TG7 cam gear and T cam gear.



2. Attaching

1. Attach the parts in the order of $\bigcirc \rightarrow \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc$.



Cam groove on the T cam gear.

Apply the grease (3 mm dia) to of cam groove (part).

T cam gear (T ratchet driving side)


PH-A : Page 5-6 PH-G : Page 5-7 (PH-T), (PH-U) : Page 5-11

5-25. MAIN CAM, TG2 SL ARM ASSEMBLY AND TENSION COIL SPRING (TG2 SL)


The two grooves on one side of the main cam drive the TG2 selection arm and the TG2 load arm assembly. Since it is difficult to attach the main cam, fix the TG2 selection arm and the TG2 load arm assembly with the main cam's phase adjusted correctly (Nearly unloading position (See 3-1. Phase Adjustment (A): page 5-6)), so that later mounting work can be performed smoothly. If fixed parts are shifted, follow "3-3. Phase Adjustment T, W: page 5-11".

①. Set the **UNLOADING** position. (Refer to page 5-3)


2. LM motor block assembly and LM worm wheel.



- 3. GL arm retainer and GL arm. (Refer to 5-22)
- 4. M slider and M slider arm. (Refer to 5-23)
- (5). TG2 spring hook.



6. Main cam, TG2 SL arm assembly and tension coil spring (TG2



2. Attaching

- 2. Adjust them according to the flowchart (START-2) on page 5-43.



Main cam (rear side)

Cam groove on the main cam.

Apply the grease (12 mm dia.) to each two of cam groove (part).


PH-A : Page 5-6 PH-F , PH-G : Page 5-7

5-26. TG3/4 ARM BLOCK ASSEMBLY (TG3/4 ARM ASSEMBLY, TG3/4 LIMITER SPRING AND TG3/4 GEAR), TG3/4 BASE BLOCK ASSEMBLY (TG3/4 BASE ASSEMBLY)


1. Removing

①. Set the **UNLOADING** position. (Refer to page 5-3)


2. LM motor block assembly and LM worm wheel.




3. Pinch arm block assembly and pinch retainer.




(4). TG3/4, TG5/6 catcher block assembly and pinch cam gear.



(5). GL arm retainer and GL arm.




(6). Rotary switch, TC gear and TG5/6 retainer.




7. Set the **LOADING** position. (Refer to page 5-3)


(8). TG3/4 arm block assembly.



TG3/4 base block assembly.



TG3/4 arm assembly, TG3/4 limiter spring, TG3/4 gear, TG3/6 roller assembly, TG3 lock spring and TG3/4 base assembly.




2. Attaching

1. Set the **LOADING** position. (Refer to page 5-3)

2. TG3/4 arm assembly, TG3/4 limiter spring, TG3/4 gear, TG3/6 roller assembly, TG3 lock spring and TG3/4 base assembly.

3. TG3/4 base block assembly.

4. TG3/4 arm block assembly.




5. Attach the parts in the order of $\textcircled{6} \rightarrow \textcircled{5} \rightarrow \textcircled{4} \rightarrow \textcircled{3} \rightarrow \textcircled{2}$.

6. Adjust them according to the flowchart (START-3) on page 5-43.


5-27. TG5/6 ARM BLOCK ASSEMBLY (TG5/6 ARM ASSEMBLY, TG5/6 LIMITER SPRING AND TG5/6 GEAR), TG5/6 BASE BLOCK ASSEMBLY (TG5/6 BASE ASSEMBLY)

1. Removing


- ①. Set the **UNLOADING** position. (Refer to page 5-3)
- 2). LM motor block assembly and LM worm wheel.




3. Pinch arm block assembly and pinch retainer.




TG3/4, TG5/6 catcher block assembly, screw of capstan motor and pinch cam gear.



(5). GL arm retainer and GL arm.




6. FL joint gear, capstan motor and TG5/6 retainer.




7. Set the **LOADING** position. (Refer to page 5-3)


3. TG3/4 arm block assembly and TG5/6 arm block assembly.



TG5/6 base block assembly.



① TG5/6 arm assembly, TG5/6 limiter spring, TG5/6 gear, TG3/6 roller assembly, TG6 lock spring and TG5/6 base assembly.




2. Attaching

1. Set the **LOADING** position. (Refer to page 5-3)

2. TG5/6 arm assembly, TG5/6 limiter spring, TG5/6 gear, TG3/6 roller assembly, TG6 lock spring and TG5/6 base assembly.

3. TG5/6 base block assembly.

4. TG3/4 arm block assembly and TG5/6 arm block assembly.




5. Attach the parts in the order of $(6) \rightarrow (5) \rightarrow (4) \rightarrow (3) \rightarrow (2)$.

6. Adjust them according to the flowchart (START-3) on page 5-43.

5-28. REEL MOTOR

· Removing/Attaching




5-29. RS ARM ASSEMBLY

1. Removing

1. Set the S/L cassette position. (Refer to page 5-2)

2. RS arm assembly.



2. Attaching

• Attach the parts in the order of \bigcirc \rightarrow \bigcirc .


5-30. RS GEAR ASSEMBLY, MIC PRESS SPRING AND MIC LEVER

1. Removing


1. Set the S/L cassette position. (Refer to page 5-2)

2. RS arm assembly. (Refer to 5-29)

3. RS gear assembly.



4. Mic press spring and Mic lever.




2. Attaching

• Attach the parts in the order of $\bigcirc \rightarrow \bigcirc \rightarrow \bigcirc \rightarrow \bigcirc$.


5-31. RACK JOINT GEAR, RACK HOLDER, MIC HOLDER, RACK (LC) AND RACK (SC)

1. Removing


1. FL joint gear, TG7 selection arm.



- 2. Set the S/L cassette position. (Refer to page 5-2)
- ③. RS arm assembly. (Refer to 5-29)
- (4). RS gear assembly, Mic press spring and Mic lever. (Refer to 5-30)
- (5). Rack joint gear, rack holder, rack (LC) and rack (SC).



- 6. FL block assembly. (Refer to page 5-2)
- 7. Push switch, Connector (4P), Mic holder.



2. Attaching


• Attach the parts in the order of $② \rightarrow ⑦ \rightarrow ⑥ \rightarrow ⑤ \rightarrow ④ \rightarrow ③ \rightarrow ①$.

5-32. PLATE LINK ASSEMBLY


1. Removing

1. Set the **L cassette** position. (Refer to page 5-2)


2. RL arm.




3. Gooseneck guard.




4. S and T band adjustor.



(5). TG1/8 base assembly. (S brake and T ratchet)




T reel plate assembly.




PH-H : Page 5-8 PH-M : Page 5-9 PH-N : Page 5-10

7. S reel plate assembly.




(B). Set the (S/L cassette) position. (Refer to page 5-2)

RS arm assembly.



10. Plate link assembly.




2. Attaching

1. Attach the parts in the order of $(8) \rightarrow (10) \rightarrow (9) \rightarrow (7) \rightarrow (6) \rightarrow (5) \rightarrow (1) \rightarrow (4) \rightarrow (3) \rightarrow (2)$.

2. Adjust them according to the flowchart (START-1) on page 5-43


5-33. ROLLER SHAFT ASSEMBLY AND ROLLER BELT

Removing/Attaching




5-34. LID OPENER

Removing/Attaching




5-35. C DOOR

Removing/Attaching



5-36. DAMPER ARM AND TENSION SPRING (DB)

· Removing/Attaching




5-37. GEAR (A), GEAR (B) AND C WORM


1. Removing

①. FL block assembly. (Refer to page 5-2)


②. Screw. $(M1.4 \times 2.5)$



3. Gear cover.



4. Gear (A), Gear (B) and C worm.




2. Attaching

• Attach the parts in the order of $\textcircled{4} \rightarrow \textcircled{3} \rightarrow \textcircled{2} \rightarrow \textcircled{1}$.


5-38. TENSION COIL SPRING (HS), TENSION SPRING (DB), SHIFT PLATE SPRING AND C SLOAT BLOCK ASSEMBLY

• Removing/Attaching (Remove the FL block assembly. (Refer to page 5-2))




5-1-6. ADJUSTMENTS AND CHECKS

6-1. ADJUSTMENT POSITION



6-2. ADJUSTMENT ORDER



6-3. ADJUSTMENT AND CHECKING METHOD

6-3-1. REEL TABLE HEIGHT ADJUSTMENT


1. Preparation before adjustment


FL block: Remove.

Position: LOADING / S cassette

Jig used: Cassette standard plate (J-7), Reel standard plate (J-8) and screwdriver (+)

2. Adjusting





6-3-2. TG2/TG7 HEIGHT ADJUSTMENT

1. Preparation before adjustment


FL block: Remove.

Position: LOADING / S cassette

Jig used: Cassette standard plate (J-7), TG2/7 preset plate (J-9) and screwdriver

(For attaching jigs, refer to page 5-5)

2. Adjusting



6-3-3. FWD/RVS POSITION ADJUSTMENT


1. Preparation before adjustment

FL block: Remove.

Position: (LOADING (The pinch roller should be stuck))/(S cassette)

Jig used: Cassette standard plate (J-7), TG2/7 preset plate (J-9) and screwdriver for tape path

2. Adjusting



6-3-4. TG2/TG7 ELECTRIC TENSION REGULATOR ADJUSTMENT

1. Preparation before adjustment


FL block: Remove.

Position: (LOADING (The pinch roller should be stuck)) / (S cassette)


Jig used: Cassette standard plate (J-7), TG2/7 preset plate (J-9) and screwdriver for tape path

(For attaching jigs, refer to page 5-5)


2. Connecting, setting, and adjusting methods



3. Adjusting



4. Adjustment order



6-3-5. FWD/RVS BACK TENSION ADJUSTMENT


1. Preparation before adjustment

Mechanism deck: Install to the unit.


Jig used : Mini DV torque cassette (J-6), pinset (For change the hooking of spring)

2. Adjusting


Note: At the FWD (TG2) side, measure the Mini DV torque cassette (J-6) in the FWD mode. At the RVS (TG7) side, measure the Mini DV torque cassette (J-6) in the RVS mode.




6-3-6. ADJUSTMENT PREPARATIONS AND RF WAVEFORM CHECK




6-3-7. TRACKING ADJUSTMENT




6-3-8. TRACKING CHECK




6-3-9. CUE AND REV CHECK



6-3-10. RISING CHECK



6-3-11. TAPE PATH CHECK



③ After adjustment and check, return the data set by the adjusting remote commander to the original value.

5-2. SERVICE MODE

5-2-1. ADJUSTING REMOTE COMMANDER

The adjusting remote commander is used for changing the calculation coefficient in signal processing, EVR data, etc. The adjusting remote commander performs bi-directional communication with the unit using the remote commander signal line (LANC). The resultant data of this bi-directional communication is written in the non-volatile memory.

1. Used Adjustment Remote Commander

- 1) With the unit set in STANDBY mode, connect the adjusting remote commander to the remote (LANC) terminal.
- Adjust the HOLD switch of the adjusting remote commander to "HOLD" (SERVICE position).
- 3) Turn on the power with the ON/STANDBY switch of the unit. If it has been properly connected, the LCD on the adjusting remote commander will display as shown in Fig. 5-2-1.




Fig. 5-2-1.

- 4) Operate the adjusting remote commander as follows.
 - Changing the page

The page increases when the EDIT SEARCH + button is pressed, and decreases when the EDIT SEARCH - button is pressed. There are altogether 16 pages, from 0 to F.

Hexadecimal notation	0 1 2 3 4 5 6 7 8 9 A B C D E F
LCD Display	0 123456789AbcdEF
Decimal notation conversion value	0 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15

Table 5-2-1.

- · Changing the address
- The address increases when the FF (►►) button is pressed, and decreases when the REW (◄◄) button is pressed. There are altogether 256 addresses, from 00 to FF.
- Changing the data (Data setting)
 The data increases when the PLAY (►) button is pressed, and decreases when the STOP (■) button is pressed. There are altogether 256 data, from 00 to FF.
- Writing the adjustment data
 The PAUSE button must be pressed to write the adjustment data (C page, D page and E page) in the nonvolatile memory.
 (The new adjustment data will not be recorded in the nonvolatile memory if this step is not performed.)

2. Precautions Upon Using The Adjusting Remote Commander

Mishandling of the adjusting remote commander may erase the correct adjustment data at times. To prevent this, it is recommended that all adjustment data be noted down before beginning adjustments and new adjustment data after each adjustment.

5-2-2. DATA PROCESSING

The calculation of the adjusting remote commander display data (hexadecimal notation) are required for obtaining the adjustment data of some adjustment items. In this case, after converting the hexadecimal notation to decimal notation, calculate and convert the result to hexadecimal notation, and use it as the adjustment data. Table 5-2-2. indicates the hexadecimal notation- the decimal notation, calculation table.

Hexa	Hexadecimal nontation-Decimal notation								②								
	The lower digits of the hexadecimal notation	0	1	2	3	4	5	6	7	8	9	А	В	С	D	E	F
	The upper digits of the hexadecimal notation											(A)	(Ь)	(c)	(d)	(E)	(F)
	0	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
1	1	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
-	. 2	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47
- 1	3	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63
	4	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79
1	5	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95
	6	96	97	98	99	100	101	102	103	104	105	106	107	108	109	110	111
	7	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126	127
	8	128	129	130	131	132	133	134	135	136	137	138	139	140	141	142	143
	9	144	145	146	147	148	149	150	151	152	153	154	155	156	157	158	159
	A (H)	160	161	162	163	164	165	166	167	168	169	170	171	172	173	174	175
①→	B (b)	176	177	178	179	180	181	182	183	184	185	186	187	188	189	190	191
	C (c)	192	193	194	195	196	197	198	199	200	201	202	203	204	205	206	207
	D (d)	208	209	210	211	212	213	214	215	216	217	218	219	220	221	222	223
	E (E)	224	225	226	227	228	229	230	231	232	233	234	235	236	237	238	239
	F (F)	240	241	242	243	244	245	246	247	248	249	250	251	252	253	254	255

Note: () indicate the adjusting remote control unit display.

(**Example**) In the case that the adjusting remote control unit display are BD (bd).

As the upper digit of the hexadecimal notation is B (b), and the lower digit is D (d), the intersection "189" of the ① and ② in the above table is the decimal notation to be calculated.

Table 5-2-2.

5-2-3. SERVICE MODE

1. Emergence Memory Address

Page C	Addresses 30 to 3B

Address	Contents			
30	EMG code when first error occurs			
32	Upper: MSW code when shift starts when first error occurs Lower: MSW code when first error occurs			
33	Lower: MSW code to be moved when first error occurs			
34	34 EMG code when second error occurs			
36	Upper: MSW code when shift starts when second error occurs Lower: MSW code when second error occurs			
37	Lower: MSW code to be moved when second error occurs			
38	EMG code when last error occurs			
3A	Upper: MSW code when shift starts when last error occurs Lower: MSW code when last error occurs			
3B	Lower: MSW code to be moved when last error occurs			

When no error occurs in the unit, data 00 is written in the above addresses (30 to 3B). When the first error occurs in the unit, the data corresponding to the error is written in the first emergency address (30 to 33). In the same way, when the second error occurs, the data corresponding to the error is written in the second emergency address (34 to 37).

Finally, when the last error occurs, the data corresponding to the error is written in the last emergency address (38 to 3B). Consequently, addresses 30 to 3B are updated each time errors occur.

Note 1: After completing adjustments, be sure to rewrite the data of addresses 30 to 3B to 00.

- 1) Select page: 0, address: 01, and set data: 01.
- Select page: C, address: 30, set data: 00, and press the PAUSE button of the adjusting remote commander.
- Select page: C, address: 31, set data: 00, and press the PAUSE button of the adjusting remote commander.
- Select page: C, address: 32, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 5) Select page: C, address: 33, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 6) Select page: C, address: 34, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 7) Select page: C, address: 35, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 8) Select page: C, address: 36, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 9) Select page: C, address: 37, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 10) Select page: C, address: 38, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 11) Select page: C, address: 39, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 12) Select page: C, address: 3A, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 13) Select page: C, address: 3B, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 14) Select page: 0, address: 01, set data: 00, and press the PAUSE button of the adjusting remote commander.

1-1. EMG Code (Emergency Code)

Codes corresponding to the errors which occur are written in addresses 30, 34, 38. The type of error indicated by the code are shown in the following table.

Code	Error Type	
00	No error (Initial state)	
· 10	Loading motor time-out during LOAD	
11	Loading motor time-out during UNLOAD	
20	Reel motor error	
22	T reel error	
23	S reel error	
24	Swing error	
32	Error during normal capstan rotation	
33	Cassette compartment LOAD error	
35	Cassette compartment UNLOAD error	
40	FG error during drum start-up	
42	FG error during normal drum rotation	
50	DEW detection	
52	Wet DEW detection	
60	Electrical tension regulator error	

5-3. VIDEO SECTION ADJUSTMENTS

When performing adjustments, refer to the layout diagrams for adjustment related parts on page 5-88.

3-1. PREPARATIONS BEFORE ADJUSTMENT

3-1-1. Equipment Used

- 1) TV monitor
- 2) Oscilloscope with 2-phenomenon, 30 MHz band, and delay mode (Unless specified otherwise, use a 10:1 probe)
- 3) Frequency counter
- 4) Digital voltmeter
- 5) Audio generator
- 6) Audio level meter
- 7) Audio distortion meter
- 8) Audio attenuator
- Pattern generator (with VIDEO OUTPUT terminal and external sync function)
- 10) Digital camera recorder NTSC: DCR-VX1000

PAL : DCR-VX1000E

11) Vectorscope

12) Alignment tape

• SW/OL reference (XH2-3)

Parts code: 8-967-997-11

• Audio operation check for NTSC (XH5-3)

Parts code: 8-967-997-51

• System operation check for NTSC (XH5-5)

Parts code: 8-967-997-61

Audio operation check for PAL (XH5-3P)

Parts code: 8-967-997-55

• System operation check for PAL (XH5-5P)

Parts code: 8-967-997-66

• BIST check for NTSC (XH5-6)

Parts code: 8-967-997-71

• BIST check for PAL (XH5-6P)

Parts code: 8-967-997-76

- 13) Adjusting remote control unit (J-6082-053-B)
- 14) Multi CPC-8 jig (J-6082-388-A). (CN775 of the RP-228 board)
- 15) Extension board
 - For extension between CN101 of the RP-228 board and CN412 of the JC-19 board.
 - For extension between CN102 of the RP-228 board and CN411 of the JC-19 board. (30P, 0.5 mm) (J-6082-270-A)
 - For extension between CN771 of the RP-228 board and drum (M901) (10P, 1 mm) (J-6082-064-A)
 - For extension between CN002 of the CM-56 board and CN501 of the VA-102 board (8P, 1 mm) (J-6082-058-A)
 - For extension between CN006 of the CM-56 board and the reel motor (M904) (15P, 1.25 mm) (J-6902-354-A)
 - For extension between CN001 of the CM-56 board and CN101 of the MD-63 board (16P, 1 mm) (J-6082-020-A)
- 16) Regulated power supply

NTAC: DSR-20MD PAL: DSR-20MDP

3-1-2. Connection of Equipment

VIDEO SELECT button on the front panel.

According to the specification for the input terminal (S VIDEO input, VIDEO input, or DV input), connect measuring equipment as shown in Fig. 5-3-1, and make adjustment.

The input terminal is specified in () of the signal column.

Any input terminal can be used unless otherwise specified.

To switch between S VIDEO INPUT and VIDEO INPUT, use the

Note 1: In adjustments specifying for the S VIDEO input to be used, using the VIDEO input would disable the product specifications of this unit from being satisfied. Always

use the input signal specified.

Note 2: If adjustments are used with the VTR with the S video output terminal as the signal source, the performance of this unit may be affected by the VTR. Use a pattern generator with a Y/C separator terminal as much as possible.




Fig. 5-3-1.

3-1-3. Adjusting Connectors (RP-228 Board CN775)

Some of the adjusting points of the video section are concentrated at CN775 of the RP-228 board. Connect the instruments via the multi CPC-8 jig (J-6082-388-A)

Pin No.	Signal Name	Pin No.	Signal Name
1	TCK	2.	TMS
3	TDI	4	GND
5	TRACK ID	6	JSWP
7	GND	8	RF MONITOR
9	VCC2	10	AGC IN
11	VCC1	12	EQ IN
. 13	LOCK	14	REF OUT
15	ENV OUT	16	GND
17	TDO	18	C1ERP
19	FLTD	20	GND

Table 5-3-1.

3-1-4. Checking the Input Signals

Because the video signal obtained from the pattern generator is used as the adjustment signal for adjustments, the video output signal must satisfy the given specifications.

1. S VIDEO Input

Connect the oscilloscope to the Y signal terminal of the S VIDEO input terminal, and check that the sync signal of the Y signal is approximately <0.286> [0.30] V and that the amplitude of the video section is approximately <0.714> [0.70] V. (When a VTR with the S VIDEO output terminal is used, also check that the chroma signal and burst signal have not remained)

Connect the oscilloscope to the chroma signal terminal of the S VIDEO input terminal, and check that the burst signal amplitude of the chroma signal is approximately <0.286> [0.30] V and flat, and that the red signal amplitude of the chroma signal is approximately <0.66> [0.67] V. The Y and chroma signals used in the adjustment are shown in Fig. 5-3-2.

< >: NTSC model

2. VIDEO Input

Connect the oscilloscope to the VIDEO input terminal, and check that the sync signal amplitude of the video signal is approximately <0.286> [0.30] V, the amplitude of the video section is approximately <0.714> [0.70] V, the amplitude of the burst signal is approximately <0.286> [0.30] V and flat, and that the red signal amplitude of the chroma signal is approximately <0.66> [0.67] V. The video signal (color bar) used for adjustments is shown in Fig. 5-3-3.

< >: NTSC model
[]: PAL model

NTAC: DSR-20MD PAL: DSR-20MDP




Fig. 5-3-2. Color Bar Signal of Pattern Generator




Fig. 5-3-3. Color Bar Signal of Pattern Generator

3-1-5. Adjustment Tapes

Use the alignment tapes shown in the following table. Use tapes specified in the signal column of each adjustment.

Name	Use
SW/OL standard (XH2-3)	Switching position adjustment
Audio operation check (XH5-3 (NTSC), XH5-3P (PAL))	Audio system adjustment
System operation check (XH5-5 (NTSC), XH5-5P (PAL))	Operation check
BIST check (XH5-6 (NTSC), XH5-6P (PAL))	BIST check

Table 5-3-2.

Fig. 5-3-4. shows the 75% color bar signals recorded on the alignment tape for Audio Operation Check (NTSC).




Fig. 5-3-4. Color Bar Signal of Alignment Tapes

3-1-6. Input/ Output Level and Impedance

LINE IN

Video input BNC connector

Input signal: 1 Vp-p

(75 ohms unbalanced)

S Video input Mini DIN 4-pin

Luminance signal: 1 Vp-p (75 ohms unbalanced)

Chrominance signal: 0.286 Vp-p

(NTSC), 0.3 Vp-p (PAL), (75 ohms unbalanced)

Audio input Phono jack (L, R)

Input level: 2 Vrms (full bit)

Input impedance: more than 47 kohms

LINE OUT

Video output BNC connector

Output signal: 1 Vp-p

(75 ohms unbalanced)

S Video output Mini DIN 4-pin

Luminance signal: 1 Vp-p

(75 ohms unbalanced) Chrominance signal: 0.286 Vp-p

(NTSC), 0.3 Vp-p (PAL),

(75 ohms unbalanced)

Audio output Phono jack (L, R)

Output level: 2 Vrms (full bit)

Output impedance: less than 10 kohms

3-2. POWER SUPPLY SYSTEM ADJUSTMENT

1. Power Supply Voltage Check Power Block (U-2 Board)

Power Block (0-2 Board)			
Mode	Playback		
Measuring Instrument	Digital voltmeter		
UNSW6V Check			
Measuring Point	Pin ① of CN11		
Specified Value	$6.0 \pm 0.5 \text{Vdc}$		
UNSW3.1V Check			
Measuring Point	Pin ② of CN11		
Specified Value	3.1 ± 0.2 Vdc		
VIDEO5V, AUDIO5V	Check		
Measuring Point	Pin ③, ⑦ of CN11		
Specified Value	5.0 ± 0.12 Vdc		
SW3.1V Check			
Measuring Point	Pin 4 of CN11		
Specified Value	3.1 ± 0.1 Vdc		
VIDEO-5V, AUDIO-5V	V Check		
Measuring Point	Pin 6, 9 of CN11		
Specified Value	-5.0 ± 0.12 Vdc		
SW5V Check			
Measuring Point	Pin ③ of CN11		
Specified Value	$5.0 \pm 0.12 \mathrm{Vdc}$		
DRUM6V Check			
Measuring Point	Pin ④ of CN10		
Specified Value	$6.0 \pm 0.5 \mathrm{Vdc}$		
MOTOR14V Check			
Measuring Point	Pin 6 of CN10		
Specified Value	14.0 ± 2 Vdc		

2. Video/Audio Block Power Supply Voltage Check Power Block (U-2 Board)

Mode	Playback
Measuring Instrument	Digital voltmeter
UNSW6V Check	
Measuring Point	Pin ② of CN12
Specified Value	$6.0 \pm 0.5 \mathrm{Vdc}$
UNSW3.1V Check	
Measuring Point	Pin ④ of CN12
Specified Value	3.1 ± 0.2 Vdc
UNSW-9V Check	
Measuring Point	Pin ⑦ of CN12
Specified Value	−9 ± 0.5 Vdc
UNSW14V Check	
Measuring Point	Pin ® of CN12
Specified Value	14 ± 2.0 Vdc
UNSW12V Check	
Measuring Point	Pin ③ of CN12
Specified Value	12.0 ± 1 Vdc

3-3. SYSTEM CONTROL SYSTEM ADJUSTMENT

1. Initializing the C, D, E Page Data

Note 1: If "Initializing the C, D, E Page Data" is performed, all data of the C page, D page and E page will be initialized.

Note 2: If the C, D, E page data has been initialized, "Modification of C, D, E page Data" and all adjustments need to be performed again.

Mode	E-E
Signal	Arbitrary
Adjustment Page	С
Adjustment Address	00 to 6F
Adjustment Page	D
Adjustment Address	00 to 4F
Adjustment Page	Е
Adjustment Address	00 to 3B

2. Input of C page Initial Data

Input method:

- 1) Select page: 0, address: 01, and set data: 01.
- Select page: 4, address: 02, set data: 01, and press the PAUSE button of the adjusting remote commander.
- 3) Select page: 4, address: 02, and confirm that the data change in the order of "01" \rightarrow "03" \rightarrow "05" \rightarrow "00".
- 4) Modify the C page data. (Refer to C page address)

3. Input of D page Initial Data

Input method:

- 1) Select page: 0, address: 01, and set data: 01.
- Select page: 2, address: 00, set data: 2D, and press the PAUSE button of the adjusting remote commander.
- 3) Select page: 2, address: 01, set data: 2D, and press the PAUSE button of the adjusting remote commander.
- 4) Select page: 2, address: 02, and confirm that the data is "01".
- 5) Modify the D page data. (Refer to D page address)

4. Input of E page Initial Data

Input method:

- 1) Select page: 0, address: 01, and set data: 01.
- 2) Select page: 5, address: 00, set data: 2D, and press the PAUSE button of the adjusting remote commander.
- 3) Select page: 5, address: 01, set data: 2D, and press the PAUSE button of the adjusting remote commander.
- 4) Select page: 5, address: 02, and confirm that the data is "01".
- 5) Modify the E page data. (Refer to E page address)

5. Modification of C, D, E, Page Data

If the C, D, E page data has been initialized, change the data of the "Fixed data-2" address shown in the following tables by manual input.

Modifying Method:

- Before changing the data, select page: 0, address: 01, and set data: 01.
- New data for changing are not shown in the tables because they are different in destination. When changing the data, copy the data built in the same model.

Note: If copy the data built in the different model, this set may not operate.

- 3) When changing the data, press the PAUSE button of the adjusting remote commander each time when setting new data to write the data in the non-volatile memory.
- 4) Check that the data of adjustment addresses is the initial value. If not, change the data to the initial value.
- After completing "Modification of C, D, E Page Data", select page: 0, address: 01, and set data: 00. Also perform all adjustments.

6. Page C Address List

Note 1: Fixed data 1: Initialized data. (Refer to 2. Input of C Page Initial Data)

Fixed data 2: Modified data. (Refer to 5. Modification of C, D, E, Page Data)

Address	Initial Value	Remark			
00	Fixed data 1 (Initial data)				
01	Fixed data 2 (Changed data, Read from same model and copy it.)				
02 to 07	Fixed data 1 (Initial data)				
08	Fixed data 2 (C copy it.)	Fixed data 2 (Changed data, Read from same model and copy it.).			
09	Fixed data 1 (Initial data)				
0A to 0C	Fixed data 1 (Ir	nitial data)			
OD.	Fixed data 2 (C copy it.)	hanged data. Read from same model and			
0E,0F	Fixed data 1 (In	nitial data)			
10 to 19	Fixed data 1 (In	nitial data)			
1A to 1F	Fixed data 1 (Ir	nitial data)			
20 to 29	Fixed data 1 (Ir	nitial data)			
2A to 2F	Fixed data 1 (Initial data)				
30 to 39	00	Emargency memory address			
3A, 3B	00	Emargency memory address			
3C, 3D	F8	PLL fo adjustment			
3E, 3F	70	Recording current adjustment			
40, 41	C0	AEQ adjustment			
42, 43	90	AEQ adjustment			
44	86	AGC Center level adjustment			
45	Fixed data 1 (In	nitial data)			
46	86	PLL Capture range adjustment			
47	C8	CLK delay adjustment			
48, 49	Fixed data 1 (In	nitial data)			
4A, 4B	Fixed data 1 (Ir	nitial data)			
4C to 4F	00	Switching position adjustment			
50	54	Capstan FG duty adjustment			
51	31	Capstan FG duty adjustment			
52 to 59	Fixed data 1 (Initial data)				
5A	00 AEQ adjustment				
5B to 5F	Fixed data 1 (Initial data)				
60 to 69	Fixed data 1 (Initial data)				
6A to 6F	Fixed data 1 (Initial data)				

Table 5-3-3.

7. Page D Address List

Note 1: Fixed data 1: Initialized data. (Refer to 3. Input of D Page Initial Data)

Fixed data 2: Modified data. (Refer to 5. Modification of C, D, E, Page Data)

Address	Initial Value	Remark				
00 to 0F						
10 to 12	Fixed data 1 (I	nitial data)				
13	Fixed data 2 (Changed data, Read from same model and copy it.)					
14	Fixed data 1 (I	nitial data)				
15 to 18	Fixed data 2 (Copy it.)	Changed data. Read from same model and				
19	Fixed data 1 (I	nitial data)				
1A to 1E	Fixed data 1 (I	nitial data)				
1F	Fixed data 2 (C	hanged data. Read from same model and				
20 to 29	Fixed data 1 (I	nitial data)				
2A, 2B	Fixed data 1 (I	nitial data)				
2C to 2F	Fixed data 2 (Copy it.)	Fixed data 2 (Changed data. Read from same model and				
30 to 32	Fixed data 1 (I	nitial data)				
33	59	IC422 27 MHz XTAL fo adjustment				
34	19	Playback CR signal level adjustment/ Encoder R-Y input level adjustment				
35	37	Playback CB signal level adjustment/ Encoder B-Y input level adjustment				
36	18	Playback Y signal level adjustment/ Y output level adjustment				
37 to 39	Fixed data 1 (I	nitial data)				
3A to 3F	Fixed data 1 (I	nitial data)				
40	Fixed data 1 (I	nitial data)				
41	00	Playback burst level adjustment				
42	Fixed data 2 (Changed data. Read from same model and copy it.)					
43	Fixed data 1 (Initial data)					
44 to 46	Fixed data 2 (Changed data. Read from same model and copy it.)					
47 to 49	Fixed data 1 (Initial data)					
1 4A	Fixed data 2 (Changed data. Read from same model and copy it.)					
4B to 4F	Fixed data 1 (I	100 - 100 -				

Table 5-3-4.

8. Page E Address List

Note 1: Fixed data 1: Initialized data. (Refer to 4. Input of E Page Initial Data)
Fixed data 2: Modified data. (Refer to 5. Modification of C, D, E, Page Data)

Address	Initial Value	Remark
00 to 1B	Fixed data 1 (1	nitial data)
ici	Fixed data 2 (Changed data. Read from same model and copy it.)	
1D 1E 1F	8D 86 80	Battery down adjustment and check
20 to 23	Fixed data 1 (1	nitial data)
24,	Fixed data 2 (Changed data, Read from same model and copy it.)	
25 to 3F	Fixed data 1 (Initial data)	

Table 5-3-5.

3-4. SERVO SYSTEM ADJUSTMENTS

1. Switching Position Adjustment (CM-56 Board)

Mode	Playback
signal	SW/OL reference tape
Measurement Point	Page: 3, address: 03 on displayed data of adjusting remote commander
Measuring Instrument	Adjusting remote commander
Adjustment Page	С
Adjustment Address	4C, 4D, 4E, 4F
Specified Value	"00"

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- 2) Select page: 3, address: 01, set data: 0E, and press the PAUSE button of the adjusting remote commander.
- Select page: 3, address: 02, and confirm that the data changes from "0E" to "00".
- 4) Select page: 3, address: 03, and confirm that the data is "00".
- 5) Turn OFF the HOLD switch on he adjusting remote commander and wait for more than 2 seconds. (The adjusted data are automatically written to page: C, address: 4C to 4F)
- 6) Turn ON the HOLD switch on the adjusting remote commander.
- 7) Select page: 0, address: 01, and set data: 00.
- Stop the tape playback.
- 9) Turn the POWER switch OFF.

2. Capstan FG Duty Adjustment (CM-56 Board)

Mode	Playback
signal	Arbitrary tape
Measurement Point	Page: 3, address: 03 on displayed data of adjusting remote commander
Measuring Instrument	Adjusting remote commander
Adjustment Page	С
Adjustment Address	50, 51
Specified Value	"00"

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- Select page: 3, address: 01, set data: 15, and press the PAUSE button of the adjusting remote commander.
- 3) Select page: 3, address: 02, and confirm that the data changes from "15" to "00".
- 4) Select page: 3, address: 03, and confirm that the data is the following value

When "00": Normal

When "01": Faulty

Perform the following adjustment only when "00" is displayed.

5) Select page: 3, address: 04 and 05, read the data, and take the values as Do4 and Do5 respectively.

(The data on page: 3, address: 05 must be 2F to 3F)

- 6) Select page: C, address: 50, set data: Do4, and press the PAUSE button of the adjusting remote commander.
- 7) Select page: C, address: 51, set data: Dos, and press the PAUSE button of the adjusting remote commander.
- 8) Select page: 3, address: 01, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 9) Select page: 0, address: 01, and set data: 00.
- 10) Stop the tape playback.
- 11) Turn the POWER switch OFF.

3-5. VIDEO SYSTEM ADJUSTMENTS

3-5-1. RP-228 Board Adjustments

1. Recording Current Adjustment (RP-228 Board)

Mode	E-E
Measurement Point	ODDch adjustment CH1: Pin (a) of CN771 (CL812) CH2: Pin (b) of CN771 (CL813) EVENch adjustment CH1: Pin (a) of CN771 (CL816) CH1: Pin (a) of CN771 (CL815)
Measuring Instrument	Oscilloscope ADD mode CH2 INV mode
Adjustment Page	С
Adjustment Address	3E, 3F
Specified Value	$A = 4.1 \pm 0.1 \text{ Vp-p}$

Connection: Disconnect CN771 and connect as follows.

- ODDch adjustment: Connect a 180 Ω resistor between Pin 6 of CN771 (CL812) and Pin 5 of CN771 (CL813).
- EVENch adjustment: Connect a 180 Ω resistor between Pin
 of CN771 (CL816) and Pin (CL815)
 of CN771 (CL815)
 of CN771 (CL815)

Adjusting method:

- Equalize the vertical range of CH1 and CH2 of the oscilloscope.
- Set the oscilloscope to the ADD mode, and set CH2 to the INV mode.
- 3) Select page: 0, address: 01, and set data: 01.
- 4) Select page: 3, address: 01, set data: 0C, and press the PAUSE button of the adjusting remote commander.
- 5) Select page: 3, address: 34, and set data: 01.
- 6) Select page: C, address: 3F (ODDch adjustment) or 3E (EVENch adjustment), change the data, and adjust the signal voltage (A) to the specified value, press the PAUSE button on the adjusting remote commander.
- 7) Select page: 3, address: 34, and set data: 04.
- 8) Select page: 3, address: 01, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 9) Select page: 0, address: 01, and set data: 00.




Fig. 5-3-5.

PLL fo Adjustment (RP-228 Board)

Mode	E-E
Measurement Point Measuring Instrument	Displayed data of page: 3, address: 04
Adjustment Page	С
Adjustment Address	3D, 3C
Specified Value	Displayed data is "FD" to "FF", "00" to "03" ("FF", "00" are center values)

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- 2) Select page: 3, address: 01, set data: 05, and press the PAUSE button of the adjusting remote commander.
- 3) Select page: 3, address: 36, and set data: 04.
- 4) Select page: 3, address: 04, and check that the average value Do4 of the displayed data is "FD" to "FF" or "00" to "03". If outside this range, select page: C, address: 3C, change the data, and check again.

[If Do4 is "80" to "FC"]

Select page: C, address: 3C, and decrease the data. (As the data is to be rewritten, press the PAUSE button of the adjusting remote commander)

[If Do4 is "04" to "7F"]

Select page: C, address: 3C, and increase the data. (As the data is to be rewritten, press the PAUSE button of the adjusting remote commander)

- 5) Select page: 3, address: 36, and set data: 05,
- 6) Select page: 3, address: 04, and check that the average value D₀₄ of displayed data is "FD" to "FF" or "00" to "03". If outside this range, select page: C, address: 3D, change the data, and check again.

[If Do4 is "80" to "FC"]

Select page: C, address: 3D, and decrease the data. (As the data is to be rewritten, press the PAUSE button of the adjusting remote commander)

[If Do4 is "04" to "7F"]

Select page: C, address: 3D, and increase the data. (As the data is to be rewritten, press the PAUSE button of the adjusting remote commander)

- Select page: 3, address: 01, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 8) Select page: 3, address: 36, and set data: 02.
- 9) Select page: 0, address: 01, and set data: 00.

3. CLK DELAY Adjustment (RP-228 Board)

Mode	Recording/playback
Mode	Recording/playback
Signal	Color bar
Measurement Point	CH1: Pin [®] of CN775 (C1ERP) CH2: Pin [©] of CN775 (JSWP)
Measuring Instrument	Oscilloscope Trigger source: CH2
Adjustment Page	С
Adjustment Address	47

Adjusting method:

- 1) Record color bar signal for two minutes on any tape.
- 2) Select page: 0, address: 01, and set data: 01.
- 3) Write the following data in page: C, address: 40 to 43, 47, 4B, 5A

To write the data, press the PAUSE button of the adjusting remote commander each time data is set.

Page: C, address: 40, data: C0

Page: C, address: 41, data: C0

Page: C, address: 42, data: 90

Page: C, address: 43, data: 90

Page: C, address: 47, data: C8

Page: C, address: 4B, data: 80

Page: C, address: 5A, data: 00

4) Playback the part recorded with the color bar.

- Select page: C, address: 47, increase the data, and read the data D1 when the CH1 pulse is set to the whole audio and video areas.
- Select page: C, address: 47, decrease the data, and read the data D2 when the CH1 pulse is set to the whole audio and video areas.
- 7) Obtain the average value of D1 and D2, and take it as D3.
- Select page: C, address: 47, set data: D3, and press the PAUSE button of the adjusting remote commander.
- Select page: C, address: 4B, set data: 0E, and press the PAUSE button of the adjusting remote commander.
- Select page: C, address: 5A, set data: 8C, and press the PAUSE button of the adjusting remote commander.
- 11) Select page: 0, address: 01, and set data: 00.
- 12) After completing the adjusting, perform "5. AEQ Adjustment".




Fig. 5-3-6.

4. AGC Center Level Adjustment (RP-228 Board)

Mode	Recording/playback
Signal	Color bar
Measurement Point	CH1: Pin [®] of CN775 (C1ERP) CH2: Pin [®] of CN775 (JSWP)
Measuring Instrument	Oscilloscope Trigger source: CH2
Adjustment Page	С
Adjustment Address	44

Adjusting method:

- 1) Record color bar signal for two minutes on any tape.
- 2) Select page: 0, address: 01, and set data: 01.
- 3) Write the following data in page: C, addresses: 40 to 44, 4B, 5A.

To write the data, press the PAUSE button of the adjusting remote commander each time data is set.

Page: C, address: 40, data: C0

Page: C, address: 41, data: C0

Page: C, address: 42, data: 90

Page: C, address: 43, data: 90

Page: C, address: 44, data: 90

Page: C, address: 4B, data: 80

Page: C, address: 5A, data: 00

- 4) Playback the part recorded with the color bar signal.
- Select page: C, address: 44, increase the data, and read the data D1 when the CH1 pulse is set to the whole audio and video areas.
- Select page: C, address: 44, decrease data, and read the data D₂ when the CH1 pulse is set to the whole audio and video areas.
- 7) Obtain the average value of D₁ and D₂, and take it as D₃.
- Select page: C, address: 44, set data: D3, and press the PAUSE button of the adjusting remote commander.
- Select page: C, address: 4B, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 10) Select page: C, address: 5A, set data: 8C, and press the PAUSE button of the adjusting remote commander.
- 11) Select page: 0, address: 01, and set data: 00.
- 12) After completing the adjusting, perform "5. AEQ Adjustment".




Fig. 5-3-7.

5. AEQ Adjustment (RP-228 Board)

Mode	Recording/playback
Signal	Arbitrary
Measurement Point	Pin ® of CN775 (RF MONITOR) (Note 1)
Measuring Instrument	Oscilloscope
Adjustment Page	С
Adjustment Address	40, 41, 42, 43, 5A

Note 1: Connect a 75 Ω resistor between Pin (3) and (7) (GND) of CN 775.

75 Ω resistor (Parts code: 1-247-804-11)

Note 2: Use the DVM60ME tape or equivalents.

Adjusting method:

1) Select page: 0, address: 01, and set data: 01.

Select page: C, address: 4B, set data: 80, and press the PAUSE button of the adjusting remote commander.

3) Write data in page: C, addresses: 40 to 43, and 5A as shown in the following table.

To write the data, press the PAUSE button of the adjusting remote commander each time data is set.

Address	Data
40	C0
41	C0
42	90
43	90
5A	00

- 4) Record color bar signal for one minute from the tape top.
- 5) Rewind the tape, and play back the part recorded.
- 6) When the RF output stabilizes, select page: 3, address: 01, and set data: 07, and press the PAUSE button of the adjusting remote commander.
- About 20 to 30 seconds after pressing the PAUSE button, select page: 3, address: 02, and check that the data changes from "07" to "00".
- Select page: 3, address: 03, and check that the data is the following value.

When "00": Normal

When "01": EVENch is faulty When "02": ODDch is faulty

When "03": EVENch and ODDch are faulty

Perform the following procedure only when "00" is displayed.

9) Select page: 3, address: 04 to 07, read the data, and take the values as Do4, Do5, Do6, and Do7.

- 10) Select page: C, address: 40, set data: Do4, and press the PAUSE button of the adjusting remote commander.
- 11) Select page: C, address: 42, set data: Dos, and press the PAUSE button of the adjusting remote commander.
- 12) Select page: C, address: 41, set data: Do6, and press the PAUSE button of the adjusting remote commander.
- 13) Select page: C, address: 43, set data: Do7, and press the PAUSE button of the adjusting remote commander.
- 14) Select page: C, address: 5A, set data: 8C, and press the PAUSE button of the adjusting remote commander.
- 15) Select page: C, address: 4B, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 16) Select page: 0, address: 01, and set data: 00.

6. PLL Capture Range Adjustment (RP-228 Board)

Mode	Recording/Playback
Signal	Color bar
Measurement Point	CH1: Pin [®] of CN775 (C1ERP) CH2: Pin [®] of CN775 (JSWP)
Measuring Instrument	Oscilloscope Trigger source: CH2
Adjustment Page	С
Adjustment Address	46

Adjusting method:

- 1) Record color bar signal for two minutes on any tape.
- 2) Select page: 0, address: 01, and set data: 01.
- 3) Write the following data in page: C, addresses: 4B and 5A.

 To write the data, press the PAUSE button of the adjusting remote commander each time data is set.

Page: C, address: 4B, data: 80

Page: C, address: 5A, data: 00

- 4) Playback the part recorded with the color bar signal.
- 5) Select page: C, address: 46, set data: 80, and press the PAUSE button of the adjusting remote commander.
- 6) Select page: C, address: 46, set the data to "60", and check that the pulse is not set at the audio area head of the ERRP waveform's ODDch of the oscilloscope (CH1).
- Select page: C, address: 46, set the data to "A0", and check that the pulse is not set at the audio area head of the C1ERP waveform's ODDch of the oscilloscope (CH1).

After confirming steps 6) and 7), select page: C, address: 46, set data: 80 again and proceed to step 12).

- 8) If the pulse is set in steps 6) and 7), select page: C, address: 46, increase the data from "80", and read the data D1 when the pulse is set at the audio area head of CH1.
- 9) Select page: C, address: 46, decrease the data from "80", and read the data D₂ when the pulse is set at the audio area head of CH1.
- 10) Obtain the average value of D1 and D2, and take it as D3.
- 11) Select page: C, address: 46, set data: D₃, and press the PAUSE button of the adjusting remote commander.
- 12) Select page: C, address: 4B, set data: 00, and press the PAUSE button of the adjusting remote commander.
- 13) Select page: C, address: 5A, set data: 8C, and press the PAUSE button of the adjusting remote commander.
- 14) Select page: 0, address: 01, and set data: 00.

When the pulse is not set at the audio area head.

CH1

CH2

When the pulse is set at the audio area head.

CH1

CH2

ODD ch

3.33 msec

Fig. 5-3-8.

7. IC774 41.85 MHz VCO Check (RP-228 Board)

Mode	E-E
Signal	Arbitrary
Measurement Point	Page: 3, address: 39 on displayed data of adjusting remote commander
Measuring Instrument	
Adjustment Value	"37" to "C9" (0.6 to 2.2 Vdc)

Checking method:

 Select page: 3, address: 39, and check that the displayed data is "37" to "C9".

3-5-2. JC-19 Board Adjustments

A/D Converter Reference Voltage Adjustment 1 (JC-19 Board)

Mode	E-E
Signal	Arbitrary
Measurement Point	Pin (5) of IC013 (CL061)
Measuring Instrument	Digital voltmeter
Adjusting Element	RV001
Specified Value	$A = 2.83 \pm 0.01 \text{ Vdc}$

Adjusting method:

1) Set the VRT voltage (A) to the specified value using RV001.

2. A/D Converter Reference Voltage Adjustment 2 (JC-19 Board)

Mode	E-E
Signal	Arbitrary
Measurement Point	Pin ③ of IC013 (CL062)
Measuring Instrument	Digital voltmeter
Adjusting Element	RV002
Specified Value	$A = 0.96 \pm 0.01 \text{ Vdc}$

Adjusting method:

1) Set the VBT voltage (A) to the specified value using RV002.

3. Y Signal Clamp Reference Voltage Adjustment (JC-19 Board)

Mode	E-E
Signal	Color bar
Measurement Point	Pin ® of IC011 (CL054)
Measuring Instrument	Digital voltmeter
Adjusting Element	RV011
Specified Value	$A = 1.150 \pm 0.005 \mathrm{Vdc}$

Connection: Connect a jumper wire between Pin ⑥ of IC018 (CL150) and GND.

Adjusting method:

1) Set the Y signal clamp reference voltage (A) to the specified value using RV011.

4. CR Signal Clamp Reference Voltage Adjustment (JC-19 Board)

Mode	E-B
Signal	Color bar
Measurement Point	Pin (8) of IC010 (CL052)
Measuring Instrument	Digital voltmeter
Adjusting Element	RV010
Specified Value	$A = 1.900 \pm 0.005 \mathrm{Vdc}$

Connection: Connect a jumper wire between Pin **(6)** of IC018 (CL150) and GND.

Adjusting method:

1) Set the CR signal clamp reference voltage (A) to the specified value using RV010.

5. CB Signal Clamp Reference Voltage Adjustment (JC-19 Board)

Mode	E-E
Signal	Color bar
Measurement Point	Pin ® of IC009 (CL053)
Measuring Instrument	Digital voltmeter
Adjusting Element	RV012
Specified Value	$A = 1.900 \pm 0.005 \text{Vdc}$

Connection: Connect a jumper wire between Pin **(6)** of IC018 (CL150) and GND.

Adjusting method:

1) Set the CB signal clamp reference voltage (A) to the specified value using RV012.

6. Playback Y Signal Level Adjustment (JC-19 Board)

Mode	Recording
Signal	DV input (Note 1)
Measurement Point	Pin (a) of CN104 or pin (a) of CN102 on VA-102 board
Measuring Instrument	Oscilloscope
Adjustment Page	D
Adjustment Address	36
Specified Value	$A = 0.43 \pm 0.04 \text{ V (NTSC)}$ $A = 0.41 \pm 0.04 \text{ V (PAL)}$

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, select page: 5, address: 02, and set data: 09. After adjustment, be sure to return the data to "00")

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- 2) Select page: D, address: 36, change data, and adjust the Y signal level (A) to the specified value.
- 3) Select page: 0, address: 01, and set data: 00.




Fig. 5-3-9.


7. Playback CR Signal Level Adjustment (JC-19 Board)

Mode	Recording
Signal	DV input (Note 1)
Measurement Point	Pin ⁽¹⁾ of CN104 or pin ⁽²⁾ of CN102 on VA-102 board
Measuring Instrument	Oscilloscope
Adjustment Page	D
Adjustment Address	34
Specified Value	$A = 540 \pm 10 \text{ mV}$

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, select page: 5, address: 02, and set data: 09. After adjustment, be sure to return the data to "00")

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- 2) Select page: D, address: 34, change data, and adjust the CR signal level (A) to the specified value.
- 3) Select page: 0, address: 01, and set data: 00.



5-72

8. Playback CB Signal Level Adjustment (JC-19 Board)


Mode	Recording
Signal	DV input (Note 1)
Measurement Point	Pin ② of CN104 or pin ③ of CN102 on VA-102 board
Measuring Instrument	Oscilloscope
Adjustment Page	D
Adjustment Address	35
Specified Value	$A = 390 \pm 10 \text{ mV}$

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, select page: 5, address: 02, and set data: 09. After adjustment, be sure to return the data to "00")

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 00.
- 2) Select page: D, address: 35, change data, and adjust the CB signal level (A) to the specified value.
- 3) Select page: 0, address: 01, and set data: 00.

For NTSC model



For PAL model




Fig. 5-3-11.

9. IC422 27MHz XTAL fo Adjustment (JC-19 Board)

Mode	Playback
Signal	Arbitrary tape
Measurement Point	Pin @ of IC442 (CL479)
Measuring Instrument	Frequency counter
Adjustment Page	D
Adjustment Address	33
Specified Value	f = 13500000 ± 100 Hz

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- Select page: D, address: 33, change data, and adjust the clock frequency (f) to the specified value.
- 3) Press the PAUSE button on the adjusting remote commander.
- 4) Select page: 0, address: 01, and set data: 00.

10. AFC Preliminary Adjustment (JC-19 Board)

	Mode	Recording
I	Signal	Color bar
	Measurement Point	Pin (9) of IC205 (CL214)
	Measuring Instrument	Digital voltmeter
	Adjusting Element	CT201
I	Specified Value	$A = 1.9 \pm 0.5 \text{ Vdc}$

Adjusting method:

1) Set the DC voltage (A) to the specified value using CT201.

11. AFC Picture Frame Adjustment (JC-19 Board)

Mode	Recording
Signal	Color bar (Video input) (Note 1)
Measurement Point	CH1: Pin ② of IC017 (CL051) CH2: Pin ② of IC205 (CL222)
Measuring Instrument	Oscilloscope
Adjusting Element	RV201
Specified Value	$T = 110 \pm 10 \text{ nsec}$

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

1) Set the time difference (T) between the H SYNC falling and AFH rising to the specified value using RV201.




Fig. 5-3-12.

12. AFC Adjustment (JC-19 Board)

Mode	Recording	
Signal	Color bar	
Measurement Point	Pin 9 of IC205 (CL214)	
Measuring Instrument	Digital voltmeter	
Adjusting Element	CT201	
Specified Value	$A = 1.80 \pm 0.05 \text{Vdc}$	

Adjusting method:

1) Set the DC voltage (A) to the specified value using CT201.

3-5-3. VA-102 Board Adjustments

1. AGC Adjustment (VA-102 Board)


Mode	E-E
Signal	Color bar (VIDEO input) (Note 1)
Measurement Point	Pin 1 of IC205 (CL220)
Measuring Instrument	Oscilloscope
Adjustment Element	RV202
Specified Value	A = 1.428 ± 0.02 V (NTSC) A = 1.400 ± 0.02 V (PAL)

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.


Adjusting method:

1) Set the Y signal level (A) to the specified value using RV202.

For NTSC model



For PAL model



2. Analog E-E VIDEO Signal Output Level Adjustment (VA-102 Board)

Mode	E-E
Signal	Color bar (S VIDEO input) (Note 1)
Measurement Point	Pin (9) of IC401 (CL434)
Measuring Instrument	Oscilloscope
Adjustment Element	RV401, RV404
Specified Value	A = 0.714 ± 0.01 V (NTSC) A = 0.700 ± 0.01 V (PAL) B = 280 ± 10 mV (NTSC) B = 300 ± 10 mV (PAL)

Note 1: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.


Note 2: Terminate the video output terminal using a 75 Ω resistor.

75 Ω resistor (Parts code: 1-247-804-11)

Adjusting method:

- 1) Set the Y signal level (A) to the specified value using RV401.
- 2) Set the burst signal level (B) to the specified value using RV404.

For NTSC model






Fig. 5-3-14.

3. Analog E-E Y Signal Output Level Check (VA-102 Board)

Mode	E-E
Signal	Color bar (S VIDEO input) (Note 1)
Measurement Point	Pin (19) of CN401 (CL436)
Measuring Instrument	Oscilloscope
Specified Value	$A = 0.714 \pm 0.02 \text{ V (NTSC)}$ $A = 0.700 \pm 0.02 \text{ V (PAL)}$

Note 1: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.


Note 2: Terminate the Y signal terminal of the S video output terminal using a 75 Ω resistor.

75 Ω resistor (Parts code: 1-247-804-11)

Checking method:

1) Check that the Y signal level (A) is the specified value.

For NTSC model



For PAL model




Fig. 5-3-15.

4. Analog E-E Chroma Signal Output Level Check (VA-102 Board)

Mode	E-E
Signal	Color bar (S VIDEO input) (Note 1)
Measurement Point	Pin (1) of IC401 (CL435)
Measuring Instrument	Oscilloscope
Specified Value	A = 286 ± 20 mV (NTSC) A = 300 ± 20 mV (PAL)


Note 1: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.

Note 2: Terminate the Chroma signal terminal of the S video output terminal using a 75 Ω resistor. 75 Ω resistor (Parts code: 1-247-804-11)

Checking method:

1) Check that the burst signal level (A) is the specified value.

For NTSC model






Fig. 5-3-16.

5. Decoder VXO Freerunning Frequency Adjustment (VA-102 Board)

Mode	E-E
Signal	Color bar (S VIDEO input) (Note 2) (Chroma signal OFF)
Measurement Point	TP201 (CL210)
Measuring Instrument	Frequency coutner
Adjustment Element	CT201
Specified Value	f = 3579545 ± 30 Hz (NTSC) f = 4433618 ± 20 Hz (PAL)

Note 1: Connect the frequency counter via high input impedance equipment such as an oscilloscope.

Note 2: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

1) Set the VXO OSC frequency (f) to the specified value using CT201.

6. Video Input Y/C Separation Adjustment (VA-102 Board)

(1) Y Signal Output Level Adjustment

Mode	E-E
Signal	Color bar (VIDEO input) (Note 1)
Measurement Point	Pin ① of IC202 (CL202)
Measuring Instrument	Oscilloscope
Adjustment Element	RV203
Specified Value	A= 0.714 ± 0.01 V (NTSC) A= 0.700 ± 0.01 V (PAL)

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.


Note 2: Terminate the video output terminal using a 75 Ω resistor.

75 Ω resistor (Parts code: 1-247-804-11)

Adjusting method:

1) Set the Y signal level (A) to the specified value using RV203.

For NTSC model






Fig. 5-3-17.

(2) Chroma Signal Output Level Adjustment

Mode	E-E
Signal	Color bar (VIDEO input) (Note 1)
Measurement Point	Emitter of Q202 (CL203)
Measuring Instrument	Oscilloscope
Adjustment Element	RV201
Specified Value	A = 286 ± 10 mV (NTSC) A = 300 ± 10 mV (PAL)

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.


Note 2: Terminate the video output terminal using a 75 Ω resistor.

75 Ω resistor (Parts code: 1-247-804-11)

Adjusting method:

1) Set the burst signal level (A) to the specified value using RV201.

For NTSC model



For PAL model




Fig. 5-3-18.

7. Decoder HUE Adjustment (VA-102 Board)


Mode	E-E
Signal	Color bar (S VIDEO input) (Note 1)
Measurement Point	Pin (1) of IC102 (CL144)
Measuring Instrument	Oscilloscope
Adjustment Element	RV207
Specified Value	A = B = C

Note 1: Set data: 00 to page: 5, address: 38

Note 2: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

1) Set the amplitude (A), (B), (C) to the same level using RV207.






Fig. 5-3-19.

8. RECY Level Adjustment (VA-102 Board)


Mode	E-E
Signal	Color bar (VIDEO input) (Note 1)
Measurement Point	Pin 1 of IC102 (CL142)
Measuring Instrument	Oscilloscope
Adjustment Element	RV205
Specified Value	$A = 1.55 \pm 0.02 \text{ V}$

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

1) Set the Y signal level (A) to the specified value using RV205.

For NTSC model



For PAL model




Fig. 5-3-20.

9. REC CR Level Adjustment (VA-102 Board)

Mode	E-E
Signal	Color bar (VIDEO input) (Note 1)
Measurement Point	Pin 1 of IC102 (CL143)
Measuring Instrument	Oscilloscope
Adjustment Element	RV204
Specified Value	$A = 1.25 \pm 0.02 \text{ V (NTSC)}$ $A = 1.20 \pm 0.02 \text{ V (PAL)}$

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

1) Set the CR signal level (A) to the specified value using RV204.




Fig. 5-3-21.

10. REC CB Level Adjustment (VA-102 Board)

Mode	E-E
Signal	Color bar (VIDEO input) (Note 1)
Measurement Point	Pin (3) of IC102 (CL144)
Measuring Instrument	Oscilloscope
Adjustment Element	RV206
Specified Value	$A = 1.20 \pm 0.02 \text{ V}$

Note 1: Set "VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

1) Set the CB signal level (A) to the specified value using RV206.




Fig. 5-3-22.

11. Encoder Freerunning Frequency Adjustment (VA-102 Board)

Mode	E-E
Signal	Color bar (S VIDEO input) (Note 2)
Measurement Point	TP401 (CL426)
Measuring Instrument	Frequency coutner
Adjustment Element	CT401
Specified Value	f = 14318182 ± 100Hz (NTSC) f = 17734475 ± 100 Hz (PAL)

Note 1: Connect the frequency counter via high input impedance equipment such as an oscilloscope.

Note 2: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

 Set the oscillation frequency (f) to the specified value using CT401.

Measurement Point




Fig. 5-3-23.

12. Playback Y Level Check (VA-102 Board)

Mode	Recording
Signal	DV input (Note 1)
Measurement Point	Pin (3) of CN401 (CL436)
Measuring Instrument	Oscilloscope
Specified Value	$A = 0.83 \pm 0.02 \text{ V (NTSC)}$ $A = 0.825 \pm 0.02 \text{ V (PAL)}$

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, and set data: 09 to page: 5, address: 02. After adjustment, be sure to return the data to "00".)

Note 2: Perform this check after confirming that the specified value in the following adjustment of the JC-19 board has been satisfied.

1. Playback Y Signal Level Adjustment.

Checking method:

1) Check that the white (75%) signal level (A) is the specified value.




Fig. 5-3-24.

13. Playback Chroma Level Adjustment (VA-102 Board)

Mode	Recording
Signal	DV input (Note 1)
Measurement Point	Pin ① of CN401 (CL435)
Measuring Instrument	Oscilloscope
Adjustment Element	RV406
Specified Value	$A = 670 \pm 10 \text{ mV (NTSC)}$ $A = 660 \pm 10 \text{ mV (PAL)}$

- Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, and set data: 09 to page: 5, address: 02. After adjustment, be sure to return the data to "00".)
- **Note 2:** Perform this adjustment after confirming that the specified value in the following adjustment of the JC-19 board has been satisfied.
 - 1. Playback CR Signal Level Adjustment.
 - 2. Playback CB Signal Level Adjustment.

Adjusting method:

1) Set the red signal level (A) to the specified value using RV406.




Fig. 5-3-25.

14. Playback Burst Level Adjustment (VA-102 Board) (PAL model only)

Mode	Recording
Signal	DV input (Note 1)
Measurement Point	Pin ① of CN401 (CL435)
Measuring Instrument	Oscilloscope
Adjustment Element	RV402
Specified Value	$A = 300 \pm 10 \text{ mV (PAL)}$

- Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, and set data: 09 to page: 5, address: 02. After adjustment, be sure to return the data to "00".)
- Note 2: Perform this adjustment after confirming that the specified value in the following adjustment of the JC-19 board has been satisfied.
 - 1. Playback CR Signal Level Adjustment.
 - 2. Playback CB Signal Level Adjustment.

Adjusting method:

1) Set the burst signal level (A) to the specified value using RV402.




Fig. 5-3-26.

3-5-4. General Adjustments

1. Y Output Level Adjustment (JC-19 Board)

Mode	E-E
Signal	Color bar (DV input) (Note 1)
Measurement Point	Video output terminal
Measuring Instrument	Oscilloscope
Adjustment Element	page: D, address: 36
Specified Value	$A = 0.83 \pm 0.01 \text{ V (NTSC)}$ $A = 0.823 \pm 0.01 \text{ V (PAL)}$

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, and set data: 09 to page: 5, address: 02. After adjustment, be sure to return the data to "00".

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- Select page: D, address: 36, change the data, adjust the Y signal level (A) to the specified value.
- 1) Select page: 0, address: 01, and set data: 00.




Fig. 5-3-27.

2. Encoder R-Y Input Level Adjustment (JC-19 Board)

Mode	E-E
Signal	Color bar (DV input) (Note 1)
Measurement Point	Video output terminal
Measuring Instrument	Vectorscope
Adjustment Element	page: D, address: 34
Specified Value	Phase: 104 ± 2° Gain: 95 ± 5%

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, and set data: 09 to page: 5, address:02. After adjustment, be sure to return the data to "00".

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- Adjust the burst luminance point to the specified position using the PHASE and GAIN knobs of the vectorscope.
- Select page: D, address: 34, change the data, adjust a red luminance point to the specified position.
- 4) Select page: 0, address: 01, and set data: 00.

For NTSC model

⊞: FOR ENCODER R-Y INPUT LEVEL ADJUSTMENT




Fig. 5-3-28.

For PAL model

⊞: FOR ENCODER R-Y INPUT LEVEL ADJUSTMENT




Fig. 5-3-29.

3. Encoder B-Y Input Level Adjustment (JC-19 Board)

Mode	E-E	
Signal	Color bar (DV input) (Note 1)	_
Measurement Point	Video output terminal	_
Measuring Instrument	Vectorscope	_
Adjustment Element	page: D, address: 35	٦
Specified Value	Phase: 348 ± 2° Gain: 66 ± 5%	

Note 1: Generate color bar signal with NTSC: DCR-VX1000/PAL: DCR-VX1000E and enter it to the DV terminal. (How to generate color bars: Connect the adjusting remote commander to the NTSC: DCR-VX1000/PAL: DCR-VX1000E, and set data: 09 to page: 5, address:02. After adjustment, be sure to return the data to "00".

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- Adjust the burst luminance point to the specified position using the PHASE and GAIN knobs of the vectorscope.
- 3) Select page: D, address: 35, change the data, adjust a blue luminance point to the specified position (inside of thick frame).
- 4) Select page: 0, address: 01, and set data: 00.

4. Decoder HUE Input Adjustment (JC-19 Board)

Mode	E-E
Signal	Color bar (S VIDEO input) (Note 1)
Measurement Point	Video output terminal
Measuring Instrument	Vectorscope
Adjustment Element	RV207
Specified Value	Phase: 104 ± 1° Gain: 95 ± 5%

Note 1: Set "S VIDEO" mode with the INPUT SELECT button on the front panel.

Adjusting method:

- 1) Select page: 0, address: 01, and set data: 01.
- 2) Select page: 2, address: 10, and set data: 20.
- Adjust the burst luminance point to the specified position using the PHASE and GAIN knobs of the vectorscope.
- 4) Adjust RV107 so that a red luminance point comes to the specified position (inside of thick frame). At this time, confirm that other color luminance points are inside each phase specified frame (± 2).
- 5) Select page: 2, address: 10, and set data: 00.
- 4) Select page: 0, address: 01, and set data: 00.

Note: When a red luminance point and a blue luminance point are not at the specified positions, adjust RV204 and RV206 so that they come to the specified positions respectively.

For NTSC model

Burst

G

FOR ENCODER

B-Y INPUT LEVEL

ADJUSTMENT

Fig. 5-3-30.

For PAL model

Fig. 5-3-31.

5. Battery Down Adjustment and Confirmation

Mode	Stop
Signal	Any
Test point	Displayed data on page 5, address 2A (LCD display of the adjusting remote commander)
Measuring Instrument	Adjusting remote commander
Adjustment page	E
Adjustment address	1D, 1E, 1F

Note 1: Make sure that the BEEP on the Menu screen is set to "ON".

Connection of Equipment:

Connect a regulated power supply and a digital voltmeter to the DC IN terminal.




Fig. 5-3-32.

Adjusting Method:

- 1) Adjust the output voltage of regulated power supply so that a digital voltmeter displays 11.0 ± 0.05 V.
- 2) Select page: 5, address: 2A, and read displayed data on the adjuting remote commander, and assume it as Do.
- 3) Select page: E, address: 1D, set data: Do, and press the PAUSE button on the adjusting remote commander.
- 4) Convert "Do" read in 2) to decimal notation, and obtain Do'.
- 5) Calculate D₁', D₂' using the following equations (decimal notation calculation).

$$D_1' = D_0' - 6$$

 $D_2' = D_0' - 12$

- 6) Convert D1' to hexadecimal notation, and obtain D1.
- 7) Select page: E, address: 1E, set data: D₁, and press the PAUSE button on the adjusting remote commander.
- 8) Convert D2' to hexadecimal notation, and obtain D2.
- 9) Select page: E, address: 1F, set data: D2, and press the PAUSE button on the adjusting remote commander.
- 10) Confirm the display and operation of the fluorescent display tube, when the voltage input from DC IN terminal is lowered from 12 V.

Input of DC 12 V : Normal operation

Input of DC 10.8 V: Beep sounds and "dc Lo" is displayed.

Input of DC 10.3 V: Beep sounds and the standby mode is activated after 2 seconds.

11) Further lower the voltage, and check the voltage when STBY indicator (red LED) turns off.

Specification: The power relay must turns off when Vdc = 9.0 V - 9.5 V.

12) On the contrary, raise the voltage, and check the voltage when STBY indicator (red LED) turns on.

Specification: The power relay must turns on when Vdc = 10.5 V - 11.0 V.

3-5-5. BIST Check

1. Playback System Check (JC-19, RP-228 Boards)

- Connect the adjusting remote commander to the LANC terminal, and turn the HOLD switch ON.
- 2) Playback the BIST check tape.

IC411(D1) Playback System Check

- Select page: 4, address: 11, set data: 04, and press the PAUSE button.
- Select page: 4, address: 11, set data: 00, and press the PAUSE button.
- Select page: 4, address: 13, set data: 03, and press the PAUSE button.
 (Data automatically returns to "00")
- 6) If IC411 (D1) → IC401 (U1) playback system is normal, the following data are displayed on page: 4, addresses: 14, 15.

Page	Address	Data
4	15	E5
4	14	. 11

7) If IC411(D1) → IC701 (IND1) playback system is normal, the following data are displayed on page: 4, addresses: 16, 17.

Page	Address	Data
4	17	C0 or BA
4	16	6E or 04

 If IC411(D1) → IC805 (A1) playback system is normal, the following data are displayed on page: 4, addresses: 18, 19.

Page	Address	Data
4	19	33 or B2
4	18	59 or 19

IC805 (A1) Playback System Check

- 9) Select page: 4, address: 11, set data: 10, and press the PAUSE button
- 10) Select page: 4, address: 11, set data: 00, and press the PAUSE button.
- 11) Select page: 4, address: 13, set data: 04, and press the PAUSE button.

(Data automatically returns to "00")

12) If IC805 (A1) playback system is normal, the following data are displayed on page: 4, addresses: 14, 15.

Page	Address	Data	_
4	15	7B	
4	14	B5	

IC401 (U1) Playback System Check

- 13) Select page: 4, address: 11, set data: 08, and press the PAUSE button.
- 14) Select page: 4, address: 42, set data: 01, and press the PAUSE button.
- 15) Select page: 4, address: 13, set data: 07, and press the PAUSE button.

(Data automatically returns to "00")

- 16) Select page: 4, address: 42, set data: 00, and press the PAUSE button.
- 17) Select page: 4, address: 11, set data: 00, and press the PAUSE button.

18) If IC401 (U1) → IC200 (S1) playback system is normal, the following data are displayed on page: 4, addresses: 14, 15.

Page	Address	Data	
4	15	1E	
4	14	F2	

19) If IC411 (D1) → IC401 (U1) playback system is normal, the following data are displayed on page: 4, addresses: 16, 17.

Page	Address	Data	
4	17	D1	
4	16	61	

20) Perform "Record System Check" successively.

Record System Check

Note: Perform "Record System Check" successively (with BIST check tape in playback status)

1) Enter the following data.

Note: Press the PAUSE button each time the data is set.

Page	Address	Data
4	41	01
4	0F	02
4	0E	01
4	40	01
4	0F	0 A
4	40	00
4	40	01
4	0F	0E
4	40	00
4	40	01
4	0F	8E
4	40	00

- With the HOLD switch on adjusting remote commander turned ON, eject the BIST check tape, and insert a record tape instead.
- 3) Set the REC mode.

IC401 (U1) Record System Check

- 4) Select page: 4, address: 11, set data: 08, and press the PAUSE
- Select page: 4, address: 42, set data: 01, and press the PAUSE
- Select page: 4, address: 13, set data: 07, and press the PAUSE button. (Data automatically returns to "00".)
- Select page: 4, address: 42, set data: 00, and press the PAUSE button.
- Select page: 4, address: 11, set data: 00, and press the PAUSE
- 9) If IC401 (U1) \rightarrow IC411 (D1) record system is normal, the following data are displayed on page: 4, addresses: 16, 17.

Page	Address	Data
4	17	. 05
4	16	80

IC411 (D1) Record System Check

- 10) Select page: 3, address: 01, set data: 0D, and press the PAUSE button.
- 11) Select page: 4, address: 1C, set data: FF, and press the PAUSE button.
- 12) Select page: 4, address: 11, set data: 04, , and press the PAUSE
- 13) Select page: 4, address: 11, set data: 00, and press the PAUSE
- 14) Select page: 4, address: 13, set data: 03, and press the PAUSE button. (Data automatically returns to "00")
- 15) If IC401 (U1) → IC411 (D1) record system is normal, the following data are displayed on page: 4, addresses: 14, 15.

Page	Address	Data
4	15	05
4	14	80

16) If IC411 (D1) → IC701 (IND1) record system is normal, the following data are displayed on page: 4, addresses: 16, 17.

Page	Address	Data
4	17	E6
4	16	BC

17) If IC805 (A1) \rightarrow IC411 (D1) record system is normal, the following data are displayed on page: 4, addresses: 18, 19.

Page	Address	Data
4	19	76
4	18	В9

18) If IC411 (D1) → IC774 (DX) record system is normal, the following data are displayed on page: 4, addresses: 1A, 1B.

Page	Address	Data
4	1B	4E
4	1A	11

3-6. AUDIO SYSTEM ADJUSTMENTS

Unless spesified otherwise, set the switches as follows. AUDIO NODE (Menu display) Fs48k INPUT SELECTVIDEO AUDIO MONITOR CH-1/2 REC LEVEL L, R Center Note 1: Set AUDIO MODE at the SET UP menu of the menu

Connection of Equipment




Fig. 5-3-33.

1. Playback Level/Indicator Check

Mode	Playback
Signal	Audio check reference tape
Measurement point	Audio output terminal (Left and Right)
Measuring Instrument	Audio level meter and frequency counter
Specified Value	32 kHz mode: 1 kHz signal should be output 48 kHz mode: 1 kHz signal level should be +6 ± 2 dBv (+8.2 ± 2 dBs) 44.1 kHz mode EMP ON: 7.35 kHz signal level is -6 ± 1 dB for 1 kHz signal level in 48 kHz mode 44.1 kHz mode EMP ON: 7.35 kHz signal level is 0 ± 1 dB for 1 kHz signal level in 48 kHz mode NS AUDIO lamp should be lit

Note: 0 dBv = 1 Vrms0 dBs = 0.775 Vrms

Checking method:

1) Check that the playback signal level satisfies the specified value.

2. E-E Level Check

Mode	E-E (LINE 1 input)					
Signal	Audio: 1 kHz -6 dBv (-3.8 dBs) Signal Audio input terminal (Left and Right)					
Signal	Video: Color bar signal Video input terminal					
Measurement point	Audio output terminal (Left and Right)					
Measuring Instrument	Audio level meter					
Specified Value	-6 ± 3 dBv (-3.8 ± 3 dBs)					

Checking method:

- 1) Check that the 1 kHz signal level satisfies the specified value.
- 2) Check that the number in the segment of the level meter (fluorescent display tube) that is lit is between 8 and 12 for both the L and R channels.

3. Recording/Playback Check (Audio Lock Mode)

Mode	Recording/Playback (LINE input)				
Signal	Audio: no signal				
	Video: Color bar Video input terminal				
Measurement point	Display data of page: 5, address: 00				
Measuring Instrument	of the adjusting remote commander				
Specified Value	① After playback pause, the changes in the data after 5 frames have been sent continuously must be in the following order. "D4"→"D6"→"D6"→"D6"→"D6"→"D6" →"D4" (NTSC) "D8"→"D8"→"D8"→"D8"→"D8" →"D8" (PAL) ② NS AUDIO lamp should be lit.				


Note 1: Check that the AUDIO MODE (menu screen) is Fs48k. Note 2: Send the frames using front panel button.

Checking method:


- 1) With no audio signal being input, record the color bar signal for about 1 minute.
- 2) Playback the recorded part, and set the playback pause mode.
- 3) Select page: 5, address: 00 using the adjusting remote commander.
- 4) Send the frames, so that the display data for page: 5, address: 00 is D4. (NTSC)
- 5) Send 5 frames continuously, and check that the display data of page: 5, address: 00 changes in the order specified.
- 6) Exit the playback pause mode, playback the recorded part, and check that the NS AUDIO lamp (front panel) is off.

3-7. ARRANGEMENT DIAGRAM FOR ADJUSTMENT PARTS


RP-228 BOARD (SIDE A)




RP-228 BOARD (SIDE B)




POWER BLOCK (U-2)





JC-19 BOARD (SIDE A)






C-19 BOARD (SIDE A) CN104 CN103 CN411 CN412 30 1 30 1 30 1 30 1 RV012 RV010 4 1 FV012 RV010 4 1 CT201 RV001 RV002 CT201 RV201 C-19 BOARD (SIDE B) CL053 RV009 14 TC009







SECTION 6 REPAIR PARTS LIST

6-1. EXPLODED VIEWS

NOTE:

- -XX and -X mean standardized parts, so they may have some difference from the original one.
- Items marked "*" are not stocked since they are seldom required for routine service. Some delay should be anticipated when ordering these items.
- The mechanical parts with no reference number in the exploded views are not supplied.
- Hardware (# mark) list and accessories are given in the last of the electrical parts list.


The components identified by mark \triangle or dotted line with mark \triangle are critical for safety

critical for safety,
Replace only with part number specified.

Les composants identifiés par une marque Δ sont critiquens pour la sécurité.


sécurité.
Ne les remplacer que par une pièce portant le numéro spécifié.

6-1-1. OVERALL ASSEMBLY




Ref. No.	Part No.	Description	<u>Remark</u>	Ref. No.	Part No.	Description	Remark
1 1 2 2 3	X-3950-022-1 X-3950-025-1	PANEL ASSY, FRONT (DSR-20MD) PANEL ASSY (P), FRONT (DSR-20MD) DOOR ASSY (DSR-20MD) DOOR ASSY (P) (DSR-20MDP) MAGNET	P)	* 7 8 * 9 10	4-886-821-01 3-987-158-01 3-987-171-01	HINGE (R) ASSY, DOOR SCREW, M3 CASE CASE, UPPER FOOT (FF-004) SUMITITE (B3), +BV	
4 5 * 6	4-942-567-01	KNOB, ROTARY EMBLEM (NO.4), SONY HINGE (L) ASSY, DOOR	•	12 13		REMOTE COMMANDER (RMT-DS20) LID, BATTERY CASE (for RMT-DSR20)	

6-1-2. CHASSIS ASSEMBLY



Ref. No.	Part No.	Description	Remark	Ref. No.	Part No.	Description	Remark
* 51	A-7074-111-A	FR-136 BOARD, COMPLETE		63	3-975-291-01	SCREW (4X6)	
52	3-970-608-41	SUMITITE (B3), +BV		* 64	3-987-157-31	PANEL. REAR	
* 53	A-7073-471-A	HP-100 BOARD, COMPLETE		65	1-958-841-11		
* 54	A-7067-251-A	VA-102 BOARD, COMPLETE (DSR-20	OMDP)	66	3-945-562-01		
* 54	A-7067-250-A	VA-102 BOARD, COMPLETE (DSR-20	OMD)	△ 67	1-958-585-11	HARNESS (AC-227)	
				*		,	
55	1-782-823-11	CABLE, FLAT (FVH-4)		68	1-958-059-11	HARNESS (VP-72)	
56		CABLE, FLAT (FVF-8)		* 69	A-7073-576-A	HG-1 BOARD, COMPLETE	
57		CABLE, FLAT (FVJ-7)		71	3-964-010-01	SCREW M2	
1 58 1 1	1-468-441-11	POWER BLOCK (U-1/U-2) (DSR-20M	D)	* 72	3-050-330-01	BRACKET, HG-1	
1 ∆ 58	1-468-442-11	POWER BLOCK (U-1/U-2) (DSR-20M	DP)	CN001	1-958-813-11	HARNESS (DH-50)	
* 59		RS-78 BOARD, COMPLETE		CN002	1-565-388-21	CONNECTOR, D-SUB 9P (REMOTE R	S-232C)
60		CABLE, FLAT (FVR-9)		J901	1-564-603-11	CONNECTOR (WITH DC SW) 4P	•
* 61	A-7073-470-A	RE-32 BOARD, COMPLETE		M905	1-698-534-31	FAN, DC	
62	1-782-826-11	CABLE, FLAT (FVR-10)					

6-1-3. MD BLOCK ASSEMBLY



Ref. No.	Part No.	Description	Remark	Ref. No.	Part No.	Description	Remark
101	1-776-148-11	CABLE, FLAT (FCM-11) 15P		* 110	A-7067-130-A	JC-19 BOARD, COMPLETE (DSR-20	MD)
102	1-776-145-11	CABLE, FLAT (FCM-8) 16P		* 111	3-987-133-01		
* 103	A-7067-127-A	CM-56 BOARD, COMPLETE (DSR-2	OMDP)	112	3-056-130-01	SCREW (M3), STEP	
* 103	A-7067-131-A	CM-56 BOARD, COMPLETE (DSR-2	20MD)	113	1-958-288-11	HARNESS (CM-130)	
* 104	3-987-138-01	FRAME, MD		.114	1-776-151-11	CABLE, FLAT (FCM-12) 14P	
105	1-764-137-11	CONNECTOR, TRANSLATION 15P		115	1-776-147-11	CABLE, FLAT (FCM-10) 15P	
* 106	A-7067-128-A	RP-228 BOARD, COMPLETE (DSR-	20MDP)	116		CABLE, FLAT (FCM-9) 9P	
* 106	A-7067-132-A	RP-228 BOARD, COMPLETE (DSR-	20MD)	117		HARNESS (CP-79)	
107	1-776-149-11	CABLE, FLEXIBLE FLAT 30P		118	1-958-061-11	HARNESS (VJ-103)	
108	1-783-376-11	CABLE, FLEXIBLE FLAT (FFC-245)		119	1-958-058-11	HARNESS (JP-55)	
109	3-732-817-01	SCREW (2X4.5), TAPPING		120	1-958-060-11	HARNESS (VJ-102)	
* 110		JC-19 BOARD, COMPLETE (DSR-20	MDP)	121		FILTER, CLAMP (FERRITE CORE)	

6-1-4. FL CASETTE COMPARTMENT ASSEMBLY

not supplied

1

2


3

4

(5)

NOTE FOR INSTALLATION

➤: Place for grease (SG-055G: 7-651-000-09) [红]: Take note of the position and specified direction.




Ref. No.	Part No.	Description	1	<u>Remark</u>	Ref. No.	Part No.	Description	Remark
701	A-7092-644-A	FL BLOCK ASSY		(5-2)	708	3-967-816-01	BELT, ROLLER	(5-38)
702	A-7092-647-A	SLOAT BLOCK ASSY, C		(5-41)	709	3-947-503-01	SCREW (M1.4X2.5)	, ,
703	3-967-604-01	SPRING (DB), TENSION	(A): 5-40/0	B: 5-41)	710	3-967-591-01	GEAR (B)	(5-40)
704	3-967-655-01	DOOR, C		(5-40)	711	3-967-590-01	GEAR (A)	(5-40)
705	3-967-613-01	SPRING (HS), TENSION C	OIL	(5-41)	712	3-967-653-01	OPENER, LID	(5-39)
706		ARM, DAMPER		(5-40)	713	3-967-592-01	WORM, C	(5-40)
707	X-3945-780-1	SHAFT ASSY, ROLLER		(5-38)	714	3-967-636-01	SPRING, SHIFT PLATE	(5-41)

6-1-5. MECHANISM CHASSIS ASSEMBLY (1) (TOP SIDE VIEW (1))

NOTE FOR INSTALLATION

PH-: Phase adjustment

: Place for grease (SG-055G: 7-651-000-09)




Ref. No.	Part No.	Description	Ref. page No.	Ref. No.	Part No.	Description	Ref. page No.
751	X-3945-801-1	BASE ASSY, TG3/4	(5-31)	765	3-967-751-01	COVER, LM	(5-15)
752	3-967-740-01	SPRING, TG3 LOCK	(5-26, 5-33)	766	3-967-767-01	WHEEL, LM WORM	(5-15)
753	X-3947-441-1	ROLLER ASSY, TG3	(5-26)	767	3-967-768-01	GEAR, PINCH DRIVING	(5-15)
754	3-703-816-74	SCRWE (M1.4X4.5), SPECIAL	HEAD	768	X-3945-803-1	BASE ASSY, TG5/6	(5-33)
755	A-7044-015-A	DRUM ASSY (DEH-08B-R)	(5-14)	769	3-967-741-01	SPRING, TG6 LOCK	(5-26, 5-33)
756	1-770-363-11	ELASTIC CONNECTOR	(5-14)	770	X-3945-810-1	ARM ASSY, PINCH	(5-16)
757	3-967-785-01	STOPPER, TAPE	(5-14)	771	3-967-645-01	SPRING (PINCH), TENSION COIL	, ,
758	3-967-728-01	SCREW (M2 X 5)		772	3-967-676-01	LIMITER, PINCH	(5-16)
759	3-967-817-01	BASE, DRUM	(5-14)	773	3-967-769-01	GEAR, PINCH CAM	(5-16)
760	A-7040-449-A	SCREW ASSY	(5-14)	774	3-967-679-01	RETAINER, PINCH	(5-16)
761	3-954-285-01	SCREW (M1.4X0.2)		775	3-967-795-03	ARM, HC	(5-16)
762	X-3945-798-1	ARM ASSY, TC	(5-15)	776	X-3945-802-1	ROLLER ASSY, TG6	(5-26)
763	3-967-675-01	HOLDER, LM	(5-15)	M901	X-3944-897-2	FPC ASSY, MOTOR	(5-14)
764	3-732-817-01	SCREW (2X4.5), TAPPING		M903		MOTOR ASSY, LM (LOADING)	(5-15)

6-1-6. MECHANISM CHASSIS ASSEMBLY (2) (TOP SIDE VIEW (2))

NOTE FOR INSTALLATION

PH:: Phase adjustment

Take note of the position and specified direction.




Ref. No.	Part No.	Description	Ref. page No.	Ref. No.	Part No.	<u>Description</u> Ref	. page No.
801	3-967-809-01	RETAINER, TG2	(5-19)	810	3-967-764-01	ARM, TG2 SELECTION	(5-18)
802	3-967-715-01	SPRING, TG2 PLATE	(5-19)	811	3-967-807-01	HOOK, TG2 SPRING	(5-18)
803	3-728-148-11	SCREW (M1.4X2.5), SPECIAL H	HEAD	812	3-967-724-01	ADJUSTOR, SPRING (5-18, @: 5-9)/®: 5-10)
804	3-967-714-01	MAGNET, ET (A: 5	5-19/®: 5-520)	813	3-967-810-01	RETAINER, TG7	(5-20)
805	X-3945-792-1	BAND ASSY, S TENSION REGUI	LATOR (5-19)	814	X-3945-806-1	ARM ASSY, TG7	(5-20)
806	X-3945-805-1	ARM ASSY, TG2	(5-19)	815	3-967-694-01	SPRING, TG7 PLATE	(5-20)
807	3-967-726-01	SPRING (TG2), TENSION COIL		816	X-3945-793-1	BAND ASSY, T TENSION REGULATO	R (5-20)
		(5-18, (): 5-9/®:5-10)	817	X-3945-783-1	ARM ASSY, TG7 LOAD	(5-20)
808	3-967-685-01		(5-18)	818	3-967-808-01	HOOK, TG7 SPRING	(5-18)
809	X-3945-782-1	ARM ASSY, TG2 LOAD	(5-19)				

6-1-7. MECHANISM CHASSIS ASSEMBLY (3) (TOP SIDE VIEW (3))

NOTE FOR INSTALLATION

PH-: Phase adjustment
: Place for grease (SG-055G: 7-651-000-09) : Take note of the position and specified direction.




Ref. No.	Part No.	Description	Ref. page No.	Ref. No.	Part No.	Description	Ref. page No.
851	1-658-990-11			861	3-967-776-01	BRAKE, S	(5-25)
852	3-318-201-11	SCREW (B) (1.4X3), TAPF	PING	862	3-967-673-01	SPRING, S BRAKE	(5-25)
853	3-967-684-01	SPRING, PLATE	(A: 5-23/B: 5-24)	863	3-967-775-01	RATCHET, T	(5-25)
854	3-967-692-01	GUARD, GOOSENECK	(5-17)	864	3-947-503-01	SCREW (M1.4X2.5)	
855	X-3945-807-1	ARM ASSY, GOOSENECK	(5-17)	865	X-3945-804-1	BASE ASSY, TG18	(5-25)
856	3-967-784-01	ARM, RL	(5-17)	866	3-967-725-01	HOLDER, T REEL	(5-22)
857	3-967-683-01	SPRING, RL PRESS	(5-17)	867	X-3945-815-1	PLATE ASSY, T REEL	(5-24)
858	X-3945-814-1	PLATE ASSY, S REEL	(5-23)	CN901		CONNECTOR 4P	(5-35)
859	3-967-680-01	LINK, RL	(5-17)	S901	1-762-551-11	SWITCH, PUSH (REC PROOF)	(5-35)
860	A-7040-441-A	TABLE BLOCK ASSY, REEL	_(A: 5-21/B: 5-22)	S902		SWITCH, PUSH (C IN SW)	(5-17)

6-1-8. MECHANISM CHASSIS ASSEMBLY (4) (BOTTOM SIDE VIEW (1))

NOTE FOR INSTALLATION

PH-: Phase adjustment

: Place for grease (SG-055G: 7-651-000-09)




Ref. No.	Part No.	Description	Ref. page No.	Ref. No.	Part No.	Description	Ref. page No.
901	3-967-678-01	GEAR, T CAM	(5-28)	914	3-967-746-01	SPRING, TG3/4 LIMITER	(5-31)
902	3-967-756-01	GEAR, TG7 CAM	(5-28)	915	X-3945-794-1	ARM ASSY, TG3/4	(5-31)
903	3-967-763-01	ARM, TG7 SELECTION	(5-28)	916	3-967-792-01	GEAR, TG5/6	(5-33)
904	3-967-677-01	SLIDER, M	(5-28)	917	3-967-748-01	SPRING, TG5/6 LIMITER	(5-33)
905	3-669-465-01	WASHER (1.5), STOPPER		918	X-3945-795-1	ARM ASSY, TG5/6	(5-33)
906	3-967-829-01	ARM, FL SELECTION	(5-28)	919	3-967-753-01	ARM, M SLIDER	(5-28)
907	3-967-828-01	GEAR, FL JOINT	(5-26)	920	3-967-819-01	CAM, MAIN	(5-29)
908	X-3945-813-1	ARM ASSY, FL JOINT	(5-27)	921	3-967-765-01	GEAR, TC	(5-27)
909	3-967-830-01	GEAR, FL RELAY	(5-27)	922	3-967-766-01	GEAR, RELAY	(5-27)
910	3-967-755-01	RETAINER, GL ARM	(5-28)	923	3-967-633-01	SPRING (TG2SL), TENSION COIL	(5-29)
911	3-947-503-01	SCREW (M1.4X2.5)		924	X-3945-781-1	ARM ASSY, TG2 SL	(5-29)
912	3-967-754-01	ARM, GL	(5-28)	M902	8-835-545-01	MOTOR, DC SCD11A/J-N (CAPSTA	
913	3-967-790-01	GEAR, TG3/4	(5-31)	S001	1-762-550-11	SWITCH, ROTARY (MODE)	(5-27)

6-1-9. MECHANISM CHASSIS ASSEMBLY (5) (BOTTOM SIDE VIEW (2))

NOTE FOR INSTALLATION

(PH-): Phase adjustment



Ref. No.	Part No.	Description	Ref. page No.	Ref. No.	Part No.	Description	Ref. page No.
* 951	A-7040-431-A	CHASSIS BLOCK ASSY, MECHA	ANICAL 63/64/65 board)	958 959	3-967-682-01 3-669-465-01	SPRING, MIC PRESS WASHER (1.5), STOPPER	(5-34)
952 953	3-967-690-01 X-3945-789-1	HOLDER, MIC ARM ASSY, RS	(5-35) (5-34)	960 961	3-967-681-01 3-947-503-01	GEAR, RACK JOINT	(5-35)
954 955		TENSION COIL SPRING LINK ASSY, PLATE	(5-34) (5-37)	962	3-967-689-01	HOLDER, RACK	(5-35)
956 957	X-3945-787-1 3-967-783-01	GEAR ASSY, RS LEVER, MIC	(5-34) (5-34)	963 964 M904	3-967-771-01 3-967-770-01 8-835-537-01	RACK (SC) RACK (LC) MOTOR, DC SRD11A/J-N (REEL)	(5-35) (5-35) (5-34)

CM-56

6-2. ELECTRICAL PARTS LIST

NOTE:

- Due to standardization, replacements in the parts list may be different from the parts specified in the diagrams or the components used on the set.
- Items marked "*" are not stocked since they are seldom required for routine service.
 Some delay should be anticipated when ordering these items.
- RESISTORS
 All resistors are in ohms.
 METAL: Metal-film resistor.
 METAL OXIDE: Metal oxide-film resistor.
 F: nonflammable
- SEMICONDUCTORS
 In each case, u: μ, for example: uA. : μA. . uPA. : μPA. . uPB. .: μPB. .: μPC. . uPD. .: μPC. .
- CAPACITORS uF: μF
- COILS uH: µH

The components identified by mark Δ or dotted line with mark Δ are critical for safety.
Replace only with part number specified.

Les composants identifiés par une marque \triangle sont critiquens pour la sécurité.

Ne les remplacer que par une pièce portant le numéro spécifié.

When indicating parts by reference number, please include the board.

Ref. No.	Part No.	Description	•	•	<u>Remark</u>	Ref. No.	Part No.	Description			Remark
*	A-7067-131-A	CM-56 BOARD, 0	OMPLETE (DSR-20	MD)	C054	1-127-530-11	ELECT	22uF	20%	20V
*		CM-56 BOARD, C				C055		CERAMIC CHIP	0.0047uF	5%	50V
	71 1001 127 71	*******			,	C056		CERAMIC CHIP	0.1uF	10%	25V
			(Re	f No. 4.0	00 Series)						
			(***	1	00 001100)	C058	1-127-530-11	ELECT	22uF	20%	20V
		< CAPACITOR >				C063		CERAMIC CHIP	0.33uF		25V
		COM MONOTON	•			C066		CERAMIC CHIP	22PF	5%	50V
0001	4 460 404 00	CEDAMIC CHID	150PF	5%	50V	C067		CERAMIC CHIP	22PF	5%	50V
C001		CERAMIC CHIP		5%	50V	C068		CERAMIC CHIP	0.01uF	0 /0	50V
C004		CERAMIC CHIP	150PF			0000	1-103-031-11	CENAIVIIC OFF	0.0141		30 V
C005		CERAMIC CHIP	0.01uF	10%	50V.	0070	4 400 040 00	OFDARMO OLUD	0.0000	100/	50V
C006	1-124-779-00		10uF	20%	16V	C073		CERAMIC CHIP	0.0068uF	10%	
C007	1-165-319-11	CERAMIC CHIP	0.1uF		50V	C075	1-163-235-11		22PF	5%	50V
						C076		CERAMIC CHIP	22PF	5%	50V
C009	1-126-206-11	ELECT CHIP	100uF	20%	6.3V	C078	1-124-779-00	ELECT CHIP	10uF	20%	16V
C011	1-163-037-11	CERAMIC CHIP	0.022uF	10%	25V	C079	1-163-009-11	CERAMIC CHIP	0.001uF	10%	50V
C017	1-164-161-11	CERAMIC CHIP	0.0022uF	10%	100V	4.				2.1	
C018		CERAMIC CHIP	0.01uF	10%	50V	C080	1-164-004-11	CERAMIC CHIP	0.1uF	10%	25V
C019		CERAMIC CHIP	0.01uF	10%	50V	C081	1-165-319-11	CERAMIC CHIP	0.1uF		50V
0010	1 100 021 01	OLI II III III	0.0.0			C082		CERAMIC CHIP	0.01uF		50V
C020	1-124-779-00	ELECT CHIP	10uF	20%	16V	C083	1-127-530-11		22uF	20%	20V
				20/0	50V	C086		CERAMIC CHIP	0.01uF	2070	50V
C021		CERAMIC CHIP	0.1uF	000/		C000	1-103-031-11	CENAIVIIC OFFI	0.01ui		307
C022	1-124-779-00		10uF	20%	16V	0007	4 407 500 44	FLFOT	00	000/	001/
C024		CERAMIC CHIP	10PF	0.5PF	50V	C087	1-127-530-11		22 u F	20%	20V
C025	1-124-779-00	ELECT CHIP	10uF	20%	16V	C088	1-126-193-11		1uF	20%	50V
					•	C090		CERAMIC CHIP	0.47uF	10%	16V
C026	1-165-319-11	CERAMIC CHIP	0.1uF		50V	C091		CERAMIC CHIP	0.01uF		50V
C027	1-163-021-91	CERAMIC CHIP	0.01uF	10%	50V	C092	1-126-205-11	ELECT CHIP	47uF	20%	6.3V
C028	1-163-227-11	CERAMIC CHIP	10PF	0.5PF	50V						
C029	1-165-319-11	CERAMIC CHIP	0.1uF		50V	C093	1-164-004-11	CERAMIC CHIP	0.1uF	10%	25V
C030		CERAMIC CHIP	0.01uF		50V	C094		CERAMIC CHIP	0.1uF	10%	25V
0000	1-100 001 11	OLI IAMIO OI III	0.0141		001	C095		CERAMIC CHIP	0.01uF		50V
0004	4 405 040 44	CERAMIC CHIP	0.1uF		50V	C096		CERAMIC CHIP	0.1uF	10%	25V
C031	1-165-319-11			000/				CERAMIC CHIP	0.0068uF	10%	50V
C032	1-126-193-11	ELECT	1uF	20%	50V	C099	1-103-019-00	GENAIVIIG GRIP	U.0000UF	10 /0	300
C034	1-163-021-91		0.01uF	10%	50V	0.101	4 400 004 44	OFDALLO OLUD	0.045		50)/
C035	1-163-021-91	CERAMIC CHIP	0.01uF	10%	50V	C101		CERAMIC CHIP	0.01uF		50V
C036	1-165-319-11	CERAMIC CHIP	0.1uF		50V	C102		CERAMIC CHIP	0.01uF		50V
						C103	1-126-204-11		47uF	20%	16V
C038	1-163-031-11	CERAMIC CHIP	0.01uF		50V	C104	1-163-009-11	CERAMIC CHIP	0.001uF	10%	50V
C039	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C105	1-164-004-11	CERAMIC CHIP	0.1uF	10%	25V
C040		CERAMIC CHIP	0.1uF	10%	25V						
C041		CERAMIC CHIP	0.047uF		50V	C106	1-163-031-11	CERAMIC CHIP	0.01uF		50V
C042		CERAMIC CHIP	180PF	5%	50V	C107		CERAMIC CHIP	0.1uF		50V
0042	1 100 201 11	OLIMANIO OIM	10071	0 / 0	001	C108		CERAMIC CHIP	0.047uF		50V
0046	1 165 210 11	CERAMIC CHIP	0.1uF		50V	C111		CERAMIC CHIP	0.01uF		50V
C046				100/		C112		CERAMIC CHIP	0.1uF		50V
C047		CERAMIC CHIP	0.1uF	10%	25V	0112	1-103-318-11	CENAIVIIO OTTIP	O. Tui		30 V
C048		CERAMIC CHIP	0.1uF	10%	25V	0440	4 400 000 00	OFDANNO OUID	0.0000	4.00/	EOV.
C049		CERAMIC CHIP	0.1uF	10%	25V	C113		CERAMIC CHIP	0.0082uF		50V
C051	1-165-319-11	CERAMIC CHIP	0.1uF		50V	C115		CERAMIC CHIP	0.0082uF		50V
						C117		CERAMIC CHIP	0.0082uF		50V
C052	1-126-193-11	ELECT	1uF	20%	50V	C118		CERAMIC CHIP	27PF	5%	50V
C053	1-126-397-11	ELECT	33uF	20%	25V	C119	1-163-009-11	CERAMIC CHIP	0.001uF	10%	50V
	1,41										

											7111 00
Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
						IC014	8-759-510-73	• • • • • • • • • • • • • • • • • • • •			jioman
C120	1-164-004-11	CERAMIC CHIP	0.1uF	10%	25V	10014	0-759-510-75	10 DA10393F-	22		
C121	1-165-319-11	CERAMIC CHIP	0.1uF		50V	IC016	8-759-510-71	IC BA10358F-	E2		
C122		CERAMIC CHIP	0.0068uF	10%	50V	IC017	8-759-011-65	IC TC74HC405	3AF (EL)		
C123	1-124-779-00		10uF	20%	16V	IC018	8-759-085-67	IC uPC339G2-	E2		
C126	1-164-004-11	CERAMIC CHIP	0.1uF	10%	25V	IC019		IC BA10358F-I			
C128	1-164-004-11	CERAMIC CHIP	0.1uF	10%	051/	IC021	8-759-335-42	IC CXA1793N-	E2		
C129		CERAMIC CHIP	0.1uF	10%	25V 25V	IC022	0 750 000 64	10 LB4007D			
C133		CERAMIC CHIP	0.01uF	10 /6	50V	IC501		IC LB1897D IC PQ30RV11			
C504		CERAMIC CHIP	0.1uF	10%	25V	IC503		IC LB1897D			
C505		CERAMIC CHIP	0.1uF	10%	25V		0.00000	10 2510075			
								< COIL >			
C506		CERAMIC CHIP	0.1uF	10%	25V						
C508	1-126-205-11		47uF	20%	6.3V	L003	1-412-282-41		470uH		
C509 C510		CERAMIC CHIP	0.47uF 0.01uF	10%	16V 50V	L004	1-414-398-11		10uH		
C511	1-126-205-11		47uF	20%	6.3V	L005 L006	1-414-398-11 1-414-398-11		10uH		
0011	7 120 200 11	ELLOT OTH	77 UI	2070	0.00	L007	1-414-402-11		10uH 47uH		
C512	1-126-193-11	ELECT	1uF	20%	50V		1 414 402 11	INDUCTOR	47 uii		
C513		CERAMIC CHIP	0.0022uF	10%	100V	L008	1-424-522-21	INDUCTOR	10uH		
C514		CERAMIC CHIP	0.01uF		50V	L010	1-424-522-21	INDUCTOR	10uH		
C515	1-126-206-11	ELECT CHIP	100uF	20%	6.3V	L011	1-409-535-41		100uH		
		0044450700				L013	1-424-524-21		47uH		
		< CONNECTOR >				L014	1-414-402-11	INDUCTOR	47uH		
CN001	1-770-690-11	CONNECTOR, FFO	/EDC 16D			L501	1 414 400 11	INDUCTOR	47		
* CN002		PIN, CONNECTOR				L501	1-414-402-11 1-414-402-11		47uH 47uH		
CN003		CONNECTOR, FFC		30P		2002	1-414-402-11	INDOCTOR	47 011		
* CN004	1-564-033-11	PIN, CONNECTOR	8P `					< IC LINK >			
CN005	1-770-692-11	CONNECTOR, FFC	FPC 9P								
. ONOOO	3 004 071 44					△ PS001	1-532-840-21	LINK, IC (1.25A)	(DSR-201	MDP)	
* CN006 * CN007		HOUSING, CONNE						TRANSIOTOR			
CN007		CONNECTOR, FFC						< TRANSISTOR	>		
011000	1 770 007 77	CONTRACTOR, 110	711 0 171			Q001	8-729-216-22	TRANSISTOR	258700	A-QRS-TX	,
		< DIODE >				Q002		TRANSISTOR	UN2213		1
						Q003		TRANSISTOR	MSD601		
D001		DIODE MA786-T				Q004	8-729-421-22	TRANSISTOR	UN2211		
D002		DIODE RD10M-T				Q008	8-729-010-25	TRANSISTOR	MSD601	I-RT1	
D004 D011		DIODE SB10-05F				0000	0.700.040.05	TDANGISTOR	Monac		
D012		DIODE MA786-T.				Q009 Q012	8-729-208-96	TRANSISTOR	MSD601		
5012	0 1 10 020 20	DIODE WITHOUT	^			Q014	8-729-421-19		UN2213	2-Y (TE16 -TY	L) .
D501		DIODE SB10-05F				Q500	8-729-216-22			A-QRS-TX	•
D502	8-719-108-24	DIODE MA151A-	TX			Q501	8-729-216-22			A-QRS-TX	
		< FUSE >				Q502	8-729-208-96			2-Y (TE16	L)
 ∆ F001	1-532-777-21	ELICE MICDO				Q503	8-729-421-19		UN2213		
221001	1-302-111-21		JDARY) (1.2)5Δ\ /D9	SR-20MD)	Q504	8-729-421-19	TRANSISTOR	UN2213	-1X	
		(02001	(I.a	.un) (DC	JI (201010)			< RESISTOR >			
		< FILTER >						\ TILOTOTOTI >			
		10 a 4 a				R001	1-216-057-00	METAL CHIP	2.2K	5%	1/10W
FL001		FILTER, BAND PAS				R002	1-216-065-91		4.7K	5%	1/10W
FL002	1-233-350-21	FILTER, BAND PAS	SS			R003	1-216-015-00		39	5%	1/10W
		10				R005	1-216-057-00		2.2K	5%	1/10W
		< IC >				R006	1-216-089-91	RES, CHIP	47K	5%	1/10W
IC001	8-759-062-66	IC TC7S66F (TE8	5B)			R009	1-216-049-91	BEC CHID	11/	E9/	1/10\4
1C002	8-759-235-19	IC TC74HC08AF	EL)			R010	1-216-049-91		1K 47K	5% 5%	1/10W 1/10W
IC003		IC CXP912032-07				R011	1-216-089-91		47K	5%	1/10W
IC005	8-759-327-00	IC CXA8044Q-T4				R012	1-216-089-91		47K	5%	1/10W
IC006	8-759-085-67	IC uPC339G2-E2				R015	1-216-295-91	•	0		
10000	0.750 100 11	10 TOT WILLS	(FL)			*					
IC008 IC009		IC TC74VHC125F IC BA6219BFP-Y-				R016	1-216-089-91		47K	5%	1/10W
IC009		IC CXA8010M-E1				R017 R018	1-216-295-91 1-216-089-91		0	E0/	4 /4 014/
IC012		IC MB3775PF-G-I				R019	1-216-089-91		47K 0	5%	1/10W
					,		. 210 200-01		U		

The components identified by mark \triangle or dotted line with mark \triangle are critical for safety. Replace only with part number specified.

Les composants identifiés par une marque \triangle sont critiques pour la sécurité. Ne les remplacer que par une pièce portant le numéro spécifié.

CM-56

Ref. No.	Part No.	Description			<u>Remark</u>	Ref. No.	Part No.	Description			<u>Remark</u>
R020	1-216-093-91	RES, CHIP	68K	5%	1/10W	R109	1-216-081-00	METAL CHIP	22K	5%	1/10W
R021	1-216-089-91	RES, CHIP	47K	5%	1/10W	R110	1-216-073-00	METAL CHIP	10K	5%	1/10W
R026	1-216-049-91	,	1K	5%	1/10W	R111	1-216-049-91				
R027	1-216-089-91		47K						1K	5%	1/10W
		•		5%	1/10W	R112	1-216-081-00		22K	5%	1/10W
R028	1-216-049-91		1K	5%	1/10W	R113	1-216-049-91		1K	5%	1/10W
R029	1-216-049-91	RES, CHIP	1K	5%	1/10W	R114	1-216-065-91	RES, CHIP	4.7K	5%	1/10W
R030	1-216-049-91	RES, CHIP	1K	5%	1/10W	R115	1-216-043-91	RES, CHIP	560	5%	1/10W
R032	1-216-049-91		1K	5%	1/10W	R116	1-216-057-00		2.2K	5%	1/10W
R033	1-216-049-91		1K	5%	1/10W	R117	1-216-043-91		560	5%	1/10W
R035	1-216-025-91	•	100	5%	1/10W	R118	1-216-057-00				
R036	1-216-025-91	•	100	5%	1/10W	R119	1-216-037-00		2.2K 10K	5% 5%	1/10W 1/10W
Dogo	1 010 075 00	METAL OLUB	4014		4.4.6111						
R039	1-216-075-00		12K	5%	1/10W	R120	1-216-073-00		10K	5%	1/10W
R040	1-216-069-00		6.8K	5%	1/10W	R121	1-219-107-91		1.5	5%	1/8W
R046	1-216-069-00		6.8K	5%	1/10W	R122	1-219-107-91	RES, CHIP	1.5	5%	1/8W
R050	1-216-077-91		15K	5%	1/10W	R125	1-219-107-91	RES, CHIP	1.5	5%	1/8W
R051	1-216-073-00	METAL CHIP	10K	5%	1/10W	R126	1-216-049-91	RES, CHIP	1K	5%	1/10W
R052	1-216-089-91	RES, CHIP	47K	5%	1/10W	R127	1-216-049-91	RES, CHIP	1K	5%	1/10W
R053	1-216-049-91	•	1K	5%	1/10W	R128	1-216-049-91	RES, CHIP			
R054	1-216-295-91	•	0	J /0	171044				1K	5%	1/10W
				F0/	4 /4 014/	R129	1-216-049-91	RES, CHIP	1K	5%	1/10W
R055	1-216-049-91		1K	5%	1/10W	R130	1-216-025-91		100	5%	1/10W
R059	1-216-043-91	RES, CHIP	560	5%	1/10W	R131	1-216-025-91	RES, CHIP	100	5%	1/10W
R060	1-216-049-91	RES, CHIP	1K	5%	1/10W	R132	1-216-025-91	BES CHID	100	5%	1/10W
R061	1-216-049-91		1K	5%	1/10W	R133	1-216-075-00		12K	5%	
R063	1-216-049-91		1K	5%	1/10W	R134	1-216-072-00				1/10W
R064	1-216-049-91		1K						9.1K	5%	1/10W
				5%	1/10W	R136	1-216-049-91		1K	5%	1/10W
R066	1-216-025-91	RES, CHIP	100	5%	1/10W	R137	1-216-049-91	RES, CHIP	1K	5%	1/10W
R067	1-216-073-00	METAL CHIP	10K	5%	1/10W	R138	1-216-049-91	RES, CHIP	1K	5%	1/10W
R069	1-216-085-00		33K	5%	1/10W	R143	1-216-073-00		10K	5%	1/10W
R070	1-216-073-00		10K	5%	1/10W	R146	1-216-295-91		0	J /0	1/1000
R071	1-216-025-91		100							F 0/	4 (4 0) 11
		,		5%	1/10W	R148	1-216-017-91		47	5%	1/10W
R075	1-216-049-91	NES, UNIP	1K	5%	1/10W	R153	1-216-295-91	SHORT	0		
R076	1-216-057-00	METAL CHIP	2.2K	5%	1/10W	R154	1-216-295-91	SHORT	0		
R077	1-216-025-91	RES, CHIP	100	5%	1/10W	R155	1-216-295-91		0		
R078	1-216-049-91		1K	5%	1/10W	R158	1-216-121-91		1M	5%	1/10W
R079	1-216-073-00		10K	5%	1/10W	R161	1-216-295-91		0	0 /0	171000
R080	1-216-025-91		100	5%	1/10W	R164				0.50/	4/4004
11000	1 210 020 01	neo, omi	100	J 76	171044	11104	1-216-672-11	WIETAL UTIF	7.5K	0.5%	1/10W
R081	1-216-049-91	RES, CHIP	1K	5%	1/10W	R165	1-216-017-91	RES, CHIP	47	5%	1/10W
R082	1-216-073-00	METAL CHIP	10K	5%	1/10W	R167	1-216-017-91	RES, CHIP	47	5%	1/10W
R083	1-216-089-91	RES, CHIP	47K	5%	1/10W	R168	1-216-055-00	METAL CHIP	1.8K	5%	1/10W
R084	1-216-025-91		100	5%	1/10W	R169	1-216-055-00		1.8K	5%	1/10W
R085	1-216-049-91		1K	5%	1/10W	R171	1-216-059-00		2.7K	5%	1/10W
		•								070	1, 1011
R086			1K	5%	1/10W	R176	1-216-033-00	METAL CHIP	220	5%	1/10W
R087	1-216-049-91	RES, CHIP	1K	5%	1/10W	R182	1-216-121-91	RES, CHIP	1M	5%	1/10W
R088	1-216-049-91	RES, CHIP	1K	5%	1/10W	R193	1-216-079-00		18K	5%	1/10W
R089	1-216-049-91	RES, CHIP	1K	5%	1/10W	R194	1-216-057-00	METAL CHIP	2.2K	5%	1/10W
R090	1-216-081-00		22K	5%	1/10W	R195	1-216-079-00		18K	5%	1/10W
D004	1 010 001 00	METAL OUR	0016	E0/	4/4031	Dica	4 040 000 00				
R091		METAL CHIP	22K	5%	1/10W	R196	1-216-057-00		2.2K	5%	1/10W
R092	1-216-089-91		47K	5%	1/10W	R201	1-216-073-00		10K	5%	1/10W
R093	1-216-049-91		1K	5%	1/10W	R203	1-216-121-91		1M	5%	1/10W
R094	1-216-671-11	METAL CHIP	6.8K	0.5%	1/10W	R206	1-216-073-00	METAL CHIP	10K	5%	1/10W
R095	1-216-645-11	METAL CHIP	560	0.5%	1/10W	R207	1-216-073-00	METAL CHIP	10K	5%	1/10W
R096	1-216-651-11	METAL CUID	1K	0.5%	1/10W	R208	1 016 045 00	METAL CLUD	600	E0/	4 /4 004/
R097						and the second second	1-216-045-00		680	5%	1/10W
	1-216-073-00		10K	5%	1/10W	R209	1-216-045-00		680	5%	1/10W
R098	1-216-121-91		1M	5%	1/10W	R211	1-216-671-11		6.8K	0.5%	1/10W
R099	1-216-105-91		220K	5%	1/10W		1-216-645-11		560	0.5%	1/10W
R102	1-216-089-91	KES, CHIP	47K	5%	1/10W	R213	1-216-651-11	METAL CHIP	1K	0.5%	1/10W
R103	1-216-089-91	RES CHIP	47K	5%	1/10W	R214	1-216-073-00	METAL CHID	101/	E0/	1/1014/
R104	1-216-295-91		0	0 /0	17 1099	R214	1-216-073-00		10K 220K	5% 5%	1/10W
R107	1-216-089-91		47K	5%	1/10W					5% 5%	1/10W
				J 70	1/1000	R216	1-216-081-00		22K	5%	1/10W
R108	1-216-295-91	onuni	0			R217	1-216-073-00	WE IAL CHIP	10K	5%	1/10W

CM-56	FP-406	FR-136

Ref. No	. Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
R218	1-216-081-00	METAL CHIP	22K	5%	1/10W	R550	1-216-073-00	METAL CHIP	10K	5%	1/10W
Dago	1 016 049 01	DEC CUID	560	5%	1/10W	R551	1-216-089-91	RES, CHIP	47K	5%	1/10W
R220											
R221			2.2K	5%	1/10W	R552	1-216-089-91		47K	5%	1/10W
R225			0			R553	1-216-073-00		10K	5%	1/10W
R227		METAL CHIP	2.2K	5%	1/10W	R554	1-216-073-00	METAL CHIP	10K	5%	1/10W
R228	3 1-216-043-91	RES, CHIP	560	5%	1/10W	R555	1-216-073-00	METAL CHIP	10K	5%	1/10W
R229	1-216-048-00	METAL CHID	910	5%	1/10W	R556	1-216-025-91	RES, CHIP	100	5%	1/10W
R230			8.2K	5%	1/10W	R560	1-216-025-91		4.7K	5%	1/10W
R231			910	5%	1/10W	R561	1-216-049-91		1K	5%	1/10W
R232			8.2K	5%	1/10W	R562	1-216-049-91		1K	5%	1/10W
R233	3 1-216-089-91	RES, CHIP	47K	5%	1/10W	R563	1-216-049-91	RES, CHIP	1K	5%	1/10W
R234	1-216-075-00	METAL CHIP	12K	5%	1/10W	R564	1-216-049-91	RES, CHIP	1K	5%	1/10W
R235			9.1K	5%	1/10W	R565	1-216-049-91		1K	5%	1/10W
R241			10K	5%	1/10W	R566	1-216-049-91		1K	5%	1/10W
R244			15K	5%	1/10W	R567	1-216-017-91	,	47	5%	1/10W
R245	1-217-671-11	METAL CHIP	1	5%	1/10W	R568	1-216-017-91	RES, CHIP	47	5%	1/10W
R247	1-216-073-00	METAL CHIP	10K	5%	1/10W	R569	1-216-017-91	RES, CHIP	47	5%	1/10W
R248			1	5%	1/10W	R570	1-216-017-91		47	5%	1/10W
R249			1	5%	1/10W	R572	1-216-295-91	•	0	570	171000
						1				5%	1/101//
R250			1	5%	1/10W	R573	1-216-049-91	•	1K		1/10W
R253	3 1-216-073-00	METAL CHIP	10K	5%	1/10W	R574	1-216-049-91	RES, CHIP	1K	5%	1/10W
R256	1-216-075-00	METAL CHIP	12K	5%	1/10W	R575	1-216-049-91	RES, CHIP	1K .	5%	1/10W
R257			18K	5%	1/10W	R901	1-216-295-91		0 (DSR-		17 1011
R259			0	0 70	17 1000	R902	1-216-295-91		0 (DSR-		
				En/	4/4/01/4	11902	1-210-233-31	אוווונו	יונים) ט	ZUMUF)	
R262			68K	5%	1/10W			. MIDDATOD			
R268	1-216-097-91	RES, CHIP	100K	5%	1/10W			< VIBRATOR >			
R270	1-216-073-00	METAL CHIP	10K	5%	1/10W	X001	1-760-655-21	VIBRATOR, CRYS	STAL (201	(Hz)	
R271			10K	5%	1/10W						
R274			10K	5%	1/10W						
R275			120K	5%	1/10W			ED ANG DOADD /	Dof No. 5	000 Caria	۵)
								FP-406 BOARD (mei No. 5,	ooo sene	s)
R283	1-216-089-91	RES, CHIP	47K	5%	1/10W			*****			
R284	1-216-025-91	RES. CHIP	100	5%	1/10W		1-658-990-11	FP-406 FLEXIBLI	E BOARD		
R285			100	5%	1/10W	1	3-318-201-11			ING	
R286			1K	5%	1/10W			HOLDER, MIC	10/1 // 1/ 1		
R290			1.8K	5%	1/10W			CLEANER, MIC			
							3-9/0-003-01	CLEANEN, MIC			
R506	1-216-049-91	neo, Unir	1K	5%	1/10W			< CONNECTOR >			
R516	1-216-295-91	SHORT	0					COMMEDICITY			
R517		METAL CHIP	3.9K	0.5%	1/10W	CN901	1-770-312-21	CONNECTOR 4P			
R518	1-216-655-11	METAL CHIP	1.5K	0.5%	1/10W						
R519	1-216-089-91	RES, CHIP	47K	5%	1/10W			< SWITCH >			
R520			47K	5%	1/10W						
						S901		SWITCH, PUSH (OF)	
R521			47K	5%	1/10W	S902	1-572-288-11	SWITCH, PUSH ((C IN)		
R524			47K	5%	1/10W						
R525	1-216-025-91	RES, CHIP	100	5%	1/10W						
R526			100	5%	1/10W	*	A-7074-111-A	FR-136 BOARD,	COMPLET	E	
R528			0					*****			
									(F	Ref.No. 5,0	000 Series)
R529			1.5	5%	1/8W	1					
R530			1.5	5%	1/8W	*	3-987-166-01	HOLDER, INDICA	TION TUB	E	
R531	1-219-107-91	RES, CHIP	1.5	5%	1/8W	.					
R535	1-216-089-91	RES, CHIP	47K	5%	1/10W			< BUZZER >			
R536			0								
						BZ101	1-529-104-11	BUZZER, PIEZOE	LECTRIC		
R537			0					OADAOITOD			
R538			0		41/200			< CAPACITOR >			
R541			10K	5%	1/10W						
R542			10K	5%	1/10W	C109		TANTALUM CHIP	33uF	20%	10V
R543	1-216-025-91	RES, CHIP	100	5%	1/10W	C110	1-164-156-11	CERAMIC CHIP	0.1uF		25V
						C111	1-164-156-11	CERAMIC CHIP	0.1uF		25V
R545	1-216-065-91	RES, CHIP	4.7K	5%	1/10W	C112		TANTALUM CHIP		20%	10V
R546			10K	5%	1/10W	C113		TANTALUM CHIP		20%	10V
R547		/ /	100	5%	1/10W						•
R549			4.7K	5%	1/10W	C114	1-164-156-11	CERAMIC CHIP	0.1uF		25V
	. 2.0 000 01	,		310					-, i ui		

FR-136

Ref.	No.	Part No.	Description			<u>Remark</u>	Ref. No.	Part No.	Description			<u>Remark</u>
C.	115 116 117 118	1-126-934-11 1-164-357-11	CERAMIC CHIP CERAMIC CHIP	1000PF 220uF 1000PF 1000PF	5% 20% 5% 5%	50V 10V 50V 50V	R103 R105 R106 R107	1-216-811-11 1-216-841-11 1-216-841-11	METAL CHIP	56K 150 47K 47K	5% 5% 5%	1/16W 1/16W 1/16W 1/16W
			< CONNECTOR >				R108	1-216-841-11	METAL CHIP	47K	5%	1/16W
CI	N104	1-774-770-11	CONNECTOR, FFO	C/FPC 30P			R109 R110 R111	1-216-841-11 1-216-841-11 1-216-841-11	METAL CHIP	47K 47K 47K	5% 5% 5%	1/16W 1/16W 1/16W
	101 108	8-719-061-58	DIODE MA151W DIODE CL-200Y	-C-TU (DUF		:	R112 R113	1-216-841-11 1-216-841-11	METAL CHIP	47K 47K	5% 5%	1/16W 1/16W
D.	109 110 111	8-719-061-58	DIODE CL-200Y DIODE CL-200Y	-C-TU (PAU			R114 R115 R116 R119	1-216-837-11 1-216-837-11 1-216-817-11 1-216-797-11	METAL CHIP	22K 22K 470 10	5% 5% 5% 5%	1/16W 1/16W 1/16W 1/16W
D.	112 113 114	8-719-027-84 8-719-106-08	DIODE CL-200Y DIODE CL-155U DIODE RD6.2M-	R/G-DT (OI T1B2		BY)	R120 R121	1-216-797-11 1-216-797-11		10	5% 5%	1/16W 1/16W
	115 116	8-719-106-08	DIODE RD6.2M- DIODE RD6.2M-	T1B2			R122 R123 R124	1-216-797-11 1-216-841-11 1-216-837-11	METAL CHIP METAL CHIP	10 47K 22K	5% 5% 5%	1/16W 1/16W 1/16W
D.	117	8-719-061-58	DIODE CL-200Y	-C-TU (REV	V)		R125	1-216-833-91	RES, CHIP	10K	5%	1/16W
FE	3101 3102 3103	1-414-445-11 1-414-445-11 1-414-445-11	FERRITE	> OuH OuH OuH			R126 R129 R130 R131 R138	1-216-833-91 1-216-837-11 1-216-833-91 1-216-833-91 1-216-825-11	METAL CHIP RES, CHIP RES, CHIP	10K 22K 10K 10K 2.2K	5% 5% 5% 5% 5%	1/16W 1/16W 1/16W 1/16W 1/16W
FE	3104 3105	1-414-445-11 1-414-445-11	FERRITE	OuH OuH		•	R139 R140	1-216-837-11 1-216-833-91	METAL CHIP	22K 10K	5% 5%	1/16W 1/16W
FE	3106 3107 3108 3109	1-414-445-11 1-414-445-11 1-414-445-11 1-414-445-11	FERRITE FERRITE	OuH OuH OuH OuH			R141 R142 R143	1-216-829-11 1-216-827-11 1-216-825-11	METAL CHIP	4.7K 3.3K 2.2K	5% 5% 5%	1/16W 1/16W 1/16W
FE FE	3110 3111 3112 3113				•		R144 R145 R146 R147 R148	1-216-845-11 1-216-797-11 1-216-797-11 1-216-797-11 1-216-797-11	METAL CHIP METAL CHIP METAL CHIP	100K 10 10 10 10	5% 5% 5% 5% 5%	1/16W 1/16W 1/16W 1/16W 1/16W
FE	3114 3115	1-414-229-11	INDUCTOR CHIP	0uH			R149 R151	1-216-821-11 1-216-811-11	METAL CHIP	1K 150	5% 5%	1/16W 1/16W
FE	3116 3117 3118		INDUCTOR CHIP INDUCTOR CHIP FERRITE				R152 R153 R154	1-216-811-11 1-216-811-11 1-216-811-11	METAL CHIP	150 150 150	5% 5% 5%	1/16W 1/16W 1/16W
			< 10 >				R155 R156	1-216-811-11 1-216-811-11		150 150	5% 5%	1/16W 1/16W
10	103 104 105		IC M66312FP-T IC uPD16311G0 IC RS-20E-T						< SWITCH >			
			< FLUORESCENT				\$101 \$102 \$104	1-572-272-11 1-572-342-11	SWITCH, TACTILI SWITCH, SLIDE (SWITCH, SLIDE (LOCAL/REN TIMER)	,	
N	D101	1-517-769-11	TUBE, FLUORESC	ENT INDIC	ATOR		S105 S106		SWITCH, SLIDE (SWITCH, TACTILI			
			< TRANSISTOR >				S107	1-579-349-11	SWITCH, SLIDE (COUNTER	SELECTY	
Q	101 103 104	8-729-424-18 8-729-421-19 8-729-421-19	TRANSISTOR	UN2113-T UN2213-T UN2213-T	X		S108 S109	1-692-838-21	SWITCH, TACTILE SWITCH, TACTILE	(RUBBER)) (EJECŤ))	R RESET)
Q.	105 106	8-729-421-19 8-729-421-19	TRANSISTOR	UN2213-T UN2213-T	Υ		S111 S112		SWITCH, TACTILI SWITCH, TACTILI	(RUBBER	(STOP)	
	107 110	8-729-421-19 8-729-421-19		UN2213-T UN2213-T			S113 S114 S115 S116	1-692-838-21 1-692-838-21 1-692-838-21	SWITCH, TACTILE SWITCH, TACTILE SWITCH, TACTILE SWITCH, TACTILE	(RUBBER) (RUBBER) (RUBBER)	(FF) (PAUSE)
							S117	1-762-333-21	SWITCH, TACTILE	(↓)		

FR-136 HG-1 HP-100

								L		1.0.	• • •	
J	Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
						***	*	A-7073-471-A	HP-100 BOARD	COMPLETE		
	S118		SWITCH, TACTIL						******			
	S119		SWITCH, TACTIL						i i	(Re	f.No. 5,0	00 Series)
	S120 S121		SWITCH, TACTIL SWITCH, TACTIL) (POWE	ER)			< CAPACITOR >			
	S122		SWITCH, TACTIL						< GAPAGITUR >			
1							C001	1-128-004-11		10uF	20%	16V
	de L	A 7070 F70 A	110 4 DOADD O	ONADI ETE			C002		CERAMIC CHIP	0.1uF		16V
4		A-7073-576-A	+**********				C003		CERAMIC CHIP ELECT CHIP	0.1uF	20%	16V
					ef.No. 8.0	000 Series)	C004		CERAMIC CHIP	10uF 47PF	20% 5%	16V 50V
	.'				,				02.0.000			
	1		< CAPACITOR >				C006		CERAMIC CHIP	0.1uF		16V
	C001	1-163-145-00	CERAMIC CHIP	0.0015uF	50/	50V	C007 C008		CERAMIC CHIP	47PF	5%	50V
	C002		CERAMIC CHIP	0.0013ul	J /0	25V	C009		ELECT CHIP CERAMIC CHIP	10uF 0.022uF	20% 10%	16V 25V
	C003	1-126-157-11		10uF	20%	16V	C010		CERAMIC CHIP	0.022uF	10%	25V
	C004	1-124-259-11		4.7uF	20%	35V						
	C005	1-126-157-11	ELECT	10uF	20%	16V	C011	1-162-966-11	CERAMIC CHIP	0.0022uF	10%	50V
			< CONNECTOR >						< CONNECTOR >	•		
	CN001 CN002		HARNESS (DH-5 PIN, CONNECTO				CN001	1-566-528-11	CONNECTOR, FF	PC (ZIF) 12P		
	ONOUZ	1 000 703-11	•	I TOF					< DIODE >		-	
			< DIODE >				D001	9-710-421-F0	DIODE MA3075	:\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\		
	D001	8-719-106-89	DIODE RD15M-	T1B2			D003	8-719-421-59	DIODE MA3075	WA- (TX)		
	D002		DIODE RD15M-				D003	8-719-421-59	DIODE MA3075	WA- (TX)		
	D003 D004		DIODE 1SS123- DIODE 1SS123-				D004		DIODE MA3075			
	D004		DIODE RD30M-				D005	8-719-421-59	DIODE MA3075	WA- (1X)		
			DIEDE TIDOOM	110			D006	8-719-421-59	DIODE MA3075	WA- (TX)		
	D006	8-719-022-76	DIODE RD30M-	T1B			D007	8-719-421-59	DIODE MA3075	WA- (TX)		
			<10>				D008	8-719-421-59	DIODE MASO75	WA- (TX)		
			(10)				D009 D010		DIODE MA3075			
	IC001		IC MM1256XF-E	3E			-,			` '		
	IC002	8-759-929-26	IC TL431CPSR				D011		DIODE MA3075			
			< TRANSISTOR >				D012	8-719-421-59	DIODE MA3075	WA- (1X)		
	Q001	8-729-120-28	TRANSISTOR	2801623-	T1-1516				< FERRITE BEAD	>		
	Q002	8-729-120-28		2SC1623-			FB001	1-500-241-22	FERRITE	0uH		
	Q003	8-729-120-28		2SC1623-			FB002	1-500-241-22		0uH		
	Q004	8-729-014-91	TRANSISTOR	2SD2185S	S-TX		FB003	1-500-241-22	FERRITE	0uH		
			< RESISTOR >						< IC >			
	R001	1-208-830-11	METAL CHIP	100K	0.5%	1/10W	IC001	8-759-369-73	IC NJM4556AM	-Δ-TF2		
	R002	1-208-830-11		100K	0.5%	1/10W	10001	3 100 000-10	NIAUCCEINIUM C.	7.16		
	R003	1-208-830-11		100K	0.5%	1/10W			< JACK >			
	R004 R005	1-208-848-11 1-208-830-11		560K 100K	0.5% 0.5%	1/10W 1/10W	J001	1-569-809-11	JACK (SMALL TY	PE) (PHONE	(S)	
	R007	1-208-830-11	METAL CHIP	100K	0.5%	1/10W			< RESISTOR >	, (
	R008	1-208-814-91		22K	0.5%	1/10W			TREGIOTOR >			
	R009	1-208-806-11		10K	0.5%	1/10W	R001	1-216-833-91				1/16W
		1-208-822-11		47K	0.5%	1/10W	R002	1-216-821-11				1/16W
	R011	1-208-822-11	WETAL CHIP	47K	0.5%	1/10W	R003 R004	1-216-821-11 1-216-833-91				1/16W
	R012	1-216-208-00	RES, CHIP	2.7K	5%	1/8W	R005	1-216-835-11				1/16W 1/16W
., .	R013	1-216-208-00	RES, CHIP	2.7K	5%	1/8W						.,
		1-208-795-11		3.6K	0.5%	1/10W		1-216-831-11				1/16W
		1-216-208-00 1-216-208-00		2.7K	5%	1/8W	R007	1-216-831-11				1/16W
	NOTO	1-210-200-00	nco, unir	2.7K	5%	1/8W		1-216-835-11 1-216-821-11				1/16W
			< RELAY >					1-216-821-11				1/16W 1/16W
	RY001	1-755-259-11	RFI AV									
	111001		TILLAI					1-216-809-11 1-216-809-11				1/16W 1/16W
									01111		- /0	., 1044

HP-100 JC-19

Ref. No.	Part No.	<u>Description</u>			Remark	Ref. No.	Part No.	<u>Description</u>			Remark	(
		< VARIABLE RES	SISTOR >			C162	1-162-974-11	CERAMIC CHIP	0.01uF		50V	
						C163		CERAMIC CHIP	0.01uF		50V	
RV001	1-238-612-11	RES, VAR, CARB	ON 20K/20	K (PHON	E LEVEL)	C164	1-162-974-11	CERAMIC CHIP	0.01uF		50V	
RV002		RES, VAR, CARE										
RV003	1-238-744-11	RES, VAR, CARB	ON, 50K (F	REC LE	VEL)	C165	1-162-974-11	CERAMIC CHIP	0.01 uF		50V	
						C166	1-162-974-11		0.01 uF		50V	
						C167	1-162-974-11		0.01uF		50V	
*		JC-19 BOARD, C				C168		CERAMIC CHIP	0.01uF		50V	
*	A-7067-126-A	JC-19 BOARD, C		(DSR-20 1	MDP)	C170	1-162-974-11	CERAMIC CHIP	0.01uF		50V	

			(H	et.No. 2,0	000 Series)	C171	1-162-974-11		0.01uF		50V	
	7-695-139-10	SCREW +P 2.6X	E TVDEO NO	TI IS MC		C172	1-162-974-11		0.01uF	000/	50V	
	1-003-132-19	30HEW #F 2.0A	J I I FEZ IV	JIN-OLI I		C173		TANTALUM CHIP		20% 20%	6.3V 6.3V	
		< CAPACITOR >				C174		TANTALUM CHIP		20%	6.3V	
		COMMONOTOR				0173	1-110-309-11	IANTALOW OTH	47 ur	20 /0	0.57	
C101	1-104-847-11	TANTALUM CHIP	22uF	20%	4V	C179	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	
C102		TANTALUM CHIP		20%	4V	C180		CERAMIC CHIP	0.01uF	2070	50V	
C103		TANTALUM CHIP		20%	4V	C181		CERAMIC CHIP	33PF	5%	50V	
C107		TANTALUM CHIP		20%	20V	C182	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	
C108	1-135-177-21	TANTALUM CHIP	1uF	20%	20V	C183	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	
C109		TANTALUM CHIP		20%	20V	C184	1-162-974-11	CERAMIC CHIP	0.01uF		50V	
C110		TANTALUM CHIP		20%	6.3V	C185		CERAMIC CHIP	0.01uF		50V	
C111		TANTALUM CHIP		20%	6.3V	C186		TANTALUM CHIP		20%	4V	
C112		TANTALUM CHIP		20%	6.3V	C187		TANTALUM CHIP		20%	4V	
C116	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C201	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	
0117	1 160 074 11	CEDAMMO CUUD	0.04		E0//	0000	1 104 057 11	OEDAMAO OUUD	10000	F0/	E01/	
C117 C118		CERAMIC CHIP	0.01uF 0.01uF		50V	C202		CERAMIC CHIP	1000PF	5%	50V	
C119		CERAMIC CHIP	10PF	0.5PF	50V 50V	C203 C204		CERAMIC CHIP CERAMIC CHIP	0.01uF 0.01uF	10% 10%	25V	
C120		CERAMIC CHIP	10PF	0.5PF	50V	C204		CERAMIC CHIP	0.01uF	10%	25V 25V	
C121		CERAMIC CHIP	10PF	0.5PF	50V	C206		CERAMIC CHIP	0.01uF	10%	25V 25V	
0.2.	1 102 010 11	OLINWIO OIM	1011	0.011	000	0200	1 102 370 11	OLITAWIO OTIII	0.0101	1070	231	
C122	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C207	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	
C123		CERAMIC CHIP	0.01uF		50V	C208		CERAMIC CHIP	0.01uF	10%	25V	
C124		CERAMIC CHIP	0.01uF		50V	C209		CERAMIC CHIP	0.01uF	10%	25V	
C125	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C210		TANTALUM CHIP		20%	6.3V	
C127	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C211	1-107-826-91	CERAMIC CHIP	0.1uF	10%	16V	
C128		CERAMIC CHIP	0.01uF		50V	C212		TANTALUM CHIP		20%	4V	
C129		TANTALUM CHIP		20%	6.3V	C214		CERAMIC CHIP	82PF	5%	50V	
C130			0.01uF		50V	C215		CERAMIC CHIP	390PF	5%	50V	
C131		CERAMIC CHIP		000/	50V	C216	No. of the Control of	CERAMIC CHIP	0.01uF		50V -	
C132	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C217	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	
C133	1 160 074 11	CERAMIC CHIP	0.01		EOV	0040	1 100 074 11	OEDAMIO OLUD	0.045		F01/	
C135		CERAMIC CHIP	0.01uF 0.01uF		50V 50V	C218 C219		CERAMIC CHIP TANTALUM CHIP		000/	50V	
C136		CERAMIC CHIP	0.01uF		50V	C219		TANTALUM CHIP		20% 20%	6.3V 6.3V	
C137		CERAMIC CHIP	0.01uF		50V	G221		CERAMIC CHIP	0.01uF	10%	25V	
C138		CERAMIC CHIP	0.01uF		50V	C223		TANTALUM CHIP		20%	20V	
• . • •			0.0.0			0220	1 100 111 21	THE THEORY OF THE	T GI	2070	201	
C139	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C224	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	
C141	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C225		TANTALUM CHIP		20%	20V	
C143	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C226		CERAMIC CHIP	0.0047uF		50V	
C148	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C227		CERAMIC CHIP	0.0033uF	10%	50V	
C149	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C229	1-104-912-11	TANTALUM CHIP	3.3uF	20%	16V	
		CERAMIC CHIP			50V	C231		CERAMIC CHIP	470PF	5%	50V	
C151		TANTALUM CHIP		20%	6.3V	C233		CERAMIC CHIP	2PF	0.25PF		
		TANTALUM CHIP		20%	6.3V	C234		CERAMIC CHIP	0.01uF		50V	
		TANTALUM CHIP		20%	6.3V	C238		CERAMIC CHIP	0.01uF		50V	
C154	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C241	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	
C155	1 105 050 14	TARITAL LINA OLUD	1005	000/	6 21/	00.40	4 464 000 44	OFDANAO OUE	0.465		1011	
		TANTALUM CHIP		20%	6.3V	C243		CERAMIC CHIP	0.1uF	E0/	16V	
		TANTALUM CHIP		20% 20%	6.3V 6.3V	C245 C246		CERAMIC CHIP	15PF	5%	50V	
		TANTALUM CHIP		20%	6.3V	C246		CERAMIC CHIP CERAMIC CHIP	0.01uF 0.01uF		50V	
		TANTALUM CHIP		20%	6.3V	C247	1-162-974-11		0.01uF		50V 50V	
0100	. 100 200-11	HAVINGOW OUR	Toul	2070	3.00	0240	1-102-014-11	OLIMNIO OTIF	v.o rur		JUV	
C160	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C249	1-110-569-11	TANTALUM CHIP	47uF	20%	4V	
C161	1-162-974-11		0.01uF		50V	C250		CERAMIC CHIP		10%	25V	
*.					,					. 0 /0		

R	ef. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
			CEDAMIC CHID	0.01.15	100/	25V	C521	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	C252	1-162-970-11	CERAMIC CHIP	0.01uF	10%			1-162-970-11	CERAMIC CHIP	0.01uF	10 /6	16V
	C253		CERAMIC CHIP	0.01uF	10%	25V	C522	1-104-360-11			20%	10V
	C254	1-162-9/0-11	CERAMIC CHIP	0.01uF	10%	25V	C523	1-104-051-11	TANTALUM CHIP	TOUF	2070	100
	C255	1-164-357-11	CERAMIC CHIP	1000PF	5%	50V	C524	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	C256	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C701	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V
	C257	1-164-357-11	CERAMIC CHIP	1000PF	5%	50V	C702	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C258		CERAMIC CHIP	0.001uF	10%	50V	C703	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C259	1-162-927-11	CERAMIC CHIP	100PF	5%	50V	C704	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V
	0000	1 104 010 11	TANTALUM CHIP	2 205	20%	16V	C705	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C260	1-104-912-11	CERAMIC CHIP	0.033uF	10%	16V	C706	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C261			0.033uF	1070	16V	C707	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C262 C401		CERAMIC CHIP	0.22ur 0.1uF	10%	16V	C707	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C402		CERAMIC CHIP	0.1uF	10 /0	16V	C709	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	0402	1-104-300-11	CENAMIO OTIII	o. rui			0700	1 104 500 11	OLINAMIO OTTO	0.747		
	C403	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C710	1-104-847-11	TANTALUM CHIP	22uF	20%	4V
	C404	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C711	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C405	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C712	1-109-982-11	CERAMIC CHIP	1uF	10%	10V
	C406		CERAMIC CHIP	0.1uF	10%	16V	C713	1-164-230-11	CERAMIC CHIP	220PF	5%	50V
	C407	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C714	1-162-919-11	CERAMIC CHIP	22PF	5%	50V
							0.717	1 100 010 11	0504440 0440	2005	Fn/	5014
	C408	1-135-259-11	TANTALUM CHIP		20%	6.3V	C715	1-162-919-11	CERAMIC CHIP	22PF	5%	50V
	C410	1-162-970-11		0.01uF	10%	25V	C801	1-162-974-11	CERAMIC CHIP	0.01uF	=0/	50V
	C412		CERAMIC CHIP	0.01uF	10%	25V	C802	1-162-923-11	CERAMIC CHIP	47PF	5%	50V
	C413		CERAMIC CHIP	0.01uF	10%	25V	C803	1-162-915-11	CERAMIC CHIP	10PF	0.5PF	50V
	C421	1-111-253-11	TANTALUM CHIP	100uF	20%	6.3V	C804	1-162-915-11	CERAMIC CHIP	10PF	0.5PF	50V
	C422	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	C807	1-164-357-11	CERAMIC CHIP	1000PF	5%	50V
	C423	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	C809	1-162-974-11	CERAMIC CHIP	0.01uF		50V
	C424	1-162-964-11		0.001uF	10%	50V	C810	1-110-569-11	TANTALUM CHIP	47uF	20%	6.3V
	C425	1-165-176-11		0.047uF	10%	16V	C811	1-162-923-11	CERAMIC CHIP	47PF	5%	50V
	C426		CERAMIC CHIP	2.2uF		16V	C812		TANTALUM CHIP	47uF	20%	6.3V
	C427	1-162-964-11		0.001uF	10%	50V	C813	1-162-974-11	CERAMIC CHIP	0.01uF		50V
	C428	1-164-360-11		0.1uF		16V	C814	1-162-974-11	CERAMIC CHIP	0.01uF	1001	50V
	C429		CERAMIC CHIP	0.1uF		16V	C815	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	C430	1-162-964-11		0.001uF	10%	50V	C816	1-162-974-11	CERAMIC CHIP	0.01uF	000/	50V
	C431	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C831	1-135-259-11	TANTALUM CHIP	TOUF	20%	6.3V
	C432	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	C832	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V
	C433		CERAMIC CHIP	0.001uF	10%	50V	C833	1-110-569-11	TANTALUM CHIP	47uF	20%	6.3V
	C434		CERAMIC CHIP	0.001uF	10%	50V	C834	1-135-151-21	TANTALUM CHIP	4.7uF	20%	4V
	C435	1-165-176-11		0.047uF	10%	16V	C835		TANTALUM CHIP		20%	4V
	C436	1-164-505-11	CERAMIC CHIP	2.2uF		16V	C837	1-162-923-11	CERAMIC CHIP	47PF	5%	50V
	0407	1 100 000 11	CERAMIC CHIP	1uF	10%	10V	C838	1-135-259-11	TANTALUM CHIP	10.15	20%	6.3V
	C437	1-109-982-11 1-109-982-11		1uF	10%	10V 10V	C839	1-162-923-11	CERAMIC CHIP	47PF	5%	50V
	C438 C439	1-162-915-11		10PF	0.5PF	50V	C840	1-135-259-11	TANTALUM CHIP		20%	6.3V
			CERAMIC CHIP	10PF	0.5PF	50V 50V	C843	1-162-923-11	CERAMIC CHIP	47PF	5%	50V
	C440 C441		CERAMIC CHIP	0.01uF	10%	25V	C844		CERAMIC CHIP	47PF	5%	50V
	0441	1-102-970-11	CENAMIC CHIP	O.OTul	10 /6	200	0044	1-102-320-11	OLITAMIO OTTI	7/11		301
	C442	1-135-259-11	TANTALUM CHIP	10uF	20%	6.3V	C845	1-104-847-11	TANTALUM CHIP	22uF	20%	4V
	C501	1-162-970-11		0.01uF	10%	25V	C846	1-104-847-11	TANTALUM CHIP	22uF	20%	4V
	C502	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C847	1-164-676-11	CERAMIC CHIP	2200PF	5%	16V
	C503		TANTALUM CHIP	10uF	20%	10V	C848	1-164-676-11	CERAMIC CHIP	2200PF	5%	16V
	C504	1-162-970-11		0.01uF	10%	25V	C849	1-135-181-21	TANTALUM CHIP	4.7uF	20%	6.3V
	0.00	4 400 070 44	OFDAMIO OUUD	0.04	400/	0514	0050	1 104 200 11	CEDAMIC CUID	20005	E0/	E01/
	C505	1-162-970-11		0.01uF	10%	25V	C850	1-164-392-11		390PF	5% 5%	50V 50V
	C506	1-162-970-11		0.01uF	10%	25V	C851 C852	1-164-392-11 1-164-360-11	CERAMIC CHIP	390PF 0.1uF	J 70	16V
	C511	1-135-259-11	TANTALUM CHIP		20%	6.3V					20%	6.3V
	C512	1-164-360-11		0.1uF	100/	16V	C853	1-135-181-21	TANTALUM CHIP TANTALUM CHIP		20%	6.3V 10V
	C513	1-164-227-11	CERAMIC CHIP	0.022uF	10%	25V	C854	1-100-149-21	TANTALUIVI UNIP	Z.ZUF	2070	100
	C514	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C855	1-163-009-11	CERAMIC CHIP	0.001uF	10%	50V
	C515	1-162-970-11		0.01uF	10%	25V	C856	1-163-009-11	CERAMIC CHIP	0.001uF	10%	50V
	C516	1-104-847-11			20%	4V	C857	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C517	1-162-970-11		0.01uF	10%	25V	C859		CERAMIC CHIP	0.1uF		16V
	C518	1-162-917-11	CERAMIC CHIP	15PF	5%	50V	C860	1-135-181-21	TANTALUM CHIP	4.7uF	20%	6.3V
	0540	4 400 070 44	OCDAMIO OLUB	0.015	100/	251	0001	1 164 260 14	CEDAMIC CLUB	0 1 ·· E		161/
	C519	1-162-970-11		0.01uF 22PF	10% 5%	25V	C861 C862	1-164-360-11	CERAMIC CHIP TANTALUM CHIP	0.1uF	20%	16V 4V
	C520	1-102-919-11	CERAMIC CHIP	2277	3%	50V	1 0002	1-100-101-21	IANIALUW CHIP	4./ UF	ZU70	4 V

E	Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description		<u>Remark</u>
	C863	1-165-176-11	CERAMIC CHIP	0.047uF	10%	16V	D501	8-719-073-01	DIODE MA11	1-TX	
	C864		CERAMIC CHIP	0.22uF	10%	10V	D503	8-719-421-27	DIODE. MA72	8-TX	
	C865	1-110-569-11	TANTALUM CHIP	47uF	20%	6.3V	D504		DIODE MA11		
							D901	8-719-073-01	DIODE MA11	1-TX	
	C901		CERAMIC CHIP	0.01uF	10%	25V	2000		DIODE MALE		
	C902		CERAMIC CHIP	0.1uF		16V	D902		DIODE KV147		
	C903		CERAMIC CHIP	0.1uF	000/	16V	D903		DIODE MA11 DIODE MA11		
	C904		TANTALUM CHIP CERAMIC CHIP	10uF 0.1uF	20%	6.3V 16V	D910	0-/19-0/3-01	DIODE MATT	1-17	
	C905	1-104-350-11	CENAIVITO CHIP	0.101		100			< FERRITE BEA	an s	
	C906	1-110-569-11	TANTALUM CHIP	47uF	20%	6.3V			V E		
	C907		CERAMIC CHIP	0.1uF		16V	FB401	1-543-955-22	FERRITE	OuH	
	C908		CERAMIC CHIP	0.1uF		16V	FB402	1-543-955-22	FERRITE	OuH	
	C909	1-135-259-11	TANTALUM CHIP		20%	6.3V					
	C910	1-115-467-11	CERAMIC CHIP	0.22uF	10%	10V			< FILTER >	•	
	0011	1 104 000 11	CERAMIC CHIP	0.1		16V	FL101	1 000.045.01	FILTER, LOW F	DACC (5 5MHz)	
	C911 C912		CERAMIC CHIP	0.1uF 0.1uF	10%	16V	FL102		FILTER, LOW F		
	C912		CERAMIC CHIP	0.0022uF	10%	50V	FL102		FILTER, LOW F		
	C915		CERAMIC CHIP	0.0022ui	10%	25V	12,00	7 200 040 21	THE LITT LOVE	7100 (0.011112)	
	C916		CERAMIC CHIP	47PF	5%	50V			< IC >		
	C917	1-135-181-21	TANTALUM CHIP		20%	6.3V	IC009		IC BA10324A		
	C919		CERAMIC CHIP	1000PF	5%	50V	IC010		IC BA10324A		
	C920		CERAMIC CHIP	1000PF	5%	50V	IC011		IC BA10324A		
	C921		CERAMIC CHIP	1000PF	5%	50V	IC012		IC BA10324A		
	C922	1-115-566-11	CERAMIC CHIP	4.7uF	10%	10V	IC013	8-759-510-71	IC BA10358F	-£2	
	C924	1-16/-260-11	CERAMIC CHIP	0.1uF		16V	10014	8-759-359-51	IC NJM431M	(TF2)	
	C926		CERAMIC CHIP	0.01uF	10%	25V	IC015		IC CXD23000		
	C927		CERAMIC CHIP	3PF	0.25PF		10016		IC CXD23000		
	C929		CERAMIC CHIP	1000PF	5%	50V	IC017	8-752-352-09	IC CXD23000	Q-T4	
	C930	1-162-927-11	CERAMIC CHIP	100PF	5%	50V	IC018	8-759-523-03	IC TC74HC40	66AFT (EL)	
	0001		OFD AND OUR	0.005	100/	4017	10040	0.750.447.75	10 0 042201	O MW T4	
	C931		CERAMIC CHIP	0.22uF	10%	10V	IC019 IC200		IC S-81322H		
	C932		CERAMIC CHIP	100PF 1000PF	5% 5%	50V 50V	IC205		IC CXD2193A		
	C933 C934		CERAMIC CHIP	1000FF	5%	50V	IC205		IC TC7S08FU		\$
	C935		CERAMIC CHIP	0.1uF	3 70	16V	IC207		IC TK11630U		
	0000	, 101 000 11	0211111110	0.747							
	C936		CERAMIC CHIP	100PF	5%	50V	IC209		IC TC74VHC1		
	C937		CERAMIC CHIP	1000PF	5%	50V	IC210		IC TC7SET08		
	C940		CERAMIC CHIP	1000PF	5%	50V		8-759-239-58			
	C941		CERAMIC CHIP	1uF	10%	10V	10212		IC TC7W00FL		
	C942	1-107-826-91	CERAMIC CHIP	0.1uF	10%	16V	IC213	8-759-523-02	IC TC74HC40	DOAFI (EL)	
			< CONNECTOR >				IC214	8-759-491-31	IC TC74VHCT	700AF (EL)	
			(00111120 1011)				IC401		IC CXD3103F		
	CN101	1-506-474-11	PIN, CONNECTOR	R 9P			IC402	8-759-328-28	IC ZA4024		
			CONNECTOR, FFO				IC403	8-759-328-28			
			CONNECTOR, FFO				IC410	8-759-433-17	IC uPD48244	5LG4-B10-9MH-E2-H	DC
			CONNECTOR, FFO				10444	0.750.505.60	10	OD 004 LVI	
	CN412	1-750-345-11	CONNECTOR, FFO	C/EPC (ZIF)	30P		IC411		IC uPD82094 IC CXP91203		
	ONTO	4 004 004 44	DIN CONNECTOR	. /1 Emm)	(CRAD) OF	,	IC421 IC422		IC CXP31203		
*			PIN, CONNECTOR CONNECTOR, BO				IC501		IC S579174P		
1)			PIN, CONNECTOR		AIID ZUI		10502		IC AK6440AN		
	CN831		PIN, CONNECTOR		(SMD) 8F		10002	0 700 110 00	70 71110 7110 111		
					. ,	•	10503		IC TC7S04FU		
			< TRIMMER >				IC504			APFV-G-BND-ER	
			040 4010005/8	(FA)			IC510		IC S-81230S		
	CT201	1-141-423-61	CAP, ADJ 20PF (A	AFC)			IC511 IC701		IC S-81350H		
			< DIODE >				10701	0.100.400.00	10 ONDE 104F	***	
							IC702		IC TSB11LV0		
	D201		DIODE KV1470T				IC703		IC HD643383		
	D421		DIODE HSM88W				IC801		IC CXD2705A		
	D422		DIODE KV1470T				IC802		IC TLV431AC		
	D423		DIODE HSM88W				IC804	ö-/59 - 465-80	IC TC74ACT0	819 (EL)	
	D424	d-119-055-86	DIODE KV1470T	L1-3			IC805	8-752-379-31	IC CXD3107F	}	
	D425	8-719-055-86	DIODE KV1470T	L1-3			10807		IC TC74LCX0		

	Ref. No.	Part No.	Description	Remark	Ref. No.	Part No.	Description			Remark
	IC831 IC832		IC NJM2115V IC NJM2115V		Q030	8-729-202-38	TRANSISTOR	2803320	6N-TE85L	-В
	IC833		IC NJM2115V		Q031	8-729-202-38	TRANSISTOR	2SC3320	6N-TE85L	-В
					Q032		TRANSISTOR		6A-T106-F	
	IC835		IC NJM2115V		Q033		TRANSISTOR		6A-T106-F	
	IC836		IC NJM2115V		Q034		TRANSISTOR		6A-T106-F	
	IC837	8-759-471-38	IC AK4520A-V	F-E2	Q035		TRANSISTOR	2SC408		•
	IC838	8-759-357-67	IC TK15125MT	L						
	IC840	8-759-358-47	IC NJM2115V	(TE2)	Q036	8-729-905-35	TRANSISTOR	2SC408	1T106R	
					Q037	8-729-905-35	TRANSISTOR	2SC408		
	IC841		IC TC75S56F (Q039	8-729-026-52	TRANSISTOR	2SA1576	6A-T106-F	₹
	IC901	8-759-523-97	IC TC74VHC12	3AFT (EL)	Q040	8-729-905-35	TRANSISTOR	2SC4081	1T106R	
1	10902	8-759-523-95	IC TC74VHC74	FT (EL))	Q041	8-729-026-52	TRANSISTOR	2SA1576	6A-T106-F	₹
	IC903		IC TC7W74FU	(TE12R)						
	IC904	8-759-429-28	IC CXD8630R		Q042		TRANSISTOR	2SC4081		
	10000		10 7074440		Q043		TRANSISTOR		5A-T106-F	₹
	IC906		IC TC74VHC12		Q044		TRANSISTOR	2SC4081		
	IC907	8-759-195-81	IC TC7S86FU (TE85R)	Q045		TRANSISTOR		6A-T106-F	₹ '
	IC908		IC TC7W08FU		Q048	8-729-402-42	TRANSISTOR	UN5213	-TX	
	IC909		IC TC74VHC74	FI (EL)			100			
	IC911	8-759-327-04	IC CXD2913Q		Q050		TRANSISTOR	2SC4081		
	1004.4	0 750 405 40	10 7111000400	21/0	Q051		TRANSISTOR	XP6501-		
	IC914		IC TLV2231CD		Q052		TRANSISTOR	2SC4081		
	IC915	0.750.050.00	IC TC4W53FU	(1E12K)	Q053		TRANSISTOR	XP6501-		
	IC916	8-759-058-62	IC TC7S08FU (1E85H)	Q200	8-729-905-35	TRANSISTOR	2SC4081	IT106R	
			< COIL >		0004	0.700.000.50	TRANSISTOR	0044576		_
			< COIL >		Q201		TRANSISTOR		6A-T106-F	1
	L011	1-414-398-11	INDUCTOR	10uH	Q501		TRANSISTOR	2SC4081		
	L012	1-414-398-11		10uH	Q502 Q504		TRANSISTOR TRANSISTOR	2SC4081		
	L013	1-414-398-11		10uH	Q505		TRANSISTOR	UN5113-		
	L014	1-414-398-11		10uH	4303	0-129-421-10	INANSISTUM	XP4401-	IVE	
	L015	1-414-398-11		10uH	Q506	8-720-101-07	TRANSISTOR	200700	T1-DLDK	
	20.0	7 771 000 11	110001011	10011	Q801		TRANSISTOR	2SC4081		
	L016	1-414-398-11	INDUCTOR	10uH	Q832		TRANSISTOR	UN5111-		
	L017	1-414-398-11		10uH	Q902		TRANSISTOR	2SC4081		
	L018	1-414-398-11		10uH	0903		TRANSISTOR	UN5213-		
	L102	1-414-398-11		10uH	(300	0-123-402-42	INANSISTON	010213-	17	
	L200	1-414-398-11		10uH	Q910	8-729-015-76	TRANSISTOR	UN5211-	TY	
					Q911		TRANSISTOR	UN5211-		
	L202	1-410-390-11	INDUCTOR CHIP	56uH		0.200.0	111111111111111111111111111111111111111	CHOLII	17.	
	L203	1-414-398-11	INDUCTOR	10uH			< RESISTOR >			
	L204	1-414-398-11	INDUCTOR	10uH						
	L205	1-411-275-21	COIL, VARIABLE		R002	1-216-864-11	METAL CHIP	0	5%	1/16W
	L206	1-410-655-31	INDUCTOR CHIP	120uH	R003	1-414-760-21	FERRITE	0uH		
					R004	1-414-760-21	FERRITE	0uH		
	L401	1-414-398-11		10uH	R005	1-414-760-21		0uH		
	L402	1-414-398-11		10uH	R009	1-414-760-21	FERRITE	0uH		
	L421		INDUCTOR CHIP							
	L422		INDUCTOR CHIP		R010	1-216-864-11	METAL CHIP	0	5%	1/16W
	L423	1-414-398-11	INDUCTOR	10uH	R011	1-216-864-11	METAL CHIP	0	5%	1/16W
					R012	1-216-864-11	METAL CHIP	0	5%	1/16W
	L424		INDUCTOR CHIP		R013	1-216-864-11		0	5%	1/16W
	L501	1-414-398-11		10uH	R014	1-216-864-11	METAL CHIP	0	5%	1/16W
	L701		INDUCTOR CHIP							
	L702	1-414-398-11		10uH	R015	1-216-864-11	METAL CHIP	0	5%	1/16W
	L703	1~410-393-11	INDUCTOR CHIP	100uH	R016		METAL CHIP	0	5%	1/16W
	1.004	4 440 000 44	MOHOTOR OUR	4.11	R017		FERRITE	OuH		
	L801		INDUCTOR CHIP		R018	1-414-760-21		0uH		
	L802		INDUCTOR CHIP	· ·	R020	1-414-760-21	FERRITE	0uH		
	L901	1-414-398-11		10uH	D004	4 44 4 700 04				
	L904 L905	1-414-398-11		10uH	R021	1-414-760-21		0uH	me.	4 (4 ***
	L903	1-411-2/3-21	COIL, VARIABLE		R023	1-216-864-11		0	5%	1/16W
	L907	1-414-398-11	INDUCTOR	10vH	R030	1-414-760-21	FERRITE	0uH	F0/	4 (4 0)
	L307	1-414-020-11	אטוטטעווו	10uH	R033	1-216-864-11		0	5%	1/16W
			< TRANSISTOR >	,	R034	1-216-864-11	METAL CHIP	0	5%	1/16W
			- THANOIOTON >		R035	1-216-864-11	METAL CHIP	0	5%	1/16W
	Q026	8-729-905-35	TRANSISTOR	2SC4081T106R	R036		METAL CHIP	0		1/16W
	Q027	8-729-905-35		2SC4081T106R	R037		METAL CHIP	0		1/16W
	Q028	8-729-905-35		2SC4081T106R	R038	The Share and the state of the	FERRITE	0 0uH	J /0	17 1 0 14
	Q029	8-729-202-38		2SC3326N-TE85L-B	R039	1-216-864-11		0	5%	1/16W
								-	5 / 0	.,

Ref. No.	Part No.	Description	: ()		Remark	Ref. No.	Part No.	Description			<u>Remark</u>
R040	1-414-760-21	FERRITE	0uH			R147	1-216-821-11	METAL CHIP	1K	5%	1/16W
R041	1-414-760-21		OuH			R148	1-216-830-11	METAL CHIP	5.6K	5%	1/16W
R042	1-414-760-21		0uH			R149	1-216-824-11	METAL CHIP			
R043	1-216-864-11		0	5%	1/16W	R150	1-216-824-11		1.8K	5%	1/16W
R044	1-414-760-21		_	376	1/1044	1		METAL CHIP	1.8K	5%	1/16W
11044	1-414-700-21	FERRITE	0uH			R151	1-216-824-11	METAL CHIP	1.8K	5%	1/16W
R045	1-216-864-11	METAL CHIP	0	5%	1/16W	D150	1 016 005 11	METAL CLUD	151/	EN	4 (4 (2) 8)
R046		METAL CHIP				R152	1-216-835-11	METAL CHIP	15K	5%	1/16W
	1-216-864-11		0	5%	1/16W	R154	1-216-809-11	METAL CHIP	100	5%	1/16W
R047	1-216-864-11		0	5%	1/16W	R155	1-216-809-11	METAL CHIP	100	5%	1/16W
R048	1-414-760-21		0uH			R156	1-216-809-11	METAL CHIP	100	5%	1/16W
R049	1-414-760-21	FERRITE	0uH			R157	1-216-864-11	METAL CHIP	0	5%	1/16W
DOFO	4 44 4 700 04	FEDDITE	0.11						·		
R050	1-414-760-21		0uH			R159	1-216-805-11	METAL CHIP	47	5%	1/16W
R052	1-414-760-21		OuH			R160	1-216-821-11	METAL CHIP	1K	5%	1/16W
R053	1-414-760-21		OuH			R161	1-216-821-11	METAL CHIP	1K	5%	1/16W
R054	1-414-760-21		0uH			R162	1-216-823-11	METAL CHIP	1.5K	5%	1/16W
R055	1-414-760-21	FERRITE	0uH			R163	1-216-864-11	METAL CHIP	0	5%	1/16W
R056	1-414-760-21		0uH			R164	1-216-816-11	METAL CHIP	390	5%	1/16W
R057	1-414-760-21		0uH			R165	1-216-864-11	METAL CHIP	0	5%	1/16W
R059	1-216-837-11	METAL CHIP	22K	5%	1/16W	R166	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
R060	1-216-818-11	METAL CHIP	560	5%	1/16W	R167	1-216-823-11	METAL CHIP	1.5K	5%	1/16W
R061	1-216-817-11	METAL CHIP	470	5%	1/16W	R168	1-216-835-11	METAL CHIP	15K	5%	1/16W
										- , -	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,
R062	1-216-821-11	METAL CHIP	1K	5%	1/16W	R169	1-216-832-11	METAL CHIP	8.2K	5%	1/16W
R063	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R170	1-216-864-11	METAL CHIP	0	5%	1/16W
R064	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R173	1-216-821-11	METAL CHIP	1K	5%	1/16W
R065	1-216-837-11	METAL CHIP	22K	5%	1/16W	R174	1-216-821-11	METAL CHIP	1K	5%	1/16W
R067	1-216-818-11	METAL CHIP	560	5%	1/16W	R175	1-216-813-11	METAL CHIP	220	5%	1/16W
				0,0		11110	1 210 010 11	WEIZE OITH	220	0 70	1/1000
R068	1-216-817-11	METAL CHIP	470	5%	1/16W	R176	1-216-821-11	METAL CHIP	1K	5%	1/16W
R069	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R177	1-216-821-11	METAL CHIP	1K	5%	1/16W
R070	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R178	1-216-817-11	METAL CHIP	470	5%	1/16W
R071	1-216-821-11	METAL CHIP	1K	5%	1/16W	R182	1-216-821-11	METAL CHIP	1K	5%	1/16W
R113	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R183	1-216-821-11	METAL CHIP			
11110	1-210-023-11	METAL OTHE	4.71	J /0	1/1000	nios	1-210-021-11	WEIAL CHIP	1K	5%	1/16W
R115	1-216-837-11	METAL CHIP	22K	5%	1/16W	R184	1-216-817-11	METAL CHIP	470	5%	1/16W
R116	1-216-837-11	METAL CHIP	22K	5%	1/16W	R185	1-216-821-11	METAL CHIP			
R117	1-216-837-11	METAL CHIP	22K	5%	1/16W				1K	5%	1/16W
R118	1-216-821-11	METAL CHIP				R186	1-216-821-11	METAL CHIP	1K	5%	1/16W
			1K	5%	1/16W	R187		METAL CHIP	470	5%	1/16W
R119	1-216-821-11	METAL CHIP	1K	5%	1/16W	R192	1-216-821-11	METAL CHIP	1K	5%	1/16W
R120	1-216-821-11	METAL CHIP	1K	5%	1/16W	R199	1 016 001 11	METAL CHID	11/	E0/	4 /4 014/
R121	1-216-864-11	METAL CHIP	0				1-216-821-11		1K	5%	1/16W
R122	1-216-864-11			5%	1/16W	R201		METAL CHIP	1K	5%	1/16W
		METAL CHIP	0	5%	1/16W	R203		METAL CHIP	0	5%	1/16W
R123		METAL CHIP	0	5%	1/16W	R204		METAL CHIP	3.3K	5%	1/16W
R124	1-216-833-91	RES, CHIP	10K	5%	1/16W	R205	1-216-864-11	METAL CHIP	0	5%	1/16W
D105	1 010 000 01	DEC CLUB	401/	F0/	4 (4 0)4(D000	4 040 004 44	METAL OLUB		==/	
R125	1-216-833-91		10K	5%	1/16W	R206		METAL CHIP	0	5%	1/16W
R126		METAL CHIP	2.2K	5%	1/16W	R207		METAL CHIP	0	5%	1/16W
R127		METAL CHIP	8.2K	5%	1/16W	R208	1-216-827-11		3.3K	5%	1/16W
R128		METAL CHIP	100	5%	1/16W	R209		METAL CHIP	0	5%	1/16W
R129	1-216-809-11	METAL CHIP	100	5%	1/16W	R210	1-216-864-11	METAL CHIP	0	5%	1/16W
D.100											
R130		METAL CHIP	100	5%	1/16W	R211	1-216-833-91	,	10K	5%	1/16W
R131		METAL CHIP	3.3K	5%	1/16W	R212	1-216-864-11	METAL CHIP	0	5%	1/16W
R133		METAL CHIP	4.7K	5%	1/16W	R213	1-216-825-11	METAL CHIP	2.2K	5%	1/16W
R134	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R214	1-216-833-91	RES, CHIP	10K	5%	1/16W
R135	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R215		METAL CHIP	39K	5%	1/16W
R136		METAL CHIP	5.6K	5%	1/16W	R216		METAL CHIP	0 .	5%	1/16W
R137		RES, CHIP	10K	5%	1/16W	R217	1-216-828-11	METAL CHIP	3.9K	5%	1/16W
R138	1-216-833-91	RES, CHIP	10K	5%	1/16W	R218	1-216-864-11	METAL CHIP	0	5%	1/16W
R139		METAL CHIP	1K	5%	1/16W	R219		METAL CHIP	0	5%	1/16W
R140		METAL CHIP	4.7K	5%	1/16W	R246		METAL CHIP	Ö	5%	1/16W
									_		
R141	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R247	1-216-821-11	METAL CHIP	1K	5%	1/16W
R142		METAL CHIP	1K	5%	1/16W	R248		METAL CHIP	470	5%	1/16W
R143		METAL CHIP	15K	5%	1/16W	R249	1-216-817-11		470	5%	1/16W
R145		METAL CHIP	1K	5%	1/16W	R250		METAL CHIP	1K	5%	1/16W
R146	1-216-821-11		1K	5%	1/16W	R256	1-216-833-91		10K	5%	1/16W
				5.5			5 000 01	, 31111	. 511	J /0	171044

Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
						D/120	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
		ALCTAL OLUB	0.01/	50 /	4.4011	R428					
R259	1-216-832-11		8.2K	5%	1/16W	R429	1-216-845-11	WETAL GRIP	100K	5%	1/16W
R260	1-216-829-11		4.7K	5%	1/16W					ma/	4 (4 0)44
R261	1-216-814-11		270	5%	1/16W	R430	1-216-805-11	METAL CHIP	47	5%	1/16W
R262	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R431	1-216-817-11		470	5%	1/16W
R264	1-216-833-91	RES, CHIP	10K	5%	1/16W	R432	1-216-845-11		100K	5%	1/16W
						R433	1-216-845-11	METAL CHIP	100K	5%	1/16W
R266	1-216-833-91	RES. CHIP	10K	5%	1/16W	R434	1-216-805-11	METAL CHIP	47	5%	1/16W
R268	1-216-829-11		4.7K	5%	1/16W						
R269	1-216-814-11		270	5%	1/16W	R435	1-216-845-11	METAL CHIP	100K	5%	1/16W
R270	1-216-829-11		4.7K	5%	1/16W	R437	1-216-864-11		0	5%	1/16W
						R438	1-216-864-11		O O	5%	1/16W
R271	1-216-864-11	WE TAL CHIP	0	5%	1/16W		1-216-864-11		0	5%	1/16W
		AAFTAL OLUD	0001/	50 /	4 (4 0)44	R439					
R272	1-216-855-11		680K	5%	1/16W	R440	1-216-833-91	RES, CHIP	10K	5%	1/16W
R273	1-216-839-11		33K	5%	1/16W				12.4		4 14 0111
R274	1-216-864-11		0	5%	1/16W	R441	1-216-833-91	•	10K	5%	1/16W
R276	1-216-833-91	RES, CHIP	10K	5%	1/16W	R442	1-216-815-11		330	5%	1/16W
R277	1-216-839-11	METAL CHIP	33K	5%	1/16W	R443	1-216-805-11	METAL CHIP	47	5%	1/16W
						R444	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
R278	1-216-864-11	METAL CHIP	0	5%	1/16W	R445	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
R279	1-216-839-11		33K	5%	1/16W						
R280	1-216-833-91		10K	5%	1/16W	R446	1-216-845-11	METAL CHIP	100K	5%	1/16W
			1K	5%	1/16W	R447	1-216-845-11		100K	5%	1/16W
R281	1-216-821-11					R448	1-216-845-11		100K	5%	1/16W
R282	1-216-839-11	WE IAL CHIP	33K	5%	1/16W						
			_			R449	1-216-821-11		1K	5%	1/16W
R285	1-216-864-11		0	5%	1/16W	R450	1-216-857-11	METAL CHIP	1M	5%	1/16W
R286	1-216-833-91		10K	5%	1/16W						
R288	1-216-833-91	RES, CHIP	10K	5%	1/16W	R451	1-216-845-11		100K	5%	1/16W
R292	1-216-833-91	RES, CHIP	10K	5%	1/16W	R452	1-216-845-11	METAL CHIP	100K	5%	1/16W
R295	1-216-864-11		0	5%	1/16W	R453	1-216-845-11	METAL CHIP	100K	5%	1/16W
						R454	1-216-845-11	METAL CHIP	100K	5%	1/16W
R297	1-216-864-11	METAL CHIP	0	5%	1/16W	R455	1-216-845-11		100K	5%	1/16W
R298	1-216-864-11		Ö	5%	1/16W		. 210 010 11			•	
				5%	1/16W	R456	1-216-845-11	METAL CHIP	100K	5%	1/16W
R299	1-216-825-11		2.2K				1-216-845-11		100K	5%	1/16W
R300	1-216-825-11		2.2K	5%	1/16W	R457					
R301	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R458	1-216-845-11		100K	5%	1/16W
						R459	1-216-845-11		100K	5%	1/16W
R302	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R460	1-216-845-11	METAL CHIP	100K	5%	1/16W
R303	1-216-821-11	METAL CHIP	1K	5%	1/16W						
R304	1-216-815-11	METAL CHIP	330	5%	1/16W	R461	1-216-864-11	METAL CHIP	0	5%	1/16W
R313	1-216-827-11		3.3K	5%	1/16W	R462	1-216-864-11	METAL CHIP	0	5%	1/16W
R314	1-216-815-11		330	5%	1/16W	R463	1-216-864-11		0	5%	1/16W
11014	1 210 010 11	WE IN LE OTHE	000	0 / 0	.,	R464	1-216-864-11		0	5%	1/16W
R315	1-216-815-11	METAL CHID	330	5%	1/16W	R465	1-216-864-11		Ö	5%	1/16W
			0	J /0	171000	11400	1 210 007 11	WEINE OIM	Ü	0 70	171011
R318	1-216-295-91			0.50/	4 /4 CM	DAGG	4 046 064 44	METAL CHID	0	5%	1/16W
R319	1-218-864-11	METAL CHIP	5.1K	0.5%	1/16W	R466	1-216-864-11				
				•	SR-20MD)	R467	1-216-864-11		0	5%	1/16W
R319	1-218-865-11	METAL CHIP	5.6K	0.5%	1/16W	R468	1-216-864-11		0	5%	1/16W
				(DS	R-20MDP)	R469	1-216-821-11		1K	5%	1/16W
R320	1-218-831-11	METAL CHIP	220	0.5%	1/16W	R470	1-216-821-11	METAL CHIP	1K	5%	1/16W
R321	1-218-851-11	METAL CHIP	1.5K	0.5%	1/16W	R471	1-216-821-11	METAL CHIP	1K	5%	1/16W
R322	1-218-831-11		220	0.5%	1/16W	R472	1-216-821-11		1K	5%	1/16W
R339	1-216-864-11		0	5%	1/16W	R502	1-216-809-11		100	5%	1/16W
	1-216-841-11		47K	5%	1/16W	R503	1-216-809-11		100	5%	1/16W
R340									0	5%	1/16W
R341	1-216-841-11	METAL CHIP	47K	5%	1/16W	R504	1-216-864-11	WE TAL UTIP	U	370	171000
			_		44011		1 010 000 11	METAL OLUB	400	50/	4 (4 0) 4 (
R342	1-216-864-11		0	5%	1/16W	R505	1-216-809-11		100	5%	1/16W
R343	1-216-840-11	METAL CHIP	39K	5%	1/16W	R506	1-216-864-11	METAL CHIP	0	5%	1/16W
R344	1-216-821-11	METAL CHIP	1K	5%	1/16W	R507	1-216-809-11		100	5%	1/16W
R345	1-216-821-11	METAL CHIP	1K	5%	1/16W	R508	1-216-864-11	METAL CHIP	0	5%	1/16W
R346	1-216-833-91		10K	5%	1/16W	R509	1-216-833-91		10K	5%	1/16W
R347	1-218-871-11	METAL CHIP	10K	0.5%	1/16W	R510	1-216-864-11	METAL CHIP	0	5%	1/16W
					1/16W	R511	1-216-833-91		10K	5%	1/16W
R401	1-216-821-11		1K	5%				•			1/16W
R422	1-216-805-11		47	5%	1/16W	R512	1-216-809-11		100	5%	
R423	1-216-805-11		47	5%	1/16W	R513	1-216-809-11	METAL CHIP	100	5%	1/16W
R424	1-216-805-11	METAL CHIP	47	5%	1/16W	R514	1-216-809-11	METAL CHIP	100	5%	1/16W
R425	1-216-833-91		10K	5%	1/16W	R515	1-216-809-11	METAL CHIP	100	5%	1/16W
R426	1-216-821-11	METAL CHIP	1K	5%	1/16W	R517	1-216-809-11	METAL CHIP	100	5%	1/16W
R427	1-216-833-91		10K	5%	1/16W	R518	1-216-809-11	METAL CHIP	100	5%	1/16W

Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
R524	1-216-841-11	METAL CHIP	47K	5%	1/16W	R701	1-216-833-91	RES, CHIP	10K	5%	1/16W
R526	1-216-841-11	METAL CHIP	47K	5%	1/16W	R702	1-216-821-11	METAL CHIP	1K	5%	1/16W
R529	1-216-827-11	METAL CHIP	3.3K	5%	1/16W	R703	1-216-857-11	METAL CHIP	1M	50/	4./4.034/
R530	1-216-827-11		3.3K	5%	1/16W	R703				5%	1/16W
R531	1-216-809-11		100	5%	1/16W	R704	1-216-833-91	RES, CHIP	10K	5%	1/16W
R532	1-216-827-11						1-216-833-91	•	10K	5%	1/16W
R533	1-216-809-11		3.3K 100	5% 5%	1/16W 1/16W	R706 R707	1-216-821-11		1K	5%	1/16W
				370	1/101/	- R/U/	1-216-845-11	METAL CHIP	100K	5%	1/16W
R534	1-216-823-11		1.5K	5%	1/16W	R708	1-216-864-11	METAL CHIP	0 .	5%	1/16W
R535	1-216-830-11		5.6K	5%	1/16W	R709	1-216-845-11	METAL CHIP	100K	5%	1/16W
R536	1-216-829-11		4.7K	5%	1/16W	R710	1-216-845-11	METAL CHIP	100K	5%	1/16W
R537	1-216-841-11		47K	5%	1/16W	R711	1-216-833-91		10K	5%	1/16W
R538	1-216-841-11	METAL CHIP	47K	5%	1/16W	R712	1-216-833-91	RES, CHIP	10K	5%	1/16W
R539	1-216-841-11		47K	5%	1/16W	R713	1-216-845-11	METAL CHIP	100K	5%	1/16W
R540	1-216-841-11		47K	5%	1/16W	R714	1-216-833-91	RES, CHIP	10K	5%	1/16W
R541	1-216-841-11	METAL CHIP	47K	5%	1/16W	R715	1-216-845-11	METAL CHIP	100K	5%	1/16W
R542	1-216-841-11	METAL CHIP	47K	5%	1/16W	R717	1-218-873-11	METAL CHIP	12K	0.5%	1/16W
R543	1-216-821-11	METAL CHIP	1K	5%	1/16W	R718	1-218-873-11	METAL CHIP	12K	0.5%	1/16W
R544	1-216-821-11	METAL CHIP	1K	5%	1/16W	R719	1-216-864-11	METAL CHIP	0	5%	1/16W
R545	1-216-821-11	METAL CHIP	1K	5%	1/16W	R720	1-218-871-11	METAL CHIP	10K	0.5%	1/16W
R546	1-216-791-11	METAL CHIP	3.3	5%	1/16W	R721	1-218-871-11	METAL CHIP	10K	0.5%	1/16W
R547	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R722			56	5%	1/16W
R548	1-216-821-11	METAL CHIP	1K	5%	1/16W	R723	1-216-806-11	RES, CHIP	56	5%	1/16W
R549	1-216-821-11	METAL CHIP	1K	5%	1/16W	R724	1-216-806-11	RES, CHIP	56	5%	1/16W
R550			47K	5%	1/16W	R725	1-216-806-11	RES, CHIP	56	5%	1/16W
R551		METAL CHIP	1K	5%	1/16W	R726	1-216-845-11	METAL CHIP	100K	5% 5%	1/16W
R553		METAL CHIP	10	5%	1/16W	R727	1-216-864-11	METAL CHIP	0	5%	1/16W
R554		METAL CHIP	10	5%	1/16W	R801	1-216-833-91	RES, CHIP	10K	5%	1/16W
R555	1 016 000 01	DEC CHID	101/	F0 /	44014	5000	4 040 045 44				
R556	1-216-833-91 1-216-833-91	RES, CHIP	10K 10K	5%	1/16W	R802		METAL CHIP	100K	5%	1/16W
R557	1-216-833-91		10K	5%	1/16W	R803	1-216-809-11	METAL CHIP	100	5%	1/16W
R558		METAL CHIP		5%	1/16W	R813	1-216-837-11	METAL CHIP	22K	5%	1/16W
R559	1-216-821-11		1K 1K	5% 5%	1/16W 1/16W	R814	1-216-142-00		4.7	5%	1/8W
11000	1-210 021-11	WILIAL OITH	IK .	J /0	1/1000	R818	1-216-845-11	METAL CHIP	100K	5%	1/16W
R560	1-216-821-11	METAL CHIP	1K	5%	1/16W	R819	1-216-837-11	METAL CHIP	22K	5%	1/16W
R561	1-216-821-11	METAL CHIP	1K	5%	1/16W	R821		METAL CHIP	4.7K	5%	1/16W
R562	1-216-821-11	METAL CHIP	1K	5%	1/16W	R822	1-216-845-11		100K	5%	1/16W
R563	1-216-821-11	METAL CHIP	1K	5%	1/16W	R823		METAL CHIP	100K	5%	1/16W
R564	1-219-570-11	RES, CHIP	10M	5%	1/16W	R824		METAL CHIP	100K	5%	1/16W
R565	1-216-864-11	METAL CHIP	0	5%	1/16W	R826	1-216-864-11	METAL CHIP	0	5%	1/16W
R566	1-216-821-11		1K	5%	1/16W	R829		METAL CHIP	0	5%	1/16W
R567	1-216-821-11		1K	5%	1/16W	R830		RES, CHIP	10K	5%	1/16W
R568	1-216-821-11		1K	5%	1/16W	R831		METAL CHIP	4.7K	5%	1/16W
R569		METAL CHIP	1K	5%	1/16W	R832		METAL CHIP	4.7K	5%	1/16W
R570	1-216-841-11	METAL CHIP	47K	5%	1/16W	R833	1-216-833-91	מבס רעום	10K	5%	4/4614
R571	1-216-841-11		47K	5%	1/16W	R834	1-216-833-91	,			1/16W
R572	1-216-845-11		100K	5%	1/16W	R835	1-216-829-11		10K	5%	1/16W
R573	1-216-821-11		1K	5%	1/16W	R836	1-216-829-11		4.7K 4.7K	5%	1/16W
R574	1-216-797-11		10	5%	1/16W	R837	1-216-830-11		5.6K	5% 5%	1/16W 1/16W
R575	1-216-821-11	METAL CHID	1K	5%	1/16/4/	Door					
					1/16W	R838	1-216-809-11		100	5%	1/16W
		METAL CHIP	10	5% 5%	1/16W	R839		METAL CHIP	5.6K	5%	1/16W
			10	5%	1/16W	R840	1-216-830-11		5.6K	5%	1/16W
		METAL CHIP	1K	5%	1/16W	R841	1-216-809-11		100	5%	1/16W
กง/ฮ	1-216-841-11	METAL CHIP	47K	5%	1/16W	R842	1-216-830-11	METAL CHIP	5.6K	5%	1/16W
	1-216-833-91		10K	5%	1/16W	R843	1-216-809-11	METAL CHIP	100	5%	1/16W
	1-216-841-11		47K	5%	1/16W	R844	1-216-833-91		10K	5%	1/16W
	1-216-841-11		47K	5%	1/16W	R845	1-216-809-11		100	5%	1/16W
	1-216-841-11		47K	5%	1/16W	R847	1-216-833-91		10K	5%	1/16W
R584	1-216-841-11	METAL CHIP	47K	5%	1/16W	R850	1-218-870-11		9.1K	0.5%	1/16W
R585	1-216-841-11	METAL CHIP	47K	5%	1/16W	R852	1-216-809-11	METAL CHIP	100	5%	1/16W
	1-216-837-11		22K	5%	1/16W	R853	1-216-833-91		10K	5%	1/16W
	1-216-837-11		22K	5%	1/16W	R854	1-216-833-91		10K	5%	1/16W
								0, 0	1011	0 /0	17 10 44

Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
R855	1-218-870-11	METAL CHIP	9.1K	0.5%	1/16W	R952	1-216-821-11	METAL CHIP	1K	5%	1/16W
R856	1-216-829-11		4.7K	5%	1/16W	R953		METAL CHIP	1K	5%	
R857	1-216-809-11	METAL CHIP	100	5%	1/16W	R954	1-216-821-11	METAL CHIP	1K	5%	1/16W
R858	1-218-707-11		4.3K	5%	1/16W	R955	1-216-821-11		1K	5%	1/16W
R859	1-216-809-11		100	5%	1/16W	R956	1-216-821-11	METAL CHIP	1K	5%	1/16W
R860	1-218-707-11	RES, CHIP	4.3K	5%	1/16W	R957	1-216-842-11	METAL CHIP	56K	5%	1/16W
R861	1-216-829-11		4.7K	5%	1/16W	R958		METAL CHIP	100K	5%	1/16W
R863	1-216-829-11		4.7K	5%	1/16W	R959	1-216-833-91	RES, CHIP	10K	5%	1/16W
R866	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R960	1-216-841-11	METAL CHIP	47K	5%	1/16W
R867	1-216-809-11		100	5%	1/16W	R962	1-216-821-11		1K	5%	1/16W
R868	1-216-829-11		·4.7K	5%	1/16W	R964	1-216-838-11		27K	5%	1/16W
R869	1-216-829-11		4.7K	5%	1/16W	R965	1-216-833-91		10K	5%	1/16W
R870	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R966	1-216-833-91	BES CHIP	10K	5%	1/16W
R871	1-216-829-11		4.7K	5%	1/16W	R967	1-216-833-91		10K		
						1				5%	1/16W
R872	1-216-829-11		4.7K	5%	1/16W	R968	1-216-841-11		47K	5%	1/16W
R873	1-216-829-11		4.7K	5%	1/16W	R969	1-216-841-11	METAL CHIP	47K	5%	1/16W
R875	1-218-839-11	METAL CHIP	470	0.5%	1/16W	R970	1-216-833-91	RES, CHIP	10K	5%	1/16W
R876	1-218-839-11		470	0.5%	1/16W	R971	1-216-813-11		220	5%	1/16W
R877	1-218-839-11	METAL CHIP	470	0.5%	1/16W	R972	1-216-841-11	METAL CHIP	47K	5%	1/16W
R878	1-218-839-11	METAL CHIP	470	0.5%	1/16W	R973	1-216-864-11	METAL CHIP	0	5%	1/16W
R879	1-216-864-11		0	5%	1/16W	R974	1-216-838-11		27K	5%	1/16W
R880			330	5%		1					
noou	1-216-815-11	WETAL UNIP	330	376	1/16W	R975	1-216-841-11	WETAL CHIP	47K	5%	1/16W
R881	1-216-815-11	METAL CHIP	330	5%	1/16W	R976	1-216-845-11	METAL CHIP	100K	5%	1/16W
R885	1-216-829-11	METAL CHIP	4.7K	5%	1/16W	R978	1-216-813-11		220	5%	1/16W
R886	1-216-829-11		4.7K	5%	1/16W	R980	1-216-841-11		47K	5%	1/16W
R887	1-216-833-91		10K	5%	1/16W	R981	1-216-821-11		1K		
						1				5%	1/16W
R888	1-216-827-11	WE IAL UNIP	3.3K	5%	1/16W	R982	1-216-833-91	RES, CHIP	10K	5%	1/16W
R889	1-216-830-11	METAL CHIP	5.6K	5%	1/16W	R983	1-216-821-11	METAL CHIP	1K	5%	1/16W
R890	1-216-849-11	METAL CHIP	220K	5%	1/16W	R984	1-216-833-91		10K	5%	1/16W
R891	1-208-813-11		20K	0.5%	1/10W	R986	1-216-864-11		0	5%	1/16W
R893											
	1-216-833-91		10K	5%	1/16W	R987	1-216-864-11		0	5%	1/16W
R894	1-216-809-11	METAL CHIP	100	5%	1/16W	R988	1-216-864-11	METAL CHIP	0	5%	1/16W
R895	1-216-833-91	RES. CHIP	10K	5%	1/16W	R989	1-216-864-11	METAL CHIP	0	5%	1/16W
R896	1-216-809-11		100	5%	1/16W	R990	1-216-813-11		220	5%	1/16W
R897	1-208-813-11										
			20K	0.5%	1/10W	R991	1-216-813-11		220	5%	1/16W
R899	1-216-864-11		0	5%	1/16W	R992	1-216-864-11		0	5%	1/16W
R903	1-216-864-11	METAL CHIP	0	5%	1/16W	RR001	1-216-864-11	METAL CHIP	0	5% (r	1/16W DSR-20MD)
R904	1-216-838-11	METAL CHIP	27K	5%	1/16W					(L	J. LOND)
R905			1K	5%	1/16W	RR002	1-216-864-11	METAL CLUD	0	En/	4 /4 () [4]
						nnouz	1-210-004-11	WE TAL CHIP	0	5%	1/16W
R906		METAL CHIP	33K	0.5%	1/16W					(DS	SR-20MDP)
R907 R908	1-216-821-11 1-216-864-11		1K 0	5% 5%	1/16W 1/16W			AVADIADI E DE	CICTOR		
11300	1-210-004-11	WE TAL OTHE	U	J 70	1/1000			< VARIABLE RE	31910K >		
R909	1-216-864-11	METAL CHIP	0	5%	1/16W	RV001	1-238-855-11	RES, ADJ, CERN	NET 4.7K		
R913	1-216-821-11	METAL CHIP	1K	5%	1/16W			,		D CONV	REF REG1)
R918	1-216-821-11		1K	5%	1/16W	RV002	1_000.055.11	RES, ADJ, CERN	ACT A ZV	D OONV.	ner near,
						RVUUZ	1-200-000-11	neo, ADJ, CENI			
R919 R920	1-216-864-11 1-218-871-11	METAL CHIP	0 10K	5% 0.5%	1/16W 1/16W	RV010	1-238-854-11	RES, ADJ, CERN		D CONV.	REF REG2)
									(0		REF REG)
R921	1-218-875-11		15K	0.5%	1/16W	RV011		RES, ADJ, CERN		LAMP RI	EF REG)
R923	1-216-178-00	RES, CHIP	150	5%	1/8W	RV012	1-238-854-11	RES, ADJ, CERN	TET 2.2K		
R924	1-216-831-11	METAL CHIP	6.8K	5%	1/16W				(0	B CLAMP	REF REG)
R926		METAL CHIP	1M	5%	1/16W				(*		
R936	1-216-845-11		100K	5%	1/16W	RV201	1-238-855-11	RES, ADJ, CERM	ET A 7V		
						117401	- 200-000-11	TIEU, ADU, UEMIV		C PICTUR	RE FRAME)
R937	1-216-845-11	METAL CHIP	100K	5%	1/16W			•	,		•
R938	1-216-847-11		150K	5%	1/16W			< VIBRATOR >			
R939	1-216-833-91		10K	5%	1/16W						
R941						V404	1 760 655 04	VIDDATOR ORY	OTAL (DOL	17-1	
	1-216-864-11		0	5%	1/16W	X421		VIBRATOR, CRY			
R943	1-216-864-11	METAL CHIP	0	5%	1/16W	X422		VIBRATOR, CRY			
						X501		VIBRATOR, CER			
R944	1-216-833-91	RES, CHIP	10K	5%	1/16W	X502		VIBRATOR, CRY			
R950	1-216-821-11	•	1K	5%	1/16W	X701		VIBRATOR, CRY			
R951	1-216-821-11		1K	5%	1/16W			The state of the other	- 11 LE (ET.U	. 0141112)	•
.1001	. 2.3 021 11	OIIII		U / U	171011						

JC-	19 MD-	63 MD-64	MD-65	POWE	ER E	BLOCK (I	J-1)			
Ref. No. X702 X801	Part No. 1-760-497-21 1-767-779-21	<u>Description</u> VIBRATOR, LITHIUM NIO VIBRATOR, CRYSTAL (49	BATE (6MHz)	Remark Re	f. No.	Part No.	Description < CONNECTOR	>		Remark
	1707 170 21	VIBILATOR, OTTOTAL (49	.152141112)		CN002	1-770-692-11	CONNECTOR, F	FC/FPC 9P		
		MD-63 BOARD (Ref No. 6	,000 Series)				< IC >		•	
C101 C102	1-126-206-11 1-163-031-11	< CAPACITOR >			C003 C004 C005 C006	8-719-820-44 8-759-510-71	IC PHOTO COL IC PHOTO COL IC BA10358F-I IC ELEMENT, I	JPLER TLP E2	907-0 (S	ONY2)
	, , , , , , , , , , , , , , , , , , , ,	< CONNECTOR >	3	0			< JUMPER RES	ISTOR >		
CN101	1-770-646-11	CONNECTOR, FFC/FPC 16	P	1	JR001 JR002	1-216-296-91 1-216-296-91		0 0	.*	
		< DIODE >					< TRANSISTOR	>		
D101	8-719-989-52	DIODE GL4600S			2001	8-729-012-46	PHOTO TRANSI	STOR PT46	OOES	
		< IC >					< RESISTOR >	010111140	,0010	
IC101 IC102 IC103 IC105	8-719-820-44 8-759-510-71	IC PHOTO COUPLER TLP IC PHOTO COUPLER TLP IC BA10358F-E2 IC ELEMENT, HALL THS1	907-0 (SONY2	2) F F	R002 R003 R004 R005 R007	1-216-031-00 1-216-069-00 1-216-081-00 1-216-069-00 1-216-081-00	METAL CHIP METAL CHIP METAL CHIP METAL CHIP	180 6.8K 22K 6.8K 22K	5% 5% 5% 5%	1/10W 1/10W 1/10W 1/10W 1/10W
		< JUMPER RESISTOR >			8008	1-216-047-91		820	5%	1/10W
JR101 JR102	1-216-296-91 1-216-296-91	SHORT 0 SHORT 0				. 2.0 0.7 01	< VARIABLE RES		3 /0	1/1049
JR103 JR104 JR105	1-216-296-91 1-216-296-91 1-216-296-91	SHORT 0 SHORT 0			IV001 IV002		RES, ADJ, CERN RES, ADJ, CERN	1ET 47K		
		< TRANSISTOR >					< SWITCH >			
Q102	8-729-012-46	PHOTO TRANSISTOR PT46	600FS	s	002	1-762-558-11	SWITCH, PUSH	(C DOWN)		
R101 R102 R103		METAL CHIP 22K	5% 1/	10W 10W			MD-65 BOARD (I		- 00 Series	3)
R107 R108	1-216-069-00 1-216-047-91	METAL CHIP 6.8K	5% 1/ ⁻ 5% 1/ ⁻	C	N201 N202	1-766-830-21 1-774-771-11	< CONNECTOR > CONNECTOR, FFI CONNECTOR, FFI	C/FPC (ZIF) C/FPC 14P	11P	
11100	7 210 001 00	< VARIABLE RESISTOR >	J /6 1/		N203 N204		PIN, CONNECTOR CONNECTOR (MI			
RV101	1-238-858-11	RES, ADJ, CERMET 47K					< JUMPER RESIS	STOR >		•
RV102	1-238-862-11	RES, ADJ, CERMET 1M		J	R201	1-216-296-91	SHORT	0	—— <u>——————————————————————————————————</u>	
S101		< SWITCH > SWITCH, PUSH (1 KEY)		<u>↑</u>		1-468-441-11 1-468-442-11	POWER BLOCK (U-1) (DSR-	20MD) 20MDP)	
		MD-64 BOARD (Ref No. 7,0	00 Series)				< CAPACITOR >		No. 10,0	00 Series)
		< CAPACITOR >					FILM	0.1uF	20%	250V
C001 C002	1-126-206-11 1-163-031-11	ELECT CHIP 100uF CERAMIC CHIP 0.01uF	20% 6.3 50\	V		1-104-706-11		0.22uF	(DS 20%	SR-20MD) 250V
			001	△ C2 △ C3 △ C4	1	1-104-705-11 1-115-383-11 1-115-383-11		0.1uF 0.001uF 0.001uF	20% 10% 10%	R-20MDP) 250V 125V 125V
				6-24	mark /	mponents identi or dotted lin are critical for e only with par crified.	e with marque safety. sécurité t num- Ne les r	mposants ic	ritiques ue par u	pour la

POWER BLOCK (U-1)

											(0 1)
Ref. No	o. Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
			•			D5	8-719-110-72		030ESB2		Homan
△ C5	1-115-383-11	CERAMIC	0.001uF	10%	125V	D6	8-719-110-72		30ESB2		
△ C6	1-104-705-11		0.1uF	20%	250V			51052 112			
C7	1-115-383-11		0.001uF	10%	125V	D8	8-719-053-20	DIODE UF	4003P		
C8	1-115-383-11		0.001uF	10%	125V	D9	8-719-109-85	DIODE RE	5.1ESB2		
 C 9	9-880-364-01	ELECT	470uF		200V	D10	8-719-510-37		LC20U		
				(E	SR-20MD)	D11	8-719-109-97		06.8ESB2		
 ∆ C9	1-117-188-11	ELECT	150uF	20%	4001	D12	8-719-110-41	DIODE RE	15ESB2		
212 03	1-117-100-11	ELECT	TOULF		400V SR-20MDP)			- TD AMOIOT	OD .		
C10	9-880-365-01	FILM	0.01uF	(DC	630V	i		< TRANSIST	UK >		
C11	9-880-366-01		680PF (D	SR-20M		Q1	8-729-037-96	TRANSISTO	B 25K2366 (F	15B-20M	ות
C11	9-880-424-01		330PF (D		,	Q1	9-880-423-01	TRANSISTO	R 2SK2483 (E	SR-20M	DP)
C12	9-880-366-01	CERAMIC	680PF (D			Q3	8-729-281-53	TRANSISTO	R 2SC1815-G	R	2. /
C12	9-880-424-01		330PF (D					< RESISTOR	>		
C13	1-107-929-11		10uF	20%	100V						
C14 C15	1-107-929-11 1-126-387-11		10uF	20%	100V	△R1	9-880-373-01	METAL OXID			DSR-20MD)
C16	9-880-367-01		2.2uF 0.1uF	20%	100V	△R1	9-880-427-01	METAL OXID	E 330K	1W	
010	3-000-001-01	I ILIVI	O. Fur			R2	9-880-374-01	METAL OXID	E 901/		SR-20MDP)
C17	9-880-368-01	FILM	470PF			R2	9-880-428-01			2W (DSR-20MD)
C18	9-880-369-01		0.15uF			11.2	3 000 720 01	MILIAL OXID	L 2/0K		SR-20MDP)
C19	9-880-370-01	FILM	0.033uF			R3	1-212-865-00	FUSEBLE	22	5%	1/4W F
C20	1-107-929-11		10uF	20%	100V					0,10	.,
C21	9-880-371-01	CERAMIC	0.001uF			R4	1-247-879-11	CARBON	100K	5%	1/4W
000	0.000.001.01					R5	1-215-884-11			5%	2W
C22	9-880-371-01		0.001uF			R6	1-215-880-11			5%	2W
C25	9-880-371-01		0.001uF (I			R7	9-880-375-01				DSR-20MD)
C26	9-880-426-01 9-880-371-01		0.001uF ([0.001uF ([R7	9-880-650-01	METAL OXID	E 150K	2W	00.001400
C26	9-880-426-01		0.001uF (I							(L	SR-20MDP)
	0 000 120 01	oc. ii iiviio	0.00141 (1	2011 2011	101)	R8	9-880-376-01	METAL OXID	E 0.22	5\M (1	OSR-20MD)
C27	1-111-066-11	ELECT	820uF	20%	25V	R8		METAL OXID		5W (1	Jon-Zulvid)
C28	1-111-066-11		820uF	20%	25V		0 000 120 01	METAL ONID	L 0.47		SR-20MDP)
C29	9-880-367-01		0.1uF			R9	1-216-377-11	METAL OXIDI	E 4.7	5%	2W
C31	9-880-367-01		0.1uF			R10	1-215-462-00	METAL	5.1K	1%	1/4W
C32	9-880-367-01	FILM	0.1uF			R11	1-249-411-11	CARBON	330	5%	1/4W
		- CONNECTOR -						19.2		• (DSR-20MD)
		< CONNECTOR >				R11	1 040 417 11	CADDON	41/	F0/	4 (4)41
* CN1	1-580-230-31	PIN, CONNECTO	R (FOR ROA	RD) 2P		KII	1-249-417-11	CARRON	1K	5%	1/4W
* CN3	9-880-386-01	BOARD IN HARN	IESS 4P	10/21		R12	1-215-383-00	METAL	27		SR-20MDP) 1/4W
						R13	1-215-385-00		33	1% 1%	1/4W
		< FUSE >				R14	1-249-430-11		12K	5%	1/4W
						R16	1-215-884-11			5%	2W
△ F1		FUSE (3.15A/12									
△F1	9-882-875-01	FUSE (T1.6AL/2	50V) (DSR-	20MDP)		R18	1-247-847-11		4.7K	5%	1/4W
 ∆ F2	9-880-385-01	FUSE (3.15A/12	5V) (DSR-20	OMD)		R19	1-247-847-11		4.7K	5%	1/4W
 △ F 2	9-002-0/0-01	FUSE (T1.6AL/2	50V) (DSR-2	ZUMDP)		R20	1-247-839-11		2.2K	5%	1/4W
		< COIL >				R21	1-249-417-11		1K	5%	1/4W
		COOIL			.	R22	1-247-843-11	CARBON	3.3K	5%	1/4W
∆L1	9-880-379-01	INDUCTOR 6mH	(DSR-20MD)	1	R23	1-249-399-11	CARRON	33	5%	1/4W
∆L1		INDUCTOR 15ml				R24	1-215-425-00		1.5K	1%	1/4 VV 1/4W
∆L2		INDUCTOR 5.6ml			i	R25	1-215-425-00		1.5K	1%	1/4W
△L2	9-880-432-01	INDUCTOR 22mH	(DSR-20MI	OP)	İ	R26	1-215-433-00		3.3K	1%	1/4W
L3	9-880-381-01	INDUCTOR			Ì	R27	1-249-417-11	CARBON	1K	5%	1/4W
1.4	0.000.000.00	INDUCTOR						•			
L4	9-880-382-01					R28	9-880-375-01				SR-20MD)
L5 L6	9-880-382-01 9-880-383-01					R28		METAL OXIDE			SR-20MD)
LU	9-000 - 000*01	אטוסטמוייו				R29	9-880-651-01	WETAL UXIDE	680	2W	
	a part a f							VADIADIED	ECICTOR -		
		< DIODE >						< VARIABLE R	ESISTUK >		
	* * * * * * * * * * * * * * * * * * * *					RV1	9-880-377-01	RES. VAAR C	ABBON 3K		
△ D1	8-719-500-58				ļ						
D4	8-719-979-63		(DSR-20MD								
D4	8-719-053-19	DIUDE UF4007	G23 (DSR-2	OMDP)	Į.						

The components identified by mark \triangle or dotted line with mark \triangle are critical for safety. Replace only with part number specified.

Les composants identifiés par une marque Δ sont critiques pour la sécurité. Ne les remplacer que par une pièce portant le numéro spécifié.

POWER BLOCK (U-1) POWER BLOCK (U-2)

Ref. No.	•										
	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
		< TRANSFORME	R >			C41	9-880-399-01		0.047uF		1101110111
∆T1	9-885-000-10	TRANSFORMER	, POWER (I	OSR-20N	1D)	C42	9-880-399-01	FILM	0.047uF		
∆T1	9-885-000-13	TRANSFORMER	, POWER (I	DSR-20N	IDP)	C44	1-117-154-11	ELECT	33uF	20%	16V
						C45	9-880-399-01	FILM	0.047uF		
		< THERMISTOR	>			C46	9-880-652-01	ELECT	100uF	20%	16V
						C47	1-117-154-11	ELECT	33uF	20%	16V
TH1	9-880-384-01	THERMISTOR	8								
						C48	9-880-652-01	ELECT	100uF	20%	16V
		< IC >				C49	9-880-399-01		0.047uF		. • •
						C50	9-880-399-01		0.047uF		
Z1	9-885-000-09	IC FA531	6P			C51	9-880-400-01		- 330uF		25V
△Z2	8-749-924-80	PHOTO COUPLE	R PS25611	L1-1-V		C54	9-880-406-01		680uF		10V
△ Z3		PHOTO COUPLE									
						C55	9-880-407-01	CERAMIC	0.0022uF		
						C56	1-115-730-11		180uF	20%	10V
Δ	1-468-441-11	POWER BLOCK	(U-2) (DSR	-20MD)		C57	9-880-403-01		0.1uF	2070	50V
\triangle						C58	9-880-399-01		0.047uF		001
		*****		zome.		C59	9-880-407-01		0.0022uF		
			(Ref	f No. 20	000 Series)	.000	0 000 401 01	OLITAVIO	0.002241		
			10		_ 50 501100)	C60	1-115-566-11	CERAMIC	4.7uF	10%	10V
		< CAPACITOR >				C61	1-115-566-11		4.7uF	10%	10V 10V
						C62	1-115-566-11		4.7uF 4.7uF	10%	10V
C1 -	1-115-781-11	ELECT	220uF	20%	25V	C63	1-115-566-11		4.7uF	10%	10V
C2	1-115-781-11	ELECT	220uF	20%	25V	C64	1-107-682-11		1uF	10%	16V
C3	9-880-399-01		0.047uF	2070	201	004	1 101-002-11	OLIMBIO	ıuı	10/0	100
C4	9-880-399-01		0.047uF			C65	1-115-566-11	CERAMIC	4.7uF	10%	10V
C5	9-880-400-01		330uF		25V	C66	1-115-566-11		4.7uF	10%	10V
00	0 000 400 01	LLLOI	oodui		201	C67	1-115-566-11		4.7uF	10%	10V
C6	9-880-401-01	FII M	0.47uF			C68	1-115-566-11		4.7uF	10%	10V
C7	1-115-787-11		820uF	20%	25V [.]	C69	1-115-566-11		4.7uF	10%	10V
C8	1-115-787-11		820uF	20%	25V 25V	009	1-110-000-11	OFTIVINIO	4./ ur	10 70	100
C10	9-880-403-01		0.1uF	£U /0	50V	C70	1-115-566-11	CEDAMIC	4 7nE	10%	101/
C11	1-115-787-11		820uF	20%	25V	C71	1-115-566-11		4.7uF 4.7uF	10%	10V 10V
311	. 110-101-11	LLLUI	OZOUF	20 /0	201	C72	1-115-566-11				
C12	9-880-404-01	FLECT	470uF		25V				1uF	10%	16V
C12	9-880-405-01		0.001uF		20 V	C73	1-107-682-11	CENAIVIIC	1uF	10%	16V
C14	1-115-737-11		0.001uF	200/	101/			COMMITTOTO			
C14 C15	1-115-737-11			20%	10V			< CONNECTOR	>		
C16			0.001uF	20%	10V	4 0814	0.000 447 04	DIM CONTROL	OD 45		
010	9-880-406-01	ELEVI	680uF		10V	* CN1	9-880-417-01	PIN, CONNECT			
C17	9-880-402-01	EH AA	0.100		EOV	* CN2	1-506-485-11	PIN, CONNECT			
	1 101 010 11		0.1uF	200/	50V	* CN10	1-506-487-11	•			
C18	1-124-942-11		180uF	20%	10V	* CN11		PIN, CONNECT			
C19	9-880-399-01		0.047uF			* CN12	1-506-48/-11	PIN, CONNECT	0K 8P		
C20	9-880-399-01		0.047uF			# ONE	4 500 401 41	0011150505	ь.		
C21	9-880-399-01	LITIM	0.047uF			* CN13	1-506-481-11	CONNECTOR 2	۲ .		
000	0 000 400 04	TH M	0.45		501/			ÉLIOE			
C22	9-880-402-01		0.1uF		50V			< FUSE >			
	9-880-399-01		0.047uF	0001	4071	A ===	0.000 445 5:	FLIDE (TEL 1055			
	1-115-730-11		180uF	20%	10V	 ♠ F2	9-880-416-01	FUSE (T5A/250	V)		
	9-880-404-01		470uF		25V						
C26	9-880-399-01	FILIVI	0.047uF					< COIL >			
	9-880-407-01		0.0022uF			L1		INDUCTOR 0.5			
	9-880-406-01		680uF		10V	L2		INDUCTOR 1ml	•		
	9-880-406-01		680uF		10V	L3		INDUCTOR 150			
	9-880-402-01		0.01uF			L4		COIL, FERRITE			
C31	9-880-406-01	ELECT	680uF		10V	L5	9-880-411-01	INDUCTOR 150	uH		
						e se la companya de l					
	1-115-785-11	ELECT	470uF		25V	L6	9-880-412-01	INDUCTOR 15u	Н		
	1-124-534-11		680uF	20%	16V	L8	9-880-413-01	INDUCTOR 330	иH		
	9-880-402-01	FILM	0.01uF			L9		INDUCTOR 15u			
005	1-115-754-11	ELECT	120uF	20%	16V	L11		INDUCTOR 150			
	9-880-403-01	FILM	0.1uF		50V	L12	9-880-415-01	INDUCTOR			
C36											
C36 C37	1-115-730-11		180uF	20%	10V	L13	9-880-412-01	INDUCTOR 15u	Н		
C36 C37 C38	9-880-399-01	FILM	0.047uF	20%	10V	L15		INDUCTOR 15u INDUCTOR 15u			
C36 C37 C38 C39		FILM FILM		20%	10V			INDUCTOR 15u	Н		

The components identified by mark \triangle or dotted line with mark \triangle are critical for safety. Replace only with part number specified.

Les composants identifiés par une marque ∆ sont critiques pour la sécurité.
Ne les remplacer que par une pièce portant le numéro spécifié.

							L				
5		5 1.1	•		D t	D (N-	Deat No.	Describetors			Demonstr
Ref. No.	Part No.	Description			<u>Remark</u>	Ref. No.	Part No.	<u>Description</u>			<u>Remark</u>
	•	< DIODE >				Z3	8-759-293-98	IC LM2576T-ADJ	LB03		
						Z4	9-880-397-01	IC uPC29L03J			
D1	8-719-500-70	DIODE D5S4M									
D2	8-719-109-89					Z5	8-759-069-28	IC PQ05RF11			
D3	8-719-500-70					Z6		IC TD62305AP			
D4	8-719-043-76					Z7		IC LM2575T-ADJ	U B03		
D5	8-719-018-83					Z8		IC LM2575T-ADJ			
D3	0-7 13-010-03	DIODE DEGAIN				Z10		IC uPC79M05HF			
DC	0.740.040.00	DIODE DOCAM				210	0-709-009-00	IC BECT SIMILISHE			
D6	8-719-018-83					744	0.750.000.04	10 00000044			
D7	8-719-107-94					Z11	0-759-090-24	IC PQ30RV11			
D8	8-719-107-94	DIODE 1SS202	2-1								
						-					
		< TRANSISTOR >				*	A-7073-470-A	RE-32 BOARD, C			

Q1		TRANSISTOR 2S				1			(He	t.No. /	,000 Series)
Q2		TRANSISTOR 2S.				1					
Q3		TRANSISTOR 2S						< CAPACITOR >			
Q4	8-729-203-76	TRANSISTOR 2S	C3328-Y								
						C101	1-163-243-11	CERAMIC CHIP	47PF	5%	50V
		< RESISTOR >				C102	1-163-243-11	CERAMIC CHIP	47PF	5%	50V
						C103	1-163-243-11	CERAMIC CHIP	47PF	5%	50V
R1	1-247-839-11	CARBON	2.2K	5%	1/4W	C104		CERAMIC CHIP	47PF	5%	50V
R2	1-215-454-00		24K	1%	1/4W	C105		CERAMIC CHIP	47PF	5%	50V
R3	1-215-427-00		1.8K	1%	1/4W	0,50				- , 4	
R4	1-215-429-00		2.2K	1%	1/4W	C106	1-163-243-11	CERAMIC CHIP	47PF	5%	50V
R5		METAL OXIDE	560	5%	1/4/V	C108		CERAMIC CHIP	0.022uF	10%	25V
กง	1-210-401-11	MIL INC OVIDE	500	J /0	IAA	C108		CERAMIC CHIP	0.022ur 0.1uF	10/0	25V 25V
De	1 015 057 14	METAL OXIDE	10	50/	1W	0110	1-109-030-81	CENAIVIIC CHIP	0.146		200
R6				5%				COMMECTOR			
R7	1-215-443-00		8.2K	1%	1/4W			< CONNECTOR >			
R8	1-215-416-00		620	1%	1/4W						
R9	1-215-429-00		2.2K	1%	1/4W	CN101		CONNECTOR, FP			
R10	1-216-448-11	METAL OXIDE	39	5%	2W	CN102	1-764-129-11	CONNECTOR, FP	C 15P		
R11	1-249-417-11		1K	5%	1/4W			< DIODE >			
R12	1-249-424-11	CARBON	3.9K	5%	1/4W						
R14	1-249-417-11	CARBON	1K	5%	1/4W	D101	8-719-421-59	DIODE MA3075	WA- (TX)		
R15	1-247-847-11	CARBON	4.7K	5%	1/4W	D102	8-719-421-59	DIODE MA3075	WA- (TX)		
R16	1-249-402-11		56	5%	1/4W	D103		DIODE MA3075			
			•	- , ,		D104		DIODE MA3075			
R17	1-249-402-11	CARRON	56	5%	1/4W	D105		DIODE MA3075			
R18	1-249-402-11		56	5%	1/4W	1 5100	0 7 10 421 00	DIODE WINDOID	1177		
R19	1-247-847-11		4.7K	5%	1/4W	D106	8-710-/21-50	DIODE MA3075	\A/A_ /TY\		
R21	1-249-417-11		1K	5%	1/4W	D107		DIODE MA3075			
								DIODE MA3075			
R22	1-215-857-11	METAL OXIDE	10	5%	1W	D108					
200		AAFTAI	••	401	4 (1) 14	D109		DIODE MA3075			
R23	1-215-387-00		39	1%	1/4W	D110	ŏ-/19-421-59	DIODE MA3075	WA- (TX)		
R24	1-215-408-00		300	1%	1/4W						
R25	1-215-405-00		220	1%	1/4W	D111		DIODE MA3075			
R26	1-215-431-00		2.7K	1%	1/4W	D112		DIODE MA3075			
R27	1-215-449-00	METAL	15K	1%	1/4W	D113	8-719-421-59	DIODE MA3075	WA- (TX)		
						D114		DIODE MA3075			
R28	1-215-430-00	METAL	2.4K	1%	1/4W	D115		DIODE MA3075			
R29	1-247-841-11		2.7K	5%	1/4W				(,		
R30	1-249-417-11		1K	5%	1/4W			< FERRITE BEAD	>		
R31	1-215-447-00		12K	1%	1/4W			TEMINIC DEAD			
						ED404	1-500 044 00	EEDDITE	Out		
R32	1-215-423-00	WEIAL	1.2K	1%	1/4W	FB101	1-500-241-22		OuH		
201	4 045 155 55	LACTA!	4.09	461	47,000	FB102	1-500-241-22		0uH		
R33	1-215-423-00		1.2K	1%	1/4W	FB103	1-500-241-22		0uH		
R43		METAL OXIDE	10	5%	1W	FB104	1-500-241-22		0uH		
R45	1-215-443-00		8.2K	1%	1/4W	FB105	1-500-241-22	FERRITE	0uH		
R46	1-215-413-00	METAL	470	1%	1/4W						
R47	1-215-429-00		2.2K	1%	1/4W	FB106	1-500-241-22	FERRITE	0uH		
						1					
R48	1-247-839-11	CARBON	2.2K	5%	1/4W			< JACK >			
R49	1-249-402-11		56	5%	1/4W						
R50	1-247-855-11		10K	5%	1/4W	J101	1-694-410-11	TERMINAL BOAR	D		
R53	1-215-857-11		10	5%	1W	"."		William BOAH		UTPLIT	(MONITOR)
	0 007 11			- /0					\ J 17 U	J., 011	
	,	< IC >	-					< JUMPER RESIS	TOR >		
Z1	8-759-164-80	IC LM2577-ADJ				JR101	1-216-296-91	SHORT	0		
Z2	8-759-520-49	IC PQ30RV21				JR102	1-216-296-91		0	4.1	

RE-32 RP-228

Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
JR103	1-216-296-91		0						0.45		
JR103	1-216-296-91		0			C772		CERAMIC CHIP	0.1uF		25V
JR104	1-216-296-91		0			C773	1-113-619-11	CERAMIC CHIP	0.47uF		10V
30100	1-210-290-91	SHUNI	U			0774	1 164 260 11	CEDAMIC CHIP	0.4		4011
JR106	1-216-296-91	SHORT	0			C774 C775		CERAMIC CHIP	0.1uF 0.01uF	100/	16V
JR107	1-216-296-91		0			C776		CERAMIC CHIP		10%	25V
JR108	1-216-296-91		0						47PF	5%	50V
JR109	1-216-296-91		0			C777		CERAMIC CHIP	47PF	5%	50V
JR110	1-216-296-91		0			C778	1-162-974-11	CERAMIC CHIP	0.01uF		50V
311110	1-210-290-91	SHUNI	U .			C779	1 160 074 11	CERAMIC CHIP	0.01uF		50V
JR111	1-216-296-91	SHORT	0		•	C780		CERAMIC CHIP	0.01uF		50V 50V
JR112	1-216-296-91		0			C780		CERAMIC CHIP	0.01uF		50V 50V
JR113	1-216-296-91		0			C782		CERAMIC CHIP	0.01uF		25V
JR114	1-216-296-91		0			C783		TANTALUM CHIP		20%	4V
JR115	1-216-296-91		0			0700	1-133-201-11	TAINTALUM GHIF	Tour	2076	40
011110	1 210 200 01	OHOIH	•			C784	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
JR116	1-216-296-91	SHORT	0			C786		CERAMIC CHIP	0.01uF	10 /0	50V
JR117	1-216-296-91		0			C788		TANTALUM CHIP		20%	10V
JR118	1-216-296-91		0			C789	1-162-974-11		0.01uF	20 /0	50V
JR119	1-216-296-91		Ö			C791		CERAMIC CHIP	0.01uF		50V
JR120	1-216-296-91		Ô			0.01	1 102 07 1 11	OLI II III III OI III	0.0141		00 V
			•			C792	1-119-750-11	TANTALUM CHIP	22uF	20%	6.3V
JR121	1-216-296-91	SHORT	0			C793		CERAMIC CHIP	0.1uF	10%	16V
JR122	1-216-296-91	SHORT	0			C794		CERAMIC CHIP	0.1uF	10%	16V
JR123	1-216-296-91		0			C795	1-128-004-11		10uF	20%	16V
						C796		CERAMIC CHIP	0.01uF		50V
		< RESISTOR >									
						C797	1-107-826-91		0.1uF	10%	16V
R105	1-216-295-91		0			C798	1-162-974-11	CERAMIC CHIP	0.01uF		50V
R106	1-216-022-00		75	5%	1/10W	C799	1-164-217-11	CERAMIC CHIP	150PF	5%	50V
R107	1-216-022-00		75	5%	1/10W	C803	1-164-217-11	CERAMIC CHIP	150PF	5%	50V
R108	1-216-022-00		75	5%	1/10W	C811	1-113-619-11	CERAMIC CHIP	0.47uF		10V
R109	1-216-295-91	SHORT	0								
						C813		CERAMIC CHIP	0.01uF		50V
R110	1-216-295-91		0			C814		CERAMIC CHIP	0.01uF	10%	25V
R111	1-216-295-91		0			C815		CERAMIC CHIP	0.001uF	10%	50V
R112	1-216-295-91		0			C816		CERAMIC CHIP	0.001uF	10%	50V
R113	1-216-295-91		0		i	C817	1-104-851-11	TANTALUM CHIP	10uF	20%	10V
R114	1-216-295-91	SHORT	0			2012	4 400 004 44				
D11E	1 010 005 01	CHODT	0		ļ	C818		CERAMIC CHIP	0.001uF	10%	50V
R115	1-216-295-91	SHUKI	U .		İ	C819		CERAMIC CHIP	0.01uF		50V
		< SWITCH >				C821		CERAMIC CHIP	0.1uF		16V
		< 24/11/u>				C822		CERAMIC CHIP	0.1uF		16V
S101	1-570-974-11	SWITCH, SLIDE (SVNC)			C823	1-104-300-11	CERAMIC CHIP	0.1uF		16V
0101	10/00/411	OWITOII, OLIDE (OTNO)			C824	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
						C825	1-164-315-11		470PF	5%	50V
*	Δ-7067-132-Δ	RP-228 BOARD,	COMPLETE	(DSR-20	IMD)	C826		CERAMIC CHIP	0.01uF	J /6	50V
*		RP-228 BOARD,				C827		CERAMIC CHIP	0.01uF		50V
		********			, ,	C828	1-162-974-11		0.01uF		50V
			(Re	f.No. 3.0	00 Series)	0020	. 102 07 111	ornamo om	0.0141		001
						C829	1-135-201-11	TANTALUM CHIP	10uF	20%	4V
	1-776-149-11	CABLE, FLEXIBLE	FLAT 30P			C830	1-164-360-11		0.1uF		16V
		CABLE, FLEXIBLE		-245)		C831		TANTALUM CHIP		20%	6.3V
				,		C832		CERAMIC CHIP	0.1uF		16V
	*	< CAPACITOR >				C833	1-162-974-11		0.01uF		50V
C146	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C834	1-162-974-11	CERAMIC CHIP	0.01uF		50V
C148	1-162-974-11	CERAMIC CHIP	0.01uF		50V	C835	1-162-974-11	CERAMIC CHIP	0.01uF		50V
C701	1-164-174-11	CERAMIC CHIP	0.0082uF	10%	25V	C836	1-162-974-11	CERAMIC CHIP	0.01uF		50V
C702	1-162-967-11	CERAMIC CHIP	0.0033uF	10%	50V	C837	1-162-970-11		0.01uF	10%	25V
C703	1-164-174-11	CERAMIC CHIP	0.0082uF	10%	25V	C838	1-162-913-11	CERAMIC CHIP	8PF	0.5PF	50V
	1-162-967-11			10%	50V	C839	1-162-913-11		8PF	0.5PF	50V
	1-164-173-11			10%	50V	C841	1-162-923-11		47PF	5%	50V
C706	1-164-173-11		0.0039uF	10%	50V	C842	1-164-360-11		0.1uF		16V
C761		CERAMIC CHIP	0.1uF		16V	C843	1-164-677-11		0.033uF	10%	16V
C762	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C844	1-164-677-11	CERAMIC CHIP	0.033uF	10%	16V
0700	4 404 000 44	OFD AND COUR	04.5		4017	0045	4.404.0==	0554440	100	ma.	
C763	1-164-360-11		0.1uF	000/	16V	C845	1-164-357-11		1000PF	5%	50V
		TANTALUM CHIP		20%	4V	C847	1-162-974-11		0.01uF	0001	50V
C771	1-162-974-11	CERAIVIIC CHIP	0.01uF		50V	C848	1-104-851-11	TANTALUM CHIP	iout	20%	10V

Ref. No	o. Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
C85	0 1-169-074-1	1 CERAMIC CHIP	0.01uF		EOV	0770	0.700 444 40		00000		
					50V	Q773		TRANSISTOR		4-T1BV4	
C85	3 1-162-9/4-1	1 CERAMIC CHIP	0.01uF		50V	Q774	8-729-141-48	TRANSISTOR	2SB62	4-T1BV4	
						Q775	8-729-037-72	TRANSISTOR	UN921	1J- (TX).	SO
C85	4 1-104-851-1	1 TANTALUM CHI	P 10uF	20%	10V	Q776		TRANSISTOR		16J-QR (
C85	5 1-104-851-1	1 TANTALUM CHI	P 10uF	20%	10V	Q777		TRANSISTOR	2SD22	16J-QR (TX) SO
C85		1 CERAMIC CHIP		5%	50V		0 . 20 00, 02		LODEL	100 011	177.00
C85		1 CERAMIC CHIP			50V	Q778	8-720-037-52	TRANSISTOR	26022	16J-QR (TV) CO
C86		1 CERAMIC CHIP		10%	25V	Q779		TRANSISTOR			
		, or white or m	0.0141	1070	201	Q784				16J-QR (
C86	2 1-169-070-11	CERAMIC CHIP	0.01uF	10%	051/	u/04	0-129-031-33	TRANSISTOR	25B14	62J-QR (1X).SU
C874		CERAMIC CHIP		10%	25V	İ					
C87				000/	50V			< RESISTOR >			
607	0 1-130-209-11	TANTALUM CHI	P 10ur	20%	6.3V						
		00111150505				R117		METAL CHIP	68	5%	1/16W
		< CONNECTOR	>			R118	1-216-833-91		10K	5%	1/16W
						R120	1-216-864-11	METAL CHIP	0	5%	1/16W
CN1		CONNECTOR, F	FC/EPC (ZIF) 30P		R121	1-216-825-11	METAL CHIP	2,2K	5%	1/16W
CN1	02 1-750-345-11	CONNECTOR, F	FC/EPC (ZIF) 30P		R122	1-216-825-11		2.2K		1/16W
CN1	03 1-750-345-11	CONNECTOR, F	FC/EPC (ZIF	30P							
CN7	71 1-770-305-11	CONNECTOR, F	FC/FPC 10P			R123	1-216-864-11	METAL CHIP	0	5%	1/16W
CN7	75 1-750-303-41	CONNECTOR, B	OARD TO B	OARD 201	Р			MEME OIM	Ū		(DSR-20MD)
						R124	1-216-864-11	METAL CHIP	0	5%	1/16W
		< DIODE >				1	1 210 004 11	WILLIAL OTHI	U		
						R137	1 016 007 11	METAL OUID	00		OSR-20MDP)
D77	8-710-072-01	DIODE MA111	TV				1-216-807-11		68	5%	1/16W
D772						R138	1-216-833-91		10K	5%	1/16W
D773		DIODE MA111				R143	1-216-833-91	RES, CHIP	10K	5%	1/16W
		DIODE KV1470									
D774		DIODE 188351				R144	1-216-831-11		6.8K	5%	1/16W
D775	8-719-052-27	DIODE 188351	-TB			R147	1-216-864-11		0	5%	1/16W
						R206	1-216-821-11	METAL CHIP	1K	5%	1/16W
D791	8-719-073-01	DIODE MA111	-TX			R308	1-216-821-11		1K	5%	1/16W
						R309	1-216-821-11		1K	5%	1/16W
		< FILTER >								0,0	,,,,,,,
						R310	1-216-821-11	METAL CHIP	1K	5%	1/16W
FL77	0 1-411-951-21	DELAY LINE, LC	(23NS)			R311	1-216-821-11		1K	5%	1/16W
FL77	1 1-233-734-21	FILTER, LOW PA	SS			R312	1-216-821-11		1K		
						R313	1-216-821-11			5%	1/16W
		< IC >				R314			1K	5%	1/16W
		(10)				NO 14	1-216-821-11	WETAL CHIP.	1K	5%	1/16W
IC770	8-759-445-03	IC AK6440AM-	Eo			D245	1 010 004 14	MATTAL OLUB			
IC77		IC MB88346LP		-D		R315	1-216-864-11		0	5%	1/16W
IC77				:R		R316	1-216-833-91	RES, CHIP	10K	5%	1/16W
		IC CXD2302Q-1				R701	1-216-825-11		2.2K	5%	1/16W
IC773		IC CXA1762Q-1				R702	1-216-829-11		4.7K	5%	1/16W
IC774	8-752-386-38	IC CXD3105R-1	16			R703	1-216-809-11	METAL CHIP	100	5%	1/16W
IC775		IC CXA2023R-T				R704	1-216-810-11	METAL CHIP	120	5%	1/16W
IC777		IC CXA2018Q-T				R705	1-216-825-11	METAL CHIP	2.2K	5%	1/16W
IC791	8-759-426-83	IC TK11228BM	CL			R706	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
						R707	1-216-809-11		100	5%	1/16W
		< COIL >				R708		METAL CHIP	120	5%	1/16W
•									120	0 70	1/1000
L105	1-414-398-11	INDUCTOR	10uH		1	R770	1-216-845-11	METAL CHIP	100K	5%	1/16W
L770	1-414-398-11		10uH			R772	1-216-296-91		0	J /0	1/1044
L773	1-414-398-11	INDUCTOR	10uH		.	R774	1-216-841-11			E0/	4 /4 014/
L774	1-414-398-11		10uH		1	R776			47K	5%	1/16W
L776	1-414-398-11		10uH		. 1		1-216-818-11		560	5%	1/16W
LIIO	1-414-030-11	INDUCTOR	TOUR			R779	1-216-847-11	METAL CHIP	150K	5%	1/16W
L779	1 410 727 21	INDUCTOR OUID	0.47.11		i					•	
		INDUCTOR CHIP				R780	1-216-837-11		22K	5%	1/16W
L781	1-412-963-11		100uH			R782	1-216-833-91		10K	5%	1/16W
L782	1-412-963-11		100uH		* *	R783	1-216-833-91	RES, CHIP	10K	5%	1/16W
L783	1-414-398-11	INDUCTOR	10uH			R786	1-216-817-11	METAL CHIP	470	5%	1/16W
L784	1-414-398-11	INDUCTOR	10uH			R787	1-202-924-11	RES, CHIP	240	5%	1/16W
					1						
L789	1-414-398-11	INDUCTOR	10uH			R788	1-202-924-11	RES. CHIP	240	5%	1/16W
						R789	1-216-824-11		1.8K	5%	1/16W
		< TRANSISTOR >				R790	1-216-841-11		47K	5%	1/16W
•						R791	1-216-815-11		330	5% 5%	
Q105	8-729-037-52	TRANSISTOR	2SD2216J	-OR (TX)	so	R792	1-216-814-11		270		1/16W
Q109	8-729-037-52		2SD2216J			11102	1 410-014-11	WEINE UNIP	Z/U	5%	1/16W
Q701	8-729-013-04		2SC4851-7		50	D700	1 010 007 11	RETAL OUR	0.017		
Q701	8-729-013-04					R793	1-216-827-11		3.3K	5%	1/16W
0772			2SC4851-1			R794	1-216-816-11		390	5%	1/16W
uliz	8-729-037-72	INANOISTUR	UN9211J-	(TX).SU		R796	1-216-809-11		100	5%	1/16W
						R797	1-216-827-11	METAL CHIP	3.3K	5%	1/16W

RP-228	RS-78	VA-102

										•	
Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
R798	1-216-815-11	METAL CHIP	330	5%	1/16W			< 10 >			
71,00	1 210 010 11	WEINE OTH	000	0 70	1/1044			(10)			
R799	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	10002	8_750_521_15	IC MAX232CW	E-TE-2		
R800	1-216-833-91		10K	5%	1/16W	10002	0-108-021-10	IO WANZOZOWE	-16-2		
R801	1-216-833-91	RES, CHIP	10K	5%	1/16W			4 COIL 5			
R802	1-216-841-11	METAL CHIP						< COIL >			
R804			47K	5%	1/16W	1.004	4 440 000 44	INDUCTOR OUR	40		
N0U4	1-216-839-11	METAL CHIP	33K	5%	1/16W	L001	1-412-029-11	INDUCTOR CHIP	TOUH		
Dooc	1 010 001 11	METAL CLUD	41/	E0/	4 (4 0)41						
R806	1-216-821-11	METAL CHIP	1K	5%	1/16W		4 7007 050 4	V4 400 DO 4 DD	001401 FTF	(DOD 0	0140)
R808	1-216-821-11		1K	5%	1/16W	*		VA-102 BOARD,			
R810	1-216-837-11		22K	5%	1/16W	*	A-7067-251-A	VA-102 BOARD,		(DSR-2	OMDP)
R812	1-216-837-11	METAL CHIP	22K	5%	1/16W			******			
R814	1-216-853-11	METAL CHIP	470K	5%	1/16W				(Re	ef.No. 1,	000 Series)
R815	1-216-853-11	METAL CHIP	470K	5%	1/16W			< CAPACITOR >			
R818	1-216-837-11		22K	5%	1/16W						
R819	1-216-839-11	METAL CHIP	33K	5%	1/16W	C051	1-113-619-11	CERAMIC CHIP	0.47uF		10V
R820	1-216-803-11	METAL CHIP	33	5%	1/16W	C052	1-164-230-11	CERAMIC CHIP	220PF	5%	50V
R822	1-216-834-11	METAL CHIP	12K	5%	1/16W	C053	1-113-619-11	CERAMIC CHIP	0.47uF		10V
						C055	1-164-360-11	CERAMIC CHIP	0.1uF		16V
R824	1-216-821-11	METAL CHIP	1K	5%	1/16W	C056	1-164-360-11	CERAMIC CHIP	0.1uF		16V
R825	1-216-841-11	METAL CHIP	47K	5%	1/16W						
R826	1-216-839-11	METAL CHIP	33K	5%	1/16W	C057	1-162-927-11	CERAMIC CHIP	100PF	5%	50V
R827	1-216-821-11	METAL CHIP	1K	5%	1/16W	C058		CERAMIC CHIP	100PF	5%	50V
R830	1-216-831-11		6.8K	5%	1/16W	C059	1-124-778-00		22uF	20%	6.3V
11000	1 210 001 11	WEI/IE OIIII	0.010		171000	C101	1-128-004-11		10uF	20%	16V
R832	1-216-807-11	METAL CHIP	68	5%	1/16W	C101		CERAMIC CHIP	0.1uF	20%	
R843	1-216-822-11		1.2K			6102	1-104-300-11	GENAIVIIG GRIP	U. IUF		16V
				5%	1/16W	0400	4 404 000 44	OFDAMAO OLUD	0.4		4014
R844	1-216-837-11	METAL CHIP	22K	5%	1/16W	C103		CERAMIC CHIP	0.1uF		16V
R849	1-218-837-11	METAL CHIP	390	0.5%	1/16W	C104		CERAMIC CHIP	0.1uF		16V
R850	1-218-835-11	METAL CHIP	330	0.5%	1/16W	C105		CERAMIC CHIP	0.1uF		16V
						C106	1-164-360-11	CERAMIC CHIP	0.1uF		16V
R851	1-218-835-11	METAL CHIP	330	0.5%	1/16W	C107	1-164-360-11	CERAMIC CHIP	0.1uF		16V
R852	1-218-837-11	METAL CHIP	390	0.5%	1/16W						
R858	1-216-816-11	METAL CHIP	390	5%	1/16W	C108	1-164-360-11	CERAMIC CHIP	0.1uF		16V
						C109	1-164-360-11	CERAMIC CHIP	0.1uF		16V
						C110		CERAMIC CHIP	0.1uF		16V
*	A-7073-472-A	RS-78 BOARD, C	OMPLETE			C111		CERAMIC CHIP	0.1uF		16V
		******				C201	1-124-778-00		22uF	20%	6.3V
			/R	ef No. 6.0	00 Series)	0201	1 124 770 00	LLLO1 OIIII	LLUI	2070	0.04
			() .	91.110. 0,0	00 001100)	C202	1-124-778-00	ELECT CHIP	22uF	20%	6.3V
		< CAPACITOR >			·	C203		CERAMIC CHIP	0.1uF	2070	16V
		V OAL AUTTOR >				C204	1-126-206-11			200/	
0001	4 464 046 44	CEDAMIC CUID	4		101				100uF	20%	6.3V
C001		CERAMIC CHIP	1uF		16V	C205		CERAMIC CHIP	0.01uF	10%	25V
C002		CERAMIC CHIP	1uF		16V	C206	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
C003		CERAMIC CHIP	1uF		16V						
C004		CERAMIC CHIP	1uF		1.6V	C207		CERAMIC CHIP	270PF	5%	50V
C005	1-164-346-11	CERAMIC CHIP	1uF		16V	C208		CERAMIC CHIP	0.1uF		16V
						C209	1-164-360-11	CERAMIC CHIP	0.1uF		16V
C006		CERAMIC CHIP	0.1uF		16V	C210		CERAMIC CHIP	47PF	5%	50V
C007	1-126-206-11	ELECT CHIP	100uF	20%	6.3V	C211	1-124-778-00	ELECT CHIP	22uF	20%	6.3V
		< CONNECTOR >				C212	1-124-778-00	ELECT CHIP	22uF	20%	6.3V
						C213	1-128-003-11	ELECT CHIP	22uF	20%	4V
CN001	1-770-693-11	CONNECTOR, FF	C/FPC 10P			C214	1-128-003-11	ELECT CHIP	22uF	20%	4V
CN002	1-565-388-21	CONNECTOR, D-	SUB 9P (R	EMOTE RS	S-232C)	C215	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
						C216	1-128-003-11	ELECT CHIP	22uF	20%	4V
		< DIODE >									
						C217	1-164-360-11	CERAMIC CHIP	0.1uF		16V
D001	8-719-062-19	DIODE MA3200	WA-TX			C218		CERAMIC CHIP	0.01uF	10%	25V
D002		DIODE MA3200			i	C219	1-162-970-11		0.01uF	10%	25V
D003		DIODE MA3200				C220		CERAMIC CHIP	0.01uF	10%	25V
D003		DIODE MA3200				C221		CERAMIC CHIP	10PF	0.5PF	50V
D004		DIODE MA3200				0221	1-102-310-11	CENAINIC CITIF	1011	0.51	30 V
נטטט	0-119-002-19	DIODE WAS200	AAW-IV			Conn	1,160,067,44	CEDAMIC CLUB	0.0000	100/	EOV.
		CEDDITE SEAS				C223	1-162-967-11			10%	50V
		< FERRITE BEAD	> .			C224	1-126-206-11		100uF	20%	6.3V
	4 800 041 =:					C225		CERAMIC CHIP	0.01uF	10%	25V
FB001	1-500-241-22		0uH			C226	1-124-778-00		22uF	20%	6.3V
FB002	1-500-241-22		0uH			C227	1-165-176-11	CERAMIC CHIP	0.047uF	10%	16V
FB003	1-500-241-22		0uH							:	
FB004	1-500-241-22		0uH			C228		CERAMIC CHIP	0.1uF		16V
FB005	1-500-241-22	FERRITE	0uH		1	C229	1-128-007-11	ELECT CHIP	2.2uF	20%	35V
1 5000											

E	Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	<u>Description</u>		. "	Remark
	C230		CERAMIC CHIP	0.01uF	000/	50V	C291	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V
	C231	1-128-006-11		4.7uF	20%	25V	2000	1 100 005 11	FUEOT OUUD	47		R-20MDP)
	C232	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C292	1-126-205-11		47uF	20%	6.3V
	0000	1-126-205-11	ELECT CHID	47.1E	20%	6.3V	C293 C294	1-126-205-11	CERAMIC CHIP	47uF 0.01uF	20% 10%	6.3V 25V
	C233 C234		CERAMIC CHIP	47uF 0.1uF	20%	16V	C294		CERAMIC CHIP	0.01uF	10%	25V 25V
	C235	1-124-778-00		22uF	20%	6.3V	0293	1-102-970-11	CENAIVIIC CHIP	U.UTUF	1076	237
	C236	1-124-776-00		4.7uF	20%	25V	C297	1-128-003-11	ELECT CHIP	22uF	20%	4V
	C237		CERAMIC CHIP	0.01uF	10%	25V	C300		CERAMIC CHIP	0.1uF	2070	16V
	0201	1 102 570 11	OLIMINO OIII	0.0741	1070	201	0000	1 104 000 11	OLITAWIO OTIII	0.141	(DS	R-20MDP)
	C238	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C301	1-162-927-11	CERAMIC CHIP	100PF	5%	50V
	C239		CERAMIC CHIP	0.1uF		16V					(DS	R-20MDP)
	C240	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C304	1-126-205-11	ELECT CHIP	47uF	20%	6.3V
	C241	1-128-004-11	ELECT CHIP	10uF	20%	16V	C305	1-126-205-11	ELECT CHIP	47uF	20%	6.3V
	C242	1-164-360-11	CERAMIC CHIP	0.1uF		16V						
							C307		CERAMIC CHIP	0.01uF	10%	25V
	C243		CERAMIC CHIP	0.1uF		16V	C308		CERAMIC CHIP	0.01uF	10%	25V
	C244		CERAMIC CHIP	0.1uF		16V	C309	1-164-360-11	CERAMIC CHIP	0.1uF		16V
	C245		CERAMIC CHIP	0.1uF	4.00/	16V	2014	4 400 004 44	EL FOT OLUD	40.5		R-20MDP)
	C246		CERAMIC CHIP	0.001uF	10%	50V	C311	1-128-004-11	ELECT CHIP	10uF	20%	16V
	C247	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	0010	4 404 000 44	OEDAMIO OUID	0.1	(บร	R-20MDP)
	C248	1 169-070-11	CERAMIC CHIP	0.01uF	10%	25V	C312	1-104-300-11	CERAMIC CHIP	0.1uF	/DC	16V R-20MDP)
	C249	1-102-970-11		4.7uF	20%	25V 25V					(טט	N-ZUIVIDE)
	C250		CERAMIC CHIP	0.1uF	20 /0	16V	C313	1-126-206-11	ELECT CHIP	100uF	20%	6.3V
	C251		CERAMIC CHIP	6PF	0.5PF		C314	1-126-206-11		100uF	20%	6.3V
	0201	1 102 011 71	oeronono omi	OI (SR-20MD)	C401	1-126-205-11		47uF	20%	6.3V
	C251	1-162-910-11	CERAMIC CHIP	5PF	0.25PF		C402		CERAMIC CHIP	0.01uF	10%	25V
						R-20MDP)	C403	1-128-003-11	ELECT CHIP	22uF	20%	4V
	C252	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C404	1-128-003-11		22uF	20%	4V
	C253		CERAMIC CHIP	100PF	5%	50V	C405	1-128-003-11		22uF	20%	4V
	C254		CERAMIC CHIP	0.01uF	10%	25V	C406		CERAMIC CHIP	0.01uF	10%	25V
	C255 C256	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V 16V	C407 C408		CERAMIC CHIP	0.01uF 0.01uF	10% 10%	25V
	0256	1-164-360-11	CERAMIC CHIP	0.1uF		100	0400	1-102-970-11	CERAMIC CHIP	U.UTUF	10%	25V
	C257	1-128-003-11	ELECT CHIP	22uF	20%	4V	C409	1-128-003-11	ELECT CHIP	22uF	20%	4V
	C258	1-128-003-11		22uF	20%	4V	C410	1-128-007-11		2.2uF	20%	35V
	C260	1-162-927-11	CERAMIC CHIP	100PF	5%	50V	C411	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	C261	1-162-919-11	CERAMIC CHIP	22PF	5%	50V	C412	1-126-205-11		47uF	20%	6.3V
	C262	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	C413	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	0004	4 404 200 44	CEDAMIC CUID	0.1		101/	0444	1.104.770.00	EL FOT OLUD	00"	0.00/	0.014
	C264		CERAMIC CHIP	0.1uF	100/	16V	C414	1-124-778-00		22uF	20%	6.3V
	C265 C266		CERAMIC CHIP	0.001uF 1uF	10%	50V 10V	C415 C416	1-128-003-11 1-128-003-11		22uF 22uF	20% 20%	4V 4V
	C268		CERAMIC CHIP	0.47uF		10V	C417		CERAMIC CHIP	330PF	5%	50V
	C269		CERAMIC CHIP	0.1uF		16V	C418		CERAMIC CHIP	22PF	5%	50V
	0200	. 101 000 11	02.17.11.11.0			101	0110	1 102 010 11	OLI II IIII OIIII	L	070	
	C270	1-115-156-11	CERAMIC CHIP	1uF		10V	C419	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	C271	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C420		CERAMIC CHIP	1uF		10V
	C272	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C421	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
	C273		CERAMIC CHIP	0.1uF		16V	C422		CERAMIC CHIP	0.056uF	10%	25V
	C274	1-165-176-11	CERAMIC CHIP	0.047uF	10%	16V	C423	1-163-139-00	CERAMIC CHIP	820PF	5%	50V
	0075	1 100 006 11	ELECT CHID	4 7F	000/	051/	0405	1 100 070 11	OFDAMIO OLUD	0.04	100/	05)/
	C275	1-128-006-11		4.7uF 100uF	20% 20%	25V 6.3V	C425		CERAMIC CHIP	0.01uF 0.022uF	10%	25V
	C276 C277		CERAMIC CHIP	0.047uF	10%	16V	C426 C427		CERAMIC CHIP	0.022ur 0.01uF	10% 10%	25V 25V
	C278		CERAMIC CHIP	0.01uF	10%	25V	C429		CERAMIC CHIP	1uF	10 /0	10V
	C279	1-128-006-11		4.7uF	20%	25V	C430		CERAMIC CHIP	1uF		10V
							3.00					
	C280	1-128-003-11	ELECT CHIP	22uF	20%	4V	C433	1-115-156-11	CERAMIC CHIP	1uF		10V
	C281		CERAMIC CHIP	6PF	0.5PF	50V	C434		CERAMIC CHIP	1uF		10V
	C282		CERAMIC CHIP	0.001uF	10%	50V	C435		CERAMIC CHIP	1uF		10V
	C283		CERAMIC CHIP	10PF	0.5PF	50V	C436		CERAMIC CHIP	0.01uF	10%	25V
	C284	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C437	1-128-003-11	ELECT CHIP	22uF	20%	4V
	C285	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	C438	1-169-070-11	CERAMIC CHIP	0.01uF	10%	25V
	C286		CERAMIC CHIP	0.00TuF	10.70	16V	C436		CERAMIC CHIP	0.01uF	10%	25V 25V
	C287	1-126-206-11		100uF	20%	6.3V	C443		CERAMIC CHIP	0.01uF	10%	25V
	C288		CERAMIC CHIP	0.001uF	10%	50V	C444	1-162-970-11		0.01uF	10%	25V
	C290	1-128-003-11		22uF	20%	4V	C445		CERAMIC CHIP	1uF		10V
				1.1								

Ref. No.	Part No.	Description			<u>Remark</u>	Ref. No.	Part No.	Description			Remark
C446	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V						
C447	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C612	1-162-920-11	CERAMIC CHIP	27PF	5%	50V
C448	1-124-778-00	ELECT CHIP	22uF	20%	6.3V						SR-20MD)
C449	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C612	1-162-921-11	CERAMIC CHIP	33PF	5%	50V
C450	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V					•	R-20MDP)
						C613	1-162-915-11	CERAMIC CHIP	10PF	0.5PF	50V
C451	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C614	1-162-915-11	CERAMIC CHIP	10PF	0.5PF	50V
C453	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C615	1-162-917-11	CERAMIC CHIP	15PF	5%	50V
C454	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V					(D	SR-20MD)
C455	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V						
C456		CERAMIC CHIP	0.01uF	10%	25V	C616	1-124-778-00	ELECT CHIP	22uF	20%	6.3V
						C618	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
C457	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C619	1-124-778-00	ELECT CHIP	22uF	20%	6.3V
C459	1-115-156-11	CERAMIC CHIP	1uF		10V	C621	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
C460	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C624	1-164-360-11	CERAMIC CHIP	0.1uF		16V
C461	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V						
C462	1-124-778-00		22uF	20%	6.3V	C628	1-164-360-11	CERAMIC CHIP	0.1uF		16V
0.00						C629	1-164-230-11	CERAMIC CHIP	220PF	5%	50V
C463	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C630	1-128-004-11	ELECT CHIP	10uF	20%	16V
C464	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C631	1-164-360-11	CERAMIC CHIP	0.1uF		16V
C465	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C651		CERAMIC CHIP	0.01uF	10%	25V
C466	1-102-970-11		22uF	20%	6.3V	0001	1 102 570 11	OLI MINIO OTTI	O.O.L	1070	201
		ELECT CHIP	47uF	20%	6.3V	C652	1-128-004-11	ELECT CHIP	10uF	20%	16V
C467	1-126-205-11	ELECT CHIP	4/UF	20%	0.37	C653	1-162-909-11	CERAMIC CHIP	4PF	0.25PF	
0400	4 404 000 44	OFDARMO OUID	0.4		101	0000	1-102-909-11	CENAIVIIC CHIP	461		R-20MDP)
C468	1-164-360-11	CERAMIC CHIP	0.1uF	000/	16V	0054	1-128-003-11	ELECT CUID	22uF	20%	4V
C469	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C654			22ur 10uF	20%	16V
C470	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C655	1-128-004-11				
C471	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C656	1-162-927-11	CERAMIC CHIP	100PF	5%	50V
C472	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	0057	4 400 000 44	05041410 01110	40000	F0/	E014
						C657	1-162-928-11		120PF	5%	50V
C473	1-126-206-11	ELECT CHIP	100uF	20%	6.3V	C658		CERAMIC CHIP	100PF	5%	50V
C474	1-126-205-11		47uF	20%	6.3V	C659	1-126-205-11		47uF	20%	6.3V
C475	1-124-778-00		22uF	20%	6.3V	C660		CERAMIC CHIP	0.01uF	10%	25V
C476	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C662	1-128-007-11	ELECT CHIP	2.2uF	20%	35V
C477	1-162-909-11	CERAMIC CHIP	4PF	0.25PF							
				(DSI	R-20MDP)	C663	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
						C664	1-126-205-11	ELECT CHIP	47uF	20%	6.3V
C478	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	C665	1-164-237-11	CERAMIC CHIP	16PF	5%	50V
C479	1-164-237-11	CERAMIC CHIP	16PF	5%	50V	C666	1-162-919-11	CERAMIC CHIP	22PF	5%	50V
C480	1-126-206-11	ELECT CHIP	100uF	20%	6.3V	C667	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
C481	1-162-919-11	CERAMIC CHIP	22PF	5%	50V						
C482	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C670	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V
						C671	1-115-156-11	CERAMIC CHIP	1uF		10V
C483	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C672	1-162-921-11	CERAMIC CHIP	33PF	5%	50V
C484		CERAMIC CHIP	33PF	5%	50V	C673	1-164-343-11	CERAMIC CHIP	0.056uF	10%	25V
C485	1-162-919-11		22PF	5%	50V	C674		CERAMIC CHIP	22PF	5%	50V
C486	1-164-360-11	CERAMIC CHIP	0.1uF		16V						
C487	1-128-004-11		10uF	20%	16V	C675	1-163-139-00	CERAMIC CHIP	820PF	5%	50V
0401	1 120 004 11	LLLOT OTT	Tour	2070	101	C676		CERAMIC CHIP	0.022uF	10%	25V
C489	1-162-920-11	CERAMIC CHIP	27PF	5%	50V	C678	1-128-003-11		22uF	20%	4V
		CERAMIC CHIP	1PF	0.25PF		C680		CERAMIC CHIP	0.01uF	10%	25V
C490						C681	1-126-206-11		100uF	20%	6.3V
C491		CERAMIC CHIP	0.01uF	10%	25V	0001	1-120-200-11	ELECT OTH	Toour	20 /0	0.54
C492		CERAMIC CHIP	10PF	0.5PF	50V	0000	1 164 960 11	CERAMIC CHIP	0.1uF		16V
C493	1-128-004-11	ELECT CHIP	10uF	20%	16V	C682				000/	
					mai /	C683	1-126-927-11		1000uF	20%	6.3V
C494		CERAMIC CHIP	0.001uF	10%	50V	C701	1-124-778-00		22uF	20%	6.3V
C496	1-128-007-11		2.2uF	20%	35V	C702	1-124-778-00		22uF	20%	6.3V
C497		CERAMIC CHIP	0.1uF		16V	C703	1-162-923-11	CERAMIC CHIP	47PF	5%	50V
C498	1-128-003-11		22uF	20%	4V	1		•			
C500	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C704		CERAMIC CHIP	0.1uF		16V
						C705		CERAMIC CHIP	47PF	5%	50V
C501	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C706		CERAMIC CHIP	0.1uF		16V
	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C707		CERAMIC CHIP	0.1uF		16V
C503	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	C708	1-128-004-11	ELECT CHIP	10uF	20%	16V
C505			10uF	20%	16V						
	1-128-004-11			20%	6.3V	C709	1-128-004-11	ELECT CHIP	10uF	20%	16V
C505 C506			22uF				4 400 040 44	OFFI LLUIS CLUIS			
C505	1-128-004-11 1-124-778-00		22 u F			C710	1-162-919-11	CERAMIC CHIP	22PF	5%	50V
C505 C506 C601	1-124-778-00	ELECT CHIP	0.01uF		25V	C710 C711		CERAMIC CHIP	22PF 22PF	5% 5%	50V 50V
C505 C506 C601	1-124-778-00 1-162-970-11	CERAMIC CHIP	0.01uF	10%			1-162-919-11				
C505 C506 C601 C602 C603	1-124-778-00 1-162-970-11 1-128-004-11	CERAMIC CHIP ELECT CHIP	0.01uF 10uF	10%	16V	C711 C713	1-162-919-11 1-164-360-11	CERAMIC CHIP	22PF		50V
C505 C506 C601	1-124-778-00 1-162-970-11 1-128-004-11	CERAMIC CHIP ELECT CHIP CERAMIC CHIP	0.01uF	10%		C711	1-162-919-11 1-164-360-11	CERAMIC CHIP CERAMIC CHIP	22PF 0.1uF		50V 16V

													_
F	Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remar	k
				0.4						0.004	100/	50V	-
	C716		CERAMIC CHIP	0.1uF	-0 /	16V	C855		CERAMIC CHIP	0.001uF	10% 10%	25V	
	C717		CERAMIC CHIP	300PF	5%	50V	C856		CERAMIC CHIP	0.01uF	10%		
	C719		CERAMIC CHIP	0.1uF		16V	C858		CERAMIC CHIP	0.1uF		16V	
	C720	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C859	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C721		CERAMIC CHIP	47PF	5%	50V	C860		CERAMIC CHIP	0.1uF		16V	
	C722		CERAMIC CHIP	0.1uF		16V	C861		CERAMIC CHIP	0.1uF		16V	
	C723	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C862	1-126-204-11		47uF	20%	16V	
	C724	1-164-389-11	CERAMIC CHIP	300PF	5%	50V	C863	1-126-400-11		22uF	20%	35V	
	C725	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C864	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C726	1-164-389-11	CERAMIC CHIP	300PF	5%	50V	C865		CERAMIC CHIP	0.1uF		16V	
	C727	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C866		CERAMIC CHIP	0.01uF	10%	25V	
	C728	1-162-923-11	CERAMIC CHIP	47PF	5%	50V	C867		CERAMIC CHIP	0.1uF		16V	
	C729	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C869		CERAMIC CHIP	0.1uF		16V	
	C730	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C870	1-128-013-11	ELECT CHIP	1uF	20%	50V	
	C731	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C871	1-162-970-11	CERAMIC CHIP	0.01uF	10%	25V	
	C732	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C872	1-162-975-11	CERAMIC CHIP	24PF	5%	50V	
	C733		CERAMIC CHIP	0.1uF		16V	C874	1-126-204-11	ELECT CHIP	47uF	20%	16V	
	C734	1-128-007-11	ELECT CHIP	2.2uF	20%	35V	C875	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C735	1-128-007-11		2.2uF	20%	35V	C876	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
											1001	0.514	
	C736		CERAMIC CHIP	0.1uF	F 0/	16V	C877		CERAMIC CHIP	0.022uF	10%	25V	
	C737		CERAMIC CHIP	300PF	5%	50V	C881	1-126-927-11		1000uF	20%	6.3V	
	C738		CERAMIC CHIP	0.1uF	mor	16V	C883		CERAMIC CHIP	0.1uF	F0/	16V	
	C739		CERAMIC CHIP	300PF	5%	50V	C886		CERAMIC CHIP	33PF	5%	50V	
	C740	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C887	1-162-921-11	CERAMIC CHIP	33PF	5%	50V	
	C741	1-124-778-00	ELECT CHIP	22uF	20%	6.3V	C888	1-128-004-11	ELECT CHIP	10uF	20%	16V	
	C742		CERAMIC CHIP	300PF	5%	50V	C889	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C743		CERAMIC CHIP	0.1uF		16V	C890	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C744	1-164-389-11	CERAMIC CHIP	300PF	5%	50V	C891	1-165-176-11	CERAMIC CHIP	0.047uF	10%	16V	
	C745	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C892	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C746	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C893	1-126-392-11	ELECT CHIP	100uF	20%	6.3V	
	C747		CERAMIC CHIP	10PF	0.5PF	50V	C894	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C748		CERAMIC CHIP	0.1uF		16V	C895	1-126-392-11		100uF	20%	6.3V	
	C749		CERAMIC CHIP	0.1uF		16V	C896	1-126-206-11		100uF	20%	6.3V	
	C750	1-164-360-11		0.1uF		16V	C897	1-164-360-11	CERAMIC CHIP	0.1uF	•	16V	
	C751	1-162-923-11	CERAMIC CHIP	47PF	5%	50V	C898	1-126-396-11	ELECT CHIP	47uF	20%	16V	
	C752		CERAMIC CHIP	0.1uF		16V	C899	1-126-396-11		47uF	20%	16V	
	C753		CERAMIC CHIP	0.1uF		16V	C900	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C754		CERAMIC CHIP	10PF	0.5PF	50V	C901		CERAMIC CHIP	0.1uF		16V	
	C755		CERAMIC CHIP	0.033uF	10%	16V	C906	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C756	1-164-360-11	CERAMIC CHIP	0.1uF		16V	C907	1-164-360-11	CERAMIC CHIP	0.1uF		16V	
	C757	1-164-677-11		0.033uF	10%	16V	C908	1-164-360-11		0.1uF		16V	
	C758		CERAMIC CHIP	0.033uF	10%	16V	C909		CERAMIC CHIP	0.1uF		16V	
	C759		CERAMIC CHIP	0.033uF	10%	16V	C910		CERAMIC CHIP	2.2uF	10%	10V	
	C760	1-162-923-11		47PF	5%	50V	C912	1-104-905-11		0.22F	1070	5.5V	
	C761	1 164 677 11	CERAMIC CHIP	0.033uF	10%	16V	C913	1-126-206-11	ELECT CHIP	100uF	20%	6.3V	
			CERAMIC CHIP	0.033uF	10%	16V	C914		CERAMIC CHIP	100PF	5%	50V	
	C762		CERAMIC CHIP	0.033aF	1070	16V	C915		CERAMIC CHIP	0.01uF	10%	25V	
	C763	1-128-004-11		10uF	20%	16V	C916	1-102-970-11		1uF	20%	50V	
	C764 C765		CERAMIC CHIP	0.1uF	20 /6	16V	C917	1-126-197-11		10uF	20%	50V	
					0001					47. 5	0001		
	C766	1-128-004-11		10uF	20%	16V	C918	1-126-204-11		47uF	20%	16V	
	C767		CERAMIC CHIP	0.033uF 0.033uF	10%	16V	C919		CERAMIC CHIP	0.1uF 0.1uF		16V 16V	
	C768		CERAMIC CHIP	0.033uF 0.1uF	10%	16V 16V	C920 C921	1-126-393-11		33uF	20%	10V	
	C771 C772		CERAMIC CHIP	47PF	5%	50V	C921		CERAMIC CHIP	0.022uF	10%	25V	
									COMMECTOR				
	C773		CERAMIC CHIP	0.1uF 47PF	5%	16V 50V			< CONNECTOR >			3 - 5	
	C774 C851		CERAMIC CHIP	47PF 0.1uF	J 70	16V	CN051	1-770-205-11	CONNECTOR, FFO	C/EDC 10D			
	C852		CERAMIC CHIP	0.1uF		16V	CN101		CONNECTOR, FFO				
	C853		CERAMIC CHIP	0.1uF		16V	CN102		CONNECTOR, FFO				
	0000	1-104-000-11	OLIMINIO OTIF	o. rui		100	CN401		CONNECTOR, FFO				
	C854	1-162-964-11	CERAMIC CHIP	0.001uF	10%	50V	* CN601		PIN, CONNECTOR				

						•	
Ref. No.	Part No.	Description	Remark	Ref. No.	Part No.	Description	Remark
CNCOO	1 770 200 11	CONNECTOR COURSE TYPE (INDI)	4D	IC104 IC105		IC TC74VHCT08FS (EL) IC TC74VHCT08FS (EL)	
CNOUZ	1-779-309-11	CONNECTOR, SQUARE TYPE (INDI)	OV IN/OUT)	IC105		IC TC74VHCT08FS (EL)	
		PIN, CONNECTOR (1.5mm) (SMD) 8	Р	IC201	8-759-337-26	IC MM1115XFBE	
		CONNECTOR, FFC/FPC 12P CONNECTOR, FFC/FPC 15P		IC202	8-759-433-44	IC MM1031XML	* *
CN851		CONNECTOR, FFC/FPC 30P		IC203	8-759-337-26	IC MM1115XFBE	
CNOSO	1 506 479 11	PIN, CONNECTOR 8P		IC204 IC205	8-759-432-78 8-759-420-62	IC MM1111XFBE	
CN852	1-300-4/3-11	FIN, CUNNECTOR OF		IC205		IC NJM2240M (TE2)	
		< TRIMMER >		10007	0.750.050.00	IO OVERGOOD (DOD COME)	
CT201	1-141-424-11	CAP, ADJ (DECODER FREERUN)		IC207 IC207		IC CXD2023Q (DSR-20MD) IC CXD2024AQ-TL (DSR-20MDP)	
CT401		CAP, ADJ (ENCODER FREERUN)		IC208	8-759-603-54	IC M51271FP-70AD	
		< D10DE >		IC210 IC213		IC TC74HC221AF (EL) (DSR-20MDF IC TC7S00FU (TE85R) (DSR-20MDF	
				102.10			′
D051		DIODE MA132WA-TX DIODE MA132WA-TX		IC401 IC402		IC MM1115XFBE IC MM1111XFBE	
D052 D201		DIODE MA111-TX		10402		IC LA7218M-TE-R	
D202	8-719-073-01	DIODE MA111-TX		IC404		IC CXA1592R-T4	
D401	8-713-101-85	DIODE 1T363-01-T8A		IC405	8-759-483-56	IC MB90089PF-G-196-BND-ER	
D403	8-719-073-01	DIODE MA111-TX		IC406	8-759-182-16	IC MM1196XFBE	
D404		DIODE MA111-TX		1C407		IC CXD2192Q-T4	
D851 D852		DIODE MA3075WA- (TX) DIODE MA3075WA- (TX)		IC602 IC651		IC MC68HC68VBIFB IC MB90089PF-G-155-BND-ER	
D853		DIODE MA3082-TX		IC652		IC LA7218M-TE-R	
Dorr	0 740 404 54	DIODE MA738-TX		10653	0 750 997.96	IC MM1115XFBE	
D855 D856		DIODE MA729- (K8).S0		IC701		IC NJM2115V (TE2)	
D858		DIODE MA111-TX		IC702	8-759-523-02	IC TC74HC4053AFT (EL)	
D859		DIODE MA111-TX		IC703		IC TC74HC4053AFT (EL)	
D861	8-719-400-56	DIODE MA3062H-TX		IC704	8 - 759-358-47	IC NJM2115V (TE2)	
D862		DIODE MA3075WA- (TX)		IC705		IC NJM2115V (TE2)	
D863		DIODE RD6.2FM-T1		IC706 IC707		IC NJM2115V (TE2) IC NJM2115V (TE2)	
D864 D866	8-719-421-27	DIODE MA728- (K8).S0 DIODE MA111-TX		IC707		IC DS1801E-014TE2	
D867		DIODE MA738-TX		10709		IC NJM2115V (TE2)	
		< DELAY LINE >		IC710	8-759-358-47	IC NJM2115V (TE2)	
				IC711	8-759-523-02	IC TC74HC4053AFT (EL)	
DL201	1-411-661-11	LINE, LC DELAY		IC712		IC NJM2115V (TE2)	
		< FERRITE BEAD >		IC713		IC NJM2115V (TE2) IC NJM2115V (TE2)	•
FB851	1-543-948-22			10715		IC NJM2115V (TE2)	
FB852 FB853	1-543-948-22 1-543-948-22			IC716		IC NJM2115V (TE2) IC NJM2115V (TE2)	
1 5000	1 040 040 22	TENNITE SUIT		IC718	8-759-358-47	IC NJM2115V (TE2)	
		< FILTER >		IC851	8-759-356-27	IC NJM2129M-TE2	
FL201	1-236-925-11	FILTER, LOW PASS		IC852	8-759-545-30	IC S579178PJ	
FL202	1-236-926-11	FILTER, BAND PASS (DSR-20MD)		10853		IC TC74HC4053AFT (EL)	
FL202		FILTER, BAND PASS (DSR-20MDP)		IC854 IC856		IC AK6440AM-E2 CONVERTER UNIT, DC/DC	
FL203 FL204		FILTER, LOW PASS FILTER, LOW PASS		1C857		IC S-3513BEFS-TB	
				10050	0 750 040 07	IC MANALOGGVE-DE	
FL205 FL401		FILTER, LOW PASS FILTER, BAND PASS (DSR-20MD)		1C858 1C861		IC MM1256XF-BE IC TL1596CPW-ELM2000	
FL401		FILTER, BAND PASS (DSR-20MDP)		IC862	8-759-822-95	IC L79M05T-FA-TL	
		10		10863		IC PQ05TZ1U	
		< IC >		IC864	o-759-929-26	IC TL431CPSR	
IC051		IC TC74HC74AF (EL)		IC865		IC TC74VHC00FT (EL)	
10052		IC HD6433837SC05H	• .	IC866	8-759-363-18	IC TC7ST04FU (TE85R)	
IC053 IC101		IC TC7WU04FU (TE12R) IC TC74VHC08FT (EL)				< JACK >	
IC102		IC TC7W125FU-TE12R			4 220 202 1:		101 C 181
IC103	8-759-594-04	IC TC74VHC125FT (EL)		J851	1-573-798-11	JACK, MINIATURE (DIA. 3.5) (CONTR	IUL S IN)
10100	0 100 024-04			•			

B	ef. No.	Part No.	Description	Remark	Ref. No.	Part No.	Description		Remark
	J852	1-573-798-11	JACK, MINIATUR		Q202		TRANSISTOR	2SC4081T106R	
	ioso	1 001 050 11	IAOK (LANO)	(CONTROL S OUT)	Q203		TRANSISTOR	2SC4081T106R	
	J853	1-091-256-11	JACK (LANC)		Q204		TRANSISTOR	2SC4081T106R	
			< COIL >		Q205	8-729-026-52	TRANSISTOR	2SA1576A-T106-R	
	1.054	4 440 000 44	. INDUSTOR OUT		Q206		TRANSISTOR	2SC4081T106R	
	L051		INDUCTOR CHIP		Q207	8-729-026-52		2SA1576A-T106-R	
	L201		INDUCTOR CHIP		Q208	8-729-427-83		XP6501-TXE	
	L202		INDUCTOR CHIP		Q209	8-729-905-35		2SC4081T106R	
	L203 L204		INDUCTOR CHIP INDUCTOR CHIP		Q210	8-729-026-52	TRANSISTOR	2SA1576A-T106-R	
	L207	1 412 020 11	INDOOTOR OTH	Tour	Q211	8-729-026-52	TRANSISTOR	2SA1576A-T106-R	
	L205	1-412-029-11	INDUCTOR CHIP	10uH	Q212	8-729-026-52		2SA1576A-T106-R	
	L206		INDUCTOR CHIP		Q213	8-729-905-35		2SC4081T106R	
	L207		INDUCTOR CHIP		Q214	8-729-905-35		2SC4081T106R	
	L208		INDUCTOR CHIP		Q215	8-729-905-35		2SC4081T106R	
	L211	1-412-808-21		470uH					
					Q216	8-729-905-35	TRANSISTOR	2SC4081T106R	
	L213		INDUCTOR CHIP		Q217	8-729-026-52	TRANSISTOR	2SA1576A-T106-R	
	L216		INDUCTOR CHIP		Q218	8-729-026-52		2SA1576A-T106-R	
	L217		INDUCTOR CHIP		Q221	8-729-905-35		2SC4081T106R	
	L218		INDUCTOR CHIP		Q222	8-729-905-35	TRANSISTOR	2SC4081T106R	
	L219	1-412-029-11	INDUCTOR CHIP	10uH					
					Q223	8-729-905-35		2SC4081T106R	
	L220			10uH (DSR-20MDP)	Q224	8-729-427-83		XP6501-TXE	
	L401		INDUCTOR CHIP		Q225	8-729-905-35		2SC4081T106R	
	L402		INDUCTOR CHIP		Q226	8-729-905-35		2SC4081T106R	
	L403 L404		INDUCTOR CHIP		0227	8-729-905-35	TRANSISTOR	2SC4081T106R	
	L404	1-412-029-11	INDUCTOR CHIP	Tour	Q228	0 700 407 00	TRANSICTOR	VDCEO1 TVE	
	L406	1-412-030-11	INDUCTOR CHIP	22uH	Q229	8-729-427-83 8-729-026-52		XP6501-TXE 2SA1576A-T106-R	
	L407		INDUCTOR CHIP		Q223	8-729-905-35		2SC4081T106R	
	L408		INDUCTOR CHIP		Q232	8-729-427-83		XP6501-TXE	
	L409		INDUCTOR CHIP		Q233	8-729-905-35		2SC4081T106R	
	L410		INDUCTOR CHIP		4250	0 720 000 00	THAIRDIOTOTT	2004001110011	
					Q403	8-729-905-35	TRANSISTOR	2SC4081T106R	
	L411	1-412-029-11	INDUCTOR CHIP	10uH	Q4Q4	8-729-905-35		2SC4081T106R	
	L412	1-412-029-11	INDUCTOR CHIP	10uH	Q405	8-729-905-35		2SC4081T106R	
	L601		INDUCTOR CHIP		Q406	8-729-905-35		2SC4081T106R	
	L602			47uH (DSR-20MD)	Q407	8-729-427-83	TRANSISTOR	XP6501-TXE	
	L602	1-410-388-31	INDUCTOR CHIP	39uH (DSR-20MDP)					
					Q408	8-729-402-42		UN5213-TX	
	L603		INDUCTOR CHIP		Q409	8-729-015-76		UN5211-TX	
	L604		INDUCTOR CHIP		Q410	8-729-905-35		2SC4081T106R	
	L651		INDUCTOR CHIP		Q411	8-729-905-35		2SC4081T106R	
	L652 L653		INDUCTOR CHIP		Q412	8-729-905-35	TRANSISTOR	2SC4081T106R	
	L033	1-410-303-11	INDUCTOR CHIP	22011	0412	9 700 005 05	TDANICICTOR	000400474000	
	L654	1_/112_020_11	INDUCTOR CHIP	1004	Q413 Q414	8-729-905-35 8-729-905-35		2SC4081T106R	
	L655		INDUCTOR CHIP		Q414	8-729-427-83		2SC4081T106R XP6501-TXE	
	L656		INDUCTOR CHIP		Q416	8-729-402-84		XN4601-TW	
	L657		INDUCTOR CHIP		Q417	8-729-427-83		XP6501-TXE	
	L851		INDUCTOR CHIP		W-111	3 120 721 00		AT OUT TAL	
					Q419	8-729-905-35	TRANSISTOR	2SC4081T106R	
	L852	1-412-026-11	INDUCTOR CHIP	1uH	Q420	8-729-905-35		2SC4081T106R	
	L853		INDUCTOR CHIP		Q601	8-729-905-35		2SC4081T106R	
	L854	1-412-029-11	INDUCTOR CHIP	10uH	Q602	8-729-905-35		2SC4081T106R	
	L855	1-412-028-11	INDUCTOR CHIP	4.7uH	Q604	8-729-905-35		2SC4081T106R	
	L856	1-412-028-11	INDUCTOR CHIP	4.7uH					
				·	Q605	8-729-905-35		2SC4081T106R	
	L857		INDUCTOR CHIP		Q606	8-729-905-35		2SC4081T106R	
	L858	1-412-028-11	INDUCTOR CHIP	4.7uH	Q607	8-729-026-52		2SA1576A-T106-R	
			TDANOICECT		Q613	8-729-026-52		2SA1576A-T106-R	
, di			< TRANSISTOR >	1 1984	Q651	8-729-905-35	TRANSISTOR	2SC4081T106R	. *
	Q101	9_720_015_70	TRANSISTOR		0000	0 700 005 05	TDANICIOTOS	000400474007	
	Q102	8-729-015-76 8-729-015-76		UN5211-TX UN5211-TX	Q653 Q655	8-729-905-35 8-729-026-52		2SC4081T106R	
	Q103	8-729-015-76		UN5211-TX	Q656			2SA1576A-T106-R	
	Q107	8-729-015-76		UN5211-TX	Q658	8-729-905-35 8-729-026-52		2SC4081T106R	
	Q108	8-729-015-76		UN5211-TX	Q701	8-729-015-76		2SA1576A-T106-R UN5211-TX	
					4.01	3 7 2 3 0 10 7 0	, individuo i Un	ONOT II-IV	
	Q201	8-729-905-35	TRANSISTOR	2SC4081T106R	0702	8-729-015-74	TRANSISTOR	UN5111-TX	

	· · · · · · · · · · · · · · · · · · ·										
Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
Q703 Q704		TRANSISTOR TRANSISTOR	2SC408 2SC408			R124	1-216-809-11	METAL CHIP	100	5%	1/16W
Q705		TRANSISTOR	UN5213			R125	1-216-809-11	METAL CHIP	100	E0/.	1/16/1/
Q706		TRANSISTOR	UN5113			R126	1-216-809-11		100	5% 5%	1/16W 1/16W
Q. 00	0 720 400 00	THAIGHTON	0110110	-17		R127	1-216-797-11		100	5% 5%	1/16W
Q707	8-720-015-76	TRANSISTOR	UN5211	_TV		R128	1-216-797-11				
Q708		TRANSISTOR	UN2225						100	5%	1/16W
Q709		TRANSISTOR				R129	1-216-797-11	WETAL CHIP	10	5%	1/16W
			UN2225			D400	4 040 000 44		400		
Q710		TRANSISTOR	UN2225			R130	1-216-809-11		100	5%	1/16W
Q711	8-729-028-70	TRANSISTOR	UN2225	I- (1X)		R131	1-216-809-11		100	5%	1/16W
0740						R132	1-216-797-11		10	5%	1/16W
Q712		TRANSISTOR	UN2225	T- (TX)		R133	1-216-809-11		100	5%	1/16W
Q713		TRANSISTOR	UN2225			R134	1-216-797-11	METAL CHIP	10	5%	1/16W
Q851		TRANSISTOR	2SC408								
Q852		TRANSISTOR	UN5213	-TX		R135	1-216-809-11	METAL CHIP	100	5%	1/16W
Q853	8-729-014-91	TRANSISTOR	2SD218	5S-TX		R136	1-216-809-11	METAL CHIP	100	5%	1/16W
						R137	1-216-809-11	METAL CHIP	100	5%	1/16W
Q854	8-729-905-35	TRANSISTOR	2SC4081	1T106R		R138	1-216-809-11		100	5%	1/16W
Q855		TRANSISTOR	2SC4081			R139	1-216-809-11		100	5%	1/16W
Q856		TRANSISTOR	UN5113-							. 070	171011
Q857		TRANSISTOR	UN5213-			R140	1-216-809-11	METAL CHIP	100	5%	1/16W
4007	0 720 102 12	111111111111111111111111111111111111111	0110210	174		R141	1-216-809-11		100		1/16W
		< RESISTOR >				R142	1-216-809-11			5%	
		C REGISTOR >							100	5%	1/16W
R051	1-216-833-91	pre cuió	101/	F0/	4 (4 0)41	R143	1-216-809-11		100	5%	1/16W
			10K	5%	1/16W	R144	1-216-809-11	METAL CHIP	100	5%	1/16W
R052	1-216-833-91		10K	5%	1/16W						
R053	1-216-833-91		10K	5%	1/16W	R145	1-216-809-11		100	5%	1/16W
R054	1-216-833-91		10K	5%	1/16W	R146	1-216-797-11		10	5%	1/16W
R056	1-216-864-11	METAL CHIP	0	5%	1/16W	R147	1-216-797-11	METAL CHIP	10	5%	1/16W
						R148	1-216-797-11	METAL CHIP	10	5%	1/16W
R059	1-216-823-11	METAL CHIP	1.5K	5%	1/16W	R152	1-216-833-91	RES, CHIP	10K	5%	1/16W
R060	1-216-801-11	METAL CHIP	22	5%	1/16W						
R061	1-216-857-11	METAL CHIP	1M	5%	1/16W	R153	1-216-833-91	RES. CHIP	10K	5%	1/16W
R062	1-216-801-11	METAL CHIP	22	5%	1/16W	R154	1-216-833-91		10K	5%	1/16W
R063	1-216-801-11		22	5%	1/16W	R155	1-216-821-11		1K	5%	1/16W
					.,	R156	1-216-809-11		100	5%	1/16W
R064	1-216-833-91	RES CHIP	10K	5%	1/16W	R157	1-216-809-11		100	5%	1/16W
R065			0	5%	1/16W	"""	1 210 000 11	WEINE OIL	100	J 70	171044
			Ū		SR-20MDP)	R158	1-216-809-11	METAL CHID	100	5%	1/16W
R066	1-216-801-11	METAL CHIP	22	5%	1/16W	R159	1-216-809-11		100		
R067	1-216-857-11		1M	5%	1/16W	5. 1				5%	1/16W
R068	1-216-801-11		22		1/16W	R160	1-216-809-11		100	5%	1/16W
nooo	1-210-001-11	WEIAL OHIP	22	5%	1/1000	R161	1-216-809-11		100	5%	1/16W
DOGO	1 010 000 11	METAL CUID	E CV	F0/	4 (4 (3))	R162	1-216-809-11	METAL CHIP	100	5%	1/16W
R069			5.6K	5%	1/16W	7400	1 010 000 11			-	
R070	1-216-833-91	RES, CHIP	10K	5%	1/16W	R163	1-216-809-11		100	5%	1/16W
R071	1-216-833-91	RES, CHIP	10K	5%	1/16W	R164			100	5%	1/16W
R101	1-216-797-11	METAL CHIP	10	5%	1/16W	R165	1-216-809-11		100	5%	1/16W
R102	1-216-797-11	METAL CHIP	10	5%	1/16W	R166	1-216-833-91		10K	5%	1/16W
						R167	1-216-809-11	METAL CHIP	100	5%	1/16W
R103	1-216-833-91	•	10K	5%	1/16W						
R104	1-216-833-91	RES, CHIP	10K	5%	1/16W	R168	1-216-809-11	METAL CHIP	100	5%	1/16W
R105	1-216-797-11	METAL CHIP	10	5%	1/16W	R169	1-216-809-11	METAL CHIP	100	5%	1/16W
R106	1-216-797-11	METAL CHIP	10	5%	1/16W	R170	1-216-833-91	RES, CHIP	10K	5%	1/16W
R107	1-216-797-11		10	5%	1/16W	R171		METAL CHIP	100	5%	1/16W
						R172	1-216-809-11		100	5%	1/16W
R108	1-216-833-91	RES CHIP	10K	5%	1/16W		1 210 000 11	MEDICE OTTO	100	0 / 0	1,71011
R109	1-216-797-11		10	5%	1/16W	R173	1-216-809-11	METAL CHID	100	5%	1/16W
R110	1-216-797-11		10	5%	1/16W	R174	1-216-809-11		100	5%	1/16W
R112	1-216-833-91		10K	5%	1/16W	R175					
R113							1-216-809-11		100	5%	1/16W
niio	1-216-864-11	METAL CHIP	0 .	5%	1/16W	R176	1-216-809-11		100	5%	1/16W
Dida	4 040 000 44	MATTAL OLUB	100	E0/	4400	R177	1-216-809-11	METAL CHIP	100	5%	1/16W
R114	1-216-809-11		100	5%	1/16W						
R115		METAL CHIP	100	5%	1/16W	R178	1-216-809-11		100	5%	1/16W
R117		METAL CHIP	100	5%	1/16W	R179	1-216-809-11		100	5%	1/16W
R118		METAL CHIP	10	5%	1/16W	R180	1-216-809-11	METAL CHIP	100	5%	1/16W
R119	1-216-797-11	METAL CHIP	10	5%	1/16W	R181	1-216-809-11	METAL CHIP	100	5%	1/16W
				(DS	R-20MDP)	R201	1-216-805-11	METAL CHIP	47	5%	1/16W
			188								
R120	1-216-809-11	METAL CHIP	100	5%	1/16W	R202	1-216-821-11	METAL CHIP	1K	5%	1/16W
R121	1-216-809-11	METAL CHIP	100	5%	1/16W	R203	1-216-810-11		120	5%	1/16W
R122	1-216-809-11	METAL CHIP	100	5%	1/16W	R204	1-216-821-11		1K	5%	1/16W
R123	1-216-809-11		100	5%	1/16W	R205		METAL CHIP	22K	5%	1/16W
			177					01111		370	17.1011

Ref. No.	Part No.	Description			<u>Remark</u>	Ref. No.	Part No.	Description			Remark
R206	1-216-833-91	RES, CHIP	10K	5%	1/16W	R272	1-216-821-11	METAL CHIP	1K	5%	1/16W
TIEGO	1 210 000 01	rico, orm	1011	070	171011	R273	1-216-819-11		680	5%	1/16W
R207	1-216-821-11	METAL CHIP	1K	5%	1/16W	R274	1-216-818-11		560	5%	1/16W
R208			2.2K	5%		R275					
	1-216-825-11				1/16W	N2/5	1-216-809-11	METAL CHIP	100	5%	1/16W
R209	1-216-821-11		1K	5%	1/16W	0070	4 040 004 44	METAL OLUB	414	F0/	4 14 0141
R210	1-216-821-11		1K	5%	1/16W	R276		METAL CHIP	1K	5%	1/16W
R211	1-216-837-11	METAL CHIP	22K	5%	1/16W	R277		METAL CHIP	0	5%	1/16W
						R279		METAL CHIP	1K	5%	1/16W
R212	1-216-834-11	METAL CHIP	12K	5%	1/16W	R280	1-216-825-11	METAL CHIP	2.2K	5%	1/16W
R213	1-216-821-11		1K	5%	1/16W	R281	1-216-817-11		470	5%	1/16W
R214	1-216-817-11		470	5%	1/16W					•	
R215	1-216-816-11		390	5%	1/16W	R282	1-218-899-11	METAL CHIP	150K	0.5%	1/16W
R216	1-216-821-11		1K	5%	1/16W	R283	1-218-871-11		10K	0.5%	1/16W
NZ10	1-210-021-11	WE TAL OTHE	IIV	370	1/1044						
D04#	4 040 007 44	METAL OLUB		=0/	4.44.00.84	R284	1-216-825-11		2.2K	5%	1/16W
R217	1-216-807-11		68	5%	1/16W	R285	1-216-864-11		0	5%	1/16W
R218	1-216-837-11		22K	5%	1/16W	R287	1-218-879-1 1	METAL CHIP	22K	0.5%	1/16W
R219	1-216-833-91	RES, CHIP	10K	5%	1/16W						
R220	1-216-821-11	METAL CHIP	1K	5%	1/16W	R288	1-216-844-11	METAL CHIP	82K	5%	1/16W
R221	1-216-817-11	METAL CHIP	470	5%	1/16W	R289	1-216-825-11		2.2K	5%	1/16W
				• , •	.,	R290	1-218-869-11		8.2K	0.5%	1/16W
R222	1-216-823-11	METAL CLID	1.5K	5%	1/16W	R291	1-216-817-11		470	5%	1/16W
R223	1-216-821-11		1K	5%	1/16W	R292	1-218-871-11	WETAL CHIP	10K	0.5%	1/16W
R224	1-216-841-11		47K	5%	1/16W						
R225	1-216-817-11	METAL CHIP	470	5%	1/16W	R293	1-216-833-91	RES, CHIP	10K	5%	1/16W
R226	1-216-809-11	METAL CHIP	100	5%	1/16W					(E	OSR-20MD)
						R294	1-216-821-11	METAL CHIP	1K	5%	1/16W
R227	1-216-821-11	METAL CHIP	1K	5%	1/16W	R295	1-216-853-11		470K	5%	1/16W
R228	1-216-809-11		100	5%	1/16W	R296	1-216-825-11		2.2K	5%	1/16W
R229	1-216-821-11		1K	5%	1/16W	R297	1-216-819-11				
						n291	1-210-019-11	METAL CHIP	680	5%	1/16W
R231	1-216-823-11		1.5K	5%	1/16W						
R232	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R298	1-216-818-11		560	5%	1/16W
						R299	1-216-829-11		4.7K	5%	1/16W
R233	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R300	1-216-805-11	METAL CHIP	47	5%	1/16W
R234	1-216-809-11	METAL CHIP	100	5%	1/16W	R301	1-216-815-11	METAL CHIP	330	5%	1/16W
R235	1-216-845-11		100K	5%	1/16W	R302	1-216-837-11	METAL CHIP	22K	5%	1/16W
R236	1-216-829-11		4.7K	5%	1/16W		1 270 007 11			0,0	.,
R237	1-216-828-11		3.9K	5%	1/16W	R303	1-216-864-11	METAL CHIP	0	5%	1/16W
N201	1-210-020-11	MILIAL OTHE	0.51	3 /0	171000	novo	1-210-004-11	METAL CHIL	Ū		
D044	4 040 004 44	METAL OLUD	•	E0/	4 (4 0) 44	D004	4 040 004 44	METAL OLUB		,	SR-20MDP)
R241	1-216-864-11	METAL CHIP	0	5%	1/16W	R304	1-216-864-11	METAL CHIP	0	5%	1/16W
					SR-20MDP)						SR-20MD)
R242	1-216-833-91		10K	5%	1/16W	R305	1-216-819-11	METAL CHIP	680	5%	1/16W
R244	1-216-823-11		1.5K	5%	1/16W	R306	1-216-811-11	METAL CHIP	150	5%	1/16W
R246	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R308	1-216-816-11	METAL CHIP	390	5%	1/16W
R248	1-216-813-11		220	5%	1/16W						
					.,	R309	1-216-845-11	METAL CHIP	100K	5%	1/16W
R249	1-216-864-11	METAL CHIP	0	5%	1/16W	R310	1-216-837-11		22K	5%	1/16W
11243	1-210-004-11	MILIAL OITH	U								
2050	1 010 000 11	LATTAL OLUB	0.017		SR-20MD)	R311	1-216-837-11		22K	5%	1/16W
R250	1-216-828-11		3.9K	5%	1/16W	R312	1-216-818-11		560	5%	1/16W
R251	1-216-813-11		220	5%	1/16W	R313	1-216-821 - 11	METAL CHIP	1K	5%	1/16W
R252	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	4					
R253	1-216-818-11	METAL CHIP	560	5%	1/16W	R314	1-216-826-11	METAL CHIP	2.7K	5%	1/16W
						R315	1-216-825-11		2.2K	5%	1/16W
R254	1-218-863-11	METAL CHIP	4.7K	0.5%	1/16W	R316	1-216-821-11				
R255	1-218-707-11		4.7K	5%	1/16W	R317			1K	5%	1/16W
							1-216-825-11		2.2K	5%	1/16W
R256	1-218-269-11	•	360	5%	1/16W	R318	1-216-837-11	METAL CHIP	22K	5%	1/16W
R257	1-216-864-11	METAL CHIP	0	5%	1/16W						
				(0	SR-20MD)	R319	1-216-840-11	METAL CHIP	39K	5%	1/16W
R258	1-218-823-11	METAL CHIP	100	0.5%	1/16W					(DS	R-20MDP)
						R320	1-216-837-11	METAL CHIP	22K	5%	1/16W
R259	1-216-864-11	METAL CHIP	0	5%	1/16W	R321	1-216-839-11		33K	5%	1/16W
R260	1-216-825-11		2.2K	5%	1/16W	R323	1-216-818-11				
									560	5%	1/16W
R261	1-216-849-11		220K	5%	1/16W	R324	1-216-809-11	METAL CHIP	100	5%	1/16W
R262	1-216-809-11		100	5%	1/16W						
R263	1-216-821-11	METAL CHIP	1K	5%	1/16W	R325	1-216-816-11	METAL CHIP	390	5%	1/16W
						R326	1-216-821-11	METAL CHIP	1K	5%	1/16W
R264	1-216-818-11	METAL CHIP	560	5%	1/16W	R327	1-216-825-11		2.2K	5%	1/16W
R265	1-216-825-11		2.2K	5%	1/16W	R328	1-216-839-11		33K	5%	1/16W
R266	1-216-825-11		2.2K	5%	1/16W	,1020	1 210 000 11	THE ITTE OTT	JUIN		
R269	1-216-837-11		2.2K	5%	1/16W	Doon	1_016_005_44	METAL CHID	אומ מ		R-20MDP)
						R329	1-216-825-11	WE TAL CHIP	2.2K	5%	1/16W
R270	1-216-864-11	WE TAL CHIP	0	5%	1/16W	Baca	4 040 055 11		0.211		
						R330	1-218-869-11		8.2K	0.5%	1/16W
R271	1-216-837-11	METAL CHIP	22K	5%	1/16W	R401	1-216-821-11	METAL CHIP	1K	5%	1/16W

ļ	Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
	R402	1-216-864-11	METAL CHIP	0	5%	1/16W	R465	1-218-847-11	METAL CHIP	1K	0.50/	
	R404	1-216-864-11		0	5%	1/16W	R466				0.5%	1/16W
	11101	1 210 004 11	WEIAL OIM	U				1-216-829-11		4.7K	5%	1/16W
	D406	1.016.064.11	METAL CUID	0		SR-20MDP)	R467	1-216-841-11		47K	5%	1/16W
	R406	1-210-004-11	METAL CHIP	0	5%	1/16W	R468	1-216-841-11	METAL CHIP	47K	5%	1/16W
					(DSR-20MD)						
							R469	1-216-839-11	METAL CHIP	33K	5%	1/16W
	R407	1-216-835-11		15K	5%	1/16W	R470	1-216-819-11	METAL CHIP	680	5%	1/16W
	R408	1-216-833-91	RES, CHIP	10K	5%	1/16W	R471	1-216-841-11	METAL CHIP	47K	5%	1/16W
	R410	1-216-817-11	METAL CHIP	470	5%	1/16W	R472	1-216-864-11		0	5%	1/16W
	R411	1-216-864-11	METAL CHIP	0	5%	1/16W	R473	1-216-821-11		1K	5%	1/16W
	R412	1-216-837-11		22K	5%	1/16W	11770	1-210-021-11	METAL OTHE	I I	376	1/10//
		1 210 001 11	WEIRE OIT	2211	J /0	17 1000	D474	4 040 000 44	METAL OLUD	0.714	=0/	
	R413	1-216-817-11	METAL CHIP	470	E0/	4/4/014/	R474	1-216-826-11		2.7K	5%	1/16W
				470	5%	1/16W	R477	1-216-825-11		2.2K	5%	1/16W
	R414	1-216-819-11	METAL CHIP	680	5%	1/16W	R480	1-216-864-11		0	5%	1/16W
	R415	1-216-864-11		0	5%	1/16W	R481	1-216-864-11	METAL CHIP	0	5%	1/16W
	R416	1-216-835-11		15K	5%	1/16W	R482	1-216-864-11	METAL CHIP	0	5%	1/16W
	R417	1-216-833-91	RES, CHIP	10K	5%	1/16W					(DS	R-20MDP)
		•							•		,	,
	R418	1-216-835-11	METAL CHIP	15K	5%	1/16W	R484	1-216-864-11	METAL CHIP	0	5%	1/16W
	R419	1-216-833-91	RES, CHIP	10K	5%	1/16W		. 210 001 11	WED IL OTH	Ū		
	R420	1-216-813-11		220	5%	1/16W	R485	1-216-864-11	METAL CLUB	•		SR-20MD)
	R421	1-216-828-11	METAL CHIP	3.9K	5%	1/16W				0	5%	1/16W
	R422	1-216-847-11					R486	1-216-821-11		1K	5%	1/16W
	N422	1-210-047-11	METAL CHIP	150K	5%	1/16W	R487	1-216-821-11		1K	5%	1/16W
	D. 400'						R488	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R423	1-216-864-11	METAL CHIP	0	5%	1/16W						
	R424	1-216-821-11	METAL CHIP	1K	5%	1/16W	R489	1-216-821-11	METAL CHIP	1K	5%	1/16W
	R425	1-216-817-11	METAL CHIP	470	5%	1/16W	R492	1-216-821-11		1K	5%	1/16W
	R426	1-216-818-11	METAL CHIP	560	5%	1/16W	R493	1-216-833-91		10K	5%	1/16W
	R427	1-216-821-11	METAL CHIP	1K	5%	1/16W	R494		METAL CHIP	75	5%	1/10W
						SR-20MDP)	R495	1-216-833-91		10K	5%	
					(50	on zolvion ,	11730	1-210-030-91	neo, one	IUK	376	1/16W
	R428	1-216-821-11	METAL CHIP	1K	5%	1/16W	D406	1 010 041 14	MACTAL OLUD	471/	F0/	4.4.004
	R429	1-216-829-11	METAL CHIP				R496	1-216-841-11	METAL CHIP	47K	5%	1/16W
				4.7K	5%	1/16W	R497	1-216-022-00		75	5%	1/10W
	R430	1-216-833-91	RES, CHIP	10K	5%	1/16W	R498	1-216-013-00		33	5%	1/10W
	R431	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R499	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R433	1-216-833-91	RES, CHIP	10K	5%	1/16W	R500	1-216-015-00	METAL CHIP	39	5%	1/10W
											• , ,	.,
	R434	1-216-821-11	METAL CHIP	1K	5%	1/16W	R501	1-216-817-11	METAL CHIP	470	5%	1/16W
						R-20MDP)	R502	1-216-845-11	METAL CHIP	100K	5%	1/16W
	R436	1-216-817-11	METAL CHIP	470	5%	1/16W	R503	1-216-821-11				
	R437	1-216-821-11	METAL CHIP	1K	5%	1/16W			METAL CHIP	1K	5%	1/16W
	R438		METAL CHIP				R504	1-216-864-11	METAL CHIP	0	5%	1/16W
				0	5%	1/16W	R505	1-216-845-11	METAL CHIP	100K	5%	1/16W
	R439	1-216-837-11	METAL CHIP	22K	5%	1/16W						
							R506	1-216-864-11	METAL CHIP	0	5%	1/16W
	R440	1-216-821-11		1K	5%	1/16W	R507	1-216-845-11	METAL CHIP	100K	5%	1/16W
	R441	1-216-819-11	METAL CHIP	680	5%	1/16W	R508	1-216-845-11	METAL CHIP	100K	5%	1/16W
	R442	1-216-837-11	METAL CHIP	22K	5%	1/16W	R509	1-216-845-11		100K	5%	1/16W
	R443		METAL CHIP	680	5%	1/16W	R510	1-216-845-11	METAL CHIP	100K	5%	1/16W
	R444		METAL CHIP	22K	5%	1/16W	,11010	1 210 010 11	WEINE OIM	10010	J /0	1/1000
					0 /0	171011	R511	1-216-845-11	METAL CUID	1001/	E0/	4/4034
	R445	1-218-838-11	METAL CHID	430	0.50/	1/16)4/			METAL CHIP	100K	5%	1/16W
i.	R446		RES, CHIP		0.5%	1/16W	R512	1-216-864-11	METAL CHIP	0	5%	1/16W
				10K	5%	1/16W	R513	1-216-829-11		4.7K	5%	1/16W
	R447		METAL CHIP	470	5%	1/16W	R514	1-216-849-11		220K	5%	1/16W
	R449	1-216-864-11	METAL CHIP	0	5%	1/16W	R517	1-216-864-11	METAL CHIP	0	5%	1/16W
					(D	SR-20MD)						
	R450	1-218-839-11	METAL CHIP	470	0.5%	1/16W	R519	1-216-864-11	METAL CHIP	0	5%	1/16W
							R520	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R451	1-218-846-11	METAL CHIP	910	0.5%	1/16W	R521		METAL CHIP	0	5%	1/16W
	R452	1-216-864-11		0	5%	1/16W	R524		METAL CHIP	1K	5%	
				·								1/16W
	R453	1-216-819-11	METAL CUID	680	5%	R-20MDP) 1/16W	R525	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R454		METAL CHIP	0 .	5%	1/16W	R526	1-216-833-91		10K	5%	1/16W
	R458	1-216-864-11	IVIETAL CHIP	0	5%	1/16W	R527	1-216-821-11	METAL CHIP	1K	5%	1/16W
											(DSR	-20MDP)
	R459		METAL CHIP	0	5%	1/16W	R528	1-216-864-11	METAL CHIP	0	5%	1/16W
	R460	1-216-864-11	METAL CHIP	0	5%	1/16W	R529		METAL CHIP	0		1/16W
	R461	1-216-864-11		0	5%	1/16W				-		R-20MD)
			= =			R-20MDP)	R531	1-216-864-11	METAL CHIP	0		
	R462	1-216-842-11	METAL CHIP	56K	5%	1/16W	11001	1-210-004-11	WILLIAL UTIL	U ,	5%	1/16W
	R463	1-216-833-91		10K			DEGG	1 010 004 11	NACTAL OUTS	•	EC.	4 (4 0) 11
	R464		•		5%	1/16W	R533		METAL CHIP	0		1/16W
	11404	1-216-864-11	METAL CHIP	0	5%	1/16W	R534		METAL CHIP	0		1/16W
					(DS	SR-20MD)	R535	1-216-864-11	METAL CHIP	0	5%	1/16W

R	ef. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description			Remark
				•	50 /			-		101/	E0/	1/16W
	R536	1-216-864-11		0	5%	1/16W	R702	1-216-833-91		10K	5%	
	R537	1-216-864-11	METAL CHIP	0	5%	1/16W	R703	1-216-833-91		10K	5%	1/16W
							R704	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R539	1-216-821-11	METAL CHIP	1K	5%	1/16W	R705	1-216-833-91	RES. CHIP	10K	5%	1/16W
	R540	1-216-864-11		0	5%	1/16W		,	,			
					5%	1/16W	R706	1-216-833-91	DEC CHID	10K	5%	1/16W
	R541	1-216-864-11		0								
	R542	1-216-864-11		0	5%	1/16W	R707	1-216-833-91		10K	5%	1/16W
	R543	1-216-864-11	METAL CHIP	0	5%	1/16W	R708	1-216-833-91		10K	5%	1/16W
					(DSR-20MD)	R709	1-216-849-11	METAL CHIP	220K	5%	1/16W
						,	R710	1-216-849-11		220K	5%	1/16W
	R543	1-211-983-11	METAL CHID	39	0.5%	1/16W						
	11040	1-211-300-11	WEINE OIM	00		SR-20MDP)	R711	1-216-839-11	METAL CHIP	33K	5%	1/16W
				40014								1/16W
	R601	1-216-841-11		47K	5%	1/16W	R712	1-216-839-11	METAL CHIP	33K	5%	
	R602	1-216-841-11	METAL CHIP	47K	5%	1/16W	R713	1-216-835-11		15K	5%	1/16W
	R603	1-216-813-11	METAL CHIP	220	5%	1/16W	R714	1-216-835-11	METAL CHIP	15K	5%	1/16W
	R604	1-216-849-11	METAL CHIP	220K	5%	1/16W	R715	1-216-839-11	METAL CHIP	33K	5%	1/16W
	11001	. 210 010 11			•	.,						
	DCCE	1 010 007 11	METAL CLUD	201/	5%	1/16W	R716	1-216-839-11	METAL CHID	33K	5%	1/16W
	R605	1-216-837-11		22K								
	R606	1-216-839-11		33K	5%	1/16W	R717	1-216-835-11		15K	5%	1/16W
	R610	1-216-817-11	METAL CHIP	470	5%	1/16W	R718	1-216-835-11		15K	5%	1/16W
	R611	1-216-816-11	METAL CHIP	390	5%	1/16W	R719	1-216-864-11	METAL CHIP	0	5%	1/16W
	R612	1-216-821-11		1K	5%	1/16W	R721	1-216-809-11	METAL CHIP	100	5%	1/16W
	11012	1 210 021 11	WE IT LE OTHE		0 / 0	., . • • •						
	DC10	1 010 017 11	METAL CHID	470	5%	1/16W	R722	1-216-841-11	METAL CHID	47K	5%	1/16W
	R613	1-216-817-11										
	R615	1-216-864-11		0	5%	1/16W	R723	1-216-837-11		22K	5%	1/16W
	R616	1-216-821-11	METAL CHIP	1K	5%	1/16W	R724	1-216-864-11		0	5%	1/16W
	R619	1-216-815-11	METAL CHIP	330	5%	1/16W	R726	1-216-809-11	METAL CHIP	100	5%	1/16W
	R621	1-216-821-11	METAL CHIP	1K	5%	1/16W	R727	1-216-845-11	METAL CHIP	100K	5%	1/16W
	HOLI	1 210 021 11			• , •	.,						
	nego	1-216-833-91	DEC CHID	10K	5%	1/16W	R728	1-216-845-11	METAL CHIP	100K	5%	1/16W
	R622		•									
	R623	1-216-853-11		470K	5%	1/16W	R729	1-216-845-11		100K	5%	1/16W
	R629	1-216-833-91	RES, CHIP	10K	5%	1/16W	R730	1-216-845-11		100K	5%	1/16W
	R630	1-216-836-11	METAL CHIP	18K	5%	1/16W	R731	1-216-837-11	METAL CHIP	22K	5%	1/16W
	R631	1-216-837-11	METAL CHIP	22K	5%	1/16W	R732	1-216-833-91	RES, CHIP	10K	5%	1/16W
	11001								• •			
	R635	1-216-864-11	METAL CHID	0	5%	1/16W	R733	1-216-833-91	BEC CHID	10K	5%	1/16W
								1-216-841-11		47K	5%	1/16W
	R636	1-216-833-91	the second second second	10K	5%	1/16W	R734					
	R639	1-216-821-11		1K	5%	1/16W	R735	1-216-833-91		10K	5%	1/16W
	R640	1-216-841-11	METAL CHIP	47K	5%	1/16W	R736	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R641	1-216-841-11	METAL CHIP	47K	5%	1/16W	R737	1-216-833-91	RES, CHIP	10K	5%	1/16W
	,,,,,,,								*			
	R642	1-216-821-11	METAL CHID	1K	5%	1/16W	R738	1-216-833-91	RES CHIP	10K	5%	1/16W
		1-216-845-11				1/16W	R739	1-216-819-11		680	5%	1/16W
	R643			100K	5%		1					
	R652	1-216-819-11		680	5%	1/16W	R740	1-216-819-11		680	5%	1/16W
	R653	1-216-864-11	METAL CHIP	0	5%	1/16W	R741	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R654	1-216-864-11	METAL CHIP	0	5%	1/16W	R742	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R657	1-216-831-11	METAL CHIP	6.8K	5%	1/16W	R743	1-216-841-11	METAL CHIP	47K	5%	1/16W
	R658	1-216-834-11		12K	5%	1/16W	R744	1-216-833-91		10K	5%	1/16W
											5%	1/16W
	R659	1-216-825-11		2.2K	5%	1/16W	R745	1-216-833-91		10K		
	R660	1-216-821-11	METAL CHIP	1K	5%	1/16W	R746	1-216-841-11		47K	5%	1/16W
	R662	1-216-821-11	METAL CHIP	1K	5%	1/16W	R747	1-216-833-91	RES, CHIP	10K	5%	1/16W
	R665	1-216-817-11	METAL CHIP	470	5%	1/16W	R748	1-216-833-91	RES, CHIP	10K	5%	1/16W
				0	5%	1/16W	R749	1-216-833-91		10K	5%	1/16W
	R666	1-216-864-11										1/16W
	R667	1-216-828-11		3.9K	5%	1/16W	R750	1-216-833-91		10K	5%	
	R668	1-216-847-11		150K	5%	1/16W	R753	1-218-332-11		130K	5%	1/16W
	R669	1-216-818-11	METAL CHIP	560	5%	1/16W	R754	1-218-332-11	RES, CHIP	130K	5%	1/16W
				:								
	R670	1-216-817-11	METAL CHIP	470	5%	1/16W	R755	1-216-833-91	RES, CHIP	10K	5%	1/16W
				100	5%	1/10W	R756	1-218-293-11		24K	5%	1/16W
	R676	1-216-025-91										
	R677	1-216-025-91		100	5%	1/10W	R757	1-218-293-11		24K	5%	1/16W
	R678	1-216-021-00		68	5%	1/10W	R758	1-216-833-91	•	10K	5%	1/16W
	R679	1-216-864-11	METAL CHIP	0	5%	1/16W	R759	1-216-841-11	METAL CHIP	47K	5%	1/16W
	R680	1-216-864-11	METAL CHID	0	5%	1/16W	R760	1-218-332-11	RES CHIP	130K	5%	1/16W
	R681	1-216-864-11	METAL CHIP	0	5%	1/16W	R761	1-218-332-11		130K	5%	1/16W
	R682	1-216-864-11		0	5%	1/16W	R762	1-216-841-11		47K	5%	1/16W
	R683	1-216-864-11		0	5%	1/16W	R763	1-216-833-91		10K	5%	1/16W
	R684	1-216-864-11	METAL CHIP	0	5%	1/16W	R764	1-216-849-11	METAL CHIP	220K	5%	1/16W
	R701	1-216-833-91	RES, CHIP	10K	5%	1/16W	R765	1-216-833-91	RES, CHIP	10K	5%	1/16W
	•				-							

											
Ref. No.	Part No.	Description			Remark	Dof No.	Dort No.	Description	•		
		•				Ref. No.	Part No.	<u>Description</u>			<u>Remark</u>
R766	1-218-882-1		30K	0.5%	1/16W	R865	1-216-821-11		1K	5%	1/16W
R767	1-218-862-1		4.3K	0.5%	1/16W	R866	1-216-864-11	METAL CHIP	0	5%	1/16W
R768	1-216-833-9		10K	5%	1/16W	R867	1-216-817-11	METAL CHIP	470	5%	1/16W
R769	1-216-833-9	RES, CHIP	10K	5%	1/16W	R868	1-216-817-11		470	5%	1/16W
											.,
R770	1-216-833-91	RES, CHIP	10K	5%	1/16W	R869	1-216-817-11	METAL CHIP	470	5%	1/16W
R771	1-216-833-91	RES, CHIP	10K	5%	1/16W	R870	1-216-817-11		470	5%	1/16W
R772	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R871	1-216-817-11		470	5%	1/16W
R773	1-216-825-11	METAL CHIP	2.2K	5%	1/16W	R872	1-216-817-11	METAL CHIP	470	5%	1/16W
R774	1-216-833-91		10K	5%	1/16W	R873	1-216-817-11		470	5%	1/16W
					., ,	1.0.0	. 210 017 11	WIETAL OTH	470	J /6	171000
R775	1-216-809-11	METAL CHIP	100	5%	1/16W	R874	1-216-817-11	METAL CHIP	470	5%	1/1014
R776	1-216-829-11		4.7K	5%	1/16W	R875	1-216-817-11				1/16W
R777	1-216-829-11		4.7K	5%	1/16W	R876			470	5%	1/16W
R778		METAL CHIP	4.7K	5%	1/16W		1-216-817-11		470	5%	1/16W
R779		METAL CHIP				R877			470	5%	1/16W
11//3	1-210-020-11	METAL UNIP	3.9K	5%	1/16W	R878	1-216-817-11	METAL CHIP	470	5%	1/16W
R780	1-216-828-11	METAL CUID	0.01/	E0/	4 /4 00 4 /						
			3.9K	5%	1/16W	R879	1-216-817-11	METAL CHIP	470	5%	1/16W
R781	1-216-849-11		220K	5%	1/16W	R880	1-216-817-11		470	5%	1/16W
R782	1-216-809-11		100	5%	1/16W	R881	1-216-817-11	METAL CHIP	470	5%	1/16W
R783	1-218-290-11		6.2K	5%	1/16W	R882	1-216-827-11	METAL CHIP	3.3K	5%	1/16W
R784	1-218-290-11	RES, CHIP	6.2K	5%	1/16W	R883	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
R785	1-216-831-11	METAL CHIP	6.8K	5%	1/16W	R884	1-216-833-91	RES, CHIP	10K	5%	1/16W
R786	1-218-882-11	METAL CHIP	30K	0.5%	1/16W	R885	1-216-817-11		470	5%	1/16W
R787	1-218-862-11	METAL CHIP	4.3K	0.5%	1/16W	R886	1-216-833-91		10K	5%	1/16W
R788	1-216-831-11		6.8K	5%	1/16W	R887	1-216-817-11	METAL CHIP	470	5%	1/16W
R789	1-216-829-11		4.7K	5%	1/16W	R888	1-216-864-11		0	5%	
		MEMILE OTH	7.710	0 /0	171000	11000	1-210-004-11	WE TAL OHIP	U	3%	1/16W
R790	1-216-809-11	METAL CHIP	100	5%	1/16W	R889	1-216-817-11	MACTAL CLUD	470	50 /	4 44 00 44
R791	1-216-809-11		100	5%	1/16W	R890			470	5%	1/16W
R792	1-216-840-11		39K				1-216-864-11		0	5%	1/16W
				5%	1/16W	R891	1-216-817-11		470	5%	1/16W
R793		METAL CHIP	100	5%	1/16W	R892	1-216-833-91		10K	5%	1/16W
R794	1-216-833-91	RES, CHIP	10K	5%	1/16W	R893	1-216-864-11	METAL CHIP	0	5%	1/16W
R795	1-216-840-11		39K	5%	1/16W	R894	1-216-833-91		10K	5%	1/16W
R796	1-216-809-11		100	5%	1/16W	R895	1-216-864-11	METAL CHIP	0	5%	1/16W
R797	1-216-821-11	METAL CHIP	1K	5%	1/16W	R896	1-216-817-11		470	5%	1/16W
R798	1-216-821-11	METAL CHIP	1K	5%	1/16W	R897		RES, CHIP	10K	5%	1/16W
R799	1-216-821-11	METAL CHIP	1K	5%	1/16W	R898	1-216-864-11		0	5%	1/16W
							. 2.0 00. 11		Ū	0 /0	171000
R800	1-216-813-11	METAL CHIP	220	5%	1/16W	R899	1-216-817-11	METAL CHIP	470	5%	1/16W
R801	1-216-813-11		220	5%	1/16W	R900	1-216-833-91		10K	5%	1/16W
R802	1-216-813-11		220	5%	1/16W	R901	1-216-833-91		10K	5%	1/16W
R803	1-216-813-11		220	5%	1/16W	R902	1-216-864-11				
R804	1-216-840-11		39K	5%	1/16W	R903	1-216-817-11		0 470	5%	1/16W
	1 270 010 11	WEINE OIII	OUIT	0 70	171000	11300	1-210-017-11	WE TAL CHIP	470	5%	1/16W
R805	1-216-840-11	METAL CHIP	39K	5%	1/16W	R904	1 016 000 01	DEC CUID	401/	F0/	4 (4 0)44
R807	1-216-864-11		0				1-216-833-91		10K	5%	1/16W
R808	1-216-864-11			5%	1/16W	R905	1-216-833-91		10K	5%	1/16W
			0	5%	1/16W	R906	1-216-864-11		0	5%	1/16W
R809	1-216-833-91		10K	5%	1/16W	R907	1-216-833-91		10K	5%	1/16W
R810	1-216-833-91	RES, CHIP	10K	5%	1/16W	R908	1-216-864-11	METAL CHIP	0	5%	1/16W
D014	1 010 000 01	DEC CUID	4016								
R811	1-216-833-91		10K	5%	1/16W	R909	1-216-817-11	METAL CHIP	470	5%	1/16W
R812	1-216-833-91		10K	5%	1/16W	R910	1-216-817-11	METAL CHIP	470	5%	1/16W
R813	1-216-833-91	RES, CHIP	10K	5%	1/16W	R911	1-216-833-91	RES, CHIP	10K	5%	1/16W
R814	1-216-833-91	RES, CHIP	10K	5%	1/16W	R912	1-216-817-11		470	5%	1/16W
R815	1-216-864-11	METAL CHIP	0	5%	1/16W	R913	1-216-833-91		10K	5%	1/16W
										•	1, 1011
R832	1-216-801-11	METAL CHIP	22	5%	1/16W	R914	1-216-841-11	METAL CHIP	47K	5%	1/16W
R855	1-216-841-11		47K	5%	1/16W	R915	1-216-817-11		470	5%	1/16W
R856	1-216-817-11		470	5%	1/16W	R916	1-216-833-91		10K	5%	1/16W
R857	1-216-841-11		47K	5%	1/16W	R917	1-216-817-11		470		
R858	1-216-833-91		10K	5%	1/16W	R918				5%	1/16W
	. 2.0 000 01	, orm	·	J /0	1/1000	11910	1-216-833-91	neo, unir	10K	5%	1/16W
R859	1-216-833-91	RES CHIP	10K	5%	1/16W	D010	1_016 0AE 44	METAL OUT	1001	En/	4 /4 03.57
R860	1-216-833-91					R919	1-216-845-11		100K	5%	1/16W
			10K	5%	1/16W	R920	1-216-841-11		47K	5%	1/16W
R861	1-216-833-91		10K	5%	1/16W	R921	1-216-817-11		470	5%	1/16W
R862	1-216-829-11		4.7K	5%	1/16W	R922		RES, CHIP	10K	5%	1/16W
R863	1-216-817-11	WEIAL CHIP	470	5%	1/16W	R923	1-216-841-11	METAL CHIP	47K	5%	1/16W
D004	4 040 04= 44	AACTA! OUT	170	Ec.'							
R864	1-216-817-11	WETAL CHIP	470	5%	1/16W	R924	1-216-833-91 F	RES, CHIP	10K	5%	1/16W

Restance Part No. Description Part No. Description Part No. Par												
R826 1-21-86-11 METAL CHIP 0 54 1/16W R930 1-21-86-11 METAL CHIP 0 55 1/16W R1002 1-21-8-11 METAL CHIP 470 55 1/16W R1002 1-21-8-8-11 METAL CHIP 0 55 1/16W R100	Ref. No.	Part No.	Description			Remark	Ref. No.	Part No.	Description	ā ·		Remark
R986 1-216-884-11 METAL CHIP 0 5% 1/16W R981 1-216-864-11 METAL CHIP 0 5% 1/16W R981 1-216-864-11 METAL CHIP 0 5% 1/16W R1001 1-216-817-11 METAL CHIP 0 5% 1/16W R1002 1-216-817-11 METAL CHIP 0 5% 1/16W R1002 1-216-817-11 METAL CHIP 0 5% 1/16W R1002 1-216-817-11 METAL CHIP 470 5% 1/16W R1002 1-216-817-11 METAL CHIP 270 5% 1/16W R1002 1-216-827-11 METAL CHIP 270 5% 1/16W R1002 1-216-828-11 M			METAL CHID	470	50/	1/16\\	B007	1-216-817-11	METAL CHIP	470	5%	1/16W
R930 1-218-88-411 METAL CHIP 0 5% 1/16W R999 1-218-89-411 METAL CHIP 0 5% 1/16W R1001 1-216-81-11 METAL CHIP 0 5% 1/16W R1002 1-216-81-11 METAL CHIP 0 5% 1/16W R1002 1-216-81-11 METAL CHIP 470 5% 1/16W R1003 1-216-82-11 METAL CHIP 470 5% 1/16W												
R331 1-216-884-11 METAL CHIP 0 5% 1/16W R332 1-216-864-11 METAL CHIP 0 5% 1/16W R332 1-216-864-11 METAL CHIP 0 5% 1/16W R334 1-216-864-11 METAL CHIP 0 5% 1/16W R1003 1-216-811-11 METAL CHIP 470 5% 1/16W R335 1-216-864-11 METAL CHIP 0 5% 1/16W R1003 1-216-811-11 METAL CHIP 270 5% 1/16W R1035 1-216-864-11 METAL CHIP 0 5% 1/16W R1035 1-216-811-11 METAL CHIP 470 5% 1/16W R1035 1-216-864-11 METAL CHIP 0 5% 1/16W R1035 1-216-811-11 METAL CHIP 470 5% 1/16W R1035 1-216-864-11 METAL CHIP 0 5% 1/16W R1028 1-216-811-11 METAL CHIP 470 5% 1/16W R1028 1-216-834-11 METAL CHIP 0 5% 1/16W R1028 1-216-824-11 METAL CHIP 0 5% 1/16W R1028 1-216-864-11 METAL CHIP 0 5% 1/16W R1028 1-216-864-11 METAL CHIP 270 5% 1/16W R1038 1-216-864-11 METAL CHIP 270 5% 1/16W R103												
R832 1-216-864-11 METAL CHIP 0 5% 1/16W R833 1-216-864-11 METAL CHIP 0 5% 1/16W R835 1-216-864-11 METAL CHIP 0 5% 1/16W R836 1-216-864-11 METAL CHIP 0 5% 1/16W R836 1-216-864-11 METAL CHIP 0 5% 1/16W R1003 1-216-817-11 METAL CHIP 22 3% 1/16W R1005 1-216-817-11 METAL CHIP 22 3% 1/16W R1005 1-216-817-11 METAL CHIP 24 3% 1/16W R1005 1-216-817-11 METAL CHIP 470 5% 1/16W R1028 1-216-864-11 METAL CHIP 0 5% 1/16W R1028 1-216-825-11 METAL CHIP 0 5% 1/16W R1028 1-216-825-11 METAL CHIP 0 5% 1/16W R1028 1-216-800 RES. CHIP 10 5% 1/16W R1028 1-216-800 RES. CHIP 1/16W R1028 1-216-800 RE												
R832 1-216-884-11 METAL CHIP 0 5% 1/16W R934 1-216-884-11 METAL CHIP 0 5% 1/16W R934 1-216-884-11 METAL CHIP 0 5% 1/16W R936 1-216-884-11 METAL CHIP 0 5% 1/16W R938 1-216-884-11 METAL CHIP 0 5% 1/16W R938 1-216-884-11 METAL CHIP 0 5% 1/16W R938 1-216-884-11 METAL CHIP 0 5% 1/16W R940 1-216-821-11 METAL CHIP 1 2 5% 1/16W R940 1-216-821-11 METAL CHIP 1 2 5% 1/16W R940 1-216-821-11 METAL CHIP 2 5% 1/16W R1031 1-216-860 RS. CHIP 47 5% 1/16W R1031 1-216-860 RS	11901	1-210-004-11	MILIAL OIII	U	0 /0	171044	1					
R833 -216-884-11 METAL CHIP 0 5% 1/16W R935 1-216-884-11 METAL CHIP 2 5% 1/16W R935 1-216-884-11 METAL CHIP 0 5% 1/16W R936 1-216-884-11 METAL CHIP 0 5% 1/16W R936 1-216-884-11 METAL CHIP 0 5% 1/16W R936 1-216-884-11 METAL CHIP 0 5% 1/16W R1020 1-216-817-11 METAL CHIP 470 5% 1/16W R1030 1-216-827-10 METAL CHIP 470 5% 1/16W R1030 1-216-827-10 METAL CHIP 2 5% 1/16W R1030 1-216-827-10 METAL CHIP 3 3 3 3 3 3 3 3 3	B932	1-216-864-11	METAL CHIP	0	5%	1/16W	"""	. 2.0 0.7 1.			• / •	
R985 1-216-88-11 METAL CHIP 0 5% 1/16W R1005 1-216-81711 METAL CHIP 27 5% 1/16W R1005 1-216-81711 METAL CHIP 470 5% 1/16W R1025 1-216-827-10 METAL CHIP 20 5% 1/16W R1025 1-216-025-11 METAL CHIP 20 5% 1/16W R1025 1-226-025-11				-			R1003	1-216-817-11	METAL CHIP	470	5%	1/16W
R395 1-216-886-11 METAL CHIP 0 5% 1/16W R1021 1-216-817-11 METAL CHIP 470 5% 1/16W R1028 1-216-817-11 METAL CHIP 270 5% 1/16W R1028 1-216-817-11 METAL CHIP 20 5% 1/16W R1028 1-216-817-11 METAL CHIP 20 5% 1/16W R1038 1-216-620-81 METAL CHIP 20 5% 1/16W R1038 1-216-680-81 METAL CHIP 20 5% 1/16W R1038 1-21				-								
R936 1-216-984-11 METAL CHIP 0 5% 1/16W R1028 1-216-817-11 METAL CHIP 470 5% 1/16W R383 1-216-884-11 METAL CHIP 0 5% 1/16W R394 1-216-884-11 METAL CHIP 0 5% 1/16W R1030 1-216-025-01 RES, CHIP 100 5% 1/16W R1030 1-216-05-01 RES, CHIP 100 5% 1/16W R1030 1-216-05-01 RES, CHIP 100 5% 1/16W R1030 1-216-05-00 RES, CHIP 470 5% 1/16W R1030 1-216-05-00 RES, CHIP 470 5% 1/16W R1031 1-216-166-00 RES, CHIP 470 5% 1/16W R1031 1-216-884-11 METAL CHIP 5% 1/16W R1031 1-216-885-11 METAL CHIP 5% 1/16W R1031 1-216-885-11 METAL CHIP 100 100 100 100 100 100 100 100												
R397 1-216-864-11 METAL CHIP 0 5% 1/16W R398 1-216-864-11 METAL CHIP 0 5% 1/16W R394 1-216-864-11 METAL CHIP 0 5% 1/16W R394 1-216-864-11 METAL CHIP 1K 5% 1/16W R1031 1-216-165-00 RES, CHIP 47 5% 1/16W R394 1-216-841-11 METAL CHIP 47K 5% 1/16W R1031 1-216-165-00 RES, CHIP 47 5% 1/16W R1031 1-216-165-00 RES, CHIP 5% 1/16W R1031 1-216-165-00 RES, CHIP 47 5% 1/16W R1031 1-216-165-00 RES, CHIP 5% 1/16W R1035 1-216-462-11 METAL CHIP 20 5% 1/16W R1035 1-228-453-11 RES, ADJ, CERNET IX (FEC OL LEVEL) R1036 1-216-482-11 METAL CHIP 20 5% 1/16W R1035 1-228-453-11 RES, ADJ, CERNET IX (FEC OL LEVEL) R1036 1-216-482-11 METAL CHIP 20 5% 1/16W R1035 1-228-453-11 RES, ADJ, CERNET IX (FEC OL LEVEL) R1036 1-216-482-11 METAL CHIP 20 5% 1/16W R1035 1-238-453-11 RES, ADJ, CERNET IX (FEC OL LEVEL) R1036 1-216-482-11 METAL CHIP 20 5% 1/16W R1035 1-238-453-11 RES, ADJ, CERNET IX (FEC OL LEVEL) R1036 1-216-482-11 METAL CHIP 470 5% 1/16W R1035 1-238-453-11 RES, ADJ, CERNET I												
R939 1-216-686-11 METAL CHIP 0 5% 1/16W R940 1-216-686-11 METAL CHIP 0 5% 1/16W R940 1-216-686-11 METAL CHIP 0 5% 1/16W R1030 1-216-680-01 RES, CHIP 100 5% 1/10W R1030 1-216-680-01 RES, CHIP 100 5% 1/10W R1031 1-216-680-01 RES, CHIP 47 5% 1/16W R1032 1-216-1660-00 RES, CHIP 47 5% 1/16W R1032 1-216-1680-01 RES, CHIP 47 5% 1/16W R1032 1-216-680-01 RES, CHIP 47 5% 1/16W R1032 1-216-886-11 METAL CHIP 5% 1/16W R1032 1-216-886-11 METAL CHIP 5% 1/16W R1035 1-216-888-11 METAL CHIP 5% 1/16W R1035 1-216-888-11 METAL CHIP 5% 1/16W R1035 1-216-881-1 METAL CHIP 100 5% 1/16W R1035 1-216-481-1 METAL CHIP 100	71000	1 210 001 11		•	0 / 0							
R988 1-216-88-4-11 METAL CHIP 0 5% 1/16W R944 1-216-88-4-11 METAL CHIP 1K 5% 1/16W R1039 1-216-16-20-01 RES, CHIP 47 5% 1/16W R1039 1-216-16-20-01 METAL CHIP 5% 1/16W R1039 1-216-16-20-01 METAL CHIP 270K 5% 1/16W R1039 1-216-16-20-01 METAL CHIP 300 5% 1/16W R1039 1-216-10-20 RES, CHIP 68 5% 1/16W R1039 1-216-10-20 RES, CHIP 1/16W R1039 1-216-10	B937	1-216-864-11	METAL CHIP	0	5%	1/16W						
R944 1-216-88-11 METAL CHIP 0 5% 1/16W R942 1-216-88-11 METAL CHIP 62K 5% 1/16W R1033 1-216-16-80-0 RES, CHIP 47 5% 1/16W R1033 1-216-86-11 METAL CHIP 47K 5% 1/16W R1033 1-216-86-11 METAL CHIP 47K 5% 1/16W R1034 1-216-86-11 METAL CHIP 47K 5% 1/16W R1034 1-216-86-11 METAL CHIP 27K 5% 1/16W R1034 1-216-86-11 METAL CHIP 27K 5% 1/16W R1034 1-216-86-11 METAL CHIP 27K 5% 1/16W R1034 1-216-86-11 METAL CHIP 10K 5% 1/16W R1034 1-216-82-11 METAL CHIP 10K 1/16W R1034 1-216-82-11 METAL CHIP 10K 1/16W R1034 1-216-82-11 METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (GG) METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (GG) METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-82-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-83-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-83-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-83-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-83-11 R15K, ADJ, CERMET 1K (FG GG) METAL CHIP 10K 1/16W R1034 1-216-83-11 R15K, ADJ, CERMET 1K (FG GG) META							R1029	1-216-027-00	METAL CHIP	120	5%	1/10W
R944 1-216-84-11 METAL CHIP 1K 5% 1/16W R1031 1-216-169-00 RES, CHIP 47 5% 1/16W R1032 1-216-169-00 RES, CHIP 68 5% 1/16W R1038 1-216-169-01 RES, CHIP 68 5% 1/16W R1038 1-216-170-00 RES, CHIP 68 63 1/16W R10				-						100		1/10W
R944 1-216-941-11 METAL CHIP 47K 5% 1/16W R946 1-216-93-91 METAL CHIP 47K 5% 1/16W R946 1-216-859-11 METAL CHIP 47K 5% 1/16W R946 1-216-859-11 METAL CHIP 27K 5% 1/16W R949 1-216-170-00 RES, CHIP 68 5% 1/16W R959 1-216-180-10 METAL CHIP 22 5% 1/16W R959 1-216-180-11 METAL CHIP 10K 0.5% 1/16W R950 1-216-180-11 METAL CHIP 0.5% 1/16W R950 1-216-180-11 METAL CHIP 0.5% 1/16W R960 1-216-180-11 METAL CHIP 0.5% 1/16W R960 1-216-180-11 METAL CHIP 0.5% 1/16W R960 1-216-180-11 METAL CHIP 47K 5% 1/16W											5%	1/8W
R948 1-216-84-11 METAL CHIP 47K 5% 1/16W R944 1-216-83-11 METAL CHIP 47K 5% 1/16W R948 1-216-828-11 METAL CHIP 47K 5% 1/16W R948 1-216-170-00 RES, CHIP 68 5% 1/16W R949 1-216-83-19 RES, CHIP 68 5% 1/16W R949 1-216-83-91 RES, CHIP 68 5% 1/16W R951 1-216-93-00 METAL CHIP 330 5% 1/16W R951 1-216-90-700 METAL CHIP 330 5% 1/16W R951 1-216-90-700 METAL CHIP 330 5% 1/16W R951 1-216-170-00 RES, CHIP 68 5% 1/16W R953 1-216-170-00 RES, CHIP 68 5% 1/16W R959 1-216-80-11 METAL CHIP 3.3K 5% 1/16W R959 1-216-80-10 METAL CHIP 3.3K 5% 1/16W R959 1-216-80-10 METAL CHIP 5.8K 0.5% 1/16W R959 1-218-80-11 METAL CHIP 6.8K 0.5% 1/16W R959 1-218-80-11 METAL CHIP 0.5% 1/16W R959 1-218-80-11 METAL CHIP 0.5% 1/16W R959 1-218-80-11 METAL CHIP 0.5% 1/16W R960 1-216-80-31 METAL CHIP 0.5% 1/16W R960 1-21										47		
R944 1-216-83-11 METAL CHIP 47K 5% 1/16W R1034 1-216-850-11 METAL CHIP 27K 5% 1/16W R1036 1-216-850-11 METAL CHIP 100K 5% 1/16W R1036 1-216-850-11 METAL CHIP 100K 5% 1/16W R1036 1-216-850-11 METAL CHIP 27K 5% 1/16W R1036 1-216-850-11 METAL CHIP 4.7K 5% 1/16W R1036 1-216-830-11 RES, ADJ, CERMET 1K (G GAIN) R1036 1-216-850-11 RES, ADJ, CERMET 1K (G GAIN) R1036 1-216-870-10 R1036 R103												1/16W
R944 12[6-85-11 METAL CHIP 470K 5% 1/16W R948 1-2[16-170-00 RES, CHIP 68 5% 1/8W R1036 1-2[6-829-11 METAL CHIP 4.7K 5% 1/16W R949 1-2[16-170-00 RES, CHIP 68 5% 1/8W R1036 1-2[6-829-11 METAL CHIP 4.7K 5% 1/16W R951 1-2[6-83-00 METAL CHIP 330 5% 1/16W R951 1-2[6-83-00 METAL CHIP 330 5% 1/16W R952 1-2[6-83-00 METAL CHIP 330 5% 1/16W R952 1-2[6-83-00 RES, CHIP 68 5% 1/16W R953 1-2[6-170-00 RES, CHIP 68 5% 1/16W R953 1-2[6-170-00 RES, CHIP 68 5% 1/16W R953 1-2[6-170-00 RES, CHIP 68 5% 1/16W R954 1-2[6-80-11 METAL CHIP 22 5% 1/16W R954 1-2[6-80-11 METAL CHIP 24 5% 1/16W R955 1-2[6-80-11 METAL CHIP 25 5% 1/16W R955 1-2[6-80-11 METAL CHIP 27 5% 1/16W R955 1-2[6-80-11 METAL CHIP 470 5% 1/16W R956 1-2[6-80-11 METAL CHIP 470 5% 1/16W R956 1-2[6-80-11 METAL CHIP 470 5% 1/16W R956 1-2[6-80-11 METAL CHIP 470 5% 1/16W R957 1-2[6-80-11 METAL CHIP 470 5% 1/16W R957 1-2[6-80-11 METAL CHIP 470 5% 1/16W R956 1-2[6-80-11 METAL CHIP 470 5% 1/16W R957 1-2[6-80-11 METAL CHIP 470 5% 1/16W R956 1-2[6-80-11 METAL CHIP 470 5% 1/16W R957 1-2[6-80-11 METAL CHIP 470 5% 1/16W R958 1-2[6-80-11 METAL CHIP 470 5% 1/16W R959 1-2[6-80-11 METAL CHIP 4	R943	1-216-841-11	METAL CHIP	47K	5%	1/16W						
R946 1-216-826-11 METAL CHIP 2.7K 5% 1/16W R949 1-216-170-00 RES, CHIP 68 5% 1/8W R949 1-216-18-33-91 RES, CHIP 68 5% 1/8W R951 1-216-83-91 RES, CHIP 68 5% 1/8W R951 1-216-83-91 RES, CHIP 68 5% 1/16W R952 1-216-170-00 RES, CHIP 68 5% 1/16W R952 1-216-170-00 RES, CHIP 68 5% 1/16W R920 1-238-853-11 RES, ADJ, GERMET 1K (G.GAIN) R952 1-216-170-00 RES, CHIP 68 5% 1/16W RV201 1-238-853-11 RES, ADJ, GERMET 1K (G.GAIN) R953 1-216-170-00 RES, CHIP 68 5% 1/16W RV201 1-238-852-11 RES, ADJ, GERMET 1K (G.GAIN) R953 1-216-810-10 REJ, CHIP 22 5% 1/16W RV206 1-238-852-11 RES, ADJ, GERMET 470 (REC G. ILEVEL) R956 1-218-871-11 METAL CHIP 22 5% 1/16W RV206 1-238-852-11 RES, ADJ, GERMET 470 (REC G. ILEVEL) R959 1-218-871-11 METAL CHIP 0 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 470 (REC G. ILEVEL) R959 1-218-871-11 METAL CHIP 0 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) R960 1-218-871-11 METAL CHIP 0 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) R960 1-218-834-11 METAL CHIP 0 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) R960 1-218-834-11 METAL CHIP 0 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) R960 1-218-834-11 METAL CHIP 470 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) R960 1-218-834-11 METAL CHIP 470 5% 1/16W RV207 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) RV207 RV208 1-238-853-11 RES, ADJ, GERMET 1K (EV LEVEL) RV208 RV20				470K	5%	1/16W	R1034	1-216-850-11	METAL CHIP	270K	5%	1/16W
R988 1-216-170-00 RES, CHIP 68 5% 1/8W R1036 1-216-829-11 METAL CHIP 4.7K 5% 1/16W R951 1-216-933-91 RES, CHIP 30 5% 1/16W R951 1-216-937-90 METAL CHIP 30 5% 1/16W R952 1-216-170-00 RES, CHIP 68 5% 1/8W R953 1-216-170-00 RES, CHIP 88 5% 1/8W R958 1-216-801-11 METAL CHIP 22 5% 1/16W R959 1-218-601-11 METAL CHIP 22 5% 1/16W R959 1-218-601-11 METAL CHIP 22 5% 1/16W R959 1-218-601-11 METAL CHIP 0 5% 1/16W R959 1-216-864-11 METAL CHIP 0 5% 1/16W R960 1-216-864-11 METAL CHIP 470 5% 1/16W				2.7K			R1035	1-216-845-11	METAL CHIP	100K	5%	1/16W
R949 1-216-18-03-91 RES, CHIP 68 5% 1/8W		1-216-170-00	RES, CHIP				R1036	1-216-829-11	METAL CHIP	4.7K	5%	1/16W
R950 1-216-833-91 RES, CHIP 10K 5% 1/16W R9201 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) R952 1-216-170-00 RES, CHIP 68 5% 1/8W RV202 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) R953 1-216-170-00 RES, CHIP 68 5% 1/8W RV203 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) R955 1-216-170-00 RES, CHIP 68 5% 1/8W RV203 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) R956 1-216-170-00 RES, CHIP 68 5% 1/8W RV204 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) RV205 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) RV206 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) RV207 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) RV208 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) RV208 1-238-853-11 RES, ADJ, CERMET IX (CBAIN) RV209 1-238-853-11 RES, ADJ, CERME					5%							
P851 1-216-070-00 METAL CHIP 330 5% 1/16W P852 1-216-170-00 RES, CHIP 68 5% 1/8W P8201 1-238-863-11 RES, ADJ, CERMET IX (G GAIN) P853 1-216-170-00 RES, CHIP 68 5% 1/8W P8201 1-238-863-11 RES, ADJ, CERMET IX (F GAIN) P859 1-216-081-00 METAL CHIP 3.3K 5% 1/8W P8201 1-238-863-11 RES, ADJ, CERMET 470 (REC OR LEVEL) P859 1-216-801-11 METAL CHIP 20 5% 1/16W P8201 1-238-852-11 RES, ADJ, CERMET 470 (REC OR LEVEL) P859 1-216-801-11 METAL CHIP 10K 0.5% 1/16W P8201 1-238-853-11 RES, ADJ, CERMET 470 (REC OR LEVEL) P859 1-216-801-11 METAL CHIP 10K 0.5% 1/16W P8201 1-238-853-11 RES, ADJ, CERMET 1X (F BURST LEVEL) P8601 1-216-804-11 METAL CHIP 0 5% 1/16W P8701 1-216-804-11 METAL CHIP 47K 5% 1/16W P8701 1-216-804-1									< VARIABLE RE	SISTOR >		
R952 1-216-170-00 RES, CHIP 68 5% 1/8W R953 1-216-170-00 RES, CHIP 68 5% 1/8W R955 1-216-170-00 RES, CHIP 68 5% 1/8W RV204 1-238-852-11 RES, ADJ, CERMET 147 (REC CA LEVEL) R957 1-216-061-00 METAL CHIP 2.5% 1/16W RV205 1-238-852-11 RES, ADJ, CERMET 470 (REC CA LEVEL) R959 1-216-887-11 METAL CHIP 2.5% 1/16W RV207 1-238-852-11 RES, ADJ, CERMET 470 (REC CA LEVEL) R959 1-218-887-11 METAL CHIP 6.8K 0.5% 1/16W RV207 1-238-853-11 RES, ADJ, CERMET 470 (REC CA LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 147 (REC CA LEVEL) RV208 1-238-853-11 RES, ADJ, CERMET 147 (REC CA LEVEL) RV209 1-238-853-11 RES	R950	1-216-833-91	RES, CHIP	10K	5%	1/16W						
R953 1-216-170-00 RES, CHIP 68 5% 1/8W RV203 1-238-883-11 RES, ADJ, CERMET 17K (PG CAL LEVEL) RV205 1-238-883-11 RES, ADJ, CERMET 470 (REC CAL LEVEL) RV205 1-238-883-11 RES, ADJ, CERMET 470 (REC CAL LEVEL) RV205 1-238-883-11 RES, ADJ, CERMET 470 (REC CAL LEVEL) RV205 1-238-883-11 RES, ADJ, CERMET 470 (REC CAL LEVEL) RV206 1-238-883-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-885-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-885-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-885-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-885-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-885-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV208 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV208 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV209 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV209 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERMET 17K (PG COR LEVEL) RV201 1-238-853-11 RES, ADJ, CERM	R951	1-216-037-00	METAL CHIP	330	5%	1/10W	RV201	1-238-853-11	RES, ADJ, CERI	MET 1K (C (GAIN)	
R955 1-216-170-00 RES, CHIP 68 5% 1/8W RV205 1-238-852-11 RES, ADJ, CERMET 470 (REC CR LEVEL) R959 1-216-610-10 METAL CHIP 22 5% 1/16W RV205 1-238-852-11 RES, ADJ, CERMET 470 (REC CR LEVEL) R959 1-218-887-11 METAL CHIP 22 5% 1/16W RV207 1-238-852-11 RES, ADJ, CERMET 470 (REC CR LEVEL) R959 1-218-887-11 METAL CHIP 22 5% 1/16W RV207 1-238-852-11 RES, ADJ, CERMET 470 (REC CR LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV207 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES, ADJ, CERMET 170 (REC CR LEVEL) (DSR-20MDP RV307 1-238-853-11 RES,	R952	1-216-170-00	RES, CHIP	68	5%	1/8W	RV202	1-238-853-11	RES, ADJ, CERI	MET 1K (AG	C)	
R056	R953	1-216-170-00	RES, CHIP	68	5%		RV203	1-238-853-11				
R956 1-216-01-00 RES, CHIP 68 5% 1/16W R957 1-216-061-00 METAL CHIP 22 5% 1/16W R959 1-216-801-11 METAL CHIP 22 5% 1/16W RV201 1-238-853-11 RES, ADJ, CERMET 470 (REC OB LEVEL) R959 1-218-867-11 METAL CHIP 0K 0.5% 1/16W RV401 1-238-853-11 RES, ADJ, CERMET 14, KIEP LEVEL) (DSR-20MDP R959 1-216-864-11 METAL CHIP 0 5% 1/16W RV402 1-238-853-11 RES, ADJ, CERMET 14, KIEP LEVEL) (DSR-20MDP R962 1-216-861-11 METAL CHIP 0 5% 1/16W RV404 1-238-853-11 RES, ADJ, CERMET 14, KIEP LEVEL) (DSR-20MDP R963 1-216-831-11 METAL CHIP 0 5% 1/16W RV406 1-238-853-11 RES, ADJ, CERMET 14, KIEP LEVEL) (DSR-20MDP R963 1-216-864-11 METAL CHIP 0 5% 1/16W R969 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-864-11 METAL CHIP 0 5% 1/16W R971 1-216-841-11 METAL CHIP 47K 5% 1/16W R971 1-216-841-11 METAL CHIP 47K 5% 1/16W R973 1-216-841-11 METAL CHIP 47K 5% 1/16W R97	R955	1-216-170-00	RES, CHIP	68	5%	1/8W	RV204					
R957 1-216-801-10 METAL CHIP 22 57 1/16W R920 1-238-853-11 RES, ADJ, CERMET 470 (REC GB LEVEL) R959 1-218-867-11 METAL CHIP 6.8K 0.5% 1/16W RV401 1-238-853-11 RES, ADJ, CERMET 47K (DECODER HUE) GOSR-20MDP R960 1-216-817-11 METAL CHIP 0 5% 1/16W RV402 1-238-853-11 RES, ADJ, CERMET 1K (FB DURST LEVEL) GOSR-20MDP R961 1-216-837-11 METAL CHIP 0 5% 1/16W RV402 1-238-853-11 RES, ADJ, CERMET 1K (FB CLEVEL) GOSR-20MDP R962 1-216-817-11 METAL CHIP 0 5% 1/16W (DSR-20MDP) R964 1-216-864-11 METAL CHIP 0 5% 1/16W (DSR-20MDP) R966 1-216-864-11 METAL CHIP 0 5% 1/16W R968 1-216-833-91 RES, CHIP 10K 5% 1/16W R969 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-841-11 METAL CHIP 47K 5% 1/16W R971 1-216-841-11 METAL CHIP 47K 5% 1/16W R980 1-216-841-11 METAL C							RV205	1-238-852-11	RES, ADJ, CERI	MET 470 (R	EC Y LEV	EL)
R958 1-216-801-11 METAL CHIP 22 5% 1/16W R959 1-218-867-11 METAL CHIP 10K 0.5% 1/16W R0401 1-238-853-11 RES, ADJ, CERMET 1K (EEV LEVEL) (DSR-20MDP R960 1-218-871-11 METAL CHIP 10K 0.5% 1/16W R0401 1-238-853-11 RES, ADJ, CERMET 1K (EEV LEVEL) (DSR-20MDP R962 1-216-864-11 METAL CHIP 470 5% 1/16W (DSR-20MDP R964 1-216-864-11 METAL CHIP 0 5% 1/16W (DSR-20MDP R964 1-216-864-11 METAL CHIP 0 5% 1/16W (DSR-20MDP R966 1-216-864-11 METAL CHIP 0 5% 1/16W R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R966 1-216-864-11 METAL CHIP 0 5% 1/16W (DSR-20MDP R966 1-216-864-11 METAL CHIP 0 5% 1/16W R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-216-864-11 METAL CHIP 0 5% 1/16W (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040 1-238-853-11 RES, ADJ, CERMET 1K (FE C LEVEL) (CSR-20MDP R040	R956											
R959 1-218-867-11 METAL CHIP 6.8K 0.5% 1/16W RV402 1-238-853-11 RES, ADJ, CERMET 1K (EF V LEVEL) CDSR-20MDP R961 1-216-864-11 METAL CHIP 470 5% 1/16W RV402 1-238-853-11 RES, ADJ, CERMET 1K (PB BURST LEVEL) CDSR-20MDP R962 1-216-817-11 METAL CHIP 470 5% 1/16W RV404 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) RV404 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) RV404 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) RV404 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) RV406 1-238-853-11 RES, ADJ, CERMET 1K (PB C LE	R957	1-216-061-00	METAL CHIP									
R960	R958											
(DSR-20MDP) R961 1-216-847-11 METAL CHIP 470 5% 1/16W R963 1-216-837-91 RES, CHIP 10K 5% 1/16W R964 1-216-864-11 METAL CHIP 470 5% 1/16W R965 1-216-817-11 METAL CHIP 470 5% 1/16W R966 1-216-864-11 METAL CHIP 0 5% 1/16W R966 1-216-864-11 METAL CHIP 0 5% 1/16W R968 1-216-839-91 RES, CHIP 10K 5% 1/16W R968 1-216-839-91 RES, CHIP 10K 5% 1/16W R970 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-864-11 METAL CHIP 0 5% 1/16W R970 1-216-841-11 METAL CHIP 470 5% 1/16W R971 1-216-841-11 METAL CHIP 47K 5% 1/16W R973 1-216-841-11 METAL CHIP 47K 5% 1/16W R975 1-216-841-11 METAL CHIP 47K 5% 1/16W R976 1-216-817-11 METAL CHIP 47C 5% 1/16W R977 1-216-841-11 METAL CHIP 47C 5% 1/16W R978 1-216-817-11 METAL CHIP 47C 5% 1/16W R979 1-216-817-11 METAL CHIP 47C 5% 1/16W R970 1-216-81-11 METAL CHIP 47C 5% 1/16W R971 1-216-841-11 METAL CHIP 47C 5% 1/16W R972 1-216-817-11 METAL CHIP 47C 5% 1/16W R973 1-216-817-11 METAL CHIP 47C 5% 1/16W R974 1-216-841-11 METAL CHIP 47C 5% 1/16W R975 1-216-817-11 METAL CHIP 47C 5% 1/16W R976 1-216-817-11 METAL CHIP 47C 5% 1/16W R977 1-216-817-11 METAL CHIP 47C 5% 1/16W R978 1-216-817-11 METAL CHIP 47C 5% 1/16W R979 1-216-841-11 METAL CHIP 47C 5% 1/16W R979 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R989 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 47C 5% 1/16W R999 1-216-841-11 METAL CHIP 4	R959											
R961 1-216-884-11 METAL CHIP 470 5% 1/16W R9404 1-238-853-11 RES, ADJ, CERMET 1K (EE C LEVEL) R962 1-216-837-11 METAL CHIP 470 5% 1/16W RV406 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) R964 1-216-864-11 METAL CHIP 470 5% 1/16W RV406 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) R965 1-216-864-11 METAL CHIP 470 5% 1/16W RV406 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) R966 1-216-864-11 METAL CHIP 470 5% 1/16W RV406 1-238-853-11 RES, ADJ, CERMET 1K (PB C LEVEL) R967 1-216-864-11 METAL CHIP 0 5% 1/16W RP401 1-535-757-11 CHIP, CHECKER R968 1-216-833-91 RES, CHIP 10K 5% 1/16W RP509 1-535-757-11 CHIP, CHECKER R969 1-216-841-11 METAL CHIP 470 5% 1/16W RP509 1-535-757-11 CHIP, CHECKER R969 1-216-841-11 METAL CHIP 47K 5% 1/16W AV509 1-579-738-21 VIBRATOR, CERAMIC (12288MHz) R976 1-216-841-11 METAL CHIP 47K 5% 1/16W AV509 1-579-738-21 VIBRATOR, CERAMIC (500KHz) R978 1-216-841-11 METAL CHIP 47K 5% 1/16W AV509 1-579-780-21 VIBRATOR, CERAMIC (500KHz) R978 1-216-841-11 METAL CHIP 47K 5% 1/16W AV509 1-579-780-21 VIBRATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R985 1-216-841-11 METAL CHIP 47K 5% 1/16W AV509 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R985 1-216-843-11 METAL CHIP 47K 5% 1/16W AV509 1-567-930-21 VIBRATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R985 1-216-843-11 METAL CHIP 47K 5% 1/16W AV509 1-567-930-21 VIBRATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R985 1-216-843-11 METAL CHIP 47K 5% 1/16W AV509 1-567-900-11 OSCILLATOR, CRYSTAL (14.33419MHz) (DSR-20MDP) R985 1-216-845-11 METAL CHIP 470 5% 1/16W AV509 1-5661-21 VIBRATOR, CRYSTAL (4.33819MHz) (DSR-20MDP) R995 1-216-845-11 METAL CHIP 470 5% 1/16W AV509 1-566-900-11 OSCILLATOR, CRYSTAL (4.33819MHz) (DSR-20MDP) R995 1-216	R960	1-218-871-11	METAL CHIP	10K	0.5%	1/16W	RV402	1-238-853-11	RES, ADJ, CERI	MET 1K (PB		
R962									DEC 404 050		•	
R963 1-216-833-91 RES, CHIP 10K 5% 1/16W 1/26-884-11 METAL CHIP 0 5% 1/16W 1/26-884-11 METAL CHIP 470 5% 1/16W 1/26-883-91 RES, CHIP 10K 5% 1/16W 1/26-883-91 RES, CHIP 470 5% 1/16W				-			RV404	1-238-853-11	RES, ADJ, CERI	VIET 1K (EE	C LEVEL)
R964 1-216-864-11 METAL CHIP 0 5% 1/16W CDSR-20MDP CREATER CDSR-20MDP CREATER CDSR-20MDP CREATER CDSR-20MDP								4 666 656 34	DEC 401 0501	AET 414 (DD	0 5 5	
CISR-20MDP CISR-20MDP							RV406	1-238-853-11	RES, ADJ, CERI	MELLIK (PB	C LEVEL	.)
R965 1-216-864-11 METAL CHIP 470 5% 1/16W TP401 1-535-757-11 CHIP, CHECKER TP401	R964	1-216-864-11	METAL CHIP	U					TECT DIN .			•
R966 1-216-864-11 METAL CHIP 0 5% 1/16W TP450 1-535-757-11 CHIP, CHECKER TP401 T-535-757-11 CHIP, CHIP TP401 T-535-757-11 CHIP, CHIP, CHIP, CHIP	DOCE	1 016 017 11	METAL CHID	470					< IEST FIN >			
R966 1-216-864-11 METAL CHIP 0 5% 1/16W TP859 1-535-757-11 CHIP, CHECKER R968 1-216-864-11 METAL CHIP 0 5% 1/16W TP859 1-535-757-11 CHIP, CHECKER R968 1-216-833-91 RES, CHIP 10K 5% 1/16W TP860 1-535-757-11 CHIP, CHECKER R969 1-216-864-11 METAL CHIP 470 5% 1/16W TP860 1-535-757-11 CHIP, CHECKER	H900	1-210-01/-11	WE TAL CHIP	470	370	1/1044	TDOOL	1.525.757.11	CHID CHECKED	1		
R967 1-216-864-11 METAL CHIP 0 5% 1/16W R968 1-216-833-91 RES, CHIP 10K 5% 1/16W R970 1-216-864-11 METAL CHIP 470 5% 1/16W R970 1-216-864-11 METAL CHIP 47 5% 1/16W X052 1-767-879-11 VIBRATOR, CERAMIC (8MHz) X052 1-767-879-11 VIBRATOR, CERAMIC (1.2288MHz) X052 1-767-879-11 VIBRATOR, CRYSTAL (14.318182MHz) X052 1-567-900-21 VIBRATOR, CRYSTAL (17.734475MHz) X402 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) X402 1-567-733-11 VIBRATOR, CRYSTAL (14.318182MHz) X402 1-567-733-11 VIBRATOR, CRYSTAL (14.318182MHz) X403 1-579-780-21 VIBRATOR, CRYSTAL (14.318182MHz) X403 1-579-780-21 VIBRATOR, CRYSTAL (14.318182MHz) X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) X601 1-579-466-11 VIBRATOR, CRYSTAL (17.734475MHz) X601 1-579-466-11 VIBRATOR, CRYSTAL (17.734475MHz) X601 1-579-466-11 VIBRATOR, CRYSTAL (14.33818MHz) X601 1-579-661-21 OSCILLATOR, CRYSTAL (14.33	DOSS	1_216_964_11	METAL CHID	0	50/	1/16///						
R968 1-216-833-91 RES, CHIP 10K 5% 1/16W R970 1-216-841-11 METAL CHIP 470 5% 1/16W 7873 1-216-841-11 METAL CHIP 47K 5% 1/16W 7874 1-216-841-11 METAL CHIP 47K 5% 1/16W 7875 1-216-841-11 METAL CHIP 47K 5% 1/16W 7876 1/16W 7876 1-216-833-91 RES, CHIP 10K 5% 1/16W 7898 1-216-837-11 METAL CHIP 470 5% 1/16W 7898 1-216-845-11 METAL CHIP 470 5% 1/16W 7898 1-216-833-91 RES, CHIP 10K							1					
R969 1-216-817-11 METAL CHIP 470 5% 1/16W												
R970 1-216-864-11 METAL CHIP 0 5% 1/16W							11 000	1 000 707 11	om, onconci	•		
R972 1-216-805-11 METAL CHIP 47 5% 1/16W X051 1-579-125-11 VIBRATOR, CERAMIC (8MHz) R973 1-216-841-11 METAL CHIP 47K 5% 1/16W X201 1-579-738-21 VIBRATOR, CERAMIC (1.2288MHz) R975 1-216-841-11 METAL CHIP 47K 5% 1/16W X201 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) (DSR-20MDP) R976 1-216-841-11 METAL CHIP 470 5% 1/16W X201 1-579-738-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R977 1-216-154-00 RES, CHIP 15 5% 1/8W R979 1-216-844-11 METAL CHIP 47K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R980 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R981 1-216-837-11 METAL CHIP 22K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (14.318182MHz) (DSR-20MDP) R986 1-216-837-11 METAL CHIP 22K 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) (DSR-20MDP) R988 1-216-837-11 METAL CHIP 470 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R992 1-216-84-11 METAL CHIP 470 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R993 1-216-84-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-567-900-11 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP									< VIBRATOR >			
R973 1-216-841-11 METAL CHIP 47K 5% 1/16W X052 1-767-879-11 VIBRATOR, CERAMIC (1.2288MHz) R974 1-216-841-11 METAL CHIP 47K 5% 1/16W X201 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R975 1-216-841-11 METAL CHIP 47K 5% 1/16W X201 1-579-738-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R976 1-216-841-11 METAL CHIP 470 5% 1/16W X201 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R977 1-216-154-00 RES, CHIP 15 5% 1/8W R978 1-216-844-11 METAL CHIP 82K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) R980 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R981 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R982 1-219-570-11 RES, CHIP 10M 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R985 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R989 1-216-864-11 METAL CHIP 470 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R992 1-216-864-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R993 1-216-845-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (14.31818MHz) R995 1-216-833-91 RES, CHIP 100K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)	11370	1-210-004-11	WEIAL OIII	U	070	1,1000	İ		(VIBIOTI OTT >			
R973 1-216-841-11 METAL CHIP 47K 5% 1/16W X052 1-767-879-11 VIBRATOR, CERAMIC (1.2288MHz) R974 1-216-841-11 METAL CHIP 47K 5% 1/16W X201 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R975 1-216-841-11 METAL CHIP 47K 5% 1/16W X201 1-579-738-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R976 1-216-841-11 METAL CHIP 470 5% 1/16W X201 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R977 1-216-154-00 RES, CHIP 15 5% 1/8W R978 1-216-844-11 METAL CHIP 82K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) R980 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R981 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R982 1-219-570-11 RES, CHIP 10M 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R985 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R989 1-216-864-11 METAL CHIP 470 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R992 1-216-864-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R993 1-216-845-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (14.31818MHz) R995 1-216-833-91 RES, CHIP 100K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)	R972	1-216-805-11	METAL CHIP	47	5%	1/16W	X051	1-579-125-11	VIBRATOR, CER	RAMIC (8MF	łz)	
R974 1-216-841-11 METAL CHIP 47K 5% 1/16W 7975 1-216-841-11 METAL CHIP 47K 5% 1/16W 7976 1-216-841-11 METAL CHIP 47C 5% 1/16W 7976 1-216-841-11 METAL CHIP 47C 5% 1/16W 7976 1-216-841-11 METAL CHIP 47C 5% 1/16W 7976 1-216-843-11 METAL CHIP 47C 5% 1/16W 7976 1-216-843-11 METAL CHIP 47C 5% 1/16W 7976 1-216-845-11 METAL CHIP 47C 5% 1/16W 7976-845-11 METAL CHIP 47C 5% 1/16W										,	•	
R975 1-216-841-11 METAL CHIP 47K 5% 1/16W R976 1-216-817-11 METAL CHIP 470 5% 1/16W X201 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R977 1-216-154-00 RES, CHIP 15 5% 1/8W R979 1-216-844-11 METAL CHIP 82K 5% 1/16W R980 1-216-841-11 METAL CHIP 47K 5% 1/16W R981 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R982 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R982 1-216-837-11 METAL CHIP 22K 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) (DSR-20MDP) R985 1-216-837-11 METAL CHIP 22K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R988 1-216-837-11 METAL CHIP 470 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MDP) R999 1-216-845-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R993 1-216-845-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R994 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)												Hz)
R976 1-216-817-11 METAL CHIP 470 5% 1/16W X201 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R977 1-216-154-00 RES, CHIP 15 5% 1/8W R979 1-216-844-11 METAL CHIP 82K 5% 1/16W R980 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz) (DSR-20MDP) R982 1-219-570-11 RES, CHIP 10M 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R985 1-216-837-11 METAL CHIP 22K 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) (DSR-20MDP) R986 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R989 1-216-864-11 METAL CHIP 470 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R990 1-216-864-11 METAL CHIP 470 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MDP) R991 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R992 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R993 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-567-900-11 OSCILLATOR, CRYSTAL (4.433818MHz)												
R977 1-216-154-00 RES, CHIP 15 5% 1/8W X401 1-577-165-11 VIBLATOR, CERAMIC (500kHz)							X201	1-579-780-21	VIBRATOR, CRY	/STAL (17.7	,	,
R978 1-216-154-00 RES, CHIP 15 5% 1/8W R979 1-216-844-11 METAL CHIP 82K 5% 1/16W R980 1-216-841-11 METAL CHIP 47K 5% 1/16W R981 1-216-841-11 METAL CHIP 47K 5% 1/16W R982 1-219-570-11 RES, CHIP 10M 5% 1/16W R985 1-216-837-11 METAL CHIP 22K 5% 1/16W R986 1-216-833-91 RES, CHIP 10K 5% 1/16W R988 1-216-817-11 METAL CHIP 470 5% 1/16W R999 1-216-845-11 METAL CHIP 470 5% 1/16W R999 1-216-833-91 RES, CHIP 10K 5% 1/16W R990 1-216-833-91 RES, CHIP 10K 5% 1/16W R991 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R996 1-216-833-91 RES, CHIP 10K 5% 1/16W R997 1-216-833-91 RES, CHIP 10K 5% 1/16W R998 1-216-833-91 RES, CHIP 10K 5% 1/16W R999 1-216-833-91 RES, CHIP 10K 5% 1/16W R990 1-216-833-91 RES, CHIP 10K 5% 1/16W									•	•		
R978 1-216-154-00 RES, CHIP 15 5% 1/8W R979 1-216-844-11 METAL CHIP 82K 5% 1/16W R980 1-216-841-11 METAL CHIP 47K 5% 1/16W R981 1-216-841-11 METAL CHIP 47K 5% 1/16W R982 1-219-570-11 RES, CHIP 10M 5% 1/16W R985 1-216-837-11 METAL CHIP 22K 5% 1/16W R986 1-216-833-91 RES, CHIP 10K 5% 1/16W R988 1-216-817-11 METAL CHIP 470 5% 1/16W R999 1-216-845-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 470 5% 1/16W R994 1-216-817-11 METAL CHIP 470 5% 1/16W R995 1-216-833-91 RES, CHIP 100K 5% 1/16W R996 1-216-845-11 METAL CHIP 470 5% 1/16W R997 1-216-845-11 METAL CHIP 470 5% 1/16W R998 1-216-845-11 METAL CHIP 470 5% 1/16W R999 1-216-845-11 METAL CHIP 470 5% 1/16W R999 1-216-845-11 METAL CHIP 470 5% 1/16W R999 1-216-845-11 METAL CHIP 470 5% 1/16W R999 1-216-833-91 RES, CHIP 100K 5% 1/16W R999 1-216-833-91 RES, CHIP 10K 5% 1/16W R990 1-216-833-91 RES, CHIP 10K 5% 1/16W R991 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R996 1-216-833-91 RES, CHIP 10K 5% 1/16W R997 1-216-833-91 RES, CHIP 10K 5% 1/16W R998 1-216-833-91 RES, CHIP 10K 5% 1/16W R999 1-216-833-91 RES, CHIP 10K 5% 1/16W R990 1-216-833-91 RES, CHIP 10K 5% 1/16W	R977	1-216-154-00	RES, CHIP	15	5%	1/8W	X401	1-577-165-11	VIBLATOR, CER	AMIC (500)	(Hz)	
R979 1-216-844-11 METAL CHIP 82K 5% 1/16W R980 1-216-841-11 METAL CHIP 47K 5% 1/16W R981 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-793-11 VIBRATOR, CRYSTAL (14.31818MHz) (DSR-20MD) (DSR-20MDP)										•		
R980 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (17.734475MHz) R981 1-219-570-11 RES, CHIP 10M 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R985 1-216-837-11 METAL CHIP 22K 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R986 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) R988 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) R992 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R993 1-216-845-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)							X402	1-567-900-11	OSCILLATOR, C	RYSTAL (14	4.31818N	AHz)
R981 1-216-841-11 METAL CHIP 47K 5% 1/16W X402 1-567-733-11 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MDP) R982 1-219-570-11 RES, CHIP 10M 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R985 1-216-837-11 METAL CHIP 22K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) R986 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) R988 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MD) R992 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R993 1-216-845-11 METAL CHIP 470 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)		1-216-841-11	METAL CHIP	47K	5%	1/16W					. ([SR-20MD)
R982 1-219-570-11 RES, CHIP 10M 5% 1/16W X403 1-579-738-21 VIBRATOR, CRYSTAL (14.318182MHz) R985 1-216-837-11 METAL CHIP 22K 5% 1/16W (DSR-20MD) R986 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) R988 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W R994 1-216-817-11 METAL CHIP 470 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MD) X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)	R981	1-216-841-11	METAL CHIP	47K	5%	1/16W	X402	1-567-733-11	VIBRATOR, CRY	/STAL (17.7	34475MI	Hz)
R985 1-216-837-11 METAL CHIP 22K 5% 1/16W R986 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) (DSR-20MD) R989 1-216-864-11 METAL CHIP 0 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MD) R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)											(DS	R-20MDP)
R986 1-216-833-91 RES, CHIP 10K 5% 1/16W X403 1-579-780-21 VIBRATOR, CRYSTAL (17.734475MHz) R988 1-216-817-11 METAL CHIP 470 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MD) R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W R994 1-216-817-11 METAL CHIP 470 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MD) X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz)							X403	1-579-738-21	VIBRATOR, CRY	/STAL (14.3		
R988 1-216-817-11 METAL CHIP 470 5% 1/16W R989 1-216-864-11 METAL CHIP 0 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MD) R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)										:	,	
R989 1-216-864-11 METAL CHIP 0 5% 1/16W X601 1-579-466-11 VIBRATOR, CRYSTAL (3.579545MHz) (DSR-20MD) R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W R994 1-216-817-11 METAL CHIP 470 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)							X403	1-579-780-21	VIBRATOR, CRY	/STAL (17.7		
R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W R994 1-216-817-11 METAL CHIP 470 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) (DSR-20MDP) (DSR-20MDP) (DSR-20MDP) (DSR-20MDP)												
R992 1-216-817-11 METAL CHIP 470 5% 1/16W R993 1-216-845-11 METAL CHIP 100K 5% 1/16W R994 1-216-817-11 METAL CHIP 470 5% 1/16W R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) (DSR-20MDP) X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)	R989	1-216-864-11	METAL CHIP	. 0	5%	1/16W	X601	1-579-466-11	VIBRATOR, CRY	STAL (3.57		
R993 1-216-845-11 METAL CHIP 100K 5% 1/16W X601 1-579-661-21 OSCILLATOR, CRYSTAL (4.433619MHz) R994 1-216-817-11 METAL CHIP 470 5% 1/16W (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)		4 040 0:=		477.0		4 14 002					(E	ISK-20MD)
R994 1-216-817-11 METAL CHIP 470 5% 1/16W (DSR-20MDP) R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)							Voc	4 270 004 04	0001114705	ייי יייייייייייייייייייייייייייייייייי	400040	*! 1~\
R995 1-216-833-91 RES, CHIP 10K 5% 1/16W X651 1-567-900-11 OSCILLATOR, CRYSTAL (14.31818MHz)							X601	1-5/9-661-21	USUILLATUR, C	nyoial (4.		•
							VCE4	1.567.000.11	OCCILIATOR O	DVCTAL /4/		
NSSO 1-210-017-11 WETAL UNIF 470 3% 1/10W (USK-20WID)							I COX	1-007-900-11	USUILLATUR, U	ntolat (14		
· ·	กษษ	1-210-01/-11	WE IAL UNIP	4/0	376	1/ 1 O VV					(L	OIT-ZUIVID)

M904 M905

S001 S901 S902 1-698-534-31 FAN, DC

1-762-550-11 SWITCH, ROTARY (MODE) 1-762-551-11 SWITCH, PUSH (REC PROOF)

1-572-288-11 SWITCH, PUSH (C IN SW)

Ref. No.	Part No.	<u>Description</u> Rema	<u>rk</u>
X651	1-567-733-11	VIBRATOR, CRYSTAL (17.734475MHz)	
		(DSR-20MI)P)
X652	1-577-165-11	VIBLATOR, CERAMIC (500kHz)	•
X851	1-767-450-11	VIBRATOR, CERAMIC (20MHz)	
X852	1-760-458-21	VIBRATOR, CRYSTAL (32.768kHz)	
X853		VIBRATOR, CRYSTAL (32.768kHz)	

		MISCELLANEOUS *************
55 56 57 ∆ 58 ∧ 58	1-782-823-11 1-782-825-11 1-782-824-11 1-468-441-11 1-468-442-11	CABLE, FLAT (FVH-4) CABLE, FLAT (FVF-8) CABLE, FLAT (FVJ-7) POWER BLOCK (U-1/U-2) (DSR-20MD) POWER BLOCK (U-1/U-2) (DSR-20MDP)
60 62 65 △ 67	1-782-822-11 1-782-826-11 1-958-841-11 1-958-585-11 1-958-059-11	CABLE, FLAT (FVR-9) CABLE, FLAT (FVR-10) HARNESS (DP-73) HARNESS (AC-227) HARNESS (VP-72)
101	1-776-148-11	CABLE, FLAT (FCM-11) 15P
102	1-776-145-11	CABLE, FLAT (FCM-8) 16P
105	1-764-137-11	CONNECTOR, TRANSLATION 15P
113	1-958-288-11	HARNESS (CM-130)
114	1-776-151-11	CABLE, FLAT (FCM-12) 14P
115	1-776-147-11	CABLE, FLAT (FCM-10) 15P
116	1-776-146-11	CABLE, FLAT (FCM-9) 9P
117	1-958-057-11	HARNESS (CP-79)
118	1-958-061-11	HARNESS (VJ-103)
119	1-958-058-11	HARNESS (JP-55)
120	1-958-060-11	HARNESS (VJ-102)
121	1-543-793-11	FILTER, CLAMP (FERRITE CORE)
755	A-7044-015-A	DRUM ASSY (DEH-08B-R)
851	1-658-990-11	FP-406 FLEXIBLE BOARD
CN901	1-770-312-21	CONNECTOR 4P
J901	1-564-603-41	CONNECTOR (WITH DC SW) 4P
M901	X-3944-897-2	FPC ASSY, MOTOR
M902	8-835-545-01	MOTOR, DC SCD11A/J-N (CAPSTAN)
M903	X-3945-784-1	MOTOR ASSY, LM (LOADING)
M904	8-835-537-01	MOTOR, DC SRD11A/J-N (REEL)

Ref. No.	Part No.	Description	Remark
,		ACCESSORIES	

	1-475-693-11	REMOTE COMMANDER (RMT-DS20)	
Δ	1-559-945-11	CORD, POWER (DSR-20MD)	
Δ	1-551-631-22	CORD, POWER (DSR-20MDP)	
	3-867-983-11	MANUAL, INSTRUCTION (ENGLISH, FF	RENCH)
,	3-867-983-21	MANUAL, INSTRUCTION (GERMAN, IT	ALIAN)
			-20MDP)

HARDWARE LIST

#1	7-685-533-19	SCREW +BTP 2.6X6 TYPE2 N-S
#2	7-682-552-09	SCREW +P 3X16
#3	7-682-547-09	SCREW +B 3X6
#4	7-685-132-19	SCREW +P 2.6X5 TYPE2 NON-SLIT
#5	7-682-147-01	SCREW +P 3X6
#6	7-628-253-20	SCREW +PS 2X6
#7	7-682-646-09	SCREW +PS 3X5
#8	7-628-253-00	SCREW +PS 2X4
#9	7-627-553-37	SCREW (M2X3), SPECIAL HEAD
#10	7-685-871-01	SCREW +BVTT 3X6 (S TIGHT)

DSR-20MD/20MDP

SONY

SERVICE MANUAL

US Model
Canadian Model
DSR-20MD
AEP Model
Australian Model
New Zealand Model

SUPPLEMENT-1

File this supplement with the service manual.

Addition of specifications.

SPECIFICATIONS

Medical Specifications

Protection against electric shock:

Class I

Protection against harmful ingress of water:

Ordinary

Degree of safety in the presence of flammable

anesthetics or oxygen:

Not suitable for use in the presence of flammable anesthetics or oxygen

Mode of operation:

Continuous